


See how you can  
guide the path  
her cancer takes



Direct access to  
same-day-results



Fully Automated Molecular Diagnostics





“We at Biocartis aim to provide direct access to personalized medicine for patients worldwide by developing fully integrated and broadly applicable molecular diagnostics. Our platform can be used in a wide variety of healthcare settings to enable rapid and high-quality care close to patients”

Rudi Pauwels, *Founder Biocartis*



# The need for improved, standardized and fast diagnostics

Cancer can hit anyone at any time and treatment remains a real challenge. Because cancer doesn't follow rules. It fights back against therapies. It adapts. It changes its path. It does whatever it can to stay ahead of us.

At the advanced edge of oncology, **rapid access** to **accurate data** about relevant cancer mutations and treatment resistance is vital and creates the opportunity for early disease interception<sup>4,5</sup>, reducing the anxiety while waiting for results and the time before starting the best possible treatment.

Current technologies in molecular oncology are complex, require a lot of hands-on time and are often difficult to implement in the local laboratory. As a consequence, most laboratories do not perform molecular tests in-house, but send them out to specialized centers, where samples are batched in order to optimize costs.<sup>1-3</sup>


This causes delay to the fast delivery of results, preventing rapid initiation of correct therapy. In the meantime the tumor grows, which is detrimental in case of aggressively growing cancers.

# The need for a rapid treatment initiation response towards patients

Fast initiation of immunotherapy or targeted therapy as first-line treatment is crucial for cancer patients, as it increases overall survival rates.<sup>9,10,11,17,22</sup> Timely detection of biomarkers therefore is very important.

Today, turnaround times of reference technologies are on average 18 days, with 14% of patients waiting longer than a month to be able to start treatment. Ninety-five percent of the patients have to wait more than a week in order to receive the biomarker results.<sup>6</sup>

Total turnaround time of reference technologies




This means that precious time is lost whereas treatment initiation could have been started and unnecessary use of chemotherapy with its side effects could have been avoided.

# Idylla™, the next level in disease interception

Idylla™, a **fully automated**, sample-to-result PCR based **molecular diagnostics** system, provides **same-day** results enabling physicians to make **timely decisions** on patients' therapy.

Idylla™, with its **compact scalable design** and **outstanding ease-of-use**, overcomes the traditional barriers of molecular diagnostics, allowing it to be used in virtually **any laboratory setting**.

Idylla™, can be used with **multiple sample types**, including **solid** and **liquid\*** biopsies. This flexibility allows use of the system for respectively **diagnosis**, and **research** or possibly future **monitoring** applications.



\* Liquid biopsy currently available only for research applications

\*\* e.g. deparaffinization for FFPE tissue samples

# Idylla™ is the first and only molecular diagnostic system that combines



## FAST RESULTS

- Less than 2 minutes hands-on time
- Short turnaround time - 40 to 150 minutes



## ACCURATE RESULTS

- High sensitivity
- Highly standardized technology
- Contamination-controlled design



## EASE-OF-USE

- Fully automated sample-to-result process
- Walk-away system (no need for any intervention during the automatic process)
- All reagents integrated in a single cartridge
- Only 1 manual step
- Storage and shipment at room temperature



## ACCESSIBLE

- Access on demand - no need for pre-processing or batching



## SAMPLE VERSATILITY

- For solid and liquid biopsy\*



## MULTIPLEXING CAPABILITY

- Detection of up to 52 relevant mutations in one cartridge
- Multiple genes and loci detection in one cartridge



## CONNECTIVITY

- Remote assistance, monitoring and upgrading
- Bi-directional LIS



Fully Automated Molecular Diagnostics




FAST ACCURATE EASY ACCESSIBLE

# The revolutionary Idylla™ workflow



The Idylla™ system in combination with the Idylla™ Molecular Oncology Assays differs from other technologies in its outstanding **ease-of-use**, leading to an unsurpassed level of **standardization**, and its **short turnaround time**, allowing immediate access to therapy.

## FFPE workflow





## Liquid biopsy\* workflow



\* Liquid biopsy currently available only for research applications

# Instruments and consumables



Instruments



Consumables



Lab infrastructure (# of rooms) 1

## Other RT-PCR

Instruments



Consumables



Lab infrastructure (# of rooms) 3

## Pyrosequencing

Instruments



Consumables



Lab infrastructure (# of rooms) 4

## Next generation sequencing

Instruments



Consumables



Lab infrastructure (# of rooms) 4





## Current oncology assays



### FFPE in - report out

#### Diagnostic products (CE IVD)

- Idylla™ BRAF Mutation Test
- Idylla™ KRAS Mutation Test
- Idylla™ NRAS-BRAF Mutation Test
- Idylla™ NRAS Mutation Test

#### Research products (RUO)

- Idylla™ BRAF Mutation Assay
- Idylla™ KRAS Mutation Assay
- Idylla™ EGFR Mutation Assay
- Idylla™ NRAS-BRAF-EGFR S492R Mutation Assay



### Plasma in - report out

#### Research products (RUO)

- Idylla™ ctBRAF Mutation Assay
- Idylla™ ctKRAS Mutation Assay
- Idylla™ ctNRAS-BRAF-EGFR S492R

## Future oncology assay targets



### FFPE in - report out

- MSI



### Plasma in - report out

- ctNRAS-BRAF
- ctEGFR

**BRAF** **ctBRAF**

# Idylla™ BRAF mutation detection on solid and liquid biopsies



Activating mutations in the *BRAF* gene are observed in about 8% of all cancers<sup>7</sup> and have been associated with sensitivity and resistance to a number of targeted anti-cancer therapeutics.

Cancers in which *BRAF* mutations are observed include: melanoma, colorectal cancer, thyroid cancer, lung cancer, hairy cell leukemia and ovarian cancer.

*BRAF* testing is recommended in all patients with metastatic melanoma and metastatic colorectal

cancer (mCRC). About 50% of all metastatic melanoma patients harbor mutations in the *BRAF* gene, making them eligible for BRAF or BRAF/MEK inhibitor therapy.<sup>8</sup> In mCRC, BRAF mutation status should be assessed alongside the assessment of tumor *RAS* mutational status for prognostic assessment (the presence of a *BRAF* mutation indicates poor prognosis). The prevalence of *BRAF* in mCRC is about 8-15%.<sup>9</sup>

**DIAGNOSTIC PRODUCT**

Idylla™ BRAF Mutation Test (CE IVD)

**BRAF**

Diagnostic use

- approx. 90 min. Sample-to-result
- < 2 min. hands-on time
- 7 mutations in codon 600
- Directly** on FFPE tissue sections (5-10µm) from **metastatic melanoma**
- Qualitative genotype call**
- Mutation detection for **baseline treatment**

**RESEARCH PRODUCT**

Idylla™ ctBRAF Mutation Assay (RUO)

**ctBRAF**

Research Use Only, not for diagnostic use

- approx. 85 min. Sample-to-result
- < 1 min. hands-on time
- 7 mutations in codon 600
- Directly** on 1 ml plasma
- Semi-quantitative genotype call + Cq values**
- Applicable in multiple cancers** harboring BRAF mutations

*Prof. B. Neyns, M.D., Ph.D*  
Medical Oncology, UZ Brussels, Belgium

"The Idylla™ system has the potential to allow the start of targeted therapy within a time window of less than 24 hours following the diagnosis of metastasis, thereby saving precious time"

## EGFR

# Idylla™ EGFR mutation detection on solid biopsy

Activating mutations in the *EGFR* gene have been associated with sensitivity and resistance to a number of targeted anti-cancer therapeutics.<sup>17</sup>

*EGFR* mutations are mainly observed in lung cancer.

*EGFR* mutation testing is recommended in all patients with advanced non-small-cell lung cancer (NSCLC) of a non-squamous subtype.

The prevalence of *EGFR* mutations in NSCLC adenocarcinomas is 10-15% of Western and up to 50% of Asian patients. Sensitizing *EGFR* mutations are predictive for response to *EGFR* tyrosine kinase inhibitors.<sup>11,17,21</sup>




### RESEARCH PRODUCT


Idylla™ EGFR Mutation Assay (RUO)

## EGFR

Research Use Only, not for diagnostic use



**Directly** on 1 FFPE tissue section (5µm) from **metastatic non-small-cell lung cancer**



**Qualitative genotype call**

*Prof Giancarlo Troncione*

University of Napoli Federico II, Naples

"Today, *EGFR* testing is a cumbersome process and it often takes several weeks before results are analyzed. This may lead to the administration of anti-*EGFR* therapy as second-line agents, which is less efficient than their use in first-line therapy. The Idylla™ *EGFR* Mutation assay technology has the potential to change that: it is a cost-effective solution, ensuring reliable and fast detection of all relevant mutations"

# Idylla™ KRAS mutation detection on solid and liquid biopsies



Activating mutations in the *RAS* genes are observed in 9-30% of all cancers and have been associated with sensitivity and resistance to a number of targeted anti-cancer therapeutics.<sup>12</sup> Cancers in which *KRAS* mutations are observed include: colorectal cancer, lung cancer and pancreatic cancer.

According to ESMO<sup>9</sup>, NCCN<sup>14</sup>, ASCO<sup>16</sup> and CAP/AMP/ASCO guidelines<sup>15</sup>, genotyping of clinically actionable mutations at a sensitivity of 5% in *RAS* genes exon 2 (codons 12 and 13), exon 3 (codons 59 and 61) and exon 4 (codons 117 and 146) is now mandatory on tumor tissue (either primary or metastasis) of all metastatic colorectal cancers, since the presence of these mutations correlate with the lack of response to certain anti-EGFR antibody therapies<sup>9</sup>. About 46% of all

metastatic colorectal tumors harbor mutations in exons 2, 3 and 4 of the *KRAS* gene.<sup>13</sup>

Several studies are ongoing to define the predictive impact of *KRAS* mutations on therapy decision for non-small-cell lung cancer (NSCLC) patients<sup>18,19,20</sup>. Currently there is evidence that *KRAS* in lung cancer has a prognostic value, indicating poor survival for patients with NSCLC, compared to the absence of *KRAS* mutations.<sup>11</sup>

Using liquid biopsies for *KRAS* testing is minimally invasive, fast and easy to perform and provides an excellent solution to study the presence of *KRAS* mutations in different cancer types.

**DIAGNOSTIC PRODUCT**

Idylla™ KRAS Mutation Test (CE IVD)



Diagnostic use

**RESEARCH PRODUCT**

Idylla™ ctKRAS Mutation Assay (RUO)



Research use only, not for diagnostic use

approx. 120 min. Sample-to-result

< 2 min. hands-on time

21 mutations in codons 12, 13, 59, 61, 117, 146

approx. 130 min. Sample-to-result

< 1 min. hands-on time

21 mutations in codons 12, 13, 59, 61, 117, 146

FFPE

Directly on FFPE tissue sections (5-10µm) from metastatic colorectal cancer

plasma

Directly on 1 ml plasma

Qualitative genotype call

Qualitative genotype call + Cq values

Mutation detection for baseline treatment

Applicable in multiple cancers harboring KRAS mutations

**Beatriz Bellosillo**

Laboratori de Biologia Molecular, Hospital del Mar, Barcelona

"Idylla™ allows very quick results with little hands-on time"

**NRAS****NRAS - BRAF****ctNRAS3**

# Idylla™ NRAS mutation detection on solid and liquid biopsies

Activating mutations in the *RAS* genes are observed in 9-30% of all cancers and have been associated with sensitivity and resistance to a number of targeted anti-cancer therapeutics.<sup>12</sup> Cancers in which *NRAS* mutations are observed include colorectal, lung, thyroid cancers and melanoma. According to ESMO<sup>9</sup>, NCCN<sup>14</sup>, ASCO<sup>16</sup> and the CAP/AMP/ASCO guidelines<sup>15</sup>, genotyping of clinically actionable mutations at a sensitivity of 5% in *RAS* genes exon 2 (codons 12 and 13), exon 3 (codons 59 and 61) and exon 4 (codons 117 and 146) is now mandatory on tumor tissue (either primary or metastasis) of all metastatic colorectal cancers, since the presence of these mutations correlate with the lack of response to certain anti-EGFR antibody therapies<sup>9</sup>. About 5% of all metastatic colorectal tumors harbor

mutations in exons 2, 3 and 4 of the *NRAS* gene.<sup>13</sup>

In metastatic colorectal cancer *BRAF* mutation status should be assessed alongside the assessment of tumor *RAS* mutational status for prognostic assessment (the presence of a *BRAF* mutation indicates poor prognosis). Using liquid biopsies for *NRAS* testing is minimally invasive, fast and easy to perform and provides an excellent solution to study these mutations in different cancer types and lesions. Recent research data<sup>23,24</sup> suggest that in about 16% of patients, mutations may develop in codon 492 of the *EGFR* gene as a mechanism of resistance, to the anti-EGFR antibody therapies such as cetuximab.



## DIAGNOSTIC PRODUCT

**NRAS - BRAF**

Idylla™ NRAS-BRAF Mutation Test (CE IVD)  
Idylla™ NRAS Mutation Test (CE IVD)

**NRAS**

Diagnostic use

## RESEARCH PRODUCT

Idylla™ ctNRAS-BRAF-EGFR S492R  
Mutation Assay (RUO)

**ctNRAS3**

Research use only, not for diagnostic use

approx. 120 min. Sample-to-result

< 2 min. hands-on time

5 in BRAF codon 600\* mutations  
\*Only available in Idylla™ NRAS-BRAF Mutation Test

18 in NRAS codons 12, 13, 59, 61, 117, 146 mutations

FFPE

Directly on FFPE tissue sections (5-10µm) from metastatic colorectal cancer

Semi-quantitative genotype call + Cq values

Mutation detection for baseline treatment

approx. 110 min. Sample-to-result

< 1 min. hands-on time

18 in NRAS codons 12, 13, 59, 61, 117, 146 mutations

5 in BRAF codon 600 mutations

2 in EGFR codon 492 mutations


Directly on 1 ml plasma

Semi-quantitative genotype call + Cq values

Applicable in multiple cancers harboring NRAS, BRAF or EGFR S492R mutations

# Idylla™ Connect


Engage in the future



 **Advanced services to ensure continuity in your laboratory workflow**

 **Automatic software updates**

New releases of assay and console software are sent to your Idylla™ console and can be installed with a single touch on the screen.

 **Immediate and remote service and support**

Idylla™ system parameters and error logs can be analyzed at anytime and anywhere to ensure swift actions and solutions.

# More insight into your data with Idylla™ Explore


Get connected and enjoy **the advantages of Idylla™ Explore**, a web-based application that allows you to analyze your data by providing

- Visualization of PCR curves from Idylla™ Test Results
- Cq values per target
- Direct Access to Console result reports



**Idylla™ Explore** can be accessed anywhere and anytime from your PC or tablet through the following link: <https://idyllaexplore.biocartis.com>

Subscribe today and **join the Idylla™ Explore community** by sending an email to [explore@biocartis.com](mailto:explore@biocartis.com)



The screenshot displays the Idylla Explore web application interface. At the top, there is a navigation bar with the Idylla logo, a settings icon, and a user profile icon labeled 'User 1'. Below this is a 'Filter' sidebar on the left and a 'Runs' table in the center. The 'Runs' table has columns for State, Sample ID, Test type, Run date, and Results. It lists six samples with their respective test types and results. A circular callout highlights a specific PCR curve for 'Sample 1' in the 'Results' column. This callout zooms into a detailed view of the PCR curves for Sample 1, showing five graphs (A-E) of fluorescence versus cycles. Below the graphs is a table of results for Sample 1, including Cq values and detection status for various targets.

State	Sample ID	Test type	Run date	Results
✓	Sample 1	KRAS	15 Jul 2016   08:15:12	MUTATION DETECTED IN KRAS CODON 12
✓	Sample 2	KRAS	18 Jul 2016   17:24:01	NO MUTATION DETECTED IN KRAS CODON 12, 13, 59, 61, 117, 146
✓	Sample 3	BRAF	22 Jul 2016   11:02:07	NO MUTATION DETECTED IN BRAF CODON 600
✓	Sample 4	KRAS	01 Aug 2016   21:07:57	MUTATION DETECTED IN KRAS CODON 146
✓	Sample 5	ctBRAF	04 Aug 2016   14:50:46	MUTATION DETECTED IN BRAF CODON 600
✓	Sample 6	KRAS	04 Aug 2016   10:08:51	MUTATION DETECTED IN KRAS CODON 146

Zoom in



## Custom-made solutions

Biocartis offers personalized solutions that fit your individual needs

- Create a network between different Idylla™ User sites and share data and knowledge
- Direct access to your data for building your own solution
- Statistical analysis on your obtained data
- Monitoring: Follow-up of your data over time
- Disease surveillance or diagnostic grid: linking of real-time molecular diagnostic test data to geo-location and sample data







## Join the investigation

There's a clear need for improved, standardized and fast diagnostics that allow faster initiation of targeted therapy for cancer patients.

Idylla™, Biocartis' fully automated molecular diagnostics system, is the first and only molecular diagnostic system that combines unsurpassed ease-of-use, speed and accuracy on multiple sample types. With its **compact, scalable design and outstanding ease-of-use**, Idylla™ overcomes the traditional barriers of molecular diagnostics, allowing it to be used in virtually any laboratory setting.

And by providing same-day-results, Idylla™ enables physicians to make timely decisions on patients' therapy.

**Idylla™, the next level in disease interception.**

# Idylla™ order information

<b>Diagnostic Products (CE-IVD)</b>		
Idylla™ BRAF Mutation Test	6 cartridges/box	Catalog# A0010/6
Idylla™ KRAS Mutation Test	6 cartridges/box	Catalog# A0020/6
Idylla™ NRAS-BRAF Mutation Test	6 cartridges/box	Catalog# A0030/6
Idylla™ NRAS Mutation Test	6 cartridges/box	Catalog# A0040/6

<b>Research Products (RUO)</b>		
Idylla™ BRAF Mutation Assay	6 cartridges/box	Catalog# A0011/6
Idylla™ KRAS Mutation Assay	6 cartridges/box	Catalog# A0021/6
Idylla™ EGFR Mutation Assay	6 cartridges/box	Catalog# A0061/6
Idylla™ NRAS-BRAF-EGFR S492R Mutation Assay	6 cartridges/box	Catalog# A0031/6
Idylla™ ctKRAS Mutation Assay	6 cartridges/box	Catalog# A0081/6
Idylla™ ctBRAF Mutation Assay	6 cartridges/box	Catalog# A0071/6
Idylla™ ctNRAS-BRAF-EGFR S492R Mutation Assay	6 cartridges/box	Catalog# A0091/6

<b>Platform (CE-IVD)</b>		
Idylla™ Instrument	1 unit	Catalog# P0010
Idylla™ Console	1 unit	Catalog# P1010

[customerservice@biocartis.com](mailto:customerservice@biocartis.com)  
[www.biocartis.com](http://www.biocartis.com)

Biocartis NV  
Generaal de Wittelaan 11 B3  
2800 Mechelen - Belgium  
T +32 15 632 888

## REFERENCES

1. Janku F et al. *Oncotarget*. 2015; 6(29): 26886–2689.
2. Sam SS et al. *Pathol Res Pract*. 2015. pii: jclinpath-2015–203345.
3. Colling R et al. *J Clin Pathol*. 2015. pii: jclinpath-2015–203345.
4. Bratzman SV et al. *Expert Rev Mol Diagn*. 2015; 15(6): 715–719.
5. Siravegna G and Bardelli A. *Genome Biol*. 2014; 15(8): 449.
6. Accès aux tests moléculaires EGFR, RAS et BRAF /Résultats d'une enquête dans 5 régions françaises, appui à la décision, INCa, janvier 2016
7. Mutations of the BRAF gene in human cancer. Helen Davies et al; *Nature* 2002, 417, 949-954
8. Clinical Practice Guidelines - Cutaneous melanoma: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up. *Annals of Oncology* 26 (Supplement 5): v126-v132, 2015.
9. ESMO consensus guidelines for the management of patients with metastatic colorectal cancer. *Annals of Oncology* 0: 1–37, 2016.
10. NCCN Clinical Practice Guidelines in Oncology – Melanoma - Version 3.2016
11. NCCN Clinical Practice Guidelines in Oncology – NSCLC – Version 4.2016
12. Adrienne D. Cox et al. Drugging the undruggable RAS: Mission Possible? *Nature Reviews Drug Discovery* Volume:13,Pages:828–851  
Year published:(2014)DOI:doi:10.1038/nrd4389
13. Jean-Yves Douillard, M.D., Ph.D., et al. Panitumumab-FOLFOX4 Treatment and RAS Mutations in Colorectal Cancer. *N Engl J Med* 2013;369:1023-34
14. NCCN Clinical Practice Guidelines in Oncology – Colon Cancer – Version 2.2016
15. [http://www.amp.org/committees/clinical\\_practice/CRCOpenComment.cfm](http://www.amp.org/committees/clinical_practice/CRCOpenComment.cfm)
16. Allegra CJ, et al. Extended RAS gene mutation testing in metastatic Colorectal Carcinoma to predict response to anti-epidermal growth factor receptor monoclonal antibody therapy: American Society of Clinical Oncology Provisional Clinical Opinion Update 2015. *Journal of Clinical Oncology* 2016; 34(2):179-85
17. M. Reck et al. Metastatic non-small-cell lung cancer (NSCLC): ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up. *Annals of Oncology*, 2014
18. ESMO @ ECC 2015: Response to EGFR Agents in Combination With Chemotherapy Demonstrated in Patients with Metastatic Colorectal Cancer of Rare KRAS Molecular Subtype. <http://www.esmo.org/Conferences/Past-Conferences/European-Cancer-Congress-2015/News/Response-to-EGFR-Agents-in-Combination-With-Chemotherapy-Demonstrated-in-Patients-with-Metastatic-Colorectal-Cancer-of-Rare-KRAS-Molecular-Subtype>. Sept 2015.
19. P A Janne et al. *BJC* 2015. Impact of KRAS codon subtypes from a randomised phase II trial of selumetinib plus docetaxel in KRAS mutant advanced non-small-cell lung cancer.
20. Alona Zer et al. *J Thor Onco* 2015. Pooled Analysis of the Prognostic and Predictive Value of KRAS Mutation Status and Mutation Subtype in Patients with NSCLC Treated with EGFR TKI's.
21. Wendy A. Cooper et al. *J Thorac Dis* 2013; 5 (S5): S479-490. Molecular Biology of lung cancer.
22. AACR 2016: 5-Year Survival Rates for Patients With Metastatic Melanoma Treated With Nivolumab Much Higher Than Historical Rates. <http://www.ascopost.com/News/39500>
23. Montagut C. et al., Identification of a mutation in the extracellular domain of the Epidermal Growth Factor Receptor conferring Cetuximab resistance in colorectal cancer. *Nature Medicine* 2012
24. Newhall K., Frequency of S492R Mutations in the Epidermal Growth Factor Receptor: Analysis of plasma DNA from Metastatic Colorectal Cancer Patients Treated with Panitumumab or Cetuximab Monotherapy. 16<sup>th</sup> World Congress on Gastrointestinal Cancer, Barcelona, Spain 2014.

## NOTICE

### **Idylla™ BRAF Mutation Test**

The MGB Probe contained in the BRAF Mutation Test is covered by applicable US patents and corresponding patents outside the US and is sold under a license from the ELITech Group. The purchase of this product includes a license to use only this amount of product solely for the purchaser's own use solely in the human in vitro diagnostic field (in accordance with applicable FDA and other regulatory requirements) and may not be used for any other commercial use, including without limitation repackaging or resale in any form (including resale by purchasers who are licensed to make and sell kits for use in the 5' Nuclease Process). No right under any other patent claim or for any other use is conveyed expressly, by implication, or by estoppel. Corresponding products conveying rights for use in all other fields may be obtained from Life Technologies under a separate catalog number. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

THE PURCHASE OF THE IDYLLA™ BRAF MUTATION TEST GRANTS THE PURCHASER RIGHTS UNDER CERTAIN ROCHE PATENTS TO USE IT SOLELY FOR PROVIDING HUMAN IN VITRO DIAGNOSTIC SERVICES. NO GENERAL PATENT OR OTHER LICENSE OF ANY KIND OTHER THAN THIS SPECIFIC RIGHT OF USE FROM PURCHASE IS GRANTED HEREBY.

### **Idylla™ BRAF Mutation Assay and Idylla™ ctBRAF Mutation Assay**

The MGB Probe contained in the Idylla™ BRAF Mutation Assay and in the Idylla™ ctBRAF Mutation Assay is covered by applicable US patents and corresponding patents outside the US and is sold under a license from the ELITech Group. The purchase of this product includes a license to use only this amount of product solely for the purchaser's own research use and may not be used for any other commercial use, including without limitation repackaging or resale in any form (including resale by purchasers who are licensed to make and sell kits for use in the 5' Nuclease Process). No right under any other patent claim or for any other use is conveyed expressly, by implication, or by estoppel. Diagnostic use rights for MGB may be obtained under a separate license from ELITech. Corresponding products conveying commercial and diagnostic use rights for MGB may be obtained from LTC only under a separate agreement. For further information contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

THE PURCHASE OF THE IDYLLA™ BRAF MUTATION ASSAY AND IDYLLA™ CTBRAF MUTATION ASSAY GRANTS THE PURCHASER RIGHTS UNDER CERTAIN ROCHE PATENTS TO USE IT SOLELY FOR THE PURCHASER'S OWN RESEARCH USE. NO GENERAL PATENT OR OTHER LICENSE OF ANY KIND OTHER THAN THIS SPECIFIC RIGHT OF USE FROM PURCHASE IS GRANTED HEREBY.

### **Idylla™ KRAS Mutation Test, Idylla™ KRAS Mutation Assay and Idylla™ ctKRAS Mutation Assay**

These assays contain PlexZyme and PlexPrime technology covered by patents granted and pending in certain jurisdictions, which are supplied under licence of SpeedX Pty Ltd. PlexZyme and Plexprime are trademarks of SpeedX Pty Ltd.

### **Idylla™ NRAS Mutation Test, Idylla™ NRAS-BRAF Mutation Test, Idylla™ NRAS-BRAF-EGFR S492R Mutation Assay and Idylla™ ctNRAS3 Mutation Assay**

The Idylla™ NRAS Mutation Test, Idylla™ NRAS-BRAF Mutation Test, Idylla™ NRAS-BRAF-EGFR S492R Mutation Assay and ctNRAS3 Mutation Assay contain PlexZyme and PlexPrime technology covered by patents granted and pending in certain jurisdictions, which are supplied under licence of SpeedX Pty Ltd. PlexZyme and Plexprime are trademarks of SpeedX Pty Ltd. The Idylla™ NRAS Mutation Test, Idylla™ NRAS-BRAF Mutation Test and the Idylla™ NRAS-BRAF-EGFR S492R Mutation Assay contain Hilyte and QXL probes. QXL and Hilyte are licensed pursuant to an agreement with Eurogentec S.A. and these licensed probes can be used solely for the purchaser's own research use. Hilyte™ is a trademark of Anaspec, Inc. QXL® is a registered trademark of Anaspec, Inc.

### **Idylla™ EGFR Mutation Assay**

The Idylla™ EGFR Mutation Assay contains PlexZyme and PlexPrime technology covered by patents granted and pending in certain jurisdictions, which are supplied under licence of SpeedX Pty Ltd. PlexZyme and Plexprime are trademarks of SpeedX Pty Ltd.

## IMPORTANT INFORMATION

Idylla™ platform, Idylla™ BRAF Mutation Test, Idylla™ KRAS Mutation Test, Idylla™ NRAS Mutation Test and Idylla™ NRAS-BRAF Mutation Test are CE-marked IVDs in Europe. Idylla™ BRAF Mutation Assay, Idylla™ ctBRAF Mutation Assay, Idylla™ KRAS Mutation Assay, Idylla™ ctKRAS Mutation Assay, Idylla™ NRAS-BRAF-EGFR S492R Mutation Assay, Idylla™ ctNRAS-BRAF-EGFR S492R Mutation Assay and Idylla™ EGFR Mutation Assay are available for Research Use Only, not for use in diagnostic procedures. Idylla™ ctBRAF Mutation Test, Idylla™ ctKRAS Mutation Test, Idylla™ EGFR Mutation Test, Idylla™ MSI Assay, Idylla™ ctNRAS-BRAF Mutation Test and Idylla™ ctEGFR Mutation Assay are under development. Idylla™ is not yet for sale in USA.

## COPYRIGHT INFORMATION

Biocartis trademark and logo are trademarks belonging to Biocartis and are used and registered in Europe. Idylla™ is a registered trademark in the United States (US) and other countries. Idylla trademark and logo are used trademarks belonging to Biocartis.



NOTES

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---



[www.biocartis.com](http://www.biocartis.com)

Biocartis NV  
Generaal de Wittelaan 11 B3  
2800 Mechelen - Belgium  
T +32 15 632 888

[customerservice@biocartis.com](mailto:customerservice@biocartis.com)  
[www.biocartis.com](http://www.biocartis.com)

Ref: ONCO\_AIO\_201704