



Urinary Cortisol

The IDS Urinary Cortisol assay is a fully automated CLIA assay for the quantitative determination of cortisol in human urine samples. Results are to be used in conjunction with other clinical and laboratory data to assist clinicians in the diagnosis and monitoring of disorders of the adrenal gland.

Cortisol is a steroid hormone synthesized in the adrenal gland regulating a variety of physiological processes^{1,2,3}.

Urinary Cortisol results are used in conjunction with a complete clinical evaluation to assist in the diagnosis and monitoring of disorders such as hypercortisolism, commonly known as Cushing's syndrome^{3,4}. Many common signs of Cushing's such as obesity, high blood pressure, and increased blood glucose are frequently observed in today's society⁵.

Cortisol present in the urine represents free serum cortisol levels which are not bound to the cortisol-binding-globulin (CBG) and are filtered into urine. Approximately 1% of unbound cortisol is excreted in the urine⁶. Urinary Cortisol has the advantage of not being affected by conditions that alter CBG levels such as oral contraceptives.

The Endocrine Society recommends the measurement of 24hr urinary cortisol as one of the screening tests for diagnosing Cushing's syndrome⁵, with a cut-off point for a normal response considered the upper limit of normal (ULN) range as established for each respective assay⁷.

Features and benefits

- Direct urine measurement, no extraction procedure required, for streamlined workflow
- No risk of biotin interference by employing direct antibody-solid phase coating principle
- Excellent functional sensitivity and precision over the clinically relevant range
- Validated reference intervals in 24-hour urine collections from healthy subjects
- Demonstrated clinical concordance in routine clinical samples


Specifications

Parameter	Description
Limit of Quantitation (LOQ)	0.36 µg/dL
Assay Measuring Range	0.36 – 52.0 µg/dL
Assay Reportable Range	0.36 – 120.0 µg/dL, via optional dilution procedure
Expected Reference Range (µg/24hr)	50.8 to 324.2 [2.5th - 97.5 percentile], n=147
Sample Volume	30 µL urine plus dead volume
Sample type	Urine samples (no preservative, or with addition of boric acid)
Time to first result	28 minutes
Calibration frequency	14 days
Reagent Stability	
On-board system:	14 days
2 – 8°C:	28 days

Traceability

The IDS Urinary Cortisol is an established non-standardized assay. In the JCTLM list, no reference measurement procedure or material are included for cortisol in human urine. The Beckman Coulter Access Cortisol (Ref 33600) was selected to guarantee the trueness and therefore the metrology and traceability of the assay.

N	Slope	95% CI	Intercept (g/dL)	95% CI	Correlation Coefficient (R ²)
205	1.00	0.96 to 1.04	-0.14	-0.46 to -0.22	0.96



Method comparison: IDS Urinary Cortisol vs. Commercially available immunoassay

205 specimens (range 0.35 – 47.12 µg/dL), including clinical routine samples were assessed in parallel with the IDS Urinary Cortisol and another commercially available automated immunoassay for measurement in direct urine samples.

The following results were obtained:

IDS Urinary Cortisol = 1.00x - 0.14 g/dL	
Slope [95 % CI]	0.96 to 1.04
Intercept [95 % CI]	-0.46 to -0.22 µg/dL
Correlation Coeff. R ²	0.96

Ordering information

Product name	Code	Description
IDS Urinary Cortisol	IS-5800	Reagent pack: 100 tests
IDS Urinary Cortisol Calibrator Set	IS-5820	Calibrator set: 6 levels
IDS Urinary Cortisol Control Set	IS-5830	Control set: 2 levels

Relevant assays

Product	Code	Description	Portfolio
IDS Cortisol	IS-4600	Reagent pack: 100 tests	Adrenal Function
IDS Salivary Cortisol	IS-4900	Reagent pack: 100 tests	Adrenal Function
IDS ACTH	IS-4500N	Reagent pack: 100 tests	Adrenal Function
IDS Aldosterone	IS-3300	Reagent pack: 100 tests	Endocrine Hypertension
IDS Direct Renin	IS-3400	Reagent pack: 100 tests	Endocrine Hypertension

References

- Ramamoorthy, S. & Cidlowski, J. A. Corticosteroids: Mechanisms of Action in Health and Disease. *Rheumatic diseases clinics of North America* 42, 15-31, [2016].
- Wild, D. *The immunoassay handbook*. [Elsevier, 2005].
- Cain, D. W. & Cidlowski, J. A. Specificity and sensitivity of glucocorticoid signalling in health and disease. *Best practice & research. Clinical endocrinology & metabolism* 29, 545-556, [2015].
- Findling, J. W. & Raff, H. Cushing's Syndrome: important issues in diagnosis and management. *The Journal of clinical endocrinology and metabolism* 91, 3746-3753, [2006].
- Nieman, L. K. et al. The diagnosis of Cushing's syndrome: an Endocrine Society Clinical Practice Guideline. *The Journal of clinical endocrinology and metabolism* 93, 1526-1540, [2008].
- Raff, H. et al. Urine free cortisol in the diagnosis of Cushing's syndrome: Is it worth doing and, if so, how? *The Journal of clinical endocrinology and metabolism* 100(2), 395-397, [2015].



Visit www.idsplc.com for an extended range of IDS assays

Connect with us

Global Headquarters

Immunodiagnostic Systems, 10 Didcot Way,
Boldon Business Park, Boldon, Tyne & Wear,
NE35 9PD, United Kingdom

+44 (0) 191 519 6155

info@idsplc.com

www.idsplc.com

Follow us

