

BOVINE cytokine Analysis simplified

The MILLIPLEX® Bovine Cytokine/Chemokine Panel 1 is the first multiplex panel designed to analyze up to 15 bovine cytokines within the same sample. Harness the power of this novel assay and save valuable time, money and sample while dramatically increasing the number of data points per plate.

MILLIPLEX® Bovine Cytokine/Chemokine Panel 1

(Cat. Nos. BCYT1-33K or BCYT1-33K-PX15/BCYT1-33K-PXBK15 for premixed beads)


Analytes: IFN γ , IL-1 α , IL-1 β , IL-4, IL-6, IL-8/CXCL8, IL-10, IL-17A, IL-36RA/IL-1F5, IP-10/CXCL10, MCP-1/CCL2, MIP-1 α /CCL3, MIP-1 β /CCL4, TNF α , VEGF-A

Samples: Serum, plasma and cell culture supernatant


Kit Formats: Select the analytes you want, choose the 15-plex premixed kit, or the custom premix option for an additional charge

Comparison of MILLIPLEX® vs. Traditional ELISAs

	MILLIPLEX® Assay	ELISA
Number of plates required	1	15
Total plate run time	1 hour	15 hours
Results per plate (duplicate samples)	570	40
Total sample used (duplicate samples)	25 μL (1:2 sample dilution)	750 μL
Lower limit of detection	~1 pg/mL	~1 pg/mL





Quantify up to 15 analytes per sample and identify trends faster



Bovine PBMCs (BioIVT, Hicksville, NY) were treated with LPS or Concanavalin A for 48 hr, after which, cell-free samples were collected and assayed with the Bovine Cytokine/Chemokine Panel 1 (n=3 mean ± SEM).

*Notes saturation on the standard curve for these sample groups.



Serum samples obtained from BioIVT (Hicksville, NY).

All with the trusted MILLIPLEX® assay specifications to guarantee consistent performance over time

Analyte	Standard Curve Range (pg/mL)	Sensitivity minDC	+2SD (pg/mL)	Precision		Accuracy (%)
				Intra-Assay (%CV)	Inter-Assay (%CV)	
IFNγ	0.1 – 2,000	0.08		<10%	<10%	94
IL-1α	0.3 – 5,000	0.36		<10%	<10%	94
IL-1β	2.6 – 40,000	4.93		<10%	<10%	91
IL-4	12.8 – 200,000	16.57		<10%	<10%	93
IL-6	2.6 – 40,000	11.23		<10%	<15%	95
IL-8/CXCL8	2.2 – 35,000	22.6		<10%	<10%	95
IL-10	0.96 – 15,000	1.05		<10%	<10%	96
IL-17A	0.6 – 10,000	0.67		<10%	<10%	103
IL-36RA/IL-1F5	1.28 – 20,000	1.23		<10%	<10%	98
IP-10/CXCL10	0.6 – 10,000	0.91		<10%	<10%	94
MCP-1/CCL2	1.92 – 30,000	2.74		<10%	<10%	95
MIP-1α/CCL3	12.8 – 200,000	98.54		<10%	<15%	96
MIP-1β/CCL4	2.2 – 35,000	2.32		<10%	<10%	97
TNFα	12.8 – 200,000	22.62		<10%	<15%	92
VEGF-A	1.28 – 20,000	1.18		<10%	<10%	97

To view MILLIPLEX® assays for protein detection across 9 species, visit SigmaAldrich.com/milliplex

To place an order or receive technical assistance

Order/Customer Service: SigmaAldrich.com/order
 Technical Service: SigmaAldrich.com/techservice
 Safety-related Information: SigmaAldrich.com/safetycenter

SigmaAldrich.com

Merck KGaA
 Frankfurter Strasse 250
 64293 Darmstadt, Germany

© 2019 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M and Millipore are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.



Lit. No. MK_FL2834EN Ver. 2.0 2019 - 23303 06/2019

Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx