ELISA Solution for Diagnosis of COVID-19

Enzyme-linked immunosorbent assays (ELISA) to detect antibodies against SARS-CoV-2, producing COVID-19, in human serum or plasma.

Vincell^Q II ON B E 2

COVID-19 ELISA IgM+IgA COVID-19 ELISA IgG

PRODUCT FEATURES

- Suitable for all automated ELISA systems.
- SARS-CoV-2 recombinant antigens.
- · Excellent sensitivity with proven PCR positive samples.
- Complete solution IgG and IgM+IgA kits available.
- Common protocol allows to perform IgG and IgM+IgA simoultaneously in the same run.
- Ready-to-use reagents.





COVID 48H 40>

VERO 24H 10X

COVID 48H 10X

VERO 24H 20X

COVID 48H 20X

VER0 24H 40X

SARS-CoV-2. Vircell culture collection

COVID-19 serology allows confirmation of infection in both symptomatic and asymptomatic patients. ELISA is a simple and easy to automate technique that allows high throughput tests at low cost.

IgM+IgA detection complements molecular tests when the viral load decreases or when false negatives are suspected due to unsuitable samples in patients with clear symptoms.

The combination of IgM+IgA significantly improves sensitivity since COVID-19 affects the mucous membrane. The early production of IgA antibodies in this kind of infections could be particularly high.

IgG detection allows monitoring and confirming recovery and possible immunity, when other molecular and serological markers become negative.

Excellent sensitivity with proven PCR positive samples: 87% were IgM+IgA positive and 83% were IgG positive after a period of 5 days with a positive PCR result

Days after PCR positive result	Number of Samples	% lgM+lgA Positive	% IgG Positive
< 1 d	84	35%	12%
1 d	142	54%	14%
2 d	67	48%	48%
3 d	45	63%	47%
4 d	39	87%	74%
5 d	23	87%	83%



Ordering information

Reference	Description	Kit contents
G1032	COVID-19 ELISA IgG	96 tests
MA1032	COVID-19 ELISA IgM+IgA	96 tests





Distributed by Abacus dx 1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com