# **ImmPRESS<sup>™</sup> Polymer Detection Reagents**

Streamline your immunohistochemistry workflow

#### Integrated antibody and enzyme technology for simplified, convenient, sensitive IHC detection

The ImmPRESS<sup>™</sup> Polymer reagents combine affinity-purified, cross-adsorbed secondary antibodies with a high-density micropolymer array of reporter enzyme. This introduces significantly more enzyme at the site of primary antibody localization to yield a stronger signal—in a single incubation step.

## Why ImmPRESS<sup>™</sup>?

- High sensitivity
- Unique micropolymer technology
- Timesaving one-step protocols
- Non-biotin based
- Ready-to-use, pre-diluted reagents
- Convenient dropper bottles
- Choice of enzyme (HRP or AP)

## When to ImmPRESS<sup>™</sup>?

- Consistent and sensitive detection
- · Easy optimization and no risk of dilution error
- Enhanced signals and diminished background staining
- High-sensitivity, crisp detection of nuclear and membrane antigens
- Simplified multiple antigen labeling
- Absence of biotin interference
- Reagents for manual or automated staining


#### **ImmPRESS<sup>™</sup>** Polymer detection

Our unique polymer conjugation chemistry enables significantly more enzyme to be introduced at the antigen target site, resulting in greater signal intensity and sensitivity.


#### Enzyme conjugated secondary antibody detection

Standard conjugation methods limit the amount of enzyme available for reactivity and hence result in lower signal intensity and substantially lower sensitivity.

#### Be ImmPRESSed!



ImmPRESS<sup>™</sup> Anti-Mouse Ig and DAB (brown) staining using cytokeratin 8/18 primary antibody.Hematoxylin counterstain (blue).



Colon: Anti-CD34 (m), ImmPRESS™ Excel Amplified Anti-Mouse IgG Staining Kit, ImmPACT™ DAB EqV (brown). Hematoxylin QS counterstain (blue).



Human tonsil (paraffin section) stained for CD3 (DAB, brown) and AE1/AE3 cytokeratin (magenta) using ImmPRESS Duet Kit MP-7714.




#### vectorlabs.com


# Ways to ImmPRESS<sup>™</sup>: ImmPRESS<sup>™</sup> Polymer Detection Kit Options

# **One-Step Detection Kits**

Species-specific, enzyme conjugated secondary antibody polymers.










## **Two-Step Detection Kits**

Anti-mouse and anti-rabbit IgG HRP conjugated secondary antibody polymers with additional amplifying reagent.











## **Double Staining Kits**

Anti-mouse and anti-rabbit IgG HRP and AP conjugated secondary polymer cocktail.







Add Enzyme Substrate

#### Select for:

- Single Antigen Staining
- Moderate to high antigen expression
- Reduce/eliminate species cross-reactivity issues

#### Select for:

- Unknown or low antigen expression
- Maximize primary antibody dilution
- Convenience of complete kit format

#### Select for:

- Double antigen staining
- Reduce, optimization requirements
- Streamline the procedure

# **Species on Species Detection**

Mouse on Mouse (M.O.M.™) ImmPRESS™ HRP Polymer Kit

Specific blocking and detection reagents for unambiguous visualization of mouse primary antibodies on mouse tissue sections.

#### Select for:

- Transgene and knock-out mice models
- Genetically modified mice
- Mouse xenograft tissue
- Normal mouse tissue

For Superior Results

Use

ImmPRESS™ Polymer Reagents

In Combination with

Vector® ImmPACT™ Enzyme Substrates

# **ImmPRESS<sup>™</sup>** Kits and Reagents

Target Species	Product	Alkaline Phosphatase (AP)	Horseradish Peroxidase (HRP)	Veterinary HRP Reagent (VR) <sup>1</sup>	Excel Amplified (HRP )²
Rabbit	ImmPRESS™ Anti-Rabbit IgG Polymer Detection Kits				
	Made in horse	<u>MP-5401</u>	<u>MP-7401</u>	<u>MP-6401</u>	<u>MP-7601</u>
	Made in goat		<u>MP-7451</u>		
Mouse	ImmPRESS™ Anti-Mouse IgG Polymer Detection Kits				
	Made in horse	<u>MP-5402</u>	<u>MP-7402</u>	<u>MP-6402</u>	<u>MP-7602</u>
	Made in horse (rat adsorbed)		<u>MP-7422</u>		
	Made in goat		<u>MP-7452</u>		
	Mouse on Mouse (M.O.M.™) ImmPRESS™ HRP Polymer Kit <sup>3</sup>		<u>MP-2400</u>		
	Mouse on Mouse (M.O.M.™) ImmPRESS™ HRP Polymer Anti-Mouse Reagent		<u>MPX-2402</u>		
Rabbit and Mouse	ImmPRESS™ Universal Antibody (Anti-Mouse IgG/Anti-Rabbit IgG) Polymer Detection Kit (made in horse)⁴		<u>MP-7500</u>		
	ImmPRESS™ Duet Double Staining HRP/AP Polymer Kit (anti-rabbit IgG-brown, anti-mouse IgG-red) <sup>5</sup>	II			
	ImmPRESS™ Duet Double Staining HRP/AP Polymer Kit (anti-mouse IgG-brown, anti-rabbit IgG-red) <sup>s</sup>	II			
Rat	ImmPRESS™ Anti-Rat IgG Polymer Detection Kit				
	Made in goat	<u>MP-5404</u>	<u>MP-7404</u>		
	Made in goat (mouse adsorbed)	<u>MP-5444</u>	<u>MP-7444</u>		
Goat	ImmPRESS™ Anti-Goat IgG Polymer Detection Kit		<u>MP-7405</u>		
Vector <sup>®</sup> ImmPA	CT™ Enzyme Substrates				
ImmPACT™ Vector® Red (Magenta)		<u>SK-5105</u>			
ImmPACT™ DAB EqV (Brown)			II		
ImmPACT™ DAB (Brown)			II		
ImmPACT™ AEC (Red)			II		
ImmPACT™ AMEC Red (Red)			II		
ImmPACT™ VIP (Purple)			II		
ImmPACT™ SG (blue/gray)			II		
ImmPACT™ NovaRed™ (Red)			I <u>SK-</u>	<u>1805</u> I	

<sup>1</sup>Eliminate off-target binding, increase signal specificity in bovine, goat, sheep, swine, horse, cat, dog, rabbit, rat, and mouse tissues.

 $^{2}$ Maximize sensitivity with a two-step protocol using an amplifier antibody; includes ImmPACT<sup>IM</sup> DAB EqV substrate; not recommended for bovine, goat, or sheep tissue (amplifier antibody made in goat, may cross-react with endogenous IgG).

<sup>3</sup>Designed for localization of mouse primary antibodies on mouse tissue samples.

<sup>4</sup>Contains both anti-rabbit and -mouse IgG, for laboratories using both types of primary antibodies (eliminates need to purchase two separate kits).

<sup>5</sup>Designed for use in dual antigen labeling studies; includes ImmPACT™ DAB and ImmPACT™ Vector® Red substrates.



#### **Distributed by Abacus ALS**