

VectaFluor™ Duet

Immunofluorescence Double Labeling Kits

Apply two colors in one step using the VectaFluor™ Duet Immunofluorescence Double Labeling Kit. These kits are designed to save time and effort in double labeling immunofluorescence (IF) protocols that can potentially be long and tedious. The kits are configured to detect a mouse and a rabbit primary antibody present on the species with green and red fluorescence in one step.

Two kit configurations are available as indicated in the Table to the right.

Selection of a VectaFluor™ Duet format is an individual preference. However, consideration should be given to certain parameters such as prevalence of the respective target antigens within a tissue section, and whether the more abundant antigen will be viewed with a green or red signal. Consider also possible overlap or col-localization of the antigens and which combination may yield the most optimal results.

Each VectaFluor™ Duet Kit Includes:

- 15 ml VectaFluor™ Duet Reagent
- 15 ml 2.5% Normal Horse Serum for blocking (ready-to-use)


Each kit is sufficient to label approximately 150 tissue sections.

Our affinity purified, extensively cross-adsorbed secondary antibodies are conjugated to DyLight® dyes in a manner that ensures a maximum degree of labeling without compromising antibody affinity or specificity. The red and green pair of DyLight® anti-mouse/anti-rabbit antibodies are then combined into a robust, stable cocktail formulation that yields sensitive, consistent dual staining. VectaFluor™ Duet Kit is compatible with fluorescence staining of cells and tissues.


Key advantages

- Two colors in one step
- Ready-To-Use
- Robust cocktail formulation of DyLight®


Product	Catalog Number
VectaFluor Duet IF Double Labeling Kit ■ DyLight 488 Anti-Rabbit (green) ■ DyLight 594 Anti-Mouse (red)	DK-8818
VectaFluor Duet IF Double Labeling Kit ■ DyLight 594 Anti-Rabbit (red) ■ DyLight 488 Anti-Mouse (green)	DK-8828



Prostate: Rabbit Anti-PSA (MAb) and Mouse Anti-Smooth Muscle Actin detected simultaneously with VectaFluor™ Duet, DyLight® 488 Anti-Rabbit (green)/DyLight® 594 Anti-Mouse (red). Mounted in VECTASHIELD® HardSet™ Mounting Medium with DAPI.



Prostate: Rabbit Anti-PSA (MAb) and Mouse Anti-Smooth Muscle Actin detected simultaneously with VectaFluor™ Duet, DyLight® 488 Anti-Rabbit (green)/DyLight 594 Anti-Mouse (red). Mounted in VECTASHIELD® HardSet™ Mounting Medium.



Colon: Mouse Anti-Cytokeratin (AE1/AE3) and Rabbit Anti-Vimentin detected simultaneously with VectaFluor™ Duet, DyLight® 488 Anti-Mouse (green)/DyLight® 594 Anti-Rabbit (red). Mounted in VECTASHIELD® HardSet™ Mounting Medium.

VectaFluor™


Ready-to-Use DyLight® Labeled Secondary Antibodies

VectaFluor™ Ready-to-Use (R.T.U.) DyLight® dye-labeled Secondary Antibodies offer maximum convenience for fluorescence staining of cells and tissues. Our affinity purified, highly cross-adsorbed secondary antibodies are conjugated to DyLight® dyes in a manner that ensures maximum degree of labeling without compromising antibody affinity or specificity. DyLight® dyes offer several advantages such as brighter fluorescence, greater photostability, and pH independence.



Each VectaFluor™ R.T.U. Secondary Antibody Kit Includes:

- 15 ml pre-diluted, stabilized, ready-to-use VectaFluor™ DyLight® Dye-labeled Secondary Antibody
- 15 ml 2.5% normal horse serum for blocking (ready-to-use).



Each kit is sufficient to label approximately 150 tissue sections.



One step detection using VectaFluor™ DyLight® conjugated secondary antibodies.



Prostate: Anti-PSA (goat polyclonal), VectaFluor™ DyLight® 488 Anti-Goat IgG. Mounted in VECTASHIELD® HardSet™ Mounting Medium with DAPI.



Prostate: Anti-PSA (goat polyclonal), VectaFluor™ DyLight® 594 Anti-Goat IgG. Mounted in VECTASHIELD® HardSet™ Mounting Medium with DAPI.

VectaFluor™ R.T.U. DyLight® dye-labeled Antibody Conjugates:

Product	DyLight® 488 (Green)	DyLight® 594 (Red)
VectaFluor™ Anti-Rabbit IgG, made in horse	DI-1788	DI-1794
VectaFluor™ Anti-Mouse IgG, made in horse	DI-2788	DI-2794
VectaFluor™ Anti-Goat IgG, made in horse	DI-3788	DI-3794

DyLight is a trademark of Thermo Fisher Scientific

Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx