



Product Catalogue

In Situ Hybridization

Arrays

EXCELLENT

PRODUCTS, SOLUTIONS AND SYSTEMS

FOR **CLINICAL CANCER DIAGNOSTICS**



Contact Details

We help you get in touch

Reach out to our international sales managers and they will bring you in contact with the right partner!

order-ISH-distribution@zytomics.com

We value your opinion and would be delighted to take your suggestions into account in future projects.

Please feel free to send us your feedback at: marketing@zytomics.com

Content

	Introduction	4
	In situ Hybridization	7
	Arrays	167

ZYTOMICS at a Glance

Founded in 2021, **ZYTOMICS** is a combination of three established and longstanding European cancer diagnostics companies: Zytomed Systems, ZytoVision and Diagomics. All original founders of these companies shared the same vision: enabling clinical pathologists to provide a precise and timely cancer diagnosis. Unified under ZYTOMICS, we proudly carry on the vision of our founders.

ZYTOMICS is a leading manufacturer and provider of products and solutions used in clinical cancer diagnostics. Our offering includes reagents, automation and analytics with a focus on immunohistochemistry, in situ hybridization and next generation sequencing diagnostics. We aim to be the partner of trust for clinical pathology laboratories by providing best-in-class product quality, regulatory compliance, clinical and scientific insights as well as supply reliability. We are certified in accordance with Regulation (EU) 2017/746 (IVDR) for a significant portion of our Class C product portfolio.

ZYTOMICS has established a leading presence in major European markets and is consistently growing beyond European borders. Operations are located across sites in Northern Germany, Switzerland, and Southern France. Our international distribution network covers more than 90 countries.

Our Mission

Be the partner of trust for clinical pathology laboratories by providing best-in-class product quality, regulatory compliance, clinical and scientific insights as well as supply reliability.

Our Vision


Aid clinical pathologists to provide a precise and timely cancer diagnosis.

Our Values

We believe in people empowerment to make the best decisions confidently, we embrace innovation to help our customers and staff achieving their goals in the most efficient and impactful manner, we foster collaborative diversity to come up with the most innovative ideas and solutions, and we promote continuous improvement to achieve lifelong professional growth as one team.



Product Portfolio

In Situ Hybridization	Arrays
	
Fluorescence In Situ Hybridization Chromogenic In Situ Hybridization	DNA Microarray-Chips

How does the catalogue work?

The products are divided into two categories based on the technology used and —as shown above —each is identified by a specific color.

In all sections, you will find guidance on product categorization to help you navigate quickly and easily. To find a specific product, we recommend consulting the table of contents to locate the appropriate category and then navigating directly to the desired product.

A digital version of the catalog is also available. You can request this from the Marketing Communications Department (marketing@zytomics.com), our Customer Support team, or your local sales representative.

Our Brand



Our Product Lines

ZytoLight®	FlexISH®	ZytoMation®
ZytoFast®	ZytoDot®	VisionArray®

IVDR Certification

We are excited to share that **more than 200 products** from our portfolio are **IVDR certified!**

As a manufacturer of CE/IVD products, we consider compliance with the In Vitro Diagnostics Regulation (EU) 2017/746 (IVDR) to be of utmost importance - and not merely a challenge. We aim to meet the highest quality standards, and we are committed to make a valuable contribution to improving patient care and safety.

By following our mission, to be your partner of trust and In line with our commitment to conformity and compliance, we remain focused on extending our IVDR certified portfolio!

On our website, you will find products that are already IVDR-marked and comply with the latest and most stringent EU requirements.



Scan the QR code or visit:
www.zytomics.com/ivdr-certified-portfolio



Product Availability and Regulatory Status

Product availability and regulatory status may vary by country and are subject to change. For detailed information on availability and the applicable regulations in your region, please visit our website or contact your local sales representative.

You can find the latest product information, portfolio expansions, and changes and new additions on our website at any time. Please note that the regulatory status of our products is also subject to change. Therefore, the most up-to-date information on our website is always authoritative.

Resources Area

We are pleased to introduce the new Resources section on our website. This section provides you with convenient, centralized access to user manuals, safety data sheets, and certificates of analysis for our ZytoVision brand products.

You can find the new Resources section here:
www.zytomics.com/de/resources

Alternatively, our customer support team is always available to assist you if you need help.

Abbreviations

CE/IVD = In Vitro Diagnostics
RUO = Research Use Only

In Situ Hybridization

ISI

Based on morphology, fast and easy. Trust what you see.

In situ hybridization, ISH, is a powerful method for visualizing nucleic acid sequences - DNA and RNA - in cells and tissues. The two main techniques used are fluorescence in situ hybridization, FISH, and chromogenic in situ hybridization, CISH.

ZYTOMICS’ extensive portfolio includes the outstanding *ZytoLight®*, *ZytoDot®*, *ZytoFast®*, *FlexISH®* and *ZytoMation®* probes from ZytoVision as well as the associated kits and detection systems.

In our product portfolio, you will find the following categories under in situ hybridisation:

	Fluorescence In Situ Hybridization	
	<i>ZytoLight®</i>	84
	Chromosome index	85
	Gene index	92
	Probes	98
	<i>FlexISH®</i>	192
	Chromosome index	193
	Gene index	194
	Probes	195
	<i>ZytoMation®</i>	202
	Chromosome index	203
	Gene index	204
	Probes	205
	Chromogenic In Situ Hybridization	
	<i>ZytoDot®</i>	212
	Chromosome index	214
	Gene index	217
	Probes	218
	<i>ZytoFast®</i>	234
	Virus Index	235
	mRNA Index	236
	Probes	237

ZytoLight®

ZytoLight® products are designed for the identification of genetic aberrations e.g. translocations, deletions, amplifications, and chromosomal aneuploidies by Fluorescence *in situ* Hybridization (FISH) in formalin-fixed, paraffin-embedded (FFPE) tissue sections or cytology specimens.

 ZytoVision GmbH

ZytoLight® SPEC and CEN Probes

ZytoLight® FISH probes are direct labeled FISH probes. ZytoLight® SPEC™ probes are designed for the detection of single copy human DNA sequences.


ZytoLight® CEN™ probes hybridize to highly repetitive human satellite DNA sequences of chromosomes producing signals specific for each individual chromosome.

ZytoLight® Kits – Convenient Solutions

For making FISH analysis reliable and user-friendly, all ZytoLight® FISH probes can be combined with the ZytoLight® FISH-Tissue Implementation Kit (Z-2028-5/-20) for FISH analyses on FFPE specimens and/or the ZytoLight® FISH-Cytology Implementation Kit (Z-2099-20) for FISH analyses on cytology specimens.

Both Implementation Kits include all necessary pretreatment solutions, wash buffers and DAPI/Dura Tect™-Solution and a detailed protocol to perform successful FISH experiments.

Method




The ZytoLight® system uses direct-labeled FISH probes [1], eliminating the need to detect the probes with fluorophore-coupled antibodies. The probes are detected by fluorescence microscopy using appropriate filter sets [2]. Due to an exciter filter [3] full-spectrum light, emitted by the microscope lamp [4] is reduced to light of a defined wavelength that specifically excites the fluorophore of the probe. This light is reflected onto the specimen by a dichroic mirror [5]. The fluorophore emits light of longer wavelengths that passes the mirror. Finally, a barrier filter [6] reduces the emitted light to a defined wavelength that can be detected.


Protocol Overview




Pretreatment




Pepsin Digestion




Apply Probe



Denature & Hybridize



Wash




Apply DAPI/DuraTect™




Analyze Slides


Chromosome Index




Chr. Band	Product Name	Reg. Status	Order No.	Page
1p36.3	ZytoLight Glioma 1p/19q Probe Set	CE/IVD	Z-2272-20	22
	ZytoLight SPEC 1p36/1q25 Dual Color Probe	CE/IVD	Z-2075-50/-200	23
	ZytoLight SPEC TNFRSF14/1q25 Dual Color Probe	RUO	Z-2323-50	113
1p36.1	ZytoLight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe	CE/IVD	Z-2019-50/-200	80
1p32.2	ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe	CE/IVD	Z-2276-50/-200	24
1p12	ZytoLight SPEC 1p12 Probe	RUO	Z-2101-200	110
	ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe	CE/IVD	Z-2102-200	32
1q21	ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe	CE/IVD	Z-2276-50/-200	24
	ZytoLight SPEC MCL1/1p12 Dual Color Probe	CE/IVD	Z-2173-200	25
1q23.1	ZytoLight SPEC NTRK1 Dual Color Break Apart Probe	CE/IVD	Z-2167-50/-200	25
1q23.3	ZytoLight SPEC NECTIN4/1p12 Dual Color Probe	CE/IVD	Z-2331-50/-200	26
	ZytoLight SPEC TCF3/PBX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2308-50	27
1q25.2	ZytoLight SPEC ABL2 Dual Color Break Apart Probe	CE/IVD	Z-2200-50	26
1q25.3	ZytoLight Glioma 1p/19q Probe Set	CE/IVD	Z-2272-20	22
	ZytoLight SPEC 1p36/1q25 Dual Color Probe	CE/IVD	Z-2075-50/-200	23
1q32.1	ZytoLight SPEC MDM4/1p12 Dual Color Probe	CE/IVD	Z-2080-200	28




Chr. Band	Product Name	Reg. Status	Order No.	Page
2p24	ZytoLight SPEC MYCN/2q11 Dual Color Probe	CE/IVD	Z-2074-50/-200	28
2p23	ZytoLight SPEC ALK/EML4 TriCheck™ Probe	CE/IVD	Z-2117-50/-200	29
	ZytoLight SPEC ALK Dual Color Break Apart Probe	CE/IVD	Z-2124-50/-200	30
	ZytoLight SPEC ALK/2q11 Dual Color Probe	CE/IVD	Z-2161-200	30
2p21	ZytoLight SPEC EML4 Dual Color Break Apart Probe	RUO	Z-2136-50	113
	ZytoLight SPEC ALK/EML4 TriCheck™ Probe	CE/IVD	Z-2117-50/-200	29
2p11.2	ZytoLight SPEC IGK Dual Color Break Apart Probe	CE/IVD	Z-2288-50	31
2q11.2	ZytoLight SPEC 2q11 Probe	RUO	Z-2049-200	110
	ZytoLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe	CE/IVD	Z-2118-200	33
2q36	ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe	CE/IVD	Z-2018-50/-200	79
	ZytoLight SPEC FOXO1/PAX3 TriCheck™ Probe	CE/IVD	Z-2185-50	79




Chr. Band	Product Name	Reg. Status	Order No.	Page
3p25	ZytoLight SPEC VHL/CEN 3 Dual Color Probe	CE/IVD	Z-2084-200	31
	ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe	CE/IVD	Z-2102-200	32
3p14.2	ZytoLight SPEC FHIT/CEN 3 Dual Color Probe	RUO	Z-2062-200	113
3p11.1-q11.1	ZytoLight Bladder Cancer Quadruple Color Probe	CE/IVD	Z-2305-50/-200	34
	ZytoLight CEN 3 Probe	RUO	Z-2001-200	110
	ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe	CE/IVD	Z-2081-50/-200	59
3q21	ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe	CE/IVD	Z-2287-50	36
3q25.1	ZytoLight SPEC WWTR1 Dual Color Break Apart Probe	CE/IVD	Z-2212-50	35
3q26.2	ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe	CE/IVD	Z-2287-50	36
	ZytoLight SPEC TERC/CEN 3 Dual Color Probe	RUO	Z-2284-200	113
3q26.3	ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe	CE/IVD	Z-2140-200	36
	ZytoLight SPEC SOX2/CEN 3 Dual Color Probe	RUO	Z-2127-200	113
	ZytoLight SPEC TP63/TBL1XR1 TriCheck™ Probe	CE/IVD	Z-2320-50	37
3q27	ZytoLight SPEC BCL6 Dual Color Break Apart Probe	CE/IVD	Z-2177-50/-200	37
3q28	ZytoLight SPEC TP63/TBL1XR1 TriCheck™ Probe	CE/IVD	Z-2320-50	37




Chr. Band	Product Name	Reg. Status	Order No.	Page
4p16.3	ZytoLight SPEC FGFR3 Dual Color Break Apart Probe	CE/IVD	Z-2170-50/-200	38
	ZytoLight SPEC FGFR3/4p11 Dual Color Probe	CE/IVD	Z-2082-200	38
	ZytoLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2282-50	39
4p11	ZytoLight SPEC 4p11 Probe	RUO	Z-2083-200	110
	ZytoLight SPEC 4p11/CEN 10/17 Triple Color Probe	CE/IVD	Z-2307-50	40
4q12	ZytoLight SPEC PDGFRA/FIP1L1 TriCheck™ Probe	CE/IVD	Z-2209-50	40




Chr. Band	Product Name	Reg. Status	Order No.	Page
5p15.3	ZytoLight SPEC TERT Dual Color Break Apart Probe	CE/IVD	Z-2273-50	41
	ZytoLight SPEC TERT/Sq31 Dual Color Probe	CE/IVD	Z-2091-50/-200	41
5p13.1	ZytoLight SPEC RICTOR/Sq31.1 Dual Color Probe	RUO	Z-2278-200	113
5q31.2	ZytoLight SPEC EGR1/D5S23,D5S721 Dual Color Probe	CE/IVD	Z-2211-50	42
5q32	ZytoLight SPEC CSF1R Dual Color Break Apart Probe	CE/IVD	Z-2202-50	42
	ZytoLight SPEC CSF1R/D5S23,D5S721 Dual Color Probe	CE/IVD	Z-2268-50	43
	ZytoLight SPEC NRG1/CD74 TriCheck™ Probe	CE/IVD	Z-2194-200	53
	ZytoLight SPEC PDGFRB Dual Color Break Apart Probe	CE/IVD	Z-2197-50	44




Chr. Band	Product Name	Reg. Status	Order No.	Page
6p25	ZytoLight SPEC IRF4,DUSP22 Dual Color Break Apart Probe	CE/IVD	Z-2210-50/-200	44
6p24	ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe	CE/IVD	Z-2152-50/-200	45
6p21.3	ZytoLight SPEC PHF1 Dual Color Break Apart Probe	CE/IVD	Z-2215-50	45
6p21.1	ZytoLight SPEC VEGFA/CEN 6 Dual Color Probe	RUO	Z-2195-200	110
6p11.1-q11	ZytoLight CEN 6 Probe	RUO	Z-2002-200	110
	ZytoLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe	CE/IVD	Z-2118-200	33
6q22.1	ZytoLight SPEC ROS1 Dual Color Break Apart Probe	CE/IVD	Z-2144-50/-200	46
	ZytoLight SPEC ROS1/CEN 6 Dual Color Probe	CE/IVD	Z-2162-200	47
6q23.3	ZytoLight SPEC MYB Dual Color Break Apart Probe	CE/IVD	Z-2143-50/-200	47
	ZytoLight SPEC MYB/CEN 6 Dual Color Probe	CE/IVD	Z-2281-50	48
	ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe	CE/IVD	Z-2152-50/-200	45
6q25.1	ZytoLight SPEC ESR1/CEN 6 Dual Color Probe	CE/IVD	Z-2069-50	48




Chr. Band	Product Name	Reg. Status	Order No.	Page
7p15.2-p15.1	ZytoLight SPEC JAZF1 Dual Color Break Apart Probe	CE/IVD	Z-2132-50	49
7p12	ZytoLight SPEC IKZF1/CEN 7 Dual Color Probe	CE/IVD	Z-2304-50	49
7p11.2	ZytoLight SPEC EGFR/CEN 7 Dual Color Probe	CE/IVD	Z-2033-50/-200	50
7p11.1-q11.1	ZytoLight Bladder Cancer Quadruple Color Probe	CE/IVD	Z-2305-50/-200	34
	ZytoLight CEN 7 Probe	RUO	Z-2003-200	110
	ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe	CE/IVD	Z-2081-50/-200	59
	ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe	CE/IVD	Z-2102-200	32
7q11.23	ZytoLight SPEC Williams-Beuren Dual Color Probe	CE/IVD	Z-2302-50	50



Chr. Band	Product Name	Reg. Status	Order No.	Page
7q22	ZytoLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe	CE/IVD	Z-2214-50	51
7q31.2	ZytoLight SPEC MET/CEN 7 Dual Color Probe	CE/IVD	Z-2087-50/-200	51
7q34	ZytoLight SPEC BRAF Dual Color Break Apart Probe	CE/IVD	Z-2189-200	52
	ZytoLight SPEC BRAF/CEN 7 Dual Color Probe	CE/IVD	Z-2191-200	52
7q36	ZytoLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe	CE/IVD	Z-2214-50	51




Chr. Band	Product Name	Reg. Status	Order No.	Page
8p12	ZytoLight SPEC NRG1 Dual Color Break Apart Probe	CE/IVD	Z-2181-200	53
	ZytoLight SPEC NRG1/CD74 TriCheck™ Probe	CE/IVD	Z-2194-200	53
8p11.2	ZytoLight SPEC FGFR1 Dual Color Break Apart Probe	CE/IVD	Z-2168-50/-200	54
	ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe	CE/IVD	Z-2072-50/-200	54
8p11.1-q11.1	ZytoLight CEN 8 Probe	CE/IVD	Z-2004-50/-200	108
8q12	ZytoLight SPEC PLAG1 Dual Color Break Apart Probe	RUO	Z-2326-50	113
8q21.3	ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe	CE/IVD	Z-2112-50/-200	55
8q24.21	ZytoLight SPEC MYC Dual Color Break Apart Probe	CE-IVD	Z-2090-50/-200	56
	ZytoLight SPEC MYC/CEN 8 Dual Color Probe	CE/IVD	Z-2092-50/-200	56
	ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2105-50/-200	57




Chr. Band	Product Name	Reg. Status	Order No.	Page
9p24	ZytoLight SPEC CD274,PDCD1LG2/CEN 9 Dual Color Probe	CE/IVD	Z-2179-50/-200	57
	ZytoLight SPEC JAK2 Dual Color Break Apart Probe	CE/IVD	Z-2294-50	58
9p21	ZytoLight Bladder Cancer Quadruple Color Probe	CE/IVD	Z-2305-50/-200	34
	ZytoLight SPEC CDKN2A/CEN 9 Dual Color Probe	CE/IVD	Z-2063-50/-200	58
	ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe	CE/IVD	Z-2081-50/-200	59
9p13	ZytoLight SPEC PAX5 Dual Color Break Apart Probe	CE/IVD	Z-2300-50	60
9q12	ZytoLight CEN 9 Probe	RUO	Z-2067-200	110
9q21.3	ZytoLight SPEC NTRK2 Dual Color Break Apart Probe	CE/IVD	Z-2205-50/-200	60
9q22.3-q31	ZytoLight SPEC NR4A3 Dual Color Break Apart Probe	CE/IVD	Z-2145-50	61
9q34.1	ZytoLight SPEC ABL1 Dual Color Break Apart Probe	CE/IVD	Z-2199-50	61
	ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe	CE/IVD	Z-2111-50/-200	63
	ZytoLight SPEC NUP214 Dual Color Break Apart Probe	CE/IVD	Z-2265-50	62

Chr. Band	Product Name	Reg. Status	Order No.	Page
10p11.2	ZytoLight SPEC KIF5B Dual Color Break Apart Probe	RUO	Z-2131-50	113
10p11.1-q11.1	ZytoLight CEN 10 Probe	RUO	Z-2079-200	110
	ZytoLight SPEC 4p11/CEN 10/17 Triple Color Probe	CE/IVD	Z-2307-50	40
10q11.2	ZytoLight SPEC RET Dual Color Break Apart Probe	CE/IVD	Z-2148-50/-200	63
10q23.3	ZytoLight SPEC PTEN/CEN 10 Dual Color Probe	CE/IVD	Z-2078-50/-200	64
10q26.1	ZytoLight SPEC FGFR2 Dual Color Break Apart Probe	CE/IVD	Z-2169-50/-200	64
	ZytoLight SPEC FGFR2/CEN 10 Dual Color Probe	CE/IVD	Z-2122-200	65




Chr. Band	Product Name	Reg. Status	Order No.	Page
11p15.4	ZytoLight SPEC NUP98 Dual Color Break Apart Probe	CE/IVD	Z-2266-50	65
11p13	ZytoLight SPEC WT1 Dual Color Break Apart Probe	CE/IVD	Z-2142-50	66
11p11.2	ZytoLight SPEC SPI1 Dual Color Break Apart Probe	RUO	Z-2291-50	67
11p11.11-q11	ZytoLight CEN 11 Probe	CE/IVD	Z-2005-200	108
11q13.3	ZytoLight SPEC CCND1 Dual Color Break Apart Probe	CE/IVD	Z-2108-50/-200	67
	ZytoLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe	CE/IVD	Z-2118-200	33
	ZytoLight SPEC CCND1/CEN 11 Dual Color Probe	CE/IVD	Z-2071-50/-200	68
	ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2125-50/-200	68
11q21	ZytoLight SPEC MAML2 Dual Color Break Apart Probe	CE/IVD	Z-2014-50/-200	70
11q22.2	ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe	CE/IVD	Z-2146-50/-200	69
11q22.3	ZytoLight SPEC ATM/CEN 11 Dual Color Probe	CE/IVD	Z-2297-50	70
	ZytoLight SPEC ATM/CEN 12 Dual Color Probe	CE/IVD	Z-2296-50	71
	ZytoLight SPEC TP53/ATM Dual Color Probe	CE/IVD	Z-2159-50/-200	71
11q23.3	ZytoLight SPEC 11q gain/loss Triple Color Probe	CE/IVD	Z-2216-50	73
	ZytoLight SPEC KMT2A Dual Color Break Apart Probe	CE/IVD	Z-2193-50/-200	74
11q24.3	ZytoLight SPEC 11q gain/loss Triple Color Probe	CE/IVD	Z-2216-50	73
	ZytoLight SPEC EWSR1/FLI1 TriCheck™ Probe	CE/IVD	Z-2183-50	103




Chr. Band	Product Name	Reg. Status	Order No.	Page
12p13.3	ZytoLight SPEC ZNF384 Dual Color Break Apart Probe	CE/IVD	Z-2275-50	74
12p13.2	ZytoLight SPEC ETV6 Dual Color Break Apart Probe	CE/IVD	Z-2176-50/-200	76
	ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2157-50/-200	75
12p12.1	ZytoLight SPEC KRAS/CEN 12 Dual Color Probe	CE/IVD	Z-2115-200	76
12p11.1-q11	ZytoLight CEN 12 Probe	CE/IVD	Z-2050-200	108
	ZytoLight SPEC ATM/CEN 12 Dual Color Probe	CE/IVD	Z-2296-50	71
	ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe	CE/IVD	Z-2160-50/-200	72
12q13.3	ZytoLight SPEC DDIT3 Dual Color Break Apart Probe	CE/IVD	Z-2100-50/-200	77
12q14	ZytoLight SPEC CDK4/CEN 12 Dual Color Probe	CE/IVD	Z-2103-50/-200	77
12q15	ZytoLight SPEC MDM2/CEN 12 Dual Color Probe	CE/IVD	Z-2013-50/-200	78



Chr. Band	Product Name	Reg. Status	Order No.	Page
13q12.1	ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe	CE/IVD	Z-2095-50/-200	108
	ZytoLight SPEC 13/21 Dual Color Probe	CE/IVD	Z-2164-200	108
	ZytoLight Aneuploidy Panel 18/X/Y and 13/21	CE/IVD	Z-2279-20	108
	ZytoLight Aneuploidy Panel X/Y and 13/18/21	CE/IVD	Z-2104-5/-20	112
13q14.1	ZytoLight SPEC FOXO1 Dual Color Break Apart Probe	CE/IVD	Z-2139-50	78
	ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe	CE/IVD	Z-2018-50/-200	79
	ZytoLight SPEC FOXO1/PAX3 TriCheck™ Probe	CE/IVD	Z-2185-50	79
	ZytoLight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe	CE/IVD	Z-2019-50/-200	80
13q14.2	ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe	CE/IVD	Z-2160-50/-200	72
	ZytoLight SPEC D13S319/13q34 Dual Color Probe	CE/IVD	Z-2280-50/-200	73
	ZytoLight SPEC RB1/13q34 Dual Color Probe	CE/IVD	Z-2324-50/-200	80




Chr. Band	Product Name	Reg. Status	Order No.	Page
14q32.3	ZytoLight SPEC IGH Dual Color Break Apart Probe	CE/IVD	Z-2110-50/-200	157
	ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2114-50/-200	172
	ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2125-50/-200	144
	ZytoLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2282-50	115
	ZytoLight SPEC MAF/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2270-50	162
	ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2271-50	163
	ZytoLight SPEC MALT1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2325-50	171
	ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2105-50/-200	133



Chr. Band	Product Name	Reg. Status	Order No.	Page
15q11.2	ZytoLight SPEC Prader-Willi Dual Color Probe	CE/IVD	Z-2318-50	82
15q11.2-q12	ZytoLight SPEC Angelman Dual Color Probe	CE/IVD	Z-2319-50	
15q14	ZytoLight SPEC NUTM1 Dual Color Break Apart Probe	CE/IVD	Z-2208-200	83
15q24	ZytoLight SPEC Angelman Dual Color Probe	CE/IVD	Z-2319-50	
	ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe	CE/IVD	Z-2113-50/-200	84
	ZytoLight SPEC Prader-Willi Dual Color Probe	CE/IVD	Z-2318-50	
15q25	ZytoLight SPEC NTRK3 Dual Color Break Apart Probe	CE/IVD	Z-2206-50/-200	83



Chr. Band	Product Name	Reg. Status	Order No.	Page
16p13.3	ZytoLight SPEC CREBBP Dual Color Break Apart Probe	CE/IVD	Z-2267-50	85
16p11.2	ZytoLight SPEC FUS Dual Color Break Apart Probe	CE/IVD	Z-2130-50	85
16q22	ZytoLight SPEC CBFβ Dual Color Break Apart Probe	CE/IVD	Z-2207-50	88
16q23	ZytoLight SPEC MAF/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2270-50	86



Chr. Band	Product Name	Reg. Status	Order No.	Page
17p13	ZytoLight SPEC TP53/17q22 Dual Color Probe	CE/IVD	Z-2198-50	88
	ZytoLight SPEC TP53/ATM Dual Color Probe	CE/IVD	Z-2159-50/-200	71
	ZytoLight SPEC TP53/CEN 17 Dual Color Probe	CE/IVD	Z-2153-50/-200	89
	ZytoLight SPEC USP6 Dual Color Break Apart Probe	CE/IVD	Z-2151-50/-200	89
	ZytoLight SPEC YWHAЕ Dual Color Break Apart Probe	CE/IVD	Z-2175-50	90
17p11.1-q11.1	ZytoLight Bladder Cancer Quadruple Color Probe	CE/IVD	Z-2305-50/-200	34
	ZytoLight CEN 17 Probe	RUO	Z-2006-200	110
	ZytoLight SPEC 4p11/CEN 10/17 Triple Color Probe	CE/IVD	Z-2307-50	40
	ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe	CE/IVD	Z-2081-50/-200	59
17q12	ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe	CE/IVD	Z-2102-200	32
	ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe	CE/IVD	Z-2015-50/-200	90
	ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit	CE/IVD	Z-2020-5/-20	90
	ZytoLight CEN 17/SPEC ERBB2 Dual Color Probe	CE/IVD	Z-2077-50/-200	91
	ZytoLight SPEC ERBB2/D17S122 Dual Color Probe	CE/IVD	Z-2190-50/-200	91
	ZytoLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe	CE/IVD	Z-2093-50/-200	92

Chr. Band	Product Name	Reg. Status	Order No.	Page
17q21.2	ZytoLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe	CE/IVD	Z-2093-50/-200	92
	ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe	CE/IVD	Z-2113-50/-200	84
17q21.3	ZytoLight SPEC COL1A1 Dual Color Break Apart Probe	CE/IVD	Z-2121-200	92
	ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe	CE/IVD	Z-2116-50/-200	93
17q22	ZytoLight SPEC TP53/17q22 Dual Color Probe	CE/IVD	Z-2198-50	88

Chr. Band	Product Name	Reg. Status	Order No.	Page
18p11.1-q11.1	ZytoLight CEN 18 Probe	RUO	Z-2007-200	110
	ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe	CE/IVD	Z-2095-50/-200	108
	ZytoLight Aneuploidy Panel 18/X/Y and 13/21	CE/IVD	Z-2279-20	108
	ZytoLight Aneuploidy Panel X/Y and 13/18/21	CE/IVD	Z-2104-5/-20	108
18q11.2	ZytoLight SPEC SS18 Dual Color Break Apart Probe	CE/IVD	Z-2097-50/-200	93
	ZytoLight SPEC SS18/SSX1 TriCheck™ Probe	CE/IVD	Z-2184-50	94
18q21.3	ZytoLight SPEC 18/CEN X/Y Triple Color Probe	CE/IVD	Z-2163-200	108
	ZytoLight SPEC BCL2 Dual Color Break Apart Probe	CE/IVD	Z-2192-50/-200	96
	ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2114-50/-200	96
	ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe	CE/IVD	Z-2146-50/-200	69
	ZytoLight SPEC MALT1 Dual Color Break Apart Probe	CE/IVD	Z-2196-50/-200	95
	ZytoLight SPEC MALT1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2325-50	95






Chr. Band	Product Name	Reg. Status	Order No.	Page
19p13.3	ZytoLight Glioma 1p/19q Probe Set	CE/IVD	Z-2272-20	22
	ZytoLight SPEC 19q13/19p13 Dual Color Probe	CE/IVD	Z-2076-50/-200	23
	ZytoLight SPEC TCF3/PBX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2308-50	27
19q13.2	ZytoLight SPEC CIC Dual Color Break Apart Probe	CE/IVD	Z-2285-50	97
19q13.3	ZytoLight Glioma 1p/19q Probe Set	CE/IVD	Z-2272-20	22
	ZytoLight SPEC 19q13/19p13 Dual Color Probe	CE/IVD	Z-2076-50/-200	23
19q13.4	ZytoLight SPEC C19MC/19p13 Dual Color Probe	CE/IVD	Z-2274-50	97


Chr. Band	Product Name	Reg. Status	Order No.	Page
20q11.2	ZytoLight SPEC BCL2L1/CEN 20 Dual Color Probe	RUO	Z-2171-200	113
20q12	ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2271-50	87
20q12-q13.1	ZytoLight SPEC PTPRT/20q11 Dual Color Probe	CE/IVD	Z-2213-50	98

Chr. Band	Product Name	Reg. Status	Order No.	Page
21q22.1	ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe	CE/IVD	Z-2112-50/-200	55
	ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2157-50/-200	75
21q22.1-q22.2	ZytoLight SPEC 21q22 Probe	RUO	Z-2086-200	110
	ZytoLight SPEC 21/CEN X/Yq12 Triple Color Probe	CE/IVD	Z-2180-200	108
	ZytoLight SPEC 13/21 Dual Color Probe	CE/IVD	Z-2164-200	108
	ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe	CE/IVD	Z-2095-50/-200	108
	ZytoLight Aneuploidy Panel 18/X/Y and 13/21	CE/IVD	Z-2279-20	111
	ZytoLight Aneuploidy Panel X/Y and 13/18/21	CE/IVD	Z-2104-5/-20	112
21q22.2	ZytoLight SPEC ERG Dual Color Break Apart Probe	CE/IVD	Z-2138-200	99
	ZytoLight SPEC ERG/TMPRSS2 TriCheck™ Probe	CE/IVD	Z-2135-200	99
21q22.3	ZytoLight SPEC ERG/TMPRSS2 TriCheck™ Probe	CE/IVD	Z-2135-200	99

Chr. Band	Product Name	Reg. Status	Order No.	Page
22q11.2	ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe	CE/IVD	Z-2111-50/-200	63
	ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe	CE/IVD	Z-2299-50	100
	ZytoLight SPEC DiGeorge Triple Color Probe	CE/IVD	Z-2289-50	101
	ZytoLight SPEC IGL Dual Color Break Apart Probe	CE/IVD	Z-2286-50	102
	ZytoLight SPEC SMARCB1/22q12 Dual Color Probe	CE/IVD	Z-2178-50	102
22q12.2	ZytoLight SPEC EWSR1 Dual Color Break Apart Probe	CE/IVD	Z-2096-50/-200	103
	ZytoLight SPEC EWSR1/FLI1 TriCheck™ Probe	CE/IVD	Z-2183-50	103
22q13.1	ZytoLight SPEC PDGFB Dual Color Break Apart Probe	CE/IVD	Z-2119-50/-200	104
	ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe	CE/IVD	Z-2116-50/-200	93
22q13.3	ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe	CE/IVD	Z-2299-50	100

Chr. Band	Product Name	Reg. Status	Order No.	Page
Xp22.33	ZytoLight SPEC CRLF2 Dual Color Break Apart Probe	CE/IVD	Z-2201-50	105
Xp11.4	ZytoLight SPEC BCOR Dual Color Break Apart Probe	CE/IVD	Z-2310-50	106
Xp11.23	ZytoLight SPEC SS18/SSX1 TriCheck™ Probe	CE/IVD	Z-2184-50	94
	ZytoLight SPEC TFE3 Dual Color Break Apart Probe	CE/IVD	Z-2109-50/-200	104
Xp11.1-q11.1	ZytoLight CEN X Probe	RUO	Z-2008-200	110
	ZytoLight CEN X/Yq12 Dual Color Probe	CE/IVD	Z-2016-50/-200	108
	ZytoLight CEN X/Y Dual Color Probe	CE/IVD	Z-2120-200	108
	ZytoLight SPEC 18/CEN X/Y Triple Color Probe	CE/IVD	Z-2163-200	108
	ZytoLight SPEC 21/CEN X/Yq12 Triple Color Probe	CE/IVD	Z-2180-200	108
	ZytoLight Aneuploidy Panel 18/X/Y and 13/21	CE/IVD	Z-2279-20	108
	ZytoLight Aneuploidy Panel X/Y and 13/18/21	CE/IVD	Z-2104-5/-20	108





Chr. Band	Product Name	Reg. Status	Order No.	Page
Yp11.32	ZytoLight SPEC CRLF2 Dual Color Break Apart Probe	CE/IVD	Z-2201-50	105
Yp11.1-q11.1	ZytoLight CEN X/Y Dual Color Probe	CE/IVD	Z-2120-200	108
	ZytoLight SPEC 18/CEN X/Y Triple Color Probe	CE/IVD	Z-2163-200	108
Yq12	ZytoLight CEN Yq12 Probe	RUO	Z-2010-200	110
	ZytoLight CEN X/Yq12 Dual Color Probe	CE/IVD	Z-2016-50/-200	108
	ZytoLight SPEC 21/CEN X/Yq12 Triple Color Probe	CE/IVD	Z-2180-200	108
	ZytoLight Aneuploidy Panel 18/X/Y and 13/21	CE/IVD	Z-2279-20	108
	ZytoLight Aneuploidy Panel X/Y and 13/18/21	CE/IVD	Z-2104-5/-20	108

Gene Index

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
ABL1	ABL, c-ABL	ZytoLight SPEC ABL1 Dual Color Break Apart Probe	CE/IVD	Z-2199-50	61
		ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe	CE/IVD	Z-2111-50/-200	63
ABL2	ARG	ZytoLight SPEC ABL2 Dual Color Break Apart Probe	CE/IVD	Z-2200-50	26
ALK	CD246	ZytoLight SPEC ALK/EML4 TriCheck™ Probe	CE/IVD	Z-2117-50/-200	29
		ZytoLight SPEC ALK Dual Color Break Apart Probe	CE/IVD	Z-2124-50/-200	30
		ZytoLight SPEC ALK/2q11 Dual Color Probe	CE/IVD	Z-2161-200	30
ARSA	ASA, MLD	ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe	CE/IVD	Z-2299-50	100
ATM	ATA, TEL1	ZytoLight SPEC ATM/CEN 11 Dual Color Probe	CE/IVD	Z-2297-50	70
		ZytoLight SPEC ATM/CEN 12 Dual Color Probe	CE/IVD	Z-2296-50	71
		ZytoLight SPEC TP53/ATM Dual Color Probe	CE/IVD	Z-2159-50/-200	71
BCL2	Bcl-2, PPP1R50	ZytoLight SPEC BCL2 Dual Color Break Apart Probe	CE/IVD	Z-2192-50/-200	96
		ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2114-50/-200	96
BCL2L1	BCLX	ZytoLight SPEC BCL2L1/CEN 20 Dual Color Probe	RUO	Z-2171-200	113
BCL6	ZNF51, LAZ3	ZytoLight SPEC BCL6 Dual Color Break Apart Probe	CE/IVD	Z-2177-50/-200	37
BCOR	KIAA1575	ZytoLight SPEC BCOR Dual Color Break Apart Probe	CE/IVD	Z-2310-50	106
BCR	ALL, BCR1	ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe	CE/IVD	Z-2111-50/-200	63
BIRC3	C-IAP, MALT2	ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe	CE/IVD	Z-2146-50/-200	69
BRAF	BRAF1, NS7	ZytoLight SPEC BRAF Dual Color Break Apart Probe	CE/IVD	Z-2189-200	52
		ZytoLight SPEC BRAF/CEN 7 Dual Color Probe	CE/IVD	Z-2191-200	52
C19MC	-	ZytoLight SPEC C19MC/19p13 Dual Color Probe	CE/IVD	Z-2274-50	97
CBFB	PEBP2B	ZytoLight SPEC CBFB Dual Color Break Apart Probe	CE/IVD	Z-2207-50	88
CCND1	BCL1, PRAD1	ZytoLight SPEC CCND1 Dual Color Break Apart Probe	CE/IVD	Z-2108-50/-200	67
		ZytoLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe	CE/IVD	Z-2118-200	33
		ZytoLight SPEC CCND1/CEN 11 Dual Color Probe	CE/IVD	Z-2071-50/-200	68
		ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2125-50/-200	68

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
CD274	PD-L1, PDL1	ZytoLight SPEC CD274,PDCD1LG2/CEN 9 Dual Color Probe	CE/IVD	Z-2179-50/-200	57
CD74	-	ZytoLight SPEC NRG1/CD74 TriCheck™ Probe	CE/IVD	Z-2194-200	53
CDK4	PSK-J3	ZytoLight SPEC CDK4/CEN 12 Dual Color Probe	CE/IVD	Z-2103-50/-200	77
CDKN2A	p16, ARF, INK4	ZytoLight Bladder Cancer Quadruple Color Probe	CE/IVD	Z-2305-50/-200	34
		ZytoLight SPEC CDKN2A/CEN 9 Dual Color Probe	CE/IVD	Z-2063-50/-200	58
		ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe	CE/IVD	Z-2081-50/-200	59
CDNK2C	-	ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe	CE/IVD	Z-2276-50/-200	24
CIC	KIAA0306	ZytoLight SPEC CIC Dual Color Break Apart Probe	CE/IVD	Z-2285-50	97
CKS1B	-	ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe	CE/IVD	Z-2276-50/-200	24
COL1A1	OI4	ZytoLight SPEC COL1A1 Dual Color Break Apart Probe	CE/IVD	Z-2121-200	92
		ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe	CE/IVD	Z-2116-50/-200	93
CREBBP	CBP, RTS	ZytoLight SPEC CREBBP Dual Color Break Apart Probe	CE/IVD	Z-2267-50	85
CRKL	-	ZytoLight SPEC DiGeorge Triple Color Probe	CE/IVD	Z-2289-50	101
CRLF2	CRL2, TSLPR	ZytoLight SPEC CRLF2 Dual Color Break Apart Probe	CE/IVD	Z-2201-50	105
CSF1R	FMS	ZytoLight SPEC CSF1R Dual Color Break Apart Probe	CE/IVD	Z-2202-50	42
		ZytoLight SPEC CSF1R/D5S23,D5S721 Dual Color Probe	CE/IVD	Z-2268-50	43
CUX1	CUT	ZytoLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe	CE/IVD	Z-2214-50	51
DDIT3	CHOP, GADD153	ZytoLight SPEC DDIT3 Dual Color Break Apart Probe	CE/IVD	Z-2100-50/-200	77
DLEU1	BCMS1, LEU1	ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe	CE/IVD	Z-2160-50/-200	72
		ZytoLight SPEC D13S319/13q34 Dual Color Probe	CE/IVD	Z-2280-50/-200	73
DUSP22	JKAP	ZytoLight SPEC IRF4,DUSP22 Dual Color Break Apart Probe	CE/IVD	Z-2210-50/-200	44
EGFR	HER1, ERBB1	ZytoLight SPEC EGFR/CEN 7 Dual Color Probe	CE/IVD	Z-2033-50/-200	50
EGR1	KROX-24	ZytoLight SPEC EGR1/D5S23,D5S721 Dual Color Probe	CE/IVD	Z-2211-50	42
ELN	WBS	ZytoLight SPEC Williams-Beuren Dual Color Probe	CE/IVD	Z-2302-50	50
EML4	ROPP120	ZytoLight SPEC EML4 Dual Color Break Apart Probe	RUO	Z-2136-50	113
		ZytoLight SPEC ALK/EML4 TriCheck™ Probe	CE/IVD	Z-2117-50/-200	29
ERBB2	HER2, HER-2, NEU	ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe	CE/IVD	Z-2015-50/-200	90
		ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit	CE/IVD	Z-2020-5/-20	90
		ZytoLight CEN 17/SPEC ERBB2 Dual Color Probe	CE/IVD	Z-2077-50/-200	91
		ZytoLight SPEC ERBB2/D17S122 Dual Color Probe	CE/IVD	Z-2190-50/-200	91
		ZytoLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe	CE/IVD	Z-2093-50/-200	92
ERG	erg-3, p55	ZytoLight SPEC ERG Dual Color Break Apart Probe	CE/IVD	Z-2138-200	99
		ZytoLight SPEC ERG/TMPRSS2 TriCheck™ Probe	CE/IVD	Z-2135-200	99
ESR1	Era, NR3A1	ZytoLight SPEC ESR1/CEN 6 Dual Color Probe	CE/IVD	Z-2069-50	48
ETV6	TEL	ZytoLight SPEC ETV6 Dual Color Break Apart Probe	CE/IVD	Z-2176-50/-200	76
		ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2157-50/-200	75
EWSR1	EWS	ZytoLight SPEC EWSR1 Dual Color Break Apart Probe	CE/IVD	Z-2096-50/-200	103
		ZytoLight SPEC EWSR1/FLI1 TriCheck™ Probe	CE/IVD	Z-2183-50	103
EZH2	KMT6A	ZytoLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe	CE/IVD	Z-2214-50	51

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
FGFR1	FLT2, BFGFR	ZytoLight SPEC FGFR1 Dual Color Break Apart Probe	CE/IVD	Z-2168-50/-200	54
		ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe	CE/IVD	Z-2072-50/-200	54
FGFR2	BEK, CD332	ZytoLight SPEC FGFR2 Dual Color Break Apart Probe	CE/IVD	Z-2169-50/-200	64
		ZytoLight SPEC FGFR2/CEN 10 Dual Color Probe	CE/IVD	Z-2122-200	65
FGFR3	CD333, JTK4	ZytoLight SPEC FGFR3 Dual Color Break Apart Probe	CE/IVD	Z-2170-50/-200	38
		ZytoLight SPEC FGFR3/4p11 Dual Color Probe	CE/IVD	Z-2082-200	38
		ZytoLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2282-50	39
FHIT	FRA3B	ZytoLight SPEC FHIT/CEN 3 Dual Color Probe	RUO	Z-2062-200	113
FIP1L1	FIP1	ZytoLight SPEC PDGFRA/FIP1L1 TriCheck™ Probe	CE/IVD	Z-2209-50	40
FLI1	EWSR2	ZytoLight SPEC EWSR1/FLI1 TriCheck™ Probe	CE/IVD	Z-2183-50	103
FOXO1	FKHR, FKH1	ZytoLight SPEC FOXO1 Dual Color Break Apart Probe	CE/IVD	Z-2139-50	78
		ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe	CE/IVD	Z-2018-50/-200	79
		ZytoLight SPEC FOXO1/PAX3 TriCheck™ Probe	CE/IVD	Z-2185-50	79
		ZytoLight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe	CE/IVD	Z-2019-50/-200	80
FUS	FUS1	ZytoLight SPEC FUS Dual Color Break Apart Probe	CE/IVD	Z-2130-50	85
GATA2	NFE1B	ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe	CE/IVD	Z-2287-50	36
HIRA	TUPLE1, TUP1	ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe	CE/IVD	Z-2299-50	100
		ZytoLight SPEC DiGeorge Triple Color Probe	CE/IVD	Z-2289-50	101
IGH	IGH@	ZytoLight SPEC IGH Dual Color Break Apart Probe	CE/IVD	Z-2110-50/-200	81
		ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2114-50/-200	96
		ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2125-50/-200	68
		ZytoLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2282-50	39
		ZytoLight SPEC MAF/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2270-50	86
		ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2271-50	87
		ZytoLight SPEC MALT1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2325-50	95
		ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2105-50/-200	57
IGK	IGK@	ZytoLight SPEC IGK Dual Color Break Apart Probe	CE/IVD	Z-2288-50	31
IGL	IGL@	ZytoLight SPEC IGL Dual Color Break Apart Probe	CE/IVD	Z-2286-50	102
IKZF1	IKAROS	ZytoLight SPEC IKZF1/CEN 7 Dual Color Probe	CE/IVD	Z-2304-50	49
IRF4	MUM1	ZytoLight SPEC IRF4,DUSP22 Dual Color Break Apart Probe	CE/IVD	Z-2210-50/-200	44
JAK2	JTK10	ZytoLight SPEC JAK2 Dual Color Break Apart Probe	CE/IVD	Z-2294-50	58
JAZF1	TIP27, ZNF802	ZytoLight SPEC JAZF1 Dual Color Break Apart Probe	CE/IVD	Z-2132-50	49
KIF5B	KNS1	ZytoLight SPEC KIF5B Dual Color Break Apart Probe	RUO	Z-2131-50	113
KMT2A	MLL	ZytoLight SPEC KMT2A Dual Color Break Apart Probe	CE/IVD	Z-2193-50/-200	74
KRAS	KRAS1	ZytoLight SPEC KRAS/CEN 12 Dual Color Probe	CE/IVD	Z-2115-200	761
LAMP1	CD107a	ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe	CE/IVD	Z-2160-50/-200	72
		ZytoLight SPEC D13S319/13q34 Dual Color Probe	CE/IVD	Z-2280-50/-200	73
MAF	-	ZytoLight SPEC MAF/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2270-50	86

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
MAFB	-	ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2271-50	87
MALT1	MLT	ZytoLight SPEC MALT1 Dual Color Break Apart Probe	CE/IVD	Z-2196-50/-200	95
		ZytoLight SPEC MALT1/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2325-50	95
		ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe	CE/IVD	Z-2146-50/-200	69
MAML2	MAM3	ZytoLight SPEC MAML2 Dual Color Break Apart Probe	CE/IVD	Z-2014-50/-200	70
MAPK1	PRKM2, ERK	ZytoLight SPEC DiGeorge Triple Color Probe	CE/IVD	Z-2289-50	101
MCL1	BCL2L3	ZytoLight SPEC MCL1/1p12 Dual Color Probe	CE/IVD	Z-2173-200	25
MDM2	HDM2	ZytoLight SPEC MDM2/CEN 12 Dual Color Probe	CE/IVD	Z-2013-50/-200	78
MDM4	MDMX	ZytoLight SPEC MDM4/1p12 Dual Color Probe	CE/IVD	Z-2080-200	28
MECOM	MDS1, EVI1	ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe	CE/IVD	Z-2287-50	36
MET	HGFR, RCCP2	ZytoLight SPEC MET/CEN 7 Dual Color Probe	CE/IVD	Z-2087-50/-200	51
MYB	c-myb	ZytoLight SPEC MYB Dual Color Break Apart Probe	CE/IVD	Z-2143-50/-200	47
		ZytoLight SPEC MYB/CEN 6 Dual Color Probe	CE/IVD	Z-2281-50	48
		ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe	CE/IVD	Z-2152-50/-200	45
MYC	CMYC, bHLHe39, c-Myc	ZytoLight SPEC MYC Dual Color Break Apart Probe	CE/IVD	Z-2090-50/-200	56
		ZytoLight SPEC MYC/CEN 8 Dual Color Probe	CE/IVD	Z-2092-50/-200	56
		ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe	CE/IVD	Z-2105-50/-200	57
MYCN	NMYC, N-myc	ZytoLight SPEC MYCN/2q11 Dual Color Probe	CE/IVD	Z-2074-50/-200	28
NECTIN4	PVRL4	ZytoLight SPEC NECTIN4/1p12 Dual Color Probe	CE/IVD	Z-2331-50/-200	26
NR4A3	CHN, CSMF	ZytoLight SPEC NR4A3 Dual Color Break Apart Probe	CE/IVD	Z-2145-50	61
NRG1	HGL, GGF	ZytoLight SPEC NRG1 Dual Color Break Apart Probe	CE/IVD	Z-2181-200	53
		ZytoLight SPEC NRG1/CD74 TriCheck™ Probe	CE/IVD	Z-2194-200	53
NTRK1	MTC, TRK	ZytoLight SPEC NTRK1 Dual Color Break Apart Probe	CE/IVD	Z-2167-50/-200	25
NTRK2	TRKB	ZytoLight SPEC NTRK2 Dual Color Break Apart Probe	CE/IVD	Z-2205-50/-200	60
NTRK3	TRKC	ZytoLight SPEC NTRK3 Dual Color Break Apart Probe	CE/IVD	Z-2206-50/-200	83
NUP98	NUP96	ZytoLight SPEC NUP98 Dual Color Break Apart Probe	CE/IVD	Z-2266-50	65
NUP214	CAN, CAIN	ZytoLight SPEC NUP214 Dual Color Break Apart Probe	CE/IVD	Z-2265-50	62
NUTM1	NUT	ZytoLight SPEC NUTM1 Dual Color Break Apart Probe	CE/IVD	Z-2208-200	83
PAX3	HUP2	ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe	CE/IVD	Z-2018-50/-200	79
		ZytoLight SPEC FOXO1/PAX3 TriCheck™ Probe	CE/IVD	Z-2185-50	79
PAX5	BSAP	ZytoLight SPEC PAX5 Dual Color Break Apart Probe	CE/IVD	Z-2300-50	60
PAX7	HUP1	ZytoLight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe	CE/IVD	Z-2019-50/-200	80
PBX1	-	ZytoLight SPEC TCF3/PBX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2308-50	27
PDCD1LG2	PD-L2, PDL2	ZytoLight SPEC CD274,PDCD1LG2/CEN 9 Dual Color Probe	CE/IVD	Z-2179-50/-200	57
PDGFB	SIS, SSV	ZytoLight SPEC PDGFB Dual Color Break Apart Probe	CE/IVD	Z-2119-50/-200	104
		ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe	CE/IVD	Z-2116-50/-200	93
PDGFRA	GAS9	ZytoLight SPEC PDGFRA/FIP1L1 TriCheck™ Probe	CE/IVD	Z-2209-50	40

ISH

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
PDGFRB	JTK12, PDGFR1	ZytoLight SPEC PDGFRB Dual Color Break Apart Probe	CE/IVD	Z-2197-50	44
PHF1	MTF2L2, PCL1	ZytoLight SPEC PHF1 Dual Color Break Apart Probe	CE/IVD	Z-2215-50	45
PIK3CA	PI3K	ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe	CE/IVD	Z-2140-200	36
PLAG1	ZNF912	ZytoLight SPEC PLAG1 Dual Color Break Apart Probe	RUO	Z-2326-50	113
PML	MYL, RNF71	ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe	CE/IVD	Z-2113-50/-200	84
PTEN	MMAC1, TEP1	ZytoLight SPEC PTEN/CEN 10 Dual Color Probe	CE/IVD	Z-2078-50/-200	64
PTPRT	KIAA0283	ZytoLight SPEC PTPRT/20q11 Dual Color Probe	CE/IVD	Z-2213-50	98
RARA	NR1B1, RAR	ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe	CE/IVD	Z-2113-50/-200	84
RB1	PPP1R130	ZytoLight SPEC RB1/13q34 Dual Color Probe	CE/IVD	Z-2324-50/-200	80
RET	HSCR1, CDHF12	ZytoLight SPEC RET Dual Color Break Apart Probe	CE/IVD	Z-2148-50/-200	63
RICTOR	AVO3, KIAA1999	ZytoLight SPEC RICTOR/5q31.1 Dual Color Probe	RUO	Z-2278-200	113
ROS1	MCF3, ROS	ZytoLight SPEC ROS1 Dual Color Break Apart Probe	CE/IVD	Z-2144-50/-200	46
		ZytoLight SPEC ROS1/CEN 6 Dual Color Probe	CE/IVD	Z-2162-200	47
RREB1	HNT	ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe	CE/IVD	Z-2152-50/-200	45
RUNX1	AML1, AMLCR1	ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe	CE/IVD	Z-2112-50/-200	55
		ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2157-50/-200	75
RUNX1T1	ETO, CDR, MTG8	ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe	CE/IVD	Z-2112-50/-200	55
SHANK3	prosap2	ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe	CE/IVD	Z-2299-50	100
SMARCB1	BAF47	ZytoLight SPEC SMARCB1/22q12 Dual Color Probe	CE/IVD	Z-2178-50	102
SNRPN	PWCR	ZytoLight SPEC Prader-Willi Dual Color Probe	CE/IVD	Z-2318-50	82
SOX2	ANOP3	ZytoLight SPEC SOX2/CEN 3 Dual Color Probe	RUO	Z-2127-200	113
SPI1	PU.1, SPI-A	ZytoLight SPEC SPI1 Dual Color Break Apart Probe	RUO	Z-2291-50	110
SS18	SYT, SSXT	ZytoLight SPEC SS18 Dual Color Break Apart Probe	CE/IVD	Z-2097-50/-200	93
		ZytoLight SPEC SS18/SSX1 TriCheck™ Probe	CE/IVD	Z-2184-50	94
SSX1	-	ZytoLight SPEC SS18/SSX1 TriCheck™ Probe	CE/IVD	Z-2184-50	94
TBL1XR1	IRA1	ZytoLight SPEC TP63/TBL1XR1 TriCheck™ Probe	CE/IVD	Z-2320-50	37
TCF3	E2A	ZytoLight SPEC TCF3/PBX1 Dual Color Dual Fusion Probe	CE/IVD	Z-2308-50	27
TERC	hTERC, TRC3	ZytoLight SPEC TERC/CEN 3 Dual Color Probe	RUO	Z-2284-200	113
TERT	EST2, TCS1	ZytoLight SPEC TERT Dual Color Break Apart Probe	CE/IVD	Z-2273-50	41
		ZytoLight SPEC TERT/5q31 Dual Color Probe	CE/IVD	Z-2091-50/-200	41
TFE3	TFEA	ZytoLight SPEC TFE3 Dual Color Break Apart Probe	CE/IVD	Z-2109-50/-200	104
TMPRSS2	PRSS10	ZytoLight SPEC ERG/TMPRSS2 TriCheck™ Probe	CE/IVD	Z-2135-200	99
TNFRSF14	HVEM	ZytoLight SPEC TNFRSF14/1q25 Dual Color Probe	RUO	Z-2323-50	113
TOP2A	TOP2	ZytoLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe	CE/IVD	Z-2093-50/-200	92
TP53	LSF1, TRP53	ZytoLight SPEC TP53/17q22 Dual Color Probe	CE/IVD	Z-2198-50	88
		ZytoLight SPEC TP53/ATM Dual Color Probe	CE/IVD	Z-2159-50/-200	71
		ZytoLight SPEC TP53/CEN 17 Dual Color Probe	CE/IVD	Z-2153-50/-200	89

ISH

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
TP63	TP73L	ZytoLight SPEC TP63/TBL1XR1 TriCheck™ Probe	CE/IVD	Z-2320-50	37
UBE3A	AS	ZytoLight SPEC Angelman Dual Color Probe	CE/IVD	Z-2319-50	81
USP6	Tre-2, TRE17	ZytoLight SPEC USP6 Dual Color Break Apart Probe	CE/IVD	Z-2151-50/-200	89
VEGFA	VEGF, VPF	ZytoLight SPEC VEGFA/CEN 6 Dual Color Probe	RUO	Z-2195-200	46
VHL	VHL1	ZytoLight SPEC VHL/CEN 3 Dual Color Probe	CE/IVD	Z-2084-200	31
		ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe	CE/IVD	Z-2102-200	32
WT1	AWT1	ZytoLight SPEC WT1 Dual Color Break Apart Probe	CE/IVD	Z-2142-50	66
WWTR1	TAZ	ZytoLight SPEC WWTR1 Dual Color Break Apart Probe	CE/IVD	Z-2212-50	35
YWHAE	14-3-3 epsilon	ZytoLight SPEC YWHAE Dual Color Break Apart Probe	CE/IVD	Z-2175-50	90
ZNF384	CIZ	ZytoLight SPEC ZNF384 Dual Color Break Apart Probe	CE/IVD	Z-2275-50	74

ZytoLight® Glioma 1p/19q Probe Set

Intended Purpose

The ZytoLight Glioma 1p/19q Probe Set is intended to be used for the qualitative detection of deletions involving the human chromosomal region 1p36.31 as well as deletions involving the human chromosomal region 19q13.32-q13.33 in formalin-fixed, paraffin-embedded specimens, such as gliomas, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of gliomas and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight Glioma 1p/19q Probe Set is a set comprising two separate probes and a quenching solution:

- ZytoLight SPEC 1p36/1q25 Dual Color Probe (Prod. No. Z-2075-200)
- ZytoLight SPEC 19q13/19p13 Dual Color Probe (Prod. No. Z-2076-200)
- ZyBlack Quenching Solution (Prod. No. BS-0002-8)

The ZytoLight SPEC 1p36/1q25 Dual Color Probe (PL34) is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1p36.31* (chr1:5,808,946-6,176,336) (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1q25.3* (chr1:184,271,714-184,986,522) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC 1p36 Probe map; Bottom: SPEC 1q25 Probe map (not to scale)

Results

Normal situation: In interphases of normal cells or cells without a deletion involving the 1p36 locus, two green and two orange signals appear. **Aberrant situation:** In a cell with deletion affecting the 1p36 locus, a reduced number of orange signals will be observed. Deletions affecting only parts of the 1p36 locus might result in a normal signal pattern with orange signals of reduced size.



Brain tissue section hybridized with the ZytoLight SPEC 1p36/1q25 Dual Color Probe without ZyBlack™ Quenching Solution.



Brain tissue section hybridized with the ZytoLight SPEC 1p36/1q25 Dual Color Probe with ZyBlack™ Quenching Solution.

The ZytoLight SPEC 19q13/19p13 Dual Color Probe (PL35) is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 19q13.32-q13.33* (chr19:47,857,776-48,374,564) .
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10.0 ng/μl), which target sequences mapping in 19p13.3* (chr19:658,555-1,144,465) .
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Fig. 2: Top: SPEC 19p13 Probe map; Bottom: SPEC 19q13 Probe map (not to scale)

Results

Normal situation: In interphases of normal cells or cells without a deletion involving the 19q13 locus, two green and two orange signals appear. **Aberrant situation:** In a cell with deletions affecting the 19q13 locus, a reduced number of orange signals will be observed. Deletions affecting only parts of the 19q13 locus might result in a normal signal pattern with orange signals of reduced size.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	1, 19	CE/IVD	20	Z-2272-20

Includes ZytoLight SPEC 1p36/1q25 Dual Color Probe, 0.2 ml; ZytoLight SPEC 19q13/19p13 Dual Color Probe, 0.2 ml; ZyBlack Quenching Solution, 8 ml

ZytoLight® SPEC 1p36/1q25 Dual Color Probe

Intended Purpose

The ZytoLight SPEC 1p36/1q25 Dual Color Probe (PL34) is intended to be used for the qualitative detection of deletions involving the human chromosomal region 1p36.31 as well as the detection of chromosome 1q25.3 specific sequences in formalin-fixed, paraffin-embedded specimens, such as glioma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of glioma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC 1p36/1q25 Dual Color Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1p36.31* (chr1:5,808,946-6,176,336) (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1q25.3* (chr1:184,271,714-184,986,522) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC 1p36 Probe map; Bottom: SPEC 1q25 Probe map (not to scale); Right: Ideogram of chromosome 1 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a deletion involving the 1p36 locus, two orange signals and two green signals appear.

Aberrant situation: In a cell with deletion affecting the 1p36 locus, a reduced number of orange signals will be observed. Deletions affecting only parts of the 1p36 locus might result in a normal signal pattern with orange signals of reduced size.



Hybridized to a glioma tissue section with 1p36 deletion as indicated by one orange signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	1	CE/IVD	5 (50 μl)	Z-2075-50
			20 (200 μl)	Z-2075-200

ZytoLight® SPEC 19q13/19p13 Dual Color Probe

Intended Purpose

The ZytoLight SPEC 19q13/19p13 Dual Color Probe (PL35) is intended to be used for the qualitative detection of deletions involving the human chromosomal region 19q13.32-q13.33 as well as the detection of chromosome 19p13.3 specific sequences in formalin-fixed, paraffin-embedded specimens, such as glioma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of glioma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC 19q13/19p13 Dual Color Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 19q13.32-q13.33* (chr19:47,857,776-48,374,564) (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 19p13.3* (chr19:658,555-1,144,465) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC 19p13 Probe map; Bottom: SPEC 19q13 Probe map (not to scale); Right: Ideogram of chromosome 19 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a deletion involving the 19q13 locus, two orange signals and two green signals appear.

Aberrant situation: In a cell with deletion affecting the 19q13 locus, a reduced number of orange signals will be observed. Deletions affecting only parts of the 19q13 locus might result in a normal signal pattern with orange signals of reduced size.



Hybridized to a glioma tissue section with 19q13 deletion as indicated by one orange signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	19	CE/IVD	5 (50 μl)	Z-2076-50
			20 (200 μl)	Z-2076-200

ZytoLight® SPEC CKS1B/CDKN2C Dual Color Probe

Intended Purpose

The ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe (PL232) is intended to be used for the qualitative detection of amplifications/gains involving the human CKS1B gene at 1q21.3-q22 and deletions involving the human CDKN2C gene at 1p32.3 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1q21.3-q22* (chr1:154,722,168-155,144,639) harboring the CKS1B gene region (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 1p32.3* (chr1:51,196,272-51,737,475) harboring the CDKN2C gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC CDKN2C Probe map; Bottom: SPEC CKS1B Probe map (not to scale); Right: Ideogram of chromosome 1 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without an amplification/gain involving the CKS1B gene region or a deletion involving the CDKN2C gene region, two orange signals and two green signals appear.

Aberrant situation: In cells with an amplification/gain of the CKS1B gene region, an increased number of orange signals will be observed. In cells with deletion affecting the CDKN2C gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the CDKN2C gene region might result in a normal signal pattern with green signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals and to metaphase chromosomes of a normal cell.



Example of an aberrant signal pattern: Bone marrow smear of a pediatric ALL case with amplification affecting the CKS1B locus as indicated by three or more orange signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	1	CE/IVD	5 (50 μl)	Z-2276-50
			20 (200 μl)	Z-2276-200

Material kindly provided by Paediatric Oncology/ Haematology, Charité – Universitätsmedizin Berlin.

ZytoLight® SPEC MCL1/1p12 Dual Color Probe

Intended Purpose

The ZytoLight SPEC MCL1/1p12 Dual Color Probe (PL129) is intended to be used for the qualitative detection of amplifications involving the human MCL1 gene as well as the detection of chromosome 1p12 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC MCL1/1p12 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1q21.3* (chr1:150,363,209-150,940,432) harboring the MCL1 gene region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1p12* (chr1:119,537,102-119,823,147) (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC MCL1 Probe map; Bottom: SPEC 1p12 Probe map (not to scale); Right: Ideogram of chromosome 1 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without an amplification involving the MCL1 gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the MCL1 gene region, an increased number of green signals or green signal clusters will be observed.



Example of an aberrant signal pattern: Paraffin-embedded H2110 cell line with interphase cells showing amplification of the MCL1 gene locus as indicated by multiple green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	1	CE/IVD	20 (200 μl)	Z-2173-200

ZytoLight® SPEC NTRK1 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC NTRK1 Dual Color Break Apart Probe (PL123) is intended to be used for the qualitative detection of translocations involving the human NTRK1 gene at 1q23.1 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC NTRK1 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1q22-q23.1* (chr1:156,245,849-156,781,745) proximal to the NTRK1 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1q23.1* (chr1:156,854,527-157,296,918) distal to the NTRK1 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC NTRK1 Probe map (not to scale) Right: Ideogram of chromosome 1 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the NTRK1 gene region, two green/orange fusion signals appear.

Aberrant situation: One NTRK1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated orange signals are the result of deletions proximal to the NTRK1 breakpoint region or are due to unbalanced translocations affecting this chromosomal region.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: Lung cancer tissue section with translocation of the NTRK1 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	1	CE/IVD	5 (50 μl)	Z-2167-50
			20 (200 μl)	Z-2167-200

Image kindly provided by Prof. Büttner, Cologne, Germany.

ZytoLight® SPEC NECTIN4/1p12 Dual Color Probe

Intended Purpose

The ZytoLight SPEC NECTIN4/1p12 Dual Color Probe (PL291) is intended to be used for the qualitative detection of amplifications involving the human NECTIN4 gene as well as the detection of chromosome 1p12 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC NECTIN4/1p12 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1q23.3* (Chr1:160,773,463-161,218,427) harboring the NECTIN4 gene region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1p12* (Chr1:119,537,102-119,823,147) (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC NECTIN4 Probe map; Bottom: SPEC 1p12 Probe map (not to scale); Right: Ideogram of chromosome 1 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without an amplification involving the NECTIN4 gene region, two green signals and two orange signals appear.

Aberrant situation: In cells with an amplification of the NECTIN4 gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: Hybridized to urothelial carcinoma tissue section with amplification of the NECTIN4 gene as indicated by green signal clusters.

Labels	Chromosome	Reg. Status	Tests	Order No.
	1	CE/IVD	5 (50 μl)	Z-2331-50
			20 (200 μl)	Z-2331-200

ZytoLight® SPEC ABL2 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC ABL2 Dual Color Break Apart Probe (PL158) is intended to be used for the qualitative detection of translocations involving the human ABL2 gene at 1q25.2 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ABL2 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 1q25.2* (chr1:179,141,608-179,659,091) distal to the ABL2 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1q25.2* (chr1:178,496,704-179,015,312) proximal to the ABL2 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ABL2 Probe map (not to scale) Right: Ideogram of chromosome 1 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the ABL2 gene region, two green/orange fusion signals appear.

Aberrant situation: One ABL2 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.



Example of an aberrant signal pattern: Blood smear with translocation of the ABL2 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	1	CE/IVD	5 (50 μl)	Z-2200-50

ZytoLight® SPEC TCF3/PBX1 Dual Color Dual Fusion Probe

Intended Purpose

The ZytoLight SPEC TCF3/PBX1 Dual Color Dual Fusion Probe is designed to detect rearrangements affecting the chromosomal region 19p13.3 harboring the TCF3 (transcription factor 3, a.k.a E2A) gene and the chromosomal region 1q23.3 harboring the PBX1 (PBX homeobox 1) gene.

TCF3 is the target of several known recurrent rearrangements in ALL that create TCF3 fusion proteins. The balanced t(1;19) (q23.3;p13.3) and the more common unbalanced der(19)t(1;19)(q23.3;p13.3), which occur in approximately 6% of pediatric B-ALL cases and in 20-25% of all pre-B-ALL cases, fuse the TCF3 gene to the PBX1 gene.

The t(17;19)(q22;p13.3) fuses TCF3 to the HLF gene in <1% of cases. TCF3-PBX1 and TCF3-HLF are chimeric transcription factors that contain the same portion of TCF3, including two transcriptional activation domains, fused to regions of PBX1 or HLF that contain unique DNA binding domains.

As a sole abnormality, t(1;19)/der(19)t(1;19) is associated with an intermediate prognosis. In the context of hyperdiploid B-ALL, this translocation is associated with a poor prognosis. TCF3-HLF fusion has an extremely poor prognosis.

In the revised 2016 WHO classification of myeloid neoplasms and acute leukemia, "B-lymphoblastic leukemia/lymphoma with t(1;19)(q23.3;p13.3);TCF3-PBX1" is classified as its own cytogenetic subgroup of ALL. Because more intensive therapy improves the outcome of patients with TCF3-PBX1 gene fusions, it is critical to identify this subset of patients so that appropriate therapy can be administered.

References

- Arber DA, *et al.* (2016) Blood 127: 2391-405.
Boomer T, *et al.* (2001) Leukemia 15: 95-102.
Mellentin JD, *et al.* (1989) Science 246: 379-82.
Rowsey RA, *et al.* (2019) Blood Cancer J 9 :81.
Shearer BM, *et al.* (2005) Br J Haematol 129: 45-52.
Tirado CA, *et al.* (2015) Biomark Res 3: 4.

Probe Description

The ZytoLight SPEC TCF3/PBX1 Dual Color Dual Fusion Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12.0 ng/μl), which target sequences mapping in 19p13.3* (chr19:1,152,432-2,233,487) harboring the TCF3 region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6.0 ng/μl), which target sequences mapping in 1q23.3* (chr1:164,223,543-164,979,228) harboring the PBX1 gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19




Fig. 1: Top: SPEC TCF3 Probe map; Bottom: SPEC PBX1 Probe map (not to scale)


Left: Ideogram of chromosome 19 indicating the hybridization locations

Right: Ideogram of chromosome 1 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.

Aberrant situation: A gene fusion is indicated by one separate green signal, one separate orange signal, and two green/orange fusion signals. A signal pattern showing one fusion signal, one separate green and two separate orange signals can occur due to an unbalanced translocation. Other relevant signal patterns may also be observed as a result of TCF3 rearrangements without the involvement of the PBX1 locus.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Bone marrow smear with rearrangement affecting the TCF3/PBX1 loci as indicated by two separate orange signals, one separate green signal, and one orange/green fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	1, 19	CE/IVD	5 (50 μl)	Z-2308-50

ZytoLight® SPEC MDM4/1p12 Dual Color Probe

Intended Purpose

The ZytoLight SPEC MDM4/1p12 Dual Color Probe (PL39) is intended to be used for the qualitative detection of amplifications involving the human MDM4 gene as well as the detection of chromosome 1p12 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight MDM4/1p12 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1q32.1* (chr1:204,126,022-204,882,307) harboring the MDM4 region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 1p12* (chr1:119,537,102-119,823,147) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC 1p12 Probe map; Bottom: SPEC MDM4 Probe map (not to scale) Right: Ideogram of chromosome 1 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without an amplification involving the MDM4 gene region, two green signals and two orange signals appear.

Aberrant situation: In cells with an amplification of the MDM4 gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	1	CE/IVD	20 (200 μl)	Z-2080-200

ZytoLight® SPEC MYCN/2q11 Dual Color Probe

Intended Purpose

The ZytoLight SPEC MYCN/2q11 Dual Color Probe (PL33) is intended to be used for the qualitative detection of amplifications involving the human MYCN gene as well as the detection of chromosome 2q11 specific sequences in formalin-fixed, paraffin-embedded specimens, such as neuroblastoma and medulloblastoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of neuroblastoma and medulloblastoma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight MYCN/2q11 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 2p24.3* (chr2:15,846,046-16,517,671) harboring the MYCN gene region.
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 2q11.2* (chr2:100,132,806-100,621,725).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC MYCN Probe map; Bottom: SPEC 2q11 Probe map (not to scale) Right: Ideogram of chromosome 2 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without an amplification involving the MYCN gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the MYCN gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	2	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2074-50 Z-2074-200

ZytoLight® SPEC ALK/EML4 TriCheck™ Probe

Intended Purpose

The ZytoLight SPEC ALK/EML4 TriCheck Probe (PL74) is intended to be used for the qualitative detection of rearrangements involving the human ALK gene at 2p23.1-p23.2 and the human EML4 gene at 2p21 in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ALK/EML4 TriCheck Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 2p23.1-p23.2* (chr2:29,460,144-30,095,822) proximal to the ALK breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 2p23.2* (chr2:29,174,204-29,383,335) distal to the ALK breakpoint region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/μl), which target sequences mapping in 2p21* (chr2:41,573,525-43,349,624) harboring the EML4 gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC ALK Probe map; Bottom: SPEC EML4 Probe map (not to scale)




Ideogram of chromosome 2 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells without rearrangement of the EML4-ALK gene region, two orange/green fusion signals and two blue signals appear.

Aberrant situation: An EML4-ALK inversion is indicated by one separate green signal, one separate orange signal, and an additional blue signal. The separate green and orange signal each co-localize with a blue signal. An ALK translocation not affecting EML4 is indicated by separated orange and green signals without an additional blue signal. An EML4-ALK inversion with deletion of 5'-ALK sequences is indicated by loss of one green signal and co-localization of the isolated orange signal with an additional blue signal.



SPEC ALK/EML4 TriCheck™ Probe on normal interphase cells with non-rearranged ALK loci (two orange/green fusion signals), and non-rearranged EML4 loci (two blue signals).



NSCLC tissue section with an EML4-ALK inversion as indicated by one green, one separated orange, and one additional blue signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	2	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2117-50 Z-2117-200

ZytoLight® SPEC ALK Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC ALK Dual Color Break Apart Probe (PL81) is intended to be used for the qualitative detection of translocations involving the human ALK gene at 2p23.1-p23.2 in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ALK Dual Color Break Apart Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 2p23.1-p23.2* (chr2:29,460,144-30,095,822) proximal to the ALK breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 2p23.2* (chr2:29,174,204-29,383,335) distal to the ALK breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without a translocation involving the ALK gene region, two green/orange fusion signals appear.

Aberrant situation: One ALK gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. An EML4-ALK inversion with deletion of 5'-ALK sequences is indicated by one or multiple isolated orange signals.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

NSCLC tissue section with translocation affecting the 2p23 locus as indicated by one orange/green fusion (non-rearranged) signal, one orange signal, and one separate green signal.

ZytoLight® SPEC ALK/2q11 Dual Color Probe

Intended Purpose

The ZytoLight SPEC ALK/2q11 Dual Color Probe (PL117) is intended to be used for the qualitative detection of amplifications involving the human ALK gene as well as the detection of chromosome 2q11 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ALK/2q11 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 2p23.1-p23.2* (chr2:29,460,144-30,095,822) harboring the ALK region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 2q11.2* (chr2:100,132,806-100,621,725) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC ALK Probe map; Bottom: SPEC 2q11 Probe map (not to scale) Right: Ideogram of chromosome 2 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without an amplification involving the ALK gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the ALK gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Example of an aberrant signal pattern: Neuroblastoma tissue section with tetrasomy of chromosome 2 as indicated by four orange (2q11) and four green (ALK) signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	2	CE/IVD	5 (50 μl)	Z-2124-50
			20 (200 μl)	Z-2124-200

ZytoLight® SPEC IGK Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC IGK Dual Color Break Apart Probe (PL243) is intended to be used for the qualitative detection of translocations involving the human IGK locus at 2p11.2 in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC IGK Dual Color Break Apart Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 2p11.2* (chr2:88,382,616-89,153,517) distal to the IGK breakpoint region (see Fig.1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 2p11.2* (chr2:89,246,977-89,609,390 and chr2:89,853,315-90,089,156) proximal to the IGK breakpoint region (see Fig. 1). Due to homologous sequence segments proximal to the IGK breakpoint region, the orange probe has two hybridization regions in close proximity.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without a translocation involving IGK locus, two green/orange fusion signals appear. Due to the two hybridization regions of the orange probe, orange signals may appear as paired signal dots.

Aberrant situation: One IGK locus affected by a translocation is indicated by one separate green signal and one separate orange signal. Due to the two hybridization regions of the orange probe, orange signals may appear as paired signal dots.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals in each nucleus and to metaphase chromosomes of a normal cell. Orange signals may appear as paired signal dots.

Example of an aberrant signal pattern: Burkitt lymphoma with an IGK translocation affecting the 2p11.2 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal (may appear as paired signal dots), and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	2	CE/IVD	20 (200 μl)	Z-2288-50

Specimen kindly provided by Dr. Brändle, Vienna, Austria.

ZytoLight® SPEC VHL/CEN 3 Dual Color Probe

Intended Purpose

The ZytoLight SPEC VHL/CEN 3 Dual Color Probe (PL43) is intended to be used for the qualitative detection of deletions involving the human VHL gene and the detection of chromosome 3 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC VHL/CEN 3 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 3p25.3* (chr3:10,051,220-10,598,496) harboring the VHL gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 3p11.1-q11.1 specific for the alpha satellite centromeric region D3Z1 of chromosome 3.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC VHL Probe map (not to scale) Right: Ideogram of chromosome 3 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a deletion involving the VHL gene region, two green and two orange signals appear.

Aberrant situation: In a cell with deletion of the VHL gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the VHL gene region might result in a normal signal pattern with green signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Example of an aberrant signal pattern: Trisomy of chromosome 3 as indicated by three orange (CEN3) and three green (VHL) signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	3	CE/IVD	20 (200 μl)	Z-2084-200

ZytoLight® SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe


Intended Purpose
The ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe (PL60) is intended to be used for the qualitative detection of deletions involving the human VHL gene as well as chromosome 1p12 specific sequences, and chromosome 7 and 17 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC VHL/1p12/7/17 Quadruple Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 3p25.3* (chr3:10,051,220-10,598,496) harboring the VHL region (see Fig. 1).
 - ZyRed (excitation 580 nm/emission 599 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 1p12* (chr1:119,537,102-119,823,147).
 - ZyGold (excitation 532 nm/emission 553 nm) labeled polynucleotides (~7 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19

Results
Normal situation: In interphases of normal cells or cells without aberration involving the respective gene regions, two green, two red, two gold, and two blue signals are expected.
Aberrant situation: In a cell with deletion affecting the VHL gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the VHL gene region might result in a normal signal pattern with green signals of reduced size. In cells with aneuploidy of chromosome 1, 7, or 17, more or less signals of the respective color will be visible.



Example of an aberrant signal pattern: Renal cell carcinoma tissue section with deletion of the VHL gene as indicated by one green signal in each nucleus.



Example of an aberrant signal pattern: Renal cell carcinoma tissue section with polysomy of the chromosome 7 and 17 as indicated by multiple gold and/or blue signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
   	1, 3, 7, 17	CE/IVD	20 (200 μl)	Z-2102-200





Fig. 1: Top: SPEC VHL Probe map; Bottom: SPEC 1p12 Probe map (not to scale)



Ideograms of chromosome 3, 1, 7 and 17 indicating the hybridization locations

ZytoLight® SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe


Intended Purpose
The ZytoLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe (PL75) is intended to be used for the qualitative detection of translocations involving the human CCND1 gene at 11q13.3 as well as for the detection of human chromosome 2q11 specific sequences as well as chromosome 6 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 11q13.2-q13.3* (chr11:68,249,010-68,705,283) proximal to CCND1 breakpoint region (see Fig. 1).
 - ZyRed (excitation 580 nm/emission 599 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 11q13.3* (chr11:69,453,301-70,031,240) distal to the CCND1 breakpoint region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37.0 ng/μl), which target sequences mapping in 2q11.2* (chr2:100,132,806-100,621,725) (see Fig. 1).
 - ZyGold (excitation 532 nm/emission 553 nm) labeled polynucleotides (~7 ng/μl), which target sequences mapping in 6p11.1-q11 specific for the alpha satellite centromeric region D6Z1 of chromosome 6.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the CCND1 gene region, two red/green fusion signals, two blue and two gold signals appear.
Aberrant situation: In a cell with translocation of the CCND1 gene region, a signal pattern consisting of one red/green fusion signal, one red, and a separate green signal indicates one normal CCND1 gene region and one CCND1 gene region affected by an 11q13.3 translocation. In cells with aneuploidy of chromosome 2 or 6, more or less signals of the respective color will be visible.



Example of an aberrant signal pattern: Renal cell carcinoma tissue section with translocation affecting the 11q13.3 locus as indicated by one non-rearranged red/green fusion signal, one red signal, and one separate green signal.



Example of an aberrant signal pattern: Renal cell carcinoma tissue section with monosomy of chromosome 2 and 6 as indicated by one blue and one gold signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
   	2, 6, 11	CE/IVD	20 (200 μl)	Z-2118-200





Fig. 1: Top: SPEC CCND1 Probe map; Bottom: SPEC 2q11 Probe map (not to scale)



Ideograms of chromosome 11, 2 and 6 indicating the hybridization locations

ZytoLight® Bladder Cancer Quadruple Color Probe

Intended Purpose
The ZytoLight Bladder Cancer Quadruple Color Probe is designed to detect CDKN2A (a.k.a. p16) deletions and aneuploidy of chromosomes 3, 7, and 17 in cytology specimens of tumors, e.g., in urine samples from patients with hematuria suspected of having bladder cancer (BC).

Moreover, it has been shown that the detection of CDKN2A deletions and/or aneusomies of chromosomes 3, 7, and/or 17 may be used for the surveillance of patients with a history of bladder cancer to early detect possible tumor recurrence.

BC represents the ninth most common cancer worldwide. About 430,000 new BC cases and 165,000 BC deaths occurred in 2012. Most of these tumors are non-invasive, well-differentiated, papillary tumors (pTa, low grade) and can be cured by endoscopic transurethral resection.

However, up to 70% of pTa and superficially invasive (pT1) tumors recur and of these, 15-30% are characterized by tumor progression. Therefore, a long-term follow-up of patients with BC is necessary.

The two standard methods used in the follow-up are either invasive (cystoscopy) or have a low sensitivity (cytology). BC cells are characterized by typical cytogenetic changes. Homozygous deletion of the CDKN2A gene at 9p21.3 and polysomy of chromosomes 3, 7, and/or 17 are common abnormalities observed in urothelial cell carcinoma, all of which can be detected by FISH. FISH on cells from urine has been shown to be highly sensitive and specific for detection of tumor cells in urine.

References
Antoni S, *et al.* (2017) Eur Urol 71: 96-108.
Dimashkieh H, *et al.* (2013) Cancer Cytopathol 121: 591-7.
Junker K, *et al.* (2006) Cytogenet Genome Res 114: 279-83.
Placer J, *et al.* (2002) Eur Urol 42: 547-52.

Probe Description
The ZytoLight Bladder Cancer Quadruple Color Probe is composed of:

- ZyGold (excitation 532 nm/emission 553 nm) labeled polynucleotides (~5,5 ng/μl), which target sequences mapping in 9p21.3* (chr9:21,742,619-22,056,863) harboring the CDKN2A gene region (see Fig. 1).
- ZyRed (excitation 580 nm/emission 599 nm) labeled polynucleotides (~0.5 ng/μl), which target sequences mapping in 3p11.1-q11.1 specific for the alpha satellite centromeric region D3Z1 of chromosome 3.
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC CDKN2A Probe map (not to scale)



Ideograms of chromosome 3, 7, 9 and 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the CDKN2A gene region, two gold, two red, two green, and two blue signals are expected.
Aberrant situation: In a cell with deletion of the CDKN2A gene region, a reduced number of gold signals will be observed. Deletions affecting only parts of the CDKN2A gene region might result in a normal signal pattern with gold signals of reduced size. In cells with aneusomy of chromosomes 3, 7, or 17, more or less signals of the respective color will be visible.



Interphase tumor cells with trisomy of chromosome 7 as indicated by three green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	3, 7, 9, 17	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2305-50 Z-2305-200

ZytoLight® SPEC WWTR1 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC WWTR1 Dual Color Break Apart Probe is designed for the detection of translocations involving the chromosomal region 3q25.1 harboring the WWTR1 (WW domain containing transcription regulator 1, a.k.a. TAZ) gene.

Epithelioid vascular tumors encompass a spectrum of diseases that includes epithelioid hemangioma (EH), a benign neoplasm, epithelioid hemangioendothelioma (EHE), a low to intermediate grade malignancy, and epithelioid angiosarcoma (EAS), a high grade malignancy. Although certain morphologic features allow to distinguish EHE from EH and EAS, the diagnosis can be challenging due to considerable morphologic overlap, particularly on small biopsies or when EAS lacks vasoformative properties. Clinical behavior and, consequently, treatment and prognosis vary significantly among vascular tumors. Therefore, it is paramount to effectively distinguish them from each other.

The recurrent translocation t(1;3)(p36.3;q25.1) was identified in approximately 90% of EHE cases, but not in other vascular tumors. t(1;3) results in the WWTR1-CAMTA1 fusion gene which encodes a putative chimeric transcription factor which is under the transcriptional control of the WWTR1 promoter. A recurrent YAP1-TFE3 gene fusion has been identified in WWTR1-CAMTA1 negative EHEs.

Thus, FISH analysis for the presence of WWTR1 translocation may serve as a useful molecular tool in the differential diagnosis of challenging cases.

References
Underson T, *et al.* (2015) Am J Surg Pathol 39: 132-9.
Errani C, *et al.* (2011) Genes Chromosomes Cancer 50: 644-53.
Mendlick MR, *et al.* (2001) Am J Surg Pathol 25: 684-7.
Puls F, *et al.* (2015) Virchows Arch 466: 473-8.
Tanas MR, *et al.* (2011) Sci Transl Med 3: 98ra82.

Probe Description
The ZytoLight SPEC WWTR1 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 3q24-q25.1* (chr3:148,533,200-149,234,601) proximal to the WWTR1 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 3q25* (chr3:149,430,325-149,933,565) distal to the WWTR1 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC WWTR1 Probe map (not to scale)




Ideogram of chromosome 3 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the WWTR1 gene region, two green/orange fusion signals appear.
Aberrant situation: One WWTR1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.



Epithelioid hemangioendothelioma interphase cells showing translocation of the WWTR1 gene as indicated by one non-rearranged orange/green fusion signal, one orange signal and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	3	CE/IVD	5 (50 μl)	Z-2212-50

ZytoLight® SPEC GATA2/MECOM Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe (PL242) is intended to be used for the qualitative detection of the translocation t(3;3)(q21.3;q26.2) as well as the inversion inv(3)(q21.3q26.2) involving the human GATA2 and MECOM genes in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6,0 ng/μl), which target sequences mapping in 3q21.3* (chr3:127,902,316-128,564,215) harboring the GATA2 gene region (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 3q26.2* (chr3:168,249,484-169,743,447) harboring the MECOM gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC GATA2 Probe map; Bottom: SPEC MECOM Probe map (not to scale); Right: Ideogram of chromosome 3 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals. Three green signals and two orange signals are the result of a translocation affecting the MECOM gene region without the involvement of the GATA2 locus.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	3	CE/IVD	5 (50 μl)	Z-2287-50

ZytoLight® SPEC PIK3CA/CEN 3 Dual Color Probe

Intended Purpose
The ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe (PL97) is intended to be used for the qualitative detection of amplifications involving the human PIK3CA gene as well as the detection of chromosome 3 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 3q26.32-q26.33* (chr3:178,535,986-179,293,464) harboring the PIK3CA gene region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 3p11.1-q11.1 specific for the alpha satellite centromeric region D3Z1 of chromosome 3.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC PIK3CA Probe map (not to scale); Right: Ideogram of chromosome 3 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the PIK3CA gene region, two green and two orange signals appear.
Aberrant situation: In cells with an amplification of the PIK3CA gene region, an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	3	CE/IVD	20 (200 μl)	Z-2140-200

ZytoLight® SPEC TP63/TBL1XR1 TriCheck™ Probe


Intended Purpose
The ZytoLight SPEC TP63/TBL1XR1 TriCheck™ Probe (PL274) is intended to be used for the qualitative detection of rearrangements involving the human TP63 gene at 3q28 and the human TBL1XR1 gene at 3q26.32 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

Interpretation of the results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

Probe Description
The ZytoLight SPEC TP63/TBL1XR1 TriCheck™ Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10.0 ng/μl), which target sequences mapping in 3q28* (chr3:189,559,557-190,097,196) distal to the TP63 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 3q28* (chr3:188,995,562-189,305,431) proximal to the TP63 breakpoint region (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37.0 ng/μl), which target sequences mapping in 3q26.32* (chr3:176,217,831-177,284,492) harboring the TBL1XR1 gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19



Top: SPEC TP63 Probe map; Bottom: SPEC TBL1XR1 Probe map (not to scale); Right: Ideogram of chromosome 3 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells without rearrangement of the TBL1XR1-TP63 gene region, two orange/green fusion signals and two blue signals appear.
Aberrant situation: A TBL1XR1-TP63 inversion is indicated by one separate green signal, one separate orange signal, and an additional blue signal. The separate green and orange signal each co-localize with a blue signal. A TP63 translocation not affecting TBL1XR1 is indicated by separated orange and green signals without an additional blue signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
  	3	CE/IVD	5 (50 μl)	Z-2320-50

ZytoLight® SPEC BCL6 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC BCL6 Dual Color Break Apart Probe (PL136) is intended to be used for the qualitative detection of translocations involving the human BCL6 gene at 3q27.3 in formalin-fixed, paraffin-embedded specimens, such as B-cell lymphoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of B-cell lymphoma and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC BCL6 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 3q27.3* (chr3:186,737,897-187,403,834) proximal to the BCL6 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 3q27.3-q28* (chr3:187,744,962-188,097,195) distal to the BCL6 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC BCL6 Probe map (not to scale); Right: Ideogram of chromosome 3 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the BCL6 gene region, two green/orange fusion signals appear.
Aberrant situation: One BCL6 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	3	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2177-50 Z-2177-200

ZytoLight® SPEC FGFR3 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC FGFR3 Dual Color Break Apart Probe (PL126) is intended to be used for the qualitative detection of translocations involving the human FGFR3 gene at 4p16.3 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC FGFR3 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 4p16.3* (chr4:1,093,149-1,729,455) distal to the FGFR3 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 4p16.3* (chr4:1,922,997-2,446,931) proximal to the FGFR3 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC FGFR3 Probe map (not to scale); Right: Ideogram of chromosome 4 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the FGFR3 gene region, two green/orange fusion signals appear.

Aberrant situation: One FGFR3 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Example of an aberrant signal pattern: Breast cancer tissue section with translocation affecting the FGFR3 gene as indicated by one non-rearranged orange/green fusion signal, one orange, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	4	CE/IVD	5 (50 μl)	Z-2170-50
			20 (200 μl)	Z-2170-200

ZytoLight® SPEC FGFR3/4p11 Dual Color Probe

Intended Purpose

The ZytoLight SPEC FGFR3/4p11 Dual Color Probe (PL41) is intended to be used for the qualitative detection of amplifications involving the human FGFR3 gene as well as the detection of chromosome 4p11 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC FGFR3/4p11 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 4p16.3* (chr4:1,531,083-2,073,649) harboring the FGFR3 region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 4p11* (chr4:48,329,914-48,762,386) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC FGFR3 Probe map; Bottom: SPEC 4p11 Probe map (not to scale); Right: Ideogram of chromosome 4 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without an amplification involving the FGFR3 gene region, two green signals and two orange signals appear.

Aberrant situation: In cells with an amplification of the FGFR3 gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Example of an aberrant signal pattern: Bladder cancer tissue section with interphase cells showing polysomy of chromosome 4 as indicated by multiple green and orange signals in the nuclei.

Labels	Chromosome	Reg. Status	Tests	Order No.
	4	CE/IVD	20 (200 μl)	Z-2082-200

ZytoLight® SPEC FGFR3/IGH Dual Color Dual Fusion Probe

Intended Purpose

The ZytoLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe is designed to detect the translocation t(4;14)(p16.3;q32.3) affecting the FGFR3 (fibroblast growth factor receptor 3, a.k.a. JTK4) gene in the chromosomal region 4p16.3 and the IGH (immunoglobulin heavy locus, a.k.a. IGH@) locus in 14q32.33.

FGFR3 encodes for a receptor tyrosine kinase, which regulates downstream signaling cascades after ligand binding. Fusion to several partner genes (including the IGH locus) can lead to a ligand-independent activation of the tyrosine kinase of the resulting FGFR3 fusion protein, frequently found in multiple myeloma (MM).

FGFR3/IGH translocations are observed in approximately 15-20% of patients with MM. The breaking points for the 4p16.3 locus are found between the FGFR3 gene and the 5' end of the NSD2 gene. The t(4;14)(p16.3;q32.3) translocation is associated with upregulation of the FGFR3 and the myeloma NSD2 (a.k.a. MMSET) domain protein. Patients with FGFR3/IGH translocation demonstrate an overall poor prognosis that is only partially mitigated by the use of the novel agents bortezomib and lenalidomide.

With conventional cytogenetics, the t(4;14)(p16.3;q32.3) translocation is difficult to identify. Thus, the detection of FGFR3/IGH translocations by fluorescence *in situ* hybridization may be of diagnostic and prognostic relevance.

References

- Bergsagel PL & Kuehl WM (2001) Oncogene 20: 5611-22.
Chesi M, *et al.* (1998) Blood 92: 3025-34.
Fabris S, *et al.* (2005) Genes Chromosomes Cancer 42: 117-27.
Fenton JA, *et al.* (2003) Oncogene 22: 1103-13.
Kalf J & Spencer A (2012) Blood Cancer J 7: e89.
Sonneveld P, *et al.* (2016) Blood 127: 2955-62.
Walker BA, *et al.* (2013) Blood 121: 3413-19.

Probe Description

The ZytoLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 4p16.3* (chr4:1,496,938-2,351,657) harboring the FGFR3 region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH locus (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC FGFR3 Probe map; Bottom: SPEC IGH Probe map (not to scale); Right: Ideograms of chromosome 4 and 14 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate orange and green signals appear.

Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	4, 14	CE/IVD	5 (50 μl)	Z-2282-50


ZytoLight® SPEC 4p11/CEN 10/17 Triple Color Probe

Intended Purpose
The ZytoLight SPEC 4p11/CEN 10/17 Triple Color Probe (PL261) is intended to be used for the qualitative detection of human chromosome 4p11 specific sequences as well as alpha satellites of chromosomes 10 and 17 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC 4p11/CEN 10/17 Triple Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 4p11* (chr4:48,329,914-48,762,386) (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 10p11.1-q11.1 specific for the alpha satellite centromeric region D10Z1 of chromosome 10.
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without an alteration of the respective chromosomes, two green, two orange, and two blue signals appear.
Aberrant situation: In a cell with aneuploidy of the chromosomes 4, 10, and/or 17, an increased or a reduced number of signals of the respective color will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
  	4, 10, 17	CE/IVD	5 (50 μl)	Z-2307-50


ZytoLight® SPEC PDGFRA/FIP1L1 TriCheck™ Probe

Intended Purpose
The ZytoLight SPEC PDGFRA/FIP1L1 TriCheck Probe (PL167) is intended to be used for the qualitative detection of rearrangements involving the human PDGFRA gene with and without participation of the human FIP1L1 gene in cytological specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC PDGFRA/FIP1L1 TriCheck Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 4q12* (chr4:53,552,536-54,238,252) proximal to the FIP1L1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 4q12* (chr4:54,351,156-54,749,671) proximal to the PDGFRA gene region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/μl), which target sequences mapping in 4q12* (chr4:55,185,968-55,915,442) distal to the PDGFRA gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without PDGFRA-FIP1L1 rearrangements, two green/orange fusion signals appear when using an appropriate dual bandpass filter set, and two blue signals appear when using an appropriate single bandpass filter set. When using an appropriate triple color bandpass filter set, two green/orange/blue fusion signals can be observed.
Aberrant situation: A PDGFRA-FIP1L1 fusion resulting from a deletion of interstitial DNA is indicated by the loss of one orange signal leading to a separate green signal co-localizing with a blue signal. A PDGFRA translocation without involvement of FIP1L1 is indicated by a green/orange fusion signal and a separate blue signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
  	4	CE/IVD	5 (50 μl)	Z-2209-50

ZytoLight® SPEC TERT Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC TERT Dual Color Break Apart Probe (PL229) is intended to be used for the qualitative detection of translocations involving the human TERT gene in 5p15.33 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC TERT Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 5p15.33* (chr5:620,184-1,161,456) distal to the TERT breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 5p15.33* (chr5:1,353,007-1,588,209) proximal to the TERT breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC TERT Probe map (not to scale); Right: Ideogram of chromosome 5 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the TERT gene region, two green/orange fusion signals appear.
Aberrant situation: One TERT gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	5	CE/IVD	5 (50 μl)	Z-2273-50

ZytoLight® SPEC TERT/5q31 Dual Color Probe

Intended Purpose
The ZytoLight SPEC TERT/5q31 Dual Color Probe (PL50) is intended to be used for the qualitative detection of amplifications involving the human TERT gene as well as the detection of chromosome 5q31 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight TERT/5q31 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 5p15.33* (chr5:974,089-1,588,209) harboring the TERT region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 5q31.2* (chr5:137,394,637-137,999,163) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC TERT Probe map (not to scale); Bottom: SPEC 5q31 Probe map (not to scale); Right: Ideogram of chromosome 5 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the TERT gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the TERT gene region, an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	5	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2091-50 Z-2091-200

ZytoLight® SPEC EGR1/D5S23,D5S721 Dual Color Probe

Intended Purpose
The ZytoLight SPEC EGR1/D5S23,D5S721 Dual Color Probe (PL169) is intended to be used for the qualitative detection of deletions involving the human EGR1 gene as well as the detection of the human D5S23,D5S721 control region at 5p15.2-p15.31 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC EGR1/D5S23,D5S721 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 5q31.2* (chr5:137,667,079-137,897,109) harboring the EGR1 region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 5p15.2-p15.31* (chr5:9,233,775-9,967,465) harboring the D5S23,D5S721 locus.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC EGR1 Probe map; Bottom: SPEC D5S23,D5S721 Probe map (not to scale); Right: Ideogram of chromosome 5 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the EGR1 gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletion affecting the EGR1 gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the EGR1 gene region might result in a normal signal pattern with orange signals of reduced size.



Hybridized to normal interphase calls as indicated by two orange and two green signals and to metaphase chromosomes of a normal cell.

Example of an aberrant signal pattern: SPEC EGR1/D5S23,D5S721 Dual Color Probe hybridized to an AML specimen with deletion of the EGR1 gene as indicated by one orange and two green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	5	CE/IVD	5 (50 μl)	Z-2211-50

ZytoLight® SPEC CSF1R Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC CSF1R Dual Color Break Apart Probe is designed to detect rearrangements involving the chromosomal region 5q32 harboring the CSF1R (colony stimulating factor 1 receptor, a.k.a. FMS) gene. The CSF1 receptor is activated by dimerization upon binding of its ligand CSF1 and is involved in macrophage development. Rearrangement of the CSF1R gene was first detected in an acute megakaryoblastic leukemia (AMKL) cell line generating the RBM6-CSF1R fusion gene. A MEF2D-CSF1R fusion gene was described in a patient with primary pre-B cell acute lymphoblastic leukemia (pre-B ALL). Both fusion proteins contain the intact kinase domain of CSF1R.

Philadelphia chromosome-like ALL (Ph-like ALL) is a subgroup of B-cell precursor ALL and is associated with a high risk of treatment failure. SSBP2-CSF1R fusions were detected in some patients with Ph-like ALL. They result from either the balanced translocation t(5;5)(q14;q32) or the duplication dup(5)(q14q32). Expression of this fusion gene results in cytokine-independent growth and enhanced STAT5 activation which are inhibited by dasatinib *in vitro*. CSF1R signaling was also shown to be suppressed by the ABL1 kinase inhibitor imatinib.

Hence, the detection of CSF1R rearrangements by FISH may help in selecting ALL patients eligible for treatment with CSF1R inhibitors.

- References**
Dewar AL, *et al.* (2005) Blood 105: 3127-32.
Gu TL, *et al.* (2007) Blood 110: 323-33.
Lilljebjörn H, *et al.* (2014) Leukemia 28: 977-9.
Roberts KG, *et al.* (2014) N Engl J Med 371: 1005-15.
Schwab C, *et al.* (2014) Blood 124: 3773.

- Probe Description**
The ZytoLight SPEC CSF1R Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 5q32-q33.1* (chr5:149,548,518-150,118,449) distal to the CSF1R breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 5q32* (chr5:149,058,515-149,421,081) proximal to the CSF1R breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC CSF1R Probe map (not to scale); Right: Ideogram of chromosome 5 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the CSF1R gene region, two green/orange fusion signals appear.
Aberrant situation: One CSF1R gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Blood smear with translocation of the CSF1R gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	5	CE/IVD	5 (50 μl)	Z-2202-50

ZytoLight® SPEC CSF1R/D5S23,D5S721 Dual Color Probe

Intended Purpose
The ZytoLight SPEC CSF1R/D5S23,D5S721 Dual Color Probe is designed for the detection of 5q deletions. The CSF1R (colony stimulating factor 1 receptor, a.k.a. C-FMS) gene is located in the chromosomal region 5q32.

The interstitial deletion of chromosome 5q is a characteristic hallmark of the myelodysplastic syndrome (MDS) with isolated del(5q). The size of the deletion as well as the breakpoints are variable but a commonly deleted region (CDR) has been narrowed to the approximately 1.5 Mb interval at 5q32-q33.1 flanked by the DNA marker D5S413 and the GLRA1 gene.

One candidate gene for the development of MDS in patients with 5q-syndrome is RPS14 (ribosomal protein 14), a tumor suppressor gene located in the chromosomal region 5q33.1. Haploinsufficiency (caused by hemizygous deletion) of RPS14 is the probable cause of the erythroid defect that characterizes the 5q-syndrome. Lenalidomide has been reported to overcome the pathogenic effect of 5q deletion in MDS.

Despite the severe phenotype of the 5q-syndrome, it has a relatively low (10%) transformation risk to acute myeloid leukemia (AML). Therefore, FISH may be a helpful tool for diagnosis and therapy decision.

- References**
Boultonwood J, *et al.* (1991) Proc Natl Acad Sci U S A 88: 6176-80.
Boultonwood J, *et al.* (2010) Blood 116: 5803-11.
Giagounidis AA, *et al.* (2004) Clin Cancer Res 12: 5-10.
Van den Berghe H & Michaux JL (1974) Nature 251: 437-8.
Swerdlow SH, *et al.* (ed.) (2008) WHO classification of tumours of haematopoietic and lymphoid tissues.

- Probe Description**
The ZytoLight SPEC CSF1R/D5S23,D5S721 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 5q32-q33.1* (chr5:149,548,518-150,118,449) distal to the CSF1R breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 5p15.2-p15.31* (chr5:9,233,775-9,967,465) harboring the D5S23,D5S721 locus (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC CSF1R Probe map (not to scale); Bottom: SPEC D5S23,D5S721 Probe map (not to scale); Right: Ideogram of chromosome 5 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the CSF1R gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletion affecting the CSF1R gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the CSF1R gene region might result in a normal signal pattern with green signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.

Bone marrow biopsy tissue section of an ALL case showing hemizygous deletion of the CSF1R gene as indicated by the loss of one green signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	5	CE/IVD	5 (50 μl)	Z-2268-50


ZytoLight® SPEC PDGFRB Dual Color Break Apart Probe

Intended Purpose
The *ZytoLight* SPEC PDGFRB Dual Color Break Apart Probe (PL155) is intended to be used for the qualitative detection of translocations involving the human PDGFRB gene at 5q32 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The *ZytoLight* SPEC PDGFRB Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 5q32-q33.1* (chr5:149,548,518-150,285,722) distal to the PDGFRB breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 5q32* (chr5:149,058,515-149,471,369) proximal to the PDGFRB breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the PDGFRB gene region, two green/orange fusion signals appear.
Aberrant situation: One PDGFRB gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.


ZytoLight® SPEC IRF4,DUSP22 Dual Color Break Apart Probe

Intended Purpose
The *ZytoLight* SPEC IRF4,DUSP22 Dual Color Break Apart Probe (PL168) is intended to be used for the qualitative detection of translocations involving the human IRF4,DUSP22 gene region at 6p25.3 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.


- Probe Description**
The *ZytoLight* SPEC IRF4,DUSP22 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 6p25.3* (chr6:557,233-940,968) proximal to the IRF4,DUSP22 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 6p25.3* (chr6:114,722-273,908) distal to the IRF4,DUSP22 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the IRF4,DUSP22 gene region, two green/orange fusion signals appear.
Aberrant situation: One IRF4,DUSP22 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: T-cell lymphoma tissue section with translocation affecting the 6p25.3 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	CE/IVD	5 (50 μl)	Z-2210-50
			20 (200 μl)	Z-2210-200

ZytoLight® SPEC RREB1/MYB/CEN 6 Triple Color Probe

Intended Purpose
The *ZytoLight* SPEC RREB1/MYB/CEN 6 Triple Color Probe (PL108) is intended to be used for the qualitative detection of amplifications involving the human RREB1 gene as well as the human MYB gene and the detection of chromosome 6 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The *ZytoLight* SPEC RREB1/MYB/CEN 6 Triple Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 6p24.3-p25.1* (chr6:6,913,938-7,406,653) harboring the RREB1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 6q23.2-q23.3* (chr6:135,141,227-135,715,246) harboring the MYB region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 6p11.1-q11 specific for the alpha satellite centromeric region D6Z1 of chromosome 6.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC RREB1 Probe map; Bottom: SPEC MYB Probe map (not to scale); Right: Ideogram of chromosome 6 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the respective gene regions, two green, two orange, and two blue signals are expected.
Aberrant situation: In cells with an amplification of the RREB1 gene region or the MYB gene region, an increased number of green or orange signals will be observed, respectively. In a cell with deletions of the RREB1 gene region or the MYB gene region, a reduced number of green or orange signals will be observed, respectively. Deletions affecting only parts of the gene regions might result in a normal signal pattern with signals of reduced size.



Hybridized to normal interphase cells as indicated by two green, two orange, and two blue signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	CE/IVD	5 (50 μl)	Z-2152-50
			20 (200 μl)	Z-2152-200


ZytoLight® SPEC PHF1 Dual Color Break Apart Probe

Intended Purpose
The *ZytoLight* SPEC PHF1 Dual Color Break Apart Probe (PL173) is intended to be used for the qualitative detection of translocations involving the human PHF1 gene at 6p21.32 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.


- Probe Description**
The *ZytoLight* SPEC PHF1 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 6p21.31-p21.32* (chr6:33,406,580-33,863,564) proximal to the PHF1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 6p21.32* (chr6:33,121,529-33,317,357) distal to the PHF1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the PHF1 gene region, two green/orange fusion signals appear.
Aberrant situation: One PHF1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.



Example of an aberrant signal pattern: Sarcoma tissue section with translocation of the PHF1 gene as indicated by one non-rearranged orange/green fusion signal, one orange, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	CE/IVD	5 (50 μl)	Z-2215-50

ZytoLight® SPEC VEGFA/CEN 6 Dual Color Probe

Intended Purpose

The ZytoLight SPEC VEGFA/CEN 6 Dual Color Probe (PL153) is intended to be used for the qualitative detection of amplifications involving the human VEGFA gene as well as the detection of chromosome 6 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

Probe Description

The ZytoLight SPEC VEGFA/CEN 6 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 6p21.1* (chr6:43,633,271-43,985,142) harboring the VEGFA gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 6p11.1-q11 specific for the alpha satellite centromeric region D6Z1 of chromosome 6.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC VEGFA Probe map (not to scale); Right: Ideogram of chromosome 6 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without an amplification involving the VEGFA gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the VEGFA gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: HCC tissue section with interphase cells showing a polysomy of chromosome 6 as indicated by multiple green (VEGFA) and orange (CEN 6) signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	RUO	20 (200 μl)	Z-2195-200

ZytoLight® SPEC ROS1 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC ROS1 Dual Color Break Apart Probe (PL101) is intended to be used for the qualitative detection of translocations involving the human ROS1 gene at 6q22.1 in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ROS1 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 6q22.1* (chr6:116,912,298-117,627,255) proximal to the ROS1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 6q22.1* (chr6:117,659,135-117,871,701) distal to the ROS1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ROS1 Probe map (not to scale); Right: Ideogram of chromosome 6 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the ROS1 gene region, two green/orange fusion signals appear.

Aberrant situation: One ROS1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated green signals are the result of deletions distal to the ROS1 breakpoint region or are due to unbalanced translocations affecting this chromosomal region.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Section of paraffin-embedded NSCLC cell line with translocation affecting the 6q22.1 locus harboring ROS1 as indicated by one orange/green fusion signal (non-rearranged), one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	CE/IVD	5 (50 μl)	Z-2144-50
			20 (200 μl)	Z-2144-200

ZytoLight® SPEC ROS1/CEN 6 Dual Color Probe

Intended Purpose

The ZytoLight SPEC ROS1/CEN 6 Dual Color Probe (PL118) is intended to be used for the qualitative detection of amplifications involving the human ROS1 gene as well as the detection of chromosome 6 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ROS1/CEN 6 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 6q22.1* (chr6:117,431,278-117,894,830) harboring the ROS1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 6p11.1-q11 specific for the alpha satellite centromeric region D6Z1 of chromosome 6.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ROS1/CEN 6 Dual Color Probe map (not to scale); Right: Ideogram of chromosome 6 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without an amplification involving the ROS1 gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the ROS1 gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: Lung cancer tissue section with interphase cells showing a polysomy of chromosome 6 as indicated by multiple orange (CEN 6) and green (ROS1) signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	CE/IVD	20 (200 μl)	Z-2162-50

ZytoLight® SPEC MYB Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC MYB Dual Color Break Apart Probe (PL100) is intended to be used for the qualitative detection of translocations involving the human MYB gene at 6q23.3 in cytologic or formalin-fixed, paraffin-embedded specimens, such as adenoid cystic carcinoma (ACC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ACC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC MYB Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 6q23.2-q23.3* (chr6:134,840,690-135,483,752) proximal to the MYB breakpoint region (see Fig.1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 6q23.3* (chr6:135,728,667-136,390,142) distal to the MYB breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC MYB Probe map (not to scale); Right: Ideogram of chromosome 6 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the MYB gene region, two green/orange fusion signals appear.

Aberrant situation: One MYB gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Adenoid cystic carcinoma tissue section with translocation affecting the 6q23.3 locus as indicated by one orange/green fusion (non-rearranged) signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	CE/IVD	5 (50 μl)	Z-2143-50
			20 (200 μl)	Z-2143-200


ZytoLight® SPEC MYB/CEN 6 Dual Color Probe

Intended Purpose
The ZytoLight SPEC MYB/CEN 6 Dual Color Probe (PL236) is intended to be used for the qualitative detection of deletions involving the human MYB gene and the detection of chromosome 6 alpha satellites in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC MYB/CEN 6 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 6q23.2-q23.3* (chr6:135,141,227-135,715,246) harboring the MYB gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 6p11.1-q11 specific for the alpha satellite centromeric region D6Z1 of chromosome 6.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a deletion involving the MYB gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletions affecting the MYB gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the MYB gene region might result in a normal signal pattern with orange signals of reduced size.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	6	CE/IVD	5 (50 μl)	Z-2281-50


ZytoLight® SPEC ESR1/CEN 6 Dual Color Probe

Intended Purpose
The ZytoLight SPEC ESR1/CEN 6 Dual Color Probe (PL27) is intended to be used for the qualitative detection of amplifications involving the human ESR1 gene as well as the detection of chromosome 6 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC ESR1/CEN 6 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 6q25.1* (chr6:152,083,365-152,478,947) harboring the ESR1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 6p11.1-q11 specific for the alpha satellite centromeric region D6Z1 of chromosome 6.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without an amplification involving the ESR1 gene region, two green and two orange signals appear.
Aberrant situation: In cells with an amplification of the ESR1 gene region, an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	6	CE/IVD	5 (50 μl)	Z-2069-50


ZytoLight® SPEC JAZF1 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC JAZF1 Dual Color Break Apart Probe (PL89) is intended to be used for the qualitative detection of translocations involving the human JAZF1 gene at 7p15.1-p15.2 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC JAZF1 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 7p15.2* (chr7:27,146,601-27,846,497) distal to the JAZF1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 7p15.1* (chr7:28,059,911-28,661,819) proximal to the JAZF1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the JAZF1 gene region, two green/orange fusion signals appear.
Aberrant situation: One JAZF1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	7	CE/IVD	5 (50 μl)	Z-2132-50

ZytoLight® SPEC IKZF1/CEN 7 Dual Color Probe

Intended Purpose
The ZytoLight SPEC IKZF1/CEN 7 Dual Color Probe is designed for the detection of deletions affecting the IKZF1 (IKAROS family zinc finger 1, a.k.a. ZNF1A1, IKAROS) gene.


The IKZF1 gene is located on 7p12.2 and encodes a zinc finger transcription factor, which is required for normal hematopoietic differentiation and proliferation, particularly in lymphoid lineages.

Genomic deletions affecting the IKZF1 gene are found in approximately 15% of pediatric and ~40% of adult B-cell precursor acute lymphoblastic leukemia (B-ALL) cases. The frequency is remarkably high in BCR-ABL1-positive (~70%) and BCR-ABL1-like (~40%) pediatric B-ALL. IKZF1 deletions were also identified in the progression of chronic myeloid leukemia to lymphoid blast crisis.

The most frequent deletions in B-ALL affect the whole gene or exons 4 to 7. Deletions affecting other exons (i.e., exons 2 to 7, exons 2 to 8, and exons 4 to 8) were also observed. IKZF1 deletions are associated with poor prognosis and high risk of relapse in cases of B-ALL. Hence, the detection of IKZF1 deletions by FISH may help in predicting the clinical outcome in patients with B-ALL.

- References**
Boer JM, *et al.* (2016) Leukemia 30: 32-8.
Hashiguchi J, *et al.* (2018) J Mol Diagn 20: 446-54.
Iacobucci I, *et al.* (2009) Blood 114: 2159-67.
Meyer C, *et al.* (2013) Am J Blood Res 3: 165-73.
Mullighan CG, *et al.* (2008) Nature 453: 110-4.

- Probe Description**
The ZytoLight SPEC IKZF1/CEN 7 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 7p12.2* (chr7:50,412,912-50,463,612) harboring the IKZF1 gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a deletion involving the IKZF1 gene region, two orange signals and two green signals appear.
Aberrant situation: In a cell with deletion affecting the IKZF1 gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the IKZF1 gene region might result in a normal signal pattern with orange signals of reduced size. In a cell with monosomy 7, a reduced number of both, the orange (IKZF1 gene region) and green (CEN 7) signals will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
 	7	CE/IVD	5 (50 μl)	Z-2304-50

ZytoLight® SPEC EGFR/CEN 7 Dual Color Probe

Intended Purpose

The *ZytoLight* SPEC EGFR/CEN 7 Dual Color Probe (PL15) is intended to be used for the qualitative detection of amplifications involving the human EGFR gene as well as the detection of chromosome 7 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC) and glioma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and glioma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC EGFR/CEN 7 Dual Color Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 7p11.2* (chr7:54,912,555-55,548,375) harboring the EGFR gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without an amplification involving the EGFR gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the EGFR gene region, multiple copies of the green signal or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Cancer cells with multiple copies of chromosome 7 and extra EGFR signals (green) in a sample from an NSCLC patient.

Labels	Chromosome	Reg. Status	Tests	Order No.
	7	CE/IVD	5 (50 μl)	Z-2033-50
			20 (200 μl)	Z-2033-200

ZytoLight® SPEC Williams-Beuren Dual Color Probe

Intended Purpose

The *ZytoLight* SPEC Williams-Beuren Dual Color Probe is designed to detect deletions affecting the chromosomal region 7q11.23 harboring the ELN (elastin, a.k.a. WBS) gene.

The Williams-Beuren syndrome (WBS) is a genetic disorder caused by a hemizygous contiguous gene deletion on chromosome 7q11.23. The estimated prevalence of the disease ranges between 1/7,500 and 1/20,000 newborns. The WBS deletion region (~1.5-1.8 Mb) consists of a single copy gene region containing app. 28 genes, including the ELN gene that is flanked by repetitive sequences known as low-copy repeats (LCRs).

The deletions arise as a consequence of misalignment of these repetitive sequences during meiosis and a following unequal crossing over due to high similarity of LCRs. Usually, WBS occurs sporadically, but some parents of WBS patients were shown to carry a paracentric inversion of the WBS locus. Presence of this inversion predisposes to chromosomal mispairing in meiosis. WBS patients clinically display a characteristic pattern of symptoms including vascular stenosis, weakness of connective tissue, a typical face, short stature, overfriendliness, and mental retardation.

FISH analysis can be performed to confirm WBS diagnosis in patients with vascular stenosis together with mental retardation.


References

- Bayés M, *et al.* (2003) Am J Hum Genet 73: 131-51.
Beuren AJ, *et al.* (1964) Am J Cardiol 13: 471-83.
Schubert C (2009) Cell Mol Life Sci 66: 1178-97.
Sugayama SM, *et al.* (2003) Arq Bras Cardiol 81: 462-73.
Williams JC, *et al.* (1961) Circulation 24: 1311-8.

Probe Description

The *ZytoLight* SPEC Williams-Beuren Dual Color Probe is composed of:


- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 7q11.23* (chr7:73,408,390-73,899,599) harboring the ELN region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Results

Normal situation: In Interphases of normal cells or cells without a deletion involving the ELN gene region, two orange signals and two green signals appear.

Aberrant situation: In a cell with deletion affecting the ELN gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the ELN gene region might result in a normal signal pattern with orange signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.



Lymphocytes and metaphase chromosomes from a Williams-Beuren syndrome case showing an ELN deletion as indicated by the loss of one orange signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	7	CE/IVD	5 (50 μl)	Z-2302-50

ZytoLight® SPEC CUX1/EZH2/CEN 7 Triple Color Probe

Intended Purpose

The *ZytoLight* SPEC CUX1/EZH2/CEN 7 Triple Color Probe (PL172) is intended to be used for the qualitative detection of deletions involving the human CUX1 gene and the human EZH2 gene as well as the detection of chromosome 7 alpha satellites in cytological specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC CUX1/EZH2/CEN 7 Triple Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10.0 ng/μl), which target sequences mapping in 7q22.1* (chr7:101,270,255-101,934,924) harboring the CUX1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 7q36.1* (chr7:148,402,839-148,647,927) harboring the EZH2 gene region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12.0 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC CUX1 Probe map; Bottom: SPEC EZH2 Probe map (not to scale); Right: Ideogram of chromosome 7 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a deletion involving the respective gene regions, two green, two orange, and two blue signals appear.

Aberrant situation: In a cell with deletion affecting the CUX1 gene region and/or the EZH2 gene region, a reduced number of green and/or orange signals will be observed. Deletions affecting only parts of the respective gene regions might result in a normal signal pattern with signals of reduced size. Monosomy 7 will result in a loss of one green, one orange, and one blue signal.



Example of an aberrant signal pattern: Bone marrow smear with deletion of the CUX1 gene as indicated by one green signal in each nucleus.



Example of an aberrant signal pattern: Bone marrow smear with deletion of the EZH2 gene as indicated by one orange signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	7	CE/IVD	5 (50 μl)	Z-2214-50

ZytoLight® SPEC MET/CEN 7 Dual Color Probe

Intended Purpose

The *ZytoLight* SPEC MET/CEN 7 Dual Color Probe (PL46) is intended to be used for the qualitative detection of amplifications involving the human MET gene as well as the detection of chromosome 7 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC MET/CEN 7 Dual Color Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 7q31.2* (chr7:115,925,700-116,718,699) harboring the MET gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Results

Normal situation: In interphases of normal cells or cells without an amplification involving the MET gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the MET gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



NSCLC specimen cells with polysomy of chromosome 7 as indicated by four orange (CEN 7) and four green (MET) signals in the nuclei.

Labels	Chromosome	Reg. Status	Tests	Order No.
	7	CE/IVD	5 (50 μl)	Z-2087-50
			20 (200 μl)	Z-2087-200

ZytoLight® SPEC BRAF Dual Color Break Apart Probe

Intended Purpose


The *ZytoLight* SPEC BRAF Dual Color Break Apart Probe (PL147) is intended to be used for the qualitative detection of translocations involving the human BRAF gene at 7q34 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description


- The *ZytoLight* SPEC BRAF Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 7q34* (chr7:140,535,100-141,233,856) distal to the BRAF breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 7q34* (chr7:139,972,721-140,469,362) proximal to the BRAF breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without a translocation involving the BRAF gene region, two green/orange fusion signals appear.

Aberrant situation: One BRAF gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. A deletion of 5' BRAF sequences resulting in a FAM131B-BRAF fusion is indicated by the loss of one green signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	7	CE/IVD	20 (200 μl)	Z-2189-200

ZytoLight® SPEC BRAF/CEN 7 Dual Color Probe

Intended Purpose


The *ZytoLight* SPEC BRAF/CEN 7 Dual Color Probe (PL149) is intended to be used for the qualitative detection of amplifications involving the human BRAF gene as well as the detection of chromosome 7 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description


- The *ZytoLight* SPEC BRAF/CEN 7 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 7q34* (chr7:140,266,210-140,792,511) harboring the BRAF gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Results

Normal situation: In interphases of normal cells or cells without an amplification involving the BRAF gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the BRAF gene region, an increased number of green signals or green signal clusters will be observed.



Normal interphase cells, BRAF (green), CEN 7 (orange).



Example of an aberrant signal pattern: NSCLC tissue section with amplification of the BRAF gene (green).

Labels	Chromosome	Reg. Status	Tests	Order No.
	7	CE/IVD	20 (200 μl)	Z-2191-200

ZytoLight® SPEC NRG1/CD74 TriCheck™ Probe

Intended Purpose

The *ZytoLight* SPEC NRG1/CD74 TriCheck Probe (PL152) is intended to be used for the qualitative detection of human NRG1 rearrangements with and without participation of the human CD74 gene in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The *ZytoLight* SPEC NRG1/CD74 TriCheck Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 8p12* (chr8:31,730,448-32,433,429) distal to the NRG1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 8p12* (chr8:32,644,505-32,985,279) proximal to the NRG1 breakpoint region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/μl), which target sequences mapping in 5q32-q33.1* (chr5:149,274,320-150,285,722) harboring the CD74 gene (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC NRG1 Probe map; Bottom: SPEC CD74 Probe map (not to scale); Right: Ideograms of chromosome 8 and 5 indicating the hybridization locations

Results

Normal situation: In Interphases of normal cells or cells without NRG1 rearrangement, two green, two orange, and two blue signals appear.

Aberrant situation: A CD74-NRG1 fusion is indicated by one separate green signal, one separate orange signal, and an additional blue signal which colocalizes with the separated orange signal. An NRG1 rearrangement not involving CD74 is indicated by separated orange and green signals without an additional blue signal. Isolated orange signals are the result of deletions distal to the NRG1 breakpoint region.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals and two blue signals per nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	5, 8	CE/IVD	20 (200 μl)	Z-2194-200

ZytoLight® SPEC NRG1 Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC NRG1 Dual Color Break Apart Probe (PL140) is intended to be used for the qualitative detection of translocations involving the human NRG1 gene at 8p12 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The *ZytoLight* SPEC NRG1 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 8p12* (chr8:31,730,448-32,433,429) distal to the NRG1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 8p12* (chr8:32,644,505-32,985,279) proximal to the NRG1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC NRG1 Probe map (not to scale); Right: Ideogram of chromosome 8 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the NRG1 gene region, two green/orange fusion signals appear.

Aberrant situation: One NRG1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated orange signals are the result of deletions distal to the NRG1 breakpoint region.



Hybridized on normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: Lung cancer tissue section with rearrangement of the NRG1 gene as indicated by extra orange signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	8	CE/IVD	20 (200 μl)	Z-2181-200

Image kindly provided by Mc Leer A, Duruisseau M, Wislez M, and colleagues, Grenoble and Paris, France.


ZytoLight® SPEC FGFR1 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC FGFR1 Dual Color Break Apart Probe (PL124) is intended to be used for the qualitative detection of translocations involving the human FGFR1 gene at 8p11.22-p11.23 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC FGFR1 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 8p11.22* (chr8:38,352,117-38,951,783) proximal to the FGFR1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 8p11.23* (chr8:37,635,912-38,239,669) distal to the FGFR1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the FGFR1 gene region, two green/orange fusion signals appear.
Aberrant situation: One FGFR1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Example of an aberrant signal pattern: 8p11 myeloproliferative syndrome (EMS) tissue section with translocation of the FGFR1 gene as indicated by one non-rearranged orange/green fusion signal, one orange, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	8	CE/IVD	5 (50 μl)	Z-2168-50
			20 (200 μl)	Z-2168-200


ZytoLight® SPEC FGFR1/CEN 8 Dual Color Probe

Intended Purpose
The ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe (PL29) is intended to be used for the qualitative detection of amplifications involving the human FGFR1 gene as well as the detection of chromosome 8 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer and squamous cell lung cancer, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and squamous cell lung cancer and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 8p11.22-p11.23* (chr8:38,063,906-38,527,745) harboring the FGFR1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 8p11.1-q11.1 specific for the alpha satellite centromeric region D8Z2 of chromosome 8.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without an amplification involving the FGFR1 gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the FGFR1 gene region an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Lung carcinoma tissue section with interphase cells showing amplification of the FGFR1 gene (green) and partly polysomy 8 (orange).

Labels	Chromosome	Reg. Status	Tests	Order No.
	8	CE/IVD	5 (50 μl)	Z-2072-50
			20 (200 μl)	Z-2072-200

ZytoLight® SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe is designed to detect the specific translocation involving the chromosomal region 21q22.12 harboring the RUNX1 (a.k.a. AML1) gene and the chromosomal region 8q21.3 harboring the RUNX1T1 (a.k.a. ETO, CBF2T1) gene.

The balanced chromosomal translocation t(8;21) is found in about 90% of acute myeloid leukemia (AML) patients. AML is a heterogeneous clonal disorder of hematopoietic progenitor cells and one of the most common malignant myeloid disorders in adults.

The runt related transcription factor 1 gene (RUNX1) and RUNX1 translocation partner 1 (RUNX1T1) gene are both involved in the transcriptional regulation of genes during normal hematopoiesis. The non-random translocation t(8;21) (q21.3;q22.1) is strongly associated with the French-American-British (FAB) phenotype M2 (AML-M2) and produces a chimeric gene consisting of the 5'-region of the RUNX1 gene fused to the 3'-region of the RUNX1T1 gene. The chimeric protein is thought to be associated with the nuclear corepressor/histone deacetylase complex to block hematopoietic differentiation.


Fluorescence *in situ* hybridization (FISH) can provide important information for the management of patients with hematologic disorders.

- References**
Dayyani F, *et al.* (2008) Blood 111: 4338-47.
Estey E & Döhner H (2006) Lancet 368: 1894-907.
Gmidene A, *et al.* (2011) Med Oncol 28 Suppl 1: 509-12.
Licht D (2001) Oncogene 20: 5560-79.
Vangala RK, *et al.* (2003) Blood 101: 270-7.

- Probe Description**
The ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 21q22.11-q22.12* (chr21:35,530,283-36,855,548) harboring the RUNX1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 8q21.3-q22.1* (chr8:92,632,490-93,746,043) harboring the RUNX1T1 region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Fig. 1: Top: SPEC RUNX1 Probe map; Bottom: SPEC RUNX1T1 Probe map (not to scale)



Ideograms of chromosome 21 and 8 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Bone marrow biopsy section with translocation affecting the RUNX1/RUNX1T1 locus as indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	8,21	CE/IVD	5 (50 μl)	Z-2112-50
			20 (200 μl)	Z-2112-200


ZytoLight® SPEC MYC Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC MYC Dual Color Break Apart Probe (PL49) is intended to be used for the qualitative detection of translocations involving the human MYC gene at 8q24.21 in cytologic or formalin-fixed, paraffin-embedded specimens, such as diffuse large B-cell lymphoma (DLBCL) or Burkitt lymphoma (BL), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of DLBCL or BL and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC MYC Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 8q24.21* (chr8:130,373,051-130,930,673) distal to the MYC breakpoint region (see Fig.1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 8q24.21* (chr8:127,888,765-128,363,281) proximal to the MYC breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In Interphases of normal cells or cells without a translocation involving the MYC gene region, two green/orange fusion signals appear.
Aberrant situation: One MYC gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
	8	CE/IVD	5 (50 μl)	Z-2090-50
			20 (200 μl)	Z-2090-200


ZytoLight® SPEC MYC/CEN 8 Dual Color Probe

Intended Purpose
The ZytoLight SPEC MYC/CEN 8 Dual Color Probe (PL51) is intended to be used for the qualitative detection of amplifications involving the human MYC gene as well as the detection of chromosome 8 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC MYC/CEN 8 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 8q24.21* (chr8:128,487,995-128,887,929) harboring the MYC gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 8p11.1-q11.1 specific for the alpha satellite centromeric region D8Z2 of chromosome 8.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without an amplification involving the MYC gene region, two green signals and two orange signals appear.
Aberrant situation: In In cells with an amplification of the MYC gene region, an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
	8	CE/IVD	5 (50 μl)	Z-2092-50
			20 (200 μl)	Z-2092-200

ZytoLight® SPEC MYC/IGH Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe (PL62) is intended to be used for the qualitative detection of the translocation t(8;14) (q24.21;q32.3) involving the human IGH and MYC genes in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH locus (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 8q24.21* (chr8:128,171,178-129,517,468) harboring the MYC gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC MYC Probe map; Bottom: SPEC IGH Probe map (not to scale); Right: Ideograms of chromosome 8 and 14 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective bands, two separate green/orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Labels	Chromosome	Reg. Status	Tests	Order No.
	8, 14	CE/IVD	5 (50 μl)	Z-2105-50
			20 (200 μl)	Z-2105-200

ZytoLight® SPEC CD274,PDCD1LG2/CEN 9 Dual Color Probe

Intended Purpose
The ZytoLight SPEC CD274,PDCD1LG2/CEN 9 Dual Color Probe (PL138) is intended to be used for the qualitative detection of amplifications involving the human CD274,PDCD1LG2 gene cluster as well as the detection of the classical satellite III region of chromosome 9 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC CD274,PDCD1LG2/CEN 9 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 9p24.1* (chr9:5,253,553-5,819,972) harboring the CD274, PDCD1LG2 gene cluster (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 9q12 specific for the classical satellite III centromeric region D9Z3 of chromosome 9.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC CD274,PDCD1LG2 Probe map (not to scale); Right: Ideogram of chromosome 9 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the CD274,PDCD1LG2 gene cluster, two green and two orange signals appear.
Aberrant situation: In cells with an amplification of the CD274,PDCD1LG2 gene cluster, an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2179-50
			20 (200 μl)	Z-2179-200

ZytoLight® SPEC JAK2 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC JAK2 Dual Color Break Apart Probe (PL248) is intended to be used for the qualitative detection of translocations involving the human JAK2 gene at 9p24.1 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC JAK2 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 9p24.1-24.2* (chr9:4,311,843-5,031,620) distal to the JAK2 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 9p24.1* (chr9:5,088,700-5,328,239) proximal to the JAK2 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC JAK2 Probe map (not to scale) Right: Ideogram of chromosome 9 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the JAK2 gene region, two green/orange fusion signals appear.

Aberrant situation: One JAK2 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.



Example of an aberrant signal pattern: Bone marrow smear with translocation of the JAK2 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2294-50

ZytoLight® SPEC CDKN2A/CEN 9 Dual Color Probe

Intended Purpose

The ZytoLight SPEC CDKN2A/CEN 9 Dual Color Probe (PL22) is intended to be used for the qualitative detection of deletions involving the human CDKN2A gene as well as the detection of the classical satellite III region of chromosome 9 in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC CDKN2A/CEN 9 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 9p21.3* (chr9:21,742,629-22,056,853) harboring the CDKN2A gene region (see Fig.1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 9q12 specific for the classical satellite III region D9Z3 of chromosome 9.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC CDKN2A Probe map (not to scale); Right: Ideogram of chromosome 9 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a deletion involving the respective gene region (CDKN2A, CDKN2B, MTAP), two green and two orange signals appear.

Aberrant situation: In a cell with deletion affecting the respective gene region, a reduced number of green signals will be observed. Please check localization of the probe (see Fig. 1) carefully. Deletions affecting only parts of the respective gene region might result in a normal signal pattern with green signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: Glioblastoma tissue section with homozygous deletion of the CDKN2A gene as indicated by the loss of both green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2063-50
			20 (200 μl)	Z-2063-200

ZytoLight® SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe

Intended Purpose

The ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe (PL40) is intended to be used for the qualitative detection of the human CDKN2A gene as well as alpha-satellites of chromosomes 3, 7, and 17 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe is composed of:

- ZyGold (excitation 532 nm/emission 553 nm) labeled polynucleotides (~5,5 ng/μl), which target sequences mapping in 9p21.3* (chr9:21,742,629-22,056,853) harboring the CDKN2A gene region (see Fig. 1).
- ZyRed (excitation 580 nm/emission 599 nm) labeled polynucleotides (~0.5 ng/μl), which target sequences mapping in 3p11.1-q11.1 specific for the alpha satellite centromeric region D3Z1 of chromosome 3.
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19




Fig. 1: SPEC CDKN2A Probe map (not to scale)




Ideograms of chromosome 9, 3, 7 and 17 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a deletion involving the CDKN2A gene region, two gold, two red, two green, and two blue signals are expected.

Aberrant situation: In a cell with deletion of the CDKN2A gene region, a reduced number of gold signals will be observed. Deletions affecting only parts of the CDKN2A gene region might result in a normal signal pattern with gold signals of reduced size. In cells with aneuploidy of chromosomes 3, 7, or 17, more or less signals of the respective color will be visible.



Normal cytological specimen hybridized with SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe as indicated by two gold (CDKN2A), two red (CEN 3), two green (CEN 7), and two blue (CEN 17) signals.



Example of an aberrant signal pattern: SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe hybridized to tumor cells showing a trisomy 9 as indicated by three CDKN2A signals (gold) in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	3, 7, 9, 17	CE/IVD	5 (50 μl)	Z-2081-50
			20 (200 μl)	Z-2081-200

ZytoLight® SPEC PAX5 Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC PAX5 Dual Color Break Apart Probe is designed to detect rearrangements involving the chromosomal region 9p13.2 harboring the PAX5 (paired box 5, a.k.a. BSAP) gene.

The transcription factor PAX5 activates crucial genes for B-cell lineage differentiation and represses genes that play an important role in other hematopoietic lineages. PAX5 is also implicated in human B-cell malignancies, as it is deregulated by chromosomal translocations in a subset of acute lymphoblastic leukemias (ALL).

B-progenitor ALL (B-ALL), a common pediatric malignancy, is characterized by the participation of PAX5 in specific chromosomal rearrangements that generate novel fusion proteins. All PAX5 fusion proteins contain the PAX5 DNA-binding domain and thus are predicted to retain the ability to bind to PAX5 transcriptional targets, but no longer provide normal transcriptional regulatory functions. The fusion proteins contribute to B-ALL *forMation* by competitively inhibiting the transcriptional activation of wildtype PAX5. PAX5 rearranged ALL patients were shown to respond well to treatment with prednisone. Hence, the identification of PAX5 rearrangements by FISH may be of therapeutic significance in ALL.


References

Busslinger M, *et al.* (1996) Proc Natl Acad Sci U S A 93: 6129-34.
Cobaleda C, *et al.* (2007) Nat Immunol 8: 463-70.
Coyaud E, *et al.* (2010) Blood 115: 3089-97.
Mullighan CG, *et al.* (2007) Nature 446: 758-64.
Nebral K, *et al.* (2009) Leukemia 23: 134-43.
Offit K, *et al.* (1992) Blood 80: 2594-9.

Probe Description

The *ZytoLight* SPEC PAX5 Dual Color Break Apart Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 9p13.2* (chr9:36,331,787-36,837,502) distal to the PAX5 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 9p13.2* (chr9:37,043,219-37,336,413) proximal to the PAX5 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without an intragenic translocation involving the PAX5 gene region, two green/orange fusion signals appear.

Aberrant situation: One PAX5 gene region affected by an intragenic translocation is indicated by one separate green signal and one separate orange signal. In a cell with an unbalanced translocation resulting in a 3' deletion, one separate orange signal and one orange/green fusion signal will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2300-50

ZytoLight® SPEC NTRK2 Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC NTRK2 Dual Color Break Apart Probe (PL163) is intended to be used for the qualitative detection of translocations involving the human NTRK2 gene at 9q21.33 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC NTRK2 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 9q21.32-q21.33* (chr9:86,569,621-87,287,312) proximal to the NTRK2 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 9q21.33* (chr9:87,589,037-88,124,082) distal to the NTRK2 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC NTRK2 Probe map (not to scale); Right: Ideogram of chromosome 9 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the NTRK2 gene region, two green/orange fusion signals appear.

Aberrant situation: One NTRK2 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated orange signals are the result of deletions proximal to the NTRK2 breakpoint region or are due to unbalanced translocations affecting this chromosomal region.



Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2205-50
			20 (200 μl)	Z-2205-200

ZytoLight® SPEC NR4A3 Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC NR4A3 Dual Color Break Apart Probe (PL102) is intended to be used for the qualitative detection of translocations involving the human NR4A3 gene at 9q22.33-q31.1 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC NR4A3 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 9q22.33* (chr9:102,070,916-102,561,593) proximal to the NR4A3 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 9q31.1* (chr9:102,636,487-103,065,504) distal to the NR4A3 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC NR4A3 Probe map (not to scale); Right: Ideogram of chromosome 9 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the NR4A3 gene region, two green/orange fusion signals appear.

Aberrant situation: One NR4A3 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2145-50

ZytoLight® SPEC ABL1 Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC ABL1 Dual Color Break Apart Probe (PL157) is intended to be used for the qualitative detection of translocations involving the human ABL1 gene at 9q34.12 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC ABL1 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 9q34.11-q34.12* (chr9:132,872,357-133,580,236) proximal to the ABL1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 9q34.12-q34.13* (chr9:133,851,960-134,273,097) distal to the ABL1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ABL1 Probe map (not to scale); Right: Ideogram of chromosome 9 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the ABL1 gene region, two green/orange fusion signals appear.

Aberrant situation: One ABL1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. One ABL1 gene region affected by a NUP214-ABL1 fusion associated with amplified episomes is indicated by one separate green signal and multiple orange signals or an orange signal cluster.



Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2199-50

ZytoLight® SPEC NUP214 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC NUP214 Dual Color Break Apart Probe is designed for the detection of translocations involving the chromosomal region 9q34.13 harboring the NUP214 (nucleoporin 214, a.k.a. CAN, CAIN) gene.

Rearrangements of the NUP214 gene have been implicated in the pathogenesis of several types of hematologic malignancies, including T-cell acute lymphoblastic leukemia (T-ALL), acute myeloid leukemia (AML), and also myelodysplastic syndrome (MDS). Several fusion partners have been identified for NUP214. The most common are the DEK, SET, and the tyrosine kinase encoding gene ABL1.

The translocation t(6;9)(p22.3;q34.1) results in a DEK-NUP214 fusion and defines a specific subcategory of AML according to the World Health Organization 2008 classification.

The SET-NUP214 fusion is associated with T-ALL, less frequently with AML, and acute undifferentiated leukemia and can result from either a translocation or a deletion. NUP214-ABL1 fusions are exclusively associated with T-ALL patients. These patients may be considered for a targeted therapy with specific tyrosine kinase inhibitors. The fusion is often located on amplified episomes and is cytogenetically cryptic but can be detected by FISH.

Malignancies with NUP214 rearrangements are associated with a poor prognosis indicating the usefulness of NUP214 also as a prognostic biomarker.

References
Fahrenkrog B (2014) New J Sci 2014: 468306.
Takeda A & Yaseen NR (2014) Semin Cancer Biol 27: 3-10.
von Lindern M, *et al.* (1992) Baillieres Clin Haematol 5: 857-79.
Zhou MH & Yang QM (2014) Oncol Lett 8: 959-62.

Probe Description
The ZytoLight SPEC NUP214 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 9q34.13* (chr9:134,153,663-134,706,700) distal to the NUP214 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 9q34.12-q34.13* (chr9:133,739,333-134,028,546) proximal to the NUP214 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC NUP214 Probe map (not to scale)



Ideogram of chromosome 9 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the NUP214 gene region, two green/orange fusion signals appear.
Aberrant situation: One NUP214 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. One locus affected by a NUP214-ABL1 fusion associated with amplified episomes is indicated by one separate green signal and multiple orange signals or an orange signal cluster. A deletion affecting the 9q34.12-q34.13 locus is indicated by one green/orange fusion signal and one separate green signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.

Labels	Chromosome	Reg. Status	Tests	Order No.
	9	CE/IVD	5 (50 μl)	Z-2265-50

ZytoLight® SPEC BCR/ABL1 Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe (PL68) is intended to be used for the qualitative detection of the translocation t(9;22)(q34.1;q11.2) involving the human BCR and ABL1 genes in cytologic specimens, such as chronic myeloid leukemia (CML), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of CML and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 22q11.22-q11.23* (chr22:23,000,029-24,431,064) harboring the BCR gene region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 9q34.11-q34.13* (chr9:133,223,081-134,103,849) harboring the ABL1 gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC BCR Probe map; Bottom: SPEC ABL1 Probe map (not to scale); Right: Ideograms of chromosome 22 and 9 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Bone marrow specimen with translocation affecting the BCR/ABL1 loci as indicated by one separate orange signal, one separate green signal and two orange/green fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	9, 22	CE/IVD	5 (50 μl)	Z-2111-50
			20 (200 μl)	Z-2111-200

ZytoLight® SPEC RET Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC RET Dual Color Break Apart Probe (PL105) is intended to be used for the qualitative detection of translocations involving the human RET gene at 10q11.21 in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC) or papillary thyroid carcinoma (PTC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC or PTC and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC RET Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 10q11.21* (chr10:43,626,274-43,902,346) distal to the RET breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 10q11.21* (chr10:43,340,888-43,510,171) proximal to the RET breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC RET Probe map (not to scale); Right: Ideogram of chromosome 10 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the RET gene region, two green/orange fusion signals appear.
Aberrant situation: One RET gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated green signals are the result of deletions proximal to the RET breakpoint region. Inversions involving genes in close proximity such as CCDC6 might not be observable due to just subtle signal separation.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Paraffin-embedded human thyroid tumor cell line (TPC-1) with translocation affecting the 10q11.21 locus as indicated by one orange/green fusion (non-rearranged) signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	10	CE/IVD	5 (50 μl)	Z-2148-50
			20 (200 μl)	Z-2148-200

ZytoLight® SPEC PTEN/CEN 10 Dual Color Probe

Intended Purpose

The *ZytoLight* SPEC PTEN/CEN 10 Dual Color Probe (PL37) is intended to be used for the qualitative detection of deletions involving the human PTEN gene as well as the detection of chromosome 10 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC PTEN/CEN 10 Dual Color Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 10q23.2-q23.31* (chr10:89,440,649-89,755,790) harboring the PTEN gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 10p11.1-q11.1 specific for the alpha satellite centromeric region D10Z1 of chromosome 10.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Results

Normal situation: In interphases of normal cells or cells without a deletion involving the PTEN gene region, two green and two orange signals appear.

Aberrant situation: In a cell with deletion affecting the PTEN gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the PTEN gene region might result in normal signal pattern with green signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: Melanoma tissue section with chromosome 10 monosomy as indicated by one orange and one green signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	10	CE/IVD	5 (50 μl)	Z-2078-50
			20 (200 μl)	Z-2078-200

ZytoLight® SPEC FGFR2 Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC FGFR2 Dual Color Break Apart Probe (PL125) is intended to be used for the qualitative detection of translocations involving the human FGFR2 gene at 10q26.13 in formalin-fixed, paraffin-embedded specimens, such as intrahepatic cholangiosarcoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of intrahepatic cholangiosarcoma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC FGFR2 Dual Color Break Apart Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 10q26.12-q26.13* (chr10:122,632,462-123,166,030) proximal to the FGFR2 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 10q26.13* (chr10:123,436,230-123,791,279) distal to the FGFR2 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Results

Normal situation: In interphases of normal cells or cells without a translocation involving the FGFR2 gene region, two green/orange fusion signals appear.

Aberrant situation: One FGFR2 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: Cholangiocellular adenocarcinoma tissue section with translocation of the FGFR2 gene as indicated by one non-rearranged orange/green fusion signal, one orange, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	10	CE/IVD	5 (50 μl)	Z-2169-50
			20 (200 μl)	Z-2169-200

Kindly provided by Prof. Dr. Büttner, Cologne, Germany.

ZytoLight® SPEC FGFR2/CEN 10 Dual Color Probe

Intended Purpose

The *ZytoLight* SPEC FGFR2/CEN 10 Dual Color Probe (PL79) is intended to be used for the qualitative detection of amplifications involving the human FGFR2 gene as well as chromosome 10 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC FGFR2/CEN 10 Dual Color Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 10q26.12-q26.13* (chr10:122,908,224-123,682,373) harboring the FGFR2 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 10p11.1-q11.1 specific for the alpha satellite centromeric region D10Z1 of chromosome 10.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Results

Normal situation: In interphases of normal cells or cells without an amplification involving the FGFR2 gene region, two green signals and two orange signals appear.

Aberrant situation: In cells with an amplification of the FGFR2 gene region an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: Breast cancer tissue section with amplification of the FGFR2 gene as indicated by green signal clusters in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	10	CE/IVD	20 (200 μl)	Z-2122-200

ZytoLight® SPEC NUP98 Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC NUP98 Dual Color Break Apart Probe (PL223) is intended to be used for the qualitative detection of translocations involving the human NUP98 gene at 11p15.4 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The *ZytoLight* SPEC NUP98 Dual Color Break Apart Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10.0 ng/μl), which target sequences mapping in 11p15.4-p15.5* (chr11:2,773,748-3,363,120) distal to the NUP98 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 11p15.4* (chr11:3,829,054-4,142,792) proximal to the NUP98 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without a translocation involving the NUP98 gene region, two green/orange fusion signals appear. In specimens from normal tissue split signals have occasionally been observed. It is therefore important to define cut-off specifications carefully in order to avoid false positive results.

Aberrant situation: One NUP98 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.

Labels	Chromosome	Reg. Status	Tests	Order No.
	11	CE/IVD	5 (50 μl)	Z-2266-50

ZytoLight® SPEC CCND1/CEN 11 Dual Color Probe

Intended Purpose
The ZytoLight SPEC CCND1/CEN 11 Dual Color Probe (PL28) is intended to be used for the qualitative detection of amplifications involving the human CCND1 gene as well as the detection of chromosome 11 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer, by fluorescence in situ hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight CCND1/CEN 11 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 11q13.3* (chr11:69,203,885-70,031,240) harboring the CCND1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 11p11.11-q11 specific for the alpha satellite centromeric region D11Z1 of chromosome 11.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPECCND1 Probe map (not to scale)



Results
Normal situation: In interphases of normal cells or cells without an amplification involving the CCND1 gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the CCND1 gene region, an increased number of green signals or green signal clusters are visible.



Labels	Chromosome	Reg. Status	Tests	Order No.
	11	CE/IVD	5 (50 μl)	Z-2071-50
			20 (200 μl)	Z-2071-200

ZytoLight® SPEC CCND1/IGH Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe (PL82) is intended to be used for the qualitative detection of the translocation t(11;14) (q13.3;q32.3) in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence in situ hybridization (FISH). The probe is intended to be used in combination with ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 11q13.3* (chr11:68,522,105-70,031,240) harboring the CCND1 gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH locus (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPECCND1 Probe map; Bottom: SPEC IGH Probe map (not to scale); Right: Ideograms of chromosome 11 and 14 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective loci, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Labels	Chromosome	Reg. Status	Tests	Order No.
	11, 14	CE/IVD	5 (50 μl)	Z-2125-50
			20 (200 μl)	Z-2125-200

ZytoLight® SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe is designed to detect translocations involving the chromosomal region 11q22.2 harboring the BIRC3 (baculoviral IAP repeat containing 3, a.k.a. API2) gene and the chromosomal region 18q21.32 harboring the MALT1 (MALT1 paracaspase, a.k.a. MLT) gene.

The recurrent translocation t(11;18)(q22.2;q21.3) is frequently found in mucosa-associated lymphoid tissue (MALT) lymphoma which represents the most common extranodal B-cell tumor and accounts for 5-10% of all non-Hodgkin lymphoma. The translocation results in the expression of chimeric fusion transcripts comprising the N-terminal end of the apoptosis inhibitor BIRC3 which is highly expressed in adult lymphoid tissue and C-terminal parts of the MALT1 protease.

The BIRC3/MALT1 fusion protein was shown to induce proteolytic cleavage of NF-kappa-B-inducing kinase (NIK) ultimately resulting in constitutive non-canonical NF-kappa-B signaling, enhanced B-cell adhesion, and apoptosis resistance. It is assumed that disruption of the BIRC3-NIK interaction and/or blocking of MALT1 protease or NIK kinase activity could represent new treatment approaches for refractory t(11;18)-positive MALT lymphoma.

- References**
Dierlamm J, *et al.* (1999) Blood 93: 3601-9.
Dierlamm J, *et al.* (2000) Blood 96: 2215-8.
Morgan JA, *et al.* (1999) Cancer Res 59: 6205-13.
Rosebeck S, *et al.* (2011) Science 331: 468-72.

- Probe Description**
The ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 11q22.1-q22.2* (chr11:101,756,072-102,581,817) harboring the BIRC3 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 18q21.31-q21.32* (chr18:56,021,766-56,724,408) harboring the MALT1 gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC BIRC3 Probe map; Bottom: SPEC MALT1 Probe map (not to scale)



Ideograms of chromosome 11 and 18 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Labels	Chromosome	Reg. Status	Tests	Order No.
	11, 18	CE/IVD	5 (50 μl)	Z-2146-50
			20 (200 μl)	Z-2146-200

ZytoLight® SPEC MAML2 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC MAML2 Dual Color Break Apart Probe (PL5) is intended to be used for the qualitative detection of translocations involving the human MAML2 gene at 11q21 in formalin-fixed, paraffin-embedded mucoepidermoid carcinoma (MEC) specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

This product is for non-automated and for laboratory professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel. This product is neither a device for near-patient testing nor a companion diagnostic.

This product is a FISH probe with ZyGreen- and ZyOrange-labeled polynucleotides intended to be used as an aid to the differential diagnosis of MEC in patients with confirmed or suspected MEC. Therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoLight SPEC MAML2 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 11q21* (chr11:96,115,829-96,797,136) distal to the MAML2 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 11q21* (chr11:95,296,828-95,668,215) proximal to the MAML2 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC MAML2 Probe map (not to scale); Right: Ideogram of chromosome 11 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the MAML2 gene region, two green/orange fusion signals appear.

Aberrant situation: One MAML2 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Mucoepidermoid carcinoma section with translocation affecting the 11q21 locus as indicated by one separate orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	11	CE/IVD	5 (50 μl)	Z-2014-50
			20 (200 μl)	Z-2014-200

ZytoLight® SPEC ATM/CEN 11 Dual Color Probe

Intended Purpose
The ZytoLight SPEC ATM/CEN 11 Dual Color Probe (PL254) is intended to be used for the qualitative detection of deletions involving the human ATM gene as well as the detection of chromosome 11 alpha satellites in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoLight SPEC ATM/CEN 11 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 11q22.3* (chr11:107,957,618-108,380,921) harboring the ATM gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 11p11.11-q11 specific for the alpha satellite centromeric region D11Z1 of chromosome 11.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ATM Probe map (not to scale); Right: Ideogram of chromosome 11 indicating the hybridization locations

Results


Normal situation: In interphases of normal cells or cells without a deletion involving the ATM gene region, two green and two orange signals appear.

Aberrant situation: In a cell with deletions affecting the ATM gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the ATM gene region might result in a normal signal pattern with orange signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.

Labels	Chromosome	Reg. Status	Tests	Order No.
	11	CE/IVD	5 (50 μl)	Z-2297-50



Example of an aberrant signal pattern: CLL with deletion affecting the ATM locus as indicated by one orange signal in each nucleus.

ZytoLight® SPEC ATM/CEN 12 Dual Color Probe

Intended Purpose
The ZytoLight SPEC ATM/CEN 12 Dual Color Probe (PL250) is intended to be used for the qualitative detection of deletions involving the human ATM gene as well as the detection of chromosome 12 alpha satellites in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoLight SPEC ATM/CEN 12 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 11q22.3* (chr11:107,957,618-108,380,921) harboring the ATM gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 12p11.1-q11 specific for the alpha satellite centromeric region D12Z3 of chromosome 12.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ATM Probe map (not to scale); Right: Ideograms of chromosome 11 and 12 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a deletion involving the ATM gene region and without aneusomy of chromosomes 12, two orange and two green signals appear.

Aberrant situation: In a cell with deletions affecting the ATM gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the ATM gene region might result in a normal signal pattern with orange signals of reduced size. In a cell with aneusomy 12, three or more copies of the green signal will be observed, respectively.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.



Example of an aberrant signal pattern: CLL with deletion of the ATM gene and amplification affecting the centromeric region of chromosome 12 as indicated by one orange signal and five or more green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	11, 12	CE/IVD	5 (50 μl)	Z-2296-50

ZytoLight® SPEC TP53/ATM Dual Color Probe

Intended Purpose
The ZytoLight SPEC TP53/ATM Dual Color Probe (PL115) is intended to be used for the qualitative detection of deletions involving the human TP53 gene as well as the human ATM gene in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoLight SPEC TP53/ATM Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 17p13.1* (chr17:7,495,749-7,663,022) harboring the TP53 gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 11q22.3* (chr11:107,957,618-108,380,921) harboring the ATM gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC TP53 Probe map; Bottom: SPEC ATM Probe map (not to scale); Right: Ideograms of chromosome 17 and 11 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a deletion involving the respective gene regions, two green and two orange signals appear.

Aberrant situation: In a cell with deletions affecting the TP53 gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the TP53 locus might result in a normal signal pattern with orange signals of reduced size. In a cell with deletions affecting the ATM gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the ATM gene region might result in a normal signal pattern with green signals of reduced size.



Example of an aberrant signal pattern: Hybridized to bone marrow biopsy with deletions of the ATM and the TP53 genes as indicated by one green and one orange signal in each nucleus.



Example of an aberrant signal pattern: SPEC TP53/ATM Dual Color Probe hybridized to bone marrow smear with deletion of the ATM gene as indicated by one green signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	11, 17	CE/IVD	5 (50 μl)	Z-2159-50
			20 (200 μl)	Z-2159-200

ZytoLight® SPEC D13S319/13q34/CEN 12 Triple Color Probe

Intended Purpose
The ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe (PL116) is intended to be used for the qualitative detection of deletions involving the human D13S319 region as well as the detection of chromosome 13q34 specific sequences and chromosome 12 alpha satellites in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 13q14.2* (chr13:50,607,438-50,781,256) harboring the D13S319 locus (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37 ng/μl), which target sequences mapping in 13q34* (chr13:113,691,216-114,323,467) harboring the LAMP1 gene region (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 12p11.1-q11 specific for the alpha satellite centromeric region D12Z3 of chromosome 12.

* according to Human Genome Assembly GRCh37/hg19

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the D13S319 locus and without aneusomy of chromosomes 12 and 13, two orange, two blue, and two green signals appear.
Aberrant situation: In a cell with deletions affecting the D13S319 locus, a reduced number of orange signals will be observed. Deletions affecting only parts of the D13S319 locus might result in a normal signal pattern with orange signals of reduced size. In a cell with monosomy 13, a reduced number of both, the orange (D13S319 locus) and blue (13q34 locus) signals will be observed. In a cell with trisomy or polysomy 12, three or more copies of the green signal will be observed, respectively.



Example of an aberrant signal pattern: SPEC D13S319/13q34/CEN 12 Triple Color Probe hybridized to bone marrow biopsy section with deletion of the D13S319 locus as indicated by one orange signal and two blue signals in each nucleus.

Example of an aberrant signal pattern: SPEC D13S319/13q34/CEN 12 Triple Color Probe hybridized to bone marrow smear with trisomy of chromosome 12 as indicated by three green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	12, 13	CE/IVD	5 (50 μl)	Z-2160-50
			20 (200 μl)	Z-2160-200





Fig. 1: Top: SPEC D13S319 Probe map; Bottom: SPEC 13q34 Probe map (not to scale)



Ideograms of chromosome 13 and 12 indicating the hybridization locations

ZytoLight® SPEC D13S319/13q34 Dual Color Probe

Intended Purpose
The ZytoLight SPEC D13S319/13q34 Dual Color Probe (PL235) is intended to be used for the qualitative detection of deletions involving the human D13S319 region and chromosome 13q34 specific sequences in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC D13S319/13q34 Dual Color Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 13q14.2* (chr13:50,607,438-50,781,256) harboring the D13S319 locus (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 13q34 (chr13:113,691,216-114,323,467) (see Fig. 1).

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC D13S319 Probe map; Bottom: SPEC 13q34 Probe map (not to scale)

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the D13S319 locus, two orange signals and two green signals appear.
Aberrant situation: In a cell with deletion affecting the D13S319 locus, a reduced number of orange signals will be observed. Deletions affecting only parts of the D13S319 locus might result in a normal signal pattern with orange signals of reduced size. In a cell with deletions affecting the 13q34 locus, a reduced number of green signals will be observed.



Hybridized to normal interphase cells as indicated by two green and two orange signals in each nucleus.

Example of an aberrant signal pattern: Bone marrow smear with deletion of the D13S319 locus as indicated by one orange and two green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	13	CE/IVD	5 (50 μl)	Z-2280-50
			20 (200 μl)	Z-2280-200

ZytoLight® SPEC 11q gain/loss Triple Color Probe

Intended Purpose
The ZytoLight SPEC 11q gain/loss Triple Color Probe (PL174) is intended to be used for the qualitative detection of human 11q alterations involving human 11q23.3 and 11q24.1-q25 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC 11q gain/loss Triple Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in the minimal gained region (MGR) at 11q23.3* (chr11:117,574,074-118,284,029) (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in the minimal loss region (MLR) at 11q24.3* (chr11:128,707,454-129,161,227) (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 11p11.1-q11 specific for the alpha satellite centromeric region D11Z1 of chromosome 11.

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC 11q Probe map (not to scale); Right: Ideogram of chromosome 11 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an 11q alteration, two green, two orange, and two blue signals appear.
Aberrant situation: In interphases of cells affected by an 11q-gain/loss alteration, one green and one orange signal is expected for the unaffected chromosome, and a gain of green signals and loss of one orange signal is expected for the chromosome affected by an 11q-gain/loss alteration. Deletions affecting only parts of the 11q loss region might result in a normal signal pattern with orange signals of reduced size.



Hybridized to normal interphase cells as indicated by two green, two orange, and two blue signals per nucleus and to metaphase chromosomes of a normal cell.

Example of an aberrant signal pattern: Burkitt-like lymphoma tissue section with 11q aberration as indicated by three green signals and one orange signal indicating the gain and loss at 11q, respectively.

Labels	Chromosome	Reg. Status	Tests	Order No.
	11	CE/IVD	5 (50 μl)	Z-2216-50

ZytoLight® SPEC KMT2A Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC KMT2A Dual Color Break Apart Probe (PL151) is intended to be used for the qualitative detection of translocations involving the human KMT2A gene at 11q23.3 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC KMT2A Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 11q23.3* (chr11:117,574,074-118,284,029) proximal to the KMT2A breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 11q23.3* (chr11:118,399,293-119,237,675) distal to the KMT2A breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC KMT2A Probe map (not to scale); Right: Ideogram of chromosome 11 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the KMT2A gene region, two green/orange fusion signals appear.
Aberrant situation: One KMT2A gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.
Example of an aberrant signal pattern: Bone marrow smear with translocation of the KMT2A gene as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	11	CE/IVD	5 (50 μl)	Z-2193-50
			20 (200 μl)	Z-2193-200

ZytoLight® SPEC ZNF384 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC ZNF384 Dual Color Break Apart Probe is designed to detect translocations involving the chromosomal region 12p13.31 harboring the ZNF384 gene.

The ZNF384 (zinc-finger protein 384, a.k.a. CIZ) gene encodes a transcription factor involved in the regulation of matrix metalloproteinases. Rearrangements of the ZNF384 gene are recurrent in acute leukemia and are most frequently found in precursor B-cell acute lymphoblastic leukemia (BCP-ALL) in children and young adults with an incidence of about 3-4%. ZNF384-related fusion genes with multiple fusion partners have been found to define a distinct subgroup of pediatric BCP-ALL with a characteristic immunophenotype. Known translocation partners are TCF3 (19p13.3), EWSR1 (22q12.2), TAF15 (17q12), EP300 (22q13.2), ARID1B (6q25.3), CREBBP (16p13.3), and BMP2K (4q21.21) with TCF3 being the most prevalent. The breakpoints are located within the ZNF384 gene. However, the balanced translocations are resulting in fusion genes including the complete protein coding *inforMation*. Since ZNF384-related fusion genes are difficult to detect by common G-banding, investigation by FISH may be of diagnostic and prognostic relevance.

- References**
Hirabayashi S, *et al.* (2017) Haematologica 102: 118-29.
Krance RA, *et al.* (1992) Leukemia 6: 251-5.
La Starza R, *et al.* (2005) Leukemia 19: 1696-9.
Shago M, *et al.* (2016) Pediatr Blood Cancer 63: 1915-21.

- Probe Description**
The ZytoLight SPEC ZNF384 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 12p13.31* (chr12:6,016,809-6,771,300) distal to the ZNF384 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 12p13.31* (chr12:6,799,546-7,175,222) proximal to the ZNF384 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ZNF384 Probe map (not to scale); Right: Ideogram of chromosome 12 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the ZNF384 gene region, two green/orange fusion signals appear.
Aberrant situation: One ZNF384 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.
Bone marrow smear of an ALL case with translocation of the ZNF384 gene as indicated by one orange/green fusion signal, one separate green, and one separate orange signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	12	CE/IVD	5 (50 μl)	Z-2275-50

ZytoLight® SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe is designed for the detection of the specific translocation involving the chromosomal region 12p13.2 harboring the ETV6 (ETS variant 6, a.k.a. TEL) gene and the chromosomal region 21q22.12 harboring the RUNX1 (runt related transcription factor 1, a.k.a. AML1) gene. The balanced chromosomal translocation t(12;21)(p13.2;q22.1), which leads to ETV6/RUNX1 fusion, represents the most frequent genetic rearrangement in initial childhood B-cell precursor (BCP) acute lymphoblastic leukemia (ALL) (19-27%) and has been associated with good prognosis.

The ETV6/RUNX1 fusion protein, comprising a putative repressor domain of ETV6, a member of the ETS family of transcription factors, fused to RUNX1, the DNA-binding subunit of the RUNX1/CBF beta transcription factor complex, acts as a trans-dominant repressor of RUNX1 regulated target genes involved in hematopoiesis.

Three secondary aberrations in ETV6/RUNX1 positive ALL have been found to negatively influence the clinical course: deletion of the second non-translocated ETV6 allele, gains of the RUNX1 gene, and duplication of the derivative chromosome 21.

Detection of t(12;21) by Fluorescence *in situ* Hybridization enables the simultaneous identification of the most common secondary changes and thus provides additional *inforMation* about the possible outcome of the disease in patients with ALL.

- References**
Fenrick R, *et al.* (1999) Mol Cell Biol 19: 6566-74.
Martínez-Ramírez A, *et al.* (2001) Haematologica 86: 1245-53.
Morrow M, *et al.* (2007) Oncogene 26: 4404-14.
Peter A, *et al.* (2009) Eur J Haematol 83: 420-32.
Shurtleff SA, *et al.* (1995) Leukemia 9: 1985-9.

- Probe Description**
The ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 12p13.2* (chr12:11,586,400-12,454,330) harboring the ETV6 region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 21q22.12* (chr21:36,106,492-36,657,941) harboring the RUNX1 gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC ETV6 Probe map; Bottom: SPEC RUNX1 Probe map (not to scale); Right: Ideograms of chromosome 12 and 21 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.
Bone marrow tissue section with translocation affecting the ETV6/RUNX1 loci as indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	12, 21	CE/IVD	5 (50 μl)	Z-2157-50
			20 (200 μl)	Z-2157-200

ZytoLight® SPEC ETV6 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC ETV6 Dual Color Break Apart Probe (PL135) is intended to be used for the qualitative detection of translocations involving the human ETV6 gene at 12p13.2 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ETV6 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 12p13.2* (chr12:12,054,737-12,634,328) proximal to the ETV6 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 12p13.2* (chr12:11,393,774-11,808,608) distal to the ETV6 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ETV6 Probe map (not to scale)




Ideogram of chromosome 12 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the ETV6 gene region, two green/orange fusion signals appear.

Aberrant situation: One ETV6 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: MASC tissue section of the salivary glands with translocation of the ETV6 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	12	CE/IVD	5 (50 μl)	Z-2176-50
			20 (200 μl)	Z-2176-200

ZytoLight® SPEC KRAS/CEN 12 Dual Color Probe

Intended Purpose

The ZytoLight SPEC KRAS/CEN 12 Dual Color Probe (PL72) is intended to be used for the qualitative detection of amplifications involving the human KRAS gene as well as chromosome 12 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC KRAS/CEN 12 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 12p12.1* (chr12:24,916,728-25,839,353) harboring the KRAS gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 12p11.1-q11 specific for the alpha satellite centromeric region D12Z3 of chromosome 12.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC KRAS Probe map (not to scale)




Ideogram of chromosome 12 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without an amplification involving the KRAS gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the KRAS gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: Lung cancer tissue section with amplification of the KRAS gene (green).

Labels	Chromosome	Reg. Status	Tests	Order No.
	12	CE/IVD	20 (200 μl)	Z-2115-200

Image kindly provided by Prof. Diebold, Lucerne, Switzerland.

ZytoLight® SPEC DDIT3 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC DDIT3 Dual Color Break Apart Probe (PL58) is intended to be used for the qualitative detection of translocations involving the human DDIT3 gene at 12q13.3 in formalin-fixed, paraffin-embedded specimens, such as myxoid liposarcomas (MLPS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of MLPS and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC DDIT3 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 12q13.3-14.1* (chr12:58,024,366-58,775,832) distal to the DDIT3 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 12q13.3* (chr12:57,386,312-57,865,800) proximal to the DDIT3 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC DDIT3 Probe map (not to scale); Right: Ideogram of chromosome 12 indicating the hybridization locations




Results

Normal situation: In interphases of normal cells or cells without a translocation involving the DDIT3 gene region, two green/orange fusion signals appear.

Aberrant situation: One DDIT3 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Myxoid liposarcoma tissue section with translocation affecting the 12q13.3-q14.1 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	12	CE/IVD	5 (50 μl)	Z-2100-50
			20 (200 μl)	Z-2100-200

ZytoLight® SPEC CDK4/CEN 12 Dual Color Probe

Intended Purpose

The ZytoLight SPEC CDK4/CEN 12 Dual Color Probe (PL61) is intended to be used for the qualitative detection of amplifications involving the human CDK4 gene as well as the detection of chromosome 12 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as atypical lipomatous tumor/well-differentiated liposarcoma (ALT/WDLPS) and dedifferentiated liposarcoma (DDLPS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ALT/WDLPS and DDLPS and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC CDK4/CEN 12 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 12q13.3-q14.1* (chr12:57,740,440-58,435,534) harboring the CDK4 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in 12p11.1-q11 specific for the alpha satellite centromeric region D12Z3 of chromosome 12.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC CDK4 Probe map (not to scale); Right: Ideogram of chromosome 12 indicating the hybridization locations




Results

Normal situation: In interphases of normal cells or cells without an amplification involving the CDK4 gene, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the CDK4 gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Liposarcoma tissue section, CDK4 signal cluster (green), CEN 12 (orange).

Labels	Chromosome	Reg. Status	Tests	Order No.
	12	CE/IVD	5 (50 μl)	Z-2103-50
			20 (200 μl)	Z-2103-200

ZytoLight® SPEC MDM2/CEN 12 Dual Color Probe

Intended Purpose
The ZytoLight SPEC MDM2/CEN 12 Dual Color Probe (PL18) is intended to be used for the qualitative detection of amplifications involving the human MDM2 gene as well as the detection of chromosome 12 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as atypical lipomatous tumor/well differentiated liposarcoma (ALT/WDLPS) and dedifferentiated liposarcoma (DDLPS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ALT/WDLPS and DDLPS and therapeutic measures should not be initiated based on the test result alone.


Probe Description
The ZytoLight SPEC MDM2/CEN 12 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 12q15* (chr12:69,071,802-69,565,421) harboring the MDM2 gene region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 12p11.1-q11 specific for the alpha satellite centromeric region D12Z3 of chromosome 12.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without an amplification involving the MDM2 gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the MDM2 gene region an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
	12	CE/IVD	5 (50 μl)	Z-2013-50
			20 (200 μl)	Z-2013-200

ZytoLight® SPEC FOXO1 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC FOXO1 Dual Color Break Apart Probe (PL96) is intended to be used for the qualitative detection of translocations involving the human FOXO1 gene at 13q14.11 in formalin-fixed, paraffin-embedded specimens, such as alveolar rhabdomyosarcoma (ARMS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ARMS and therapeutic measures should not be initiated based on the test result alone.


Probe Description
The ZytoLight SPEC FOXO1 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 13q14.11* (chr13:40,285,558-41,132,595) proximal to the FOXO1 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 13q14.11* (chr13:41,246,917-42,054,781) distal to the FOXO1 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the FOXO1 gene region, two green/orange fusion signals appear.
Aberrant situation: One FOXO1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
	13	CE/IVD	5 (50 μl)	Z-2139-50

ZytoLight® SPEC FOXO1/PAX3 Dual Color Single Fusion Probe

Intended Purpose
The ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe (PL16) is intended to be used for the qualitative detection of translocation t(2;13) (q36;q14.1) involving the human FOXO1 and PAX3 genes in formalin-fixed, paraffin-embedded specimens, such as alveolar rhabdomyosarcoma (ARMS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ARMS and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 13q14.11* (chr13:40,816,168-41,132,595) proximal to the FOXO1 breakpoint region (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 2q36.1* (chr2:223,196,078-223,539,352) distal to the PAX3 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC FOXO1 Probe map; Bottom: SPEC PAX3 Probe map (not to scale); Right: Ideograms of chromosome 13 and 2 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and one orange/green fusion signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
	2, 13	CE/IVD	5 (50 μl)	Z-2018-50
			20 (200 μl)	Z-2018-200

ZytoLight® SPEC FOXO1/PAX3 TriCheck™ Probe

Intended Purpose
The ZytoLight SPEC FOXO1/PAX3 TriCheck Probe (PL143) is intended to be used for the qualitative detection of rearrangements involving the human FOXO1 gene at 13q14.11 and the human PAX3 gene at 2q36.1 in formalin-fixed, paraffin-embedded specimens, such as alveolar rhabdomyosarcoma (ARMS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ARMS and therapeutic measures should not be initiated based on the test result alone.


Probe Description
The ZytoLight SPEC FOXO1/PAX3 TriCheck Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 13q14.11* (chr13:41,246,917-42,054,781) distal to the FOXO1 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 13q14.11* (chr13:40,578,036-41,132,595) proximal to the FOXO1 breakpoint region (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/μl), which target sequences mapping in 2q36.1* (chr2:223,196,078-223,936,825) distal to the PAX3 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19



Results
Normal situation: In interphases of normal cells or cells without PAX3-FOXO1 rearrangement, green/orange fusion signals and two blue signals appear.
Aberrant situation: A PAX3-FOXO1 fusion is indicated by one separate orange signal co-localizing with one blue signal and one separate green signal. A FOXO1 translocation without involvement of PAX3 is indicated by the split of one green/orange fusion signal without co-localization of the separated orange signal with one blue signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
	2, 13	CE/IVD	5 (50 μl)	Z-2185-50

ZytoLight® SPEC FOXO1/PAX7 Dual Color Single Fusion Probe

Intended Purpose
The *ZytoLight* SPEC FOXO1/PAX7 Dual Color Single Fusion Probe (PL17) is intended to be used for the qualitative detection of translocation t(1;13) (p36.1;q14.1) involving the human FOXO1 and PAX7 genes in formalin-fixed, paraffin-embedded specimens, such as alveolar rhabdomyosarcoma (ARMS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ARMS and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The *ZytoLight* SPEC FOXO1/PAX7 Dual Color Single Fusion Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 13q14.11* (chr13:40,816,168-41,132,595) proximal to the FOXO1 breakpoint region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1p36.13* (chr1:18,139,970-18,956,785) distal to the PAX7 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Fig. 1: Top: SPEC FOXO1 Probe map; Bottom: SPEC PAX7 Probe map (not to scale); Right: Ideograms of chromosome 13 and 1 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and one orange/green fusion signal.



Hybridized to normal interphase cells as indicated by two orange and two green signals.

ZytoLight® SPEC RB1/13q34 Dual Color Probe

Intended Purpose
The *ZytoLight* SPEC RB1/13q34 Dual Color Probe is designed for the detection of deletions affecting the RB1 gene. The RB1 (RB transcriptional corepressor 1, a.k.a. pRb) gene is located on 13q14.2 and encodes a protein which acts as a tumor suppressor playing a crucial role in cell cycle regulation and genome stability.

Deletions of RB1 are frequently found in retinoblastoma. However, either monoallelic or biallelic deletions of RB1 are also common in a wide variety of solid tumors and hematologic malignancies such as multiple myeloma (MM) and chronic lymphocytic leukemia (CLL). While 13q14 deletions exclusive of RB1 confer a more favorable prognosis in CLL patients, 13q14 deletions that encompass the RB1 locus (present in approx. 20% of all CLL cases) are associated with shortened survival.

Hence, FISH is a valuable tool for the detection of RB1 gene deletions and can be used in combination with further biological markers, morphology and clinical information for the prediction of disease progression and overall survival.

- References**
Dal Bo M, *et al.* (2011) Genes Chromosomes Cancer 50: 633-43.
Dao DD, *et al.* (1994) Leukemia 8: 1280-4.
Di Fiore R, *et al.* (2013) J Cell Physiol 228: 1676-87.
Juge-Morineau N, *et al.* (1997) Leuk Lymphoma 24: 229-37.
Orlandi EM, *et al.* (2013) Hematol Oncol 31: 136-42.
Ouillette P, *et al.* (2011) Clin Cancer Res 17: 6778-90.

- Probe Description**
The *ZytoLight* SPEC RB1/13q34 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 13q14.2* (chr13:48,776,918-49,092,570) harboring the RB1 gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10.0 ng/μl), which target sequences mapping in 13q34* (chr13:113,691,216-114,323,467) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Fig. 1: Top: SPEC RB1 Probe map; Bottom: SPEC 13q34 Probe map

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the RB1 gene region, two orange signals and two green signals appear.
Aberrant situation: In a cell with deletion affecting the RB1 gene region, a reduced number of the orange signals will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Hybridized to lipoma tissue section with deletion of the RB1 gene as indicated by one orange signal and two green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	13	CE/IVD	5 (50 μl)	Z-2324-50
			20 (200 μl)	Z-2324-200

ZytoLight® SPEC IGH Dual Color Break Apart Probe

Intended Purpose
The *ZytoLight* SPEC IGH Dual Color Break Apart Probe (PL67) is intended to be used for the qualitative detection of translocations involving the human IGH locus at 14q32.33 in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with *ZytoLight* FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).


The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.


- Probe Description**
The *ZytoLight* SPEC IGH Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 14q32.33* (chr14:106,690,778-107,268,412) distal to the IGH breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,296,741-105,909,611) proximal to the IGH breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC IGH Probe map (not to scale)



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the IGH locus, two green/orange fusion signals appear.
Aberrant situation: One IGH locus affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Example of an aberrant signal pattern: Burkitt lymphoma tissue section with translocation affecting the 14q32.33 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	14	CE/IVD	5 (50 μl)	Z-2110-50
			20 (200 μl)	Z-2110-200

ZytoLight® SPEC Angelman Dual Color Probe

Intended Purpose
The *ZytoLight* SPEC Angelman Dual Color Probe (PL273) is intended to be used for the qualitative detection of human UBE3A gene deletions as well as the detection of the human PML gene in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the above mentioned product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

Interpretation of the results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist/human geneticist.

- Probe Description**
The *ZytoLight* SPEC Angelmann Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 15q11.2-q12* (chr15:25,550,004-25,748,900) harboring the UBE3A gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10.0 ng/μl), which target sequences mapping in 15q24.1* (chr15:73,910,690-74,699,298) harboring the PML gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC UBE3A Probe map; Bottom: SPEC PML Probe map (not to scale); Right: Ideogram of chromosome 15 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the UBE3A gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletion affecting the UBE3A gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the UBE3A gene region might result in normal signal pattern with orange signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.

Lymphocytes and metaphase chromosomes from an Angelman syndrome case showing deletion affecting the chromosomal region 15q11.2-q12 as indicated by the loss of one orange signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	15	CE/IVD	5 (50 μl)	Z-2319-50

ZytoLight® SPEC Prader-Willi Dual Color Probe

Intended Purpose
The ZytoLight SPEC Prader-Willi Dual Color Probe is designed to detect deletions affecting the chromosomal region 15q11.2 harboring the SNRPN (small nuclear ribonucleoprotein polypeptide N, a.k.a. PWCR) gene.

The Prader-Willi syndrome (PWS) is a sporadic genetic disorder caused by genomic errors that inactivate paternally-inherited genes in the PWS critical region on chromosome 15q11-q13. The absence of expression of one or more of these genes contributes to different phenotypes of PWS. There are three main genetic causes: paternal 5-7 Mb deletion of the 15q11-q13 region, maternal uniparental disomy 15, or imprinting defects in the PWS critical region.

The SNRPN gene is located within the PWS region and has an important regulatory role over the imprinted genes located in chromosome 15. The estimated prevalence of the disease ranges between 1/15,000 and 1/30,000 newborns. PWS patients clinically display a characteristic pattern of symptoms including hypotonia with poor suck and poor weight gain in infancy, mild mental retardation, hypogonadism, growth hormone insufficiency causing short stature, early childhood-onset hyperphagia and obesity, characteristic appearance, and behavioral and sometimes psychiatric disturbance.

Early diagnosis offers the opportunity to significantly improve health and quality of life of people with PWS. FISH analysis can be performed to detect deletions within the PWS critical region and can help to confirm PWS diagnosis in patients with clinical features characteristic for PWS.

References
Cassidy AB & Driscoll DJ (2009) Eur J Hum Genet 17: 3-13.
Costa RA, *et al.* (2019) Front Endocrinol (Lausanne) 10: 864.

- Probe Description**
The ZytoLight SPEC Prader-Willi Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 15q11.2* (chr15:25,097,811-25,270,969) harboring the SNRPN gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10.0 ng/μl), which target sequences mapping in 15q24.1* (chr15:73,910,690-74,699,298) harboring the PML gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC SNRPN Probe map; Bottom: SPEC PML Probe map (not to scale)




Ideogram of chromosome 15 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the SNRPN gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletion affecting the SNRPN gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the SNRPN gene region might result in normal signal pattern with orange signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.



Lymphocytes and metaphase chromosomes from a Prader-Willi syndrome case showing deletion affecting the chromosomal region 15q11.2 as indicated by the loss of one orange signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	15	CE/IVD	5 (50 μl)	Z-2318-50

ZytoLight® SPEC NUTM1 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC NUTM1 Dual Color Break Apart Probe (PL166) is intended to be used for the qualitative detection of translocations involving the human NUTM1 gene at 15q14 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC NUTM1 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 15q14* (chr15:33,999,791-34,587,649) proximal to the NUTM1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 15q14* (chr15:34,873,659-35,326,986) distal to the NUTM1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Fig. 1: SPEC NUTM1 Probe map (not to scale); Right: Ideogram of chromosome 15 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the NUTM1 gene region, two green/orange fusion signals appear.
Aberrant situation: One NUTM1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: NMC tissue section with translocation of the NUTM1 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	15	CE/IVD	20 (200 μl)	Z-2208-200

ZytoLight® SPEC NTRK3 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC NTRK3 Dual Color Break Apart Probe (PL164) is intended to be used for the qualitative detection of translocations involving the human NTRK3 gene at 15q25.3 in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC NTRK3 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 15q25.3-q26.1* (chr15:88,825,346-89,475,889) distal to the NTRK3 breakpoint region (see Fig.1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 15q25.3* (chr15:87,976,717-88,471,002) proximal to the NTRK3 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC NTRK3 Probe map (not to scale); Right: Ideogram of chromosome 15 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the NTRK3 gene region, two green/orange fusion signals appear.
Aberrant situation: One NTRK3 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated orange signals are the result of deletions distal to the NTRK3 breakpoint region or are due to unbalanced translocations affecting this chromosomal region.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: Secretory breast carcinoma tissue section with translocation affecting the 15q25.3-q26.1 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	15	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2206-50 Z-2206-200

ZytoLight® SPEC PML/RARA Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe is designed to detect the translocation t(15;17)(q24;q21.2) affecting the PML gene in the chromosomal region 15q24.1 and the RARA locus in 17q21.2.

Translocations involving the PML (promyelocytic leukemia, a.k.a. MYL) gene and the RARA (retinoic acid receptor alpha, a.k.a. RARα) gene are considered to be characteristic for acute promyelocytic leukemia (APL), a subtype of acute myeloid leukemia.

Various fusion partners of RARA have been identified, however, in 95% of all APL cases, rearrangements involving the PML gene are detectable. This translocation t(15;17)(q24;q21) leads to a gene fusion of the PML and the RARA gene. The fusion is supposed to play a fundamental role in induction, development, and progression of APL.

Since the PML/RARA fusion accounts for the response of these neoplasms to all-trans retinoic acid (ATRA) therapy and other conventional chemotherapy it is important to accurately distinguish between t(15;17) translocations and translocations involving other partners of RARA.

References
Abe S, *et al.* (2008) Cancer Genet and Cytogenet 184: 44-7.
Brockmann SR, *et al.* (2003) Cancer Genet and Cytogenet 145: 144-51.
Reiter A, *et al.* (2004) Acta Hematol 112: 55-67.
Sanz MA, *et al.* (2009) Blood 113: 1875-91.

Probe Description
The ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 17q12-q21.2* (chr17:37,953,503-38,818,030) (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 15q24.1* (chr15:73,910,690-74,699,298) (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC PML Probe map; Bottom: SPEC RARA Probe map (not to scale)



Ideograms of chromosome 15 and 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Bone marrow biopsy section with translocation affecting the PML/RARA loci as indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	15, 17	CE/IVD	5 (50 μl)	Z-2113-50
			20 (200 μl)	Z-2113-200

ZytoLight® SPEC CREBBP Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC CREBBP Dual Color Break Apart Probe is designed for the detection of translocations involving the chromosomal region 16p13.3 harboring the CREBBP (CREB binding protein, a.k.a. CBP, RTS) gene. The CREBBP protein regulates transcription by means of histone acetyltransferase activity and by binding to several proteins with key cell cycle functions, such as p53 and NFκB.

Rearrangements of the CREBBP gene have been observed in several hematologic malignancies. Three different fusion partners have been described so far.

KMT2A (a.k.a. MLL) is fused to CREBBP in therapy-related acute myeloid (AML) or lymphoid leukemia (ALL) and myelodysplastic syndrome (MDS) with t(11;16)(q23.3;p13.3). The translocation t(10;16)(q22.2;p13.3) was reported in some AML cases and fuses KAT6B (a.k.a. MORF) to CREBBP. CREBBP is also rearranged with KAT6A (a.k.a. MOZ) in de novo and therapy-related AML with t(8;16)(p11.2;p13.3) after treatment with topoisomerase II inhibitors. This rearrangement is associated with an infrequent but well-defined type of AML that has characteristic morphocytochemical features. The prognosis is usually extremely poor, with a median survival of two months. The KAT6A/CREBBP AML tends to develop within two years of adjuvant chemotherapy, especially in former breast cancer patients.

Thus, FISH analysis for the detection of CREBBP translocation may serve as a diagnostic tool to identify cases with hematologic malignancies with an aggressive presentation.

References
Borrow J, *et al.* (1996) Nat Genet 14: 33-41.
Camós M, *et al.* (2006) Cancer Res 66: 6947-54.
Gupta A, *et al.* (2014) Case Rep Oncol Med 2014: 361748.
Rozman M, *et al.* (2004) Genes Chromosomes Cancer 40: 140-5.
Taki T, *et al.* (1997) Blood 89: 3945-50.
Vizmanos JL, *et al.* (2003) Genes Chromosomes Cancer 36: 402-5.

Probe Description
The ZytoLight SPEC CREBBP Dual Color Break Apart Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 16p13.3* (chr16:3,978,217-4,521,684) proximal to the CREBBP breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 16p13.3* (chr16:3,287,067-3,762,188) distal to the CREBBP breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC CREBBP Probe map (not to scale); Right: Ideogram of chromosome 16 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the CREBBP gene region, two green/orange fusion signals appear.
Aberrant situation: One CREBBP gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
	16	CE/IVD	5 (50 μl)	Z-2267-50

ZytoLight® SPEC FUS Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC FUS Dual Color Break Apart Probe (PL87) is intended to be used for the qualitative detection of translocations involving the human FUS gene at 16p11.2 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.


Probe Description
The ZytoLight SPEC FUS Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 16p11.2* (chr16:30,383,304-31,007,836) distal to the FUS breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 16p11.2* (chr16:31,213,259-31,927,155) proximal to the FUS breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC FUS Probe map (not to scale)



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the FUS gene region, two green/orange fusion signals appear.
Aberrant situation: One FUS gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Example of an aberrant signal pattern: Myxoid liposarcoma tissue section with translocation affecting the 16p11.2 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	16	CE/IVD	5 (50 μl)	Z-2130-50

ZytoLight® SPEC MAF/IGH Dual Color Dual Fusion Probe

Intended Purpose

The ZytoLight SPEC MAF/IGH Dual Color Dual Fusion Probe is designed to detect the translocations affecting the MAF gene in the chromosomal region 16q23.2 and the IGH locus in 14q32.33. The translocation t(14;16)(q32.3;q23) is frequently found in multiple myeloma (MM).

MM is a malignant post-germinal center tumor of somatically-mutated, isotype-switched plasma cells that accumulate in the bone marrow. It is often preceded by a premalignant state known as monoclonal gammopathy of undetermined significance (MGUS).

Five recurrent primary translocations involving the immunoglobulin heavy locus (IGH) have been identified in 40% of MGUS and MM tumors. They include t(11;14)(q13.3;q32.3), t(6;14)(p21.1;q32.3), t(4;14)(p16.3;q32.3), t(14;16)(q32.3;q23), and t(14;20)(q32.3;q12), which involve the genes CCND1, CCND3, FGFR3 and NSD2, MAF, and MAFB, respectively. All of these translocations lead to the dysregulation and overexpression of the target genes as a consequence of their juxtaposition to regulatory sequences of the IGH locus. t(14;16) occurs in approximately 5% of MM patients and is associated with a more aggressive clinical outcome. The 16q23 breakpoints have been found to be scattered 550-1280kb centromerically to the MAF gene within the WWOX gene.

Hence, detection of t(14;16) by FISH represents a useful prognostic tool and may aid in therapeutic decision making in MM.

References

Chesi M, *et al.* (1998) Blood 91: 4457-63.
Fabris S, *et al.* (2005) Genes Chromosomes Cancer 42: 117-27.
Fonseca R, *et al.* (2009) Leukemia 23: 2210-21.
Gabrea A, *et al.* (2006) DNA Repair (Amst) 5: 1225-33.

Probe Description

- The ZytoLight SPEC MAF/IGH Dual Color Dual Fusion Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 16q23.1-q23.2* (chr16:78,089,697-79,657,277) harboring the MAF gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH locus (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC MAF Probe map; Bottom: SPEC IGH Probe map (not to scale); Right: Ideograms of chromosome 16 and 14 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the respective bands, two separate green and orange signals appear.

Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Bone marrow CD138+ cells with translocation affecting the MAF/IGH loci as indicated by two orange/green fusion signals, a single orange, and a separate green signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	14, 16	CE/IVD	5 (50 μl)	Z-2270-50

Kindly provided by Prof. Dr. Oskar A. Haas, Vienna, Austria.

ZytoLight® SPEC MAFB/IGH Dual Color Dual Fusion Probe

Intended Purpose

The ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe is designed to detect the translocations affecting the MAFB gene in the chromosomal region 20q12 and the IGH locus in 14q32.33. The translocation t(14;20)(q32.3;q12) is frequently found in multiple myeloma (MM).

MM is a low proliferative, malignant post-germinal center tumor of somatically mutated, isotype-switched plasma cells that accumulate in the bone marrow. It is often preceded by a premalignant state known as monoclonal gammopathy of undetermined significance (MGUS).

Five recurrent primary translocations involving the immunoglobulin heavy locus (IGH) have been identified in 40% of MGUS and MM tumors. They include t(11;14)(q13.3;q32.3), t(6;14)(p21.1;q32.3), t(4;14)(p16.3;q32.3), t(14;16)(q32.3;q23), and t(14;20)(q32.3;q12), which involve the genes CCND1, CCND3, FGFR3 and NSD2, MAF, and MAFB, respectively. All of these translocations lead to the dysregulation and overexpression of the target genes as a consequence of their juxtaposition to regulatory sequences of the IGH locus.

The t(14;20) occurs in approximately 1-2% of MM patients and is associated with an adverse prognosis. Thus, currently, detection of t(14;20) by FISH is a reliable prognostic tool and may sustain therapeutic decision making in MM.

References

Boersma-Vreugdenhil GR, *et al.* (2004) Br J Haematol 126: 355-63.
Chesi M, *et al.* (1998) Blood 92: 4457-63.
Fabris S, *et al.* (2005) Genes Chromosomes Cancer 42: 117-27.
Fonseca R, *et al.* (2009) Leukemia 23: 2210-21.
Gabrea A, *et al.* (2006) DNA Repair (Amst) 5: 1225-33.
Hanamura I, *et al.* (2001) Jpn N Cancer Res 92: 638-44.

Probe Description

- The ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 20q12* (chr20:37,782,012-39,385,613) harboring the MAFB gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH locus (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC MAFB Probe map; Bottom: SPEC IGH Probe map (not to scale); Right: Ideograms of chromosome 20 and 14 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the respective loci, two separate green and orange signals appear.

Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.



Bone marrow CD138+ cells with translocation affecting the MAFB/IGH loci as indicated by two orange/green fusion signals, a single orange, and a separate green signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	14, 20	CE/IVD	5 (50 μl)	Z-2271-50

Kindly provided by Prof. Dr. Oskar A. Haas, Vienna, Austria.

ZytoLight® SPEC CBFB Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC CBFB Dual Color Break Apart Probe (PL165) is intended to be used for the qualitative detection of translocations involving the human CBFB gene at 16q22.1 in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC CBFB Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 16q22.1* (chr16:67,161,347-67,605,304) distal to the CBFB breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 16q22.1* (chr16:66,882,262-67,102,895) proximal to the CBFB breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC CBFB Probe map (not to scale); Right: Ideogram of chromosome 16 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the CBFB gene region, two green/orange fusion signals appear.
Aberrant situation: One CBFB gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	16	CE/IVD	5 (50 μl)	Z-2207-50

ZytoLight® SPEC TP53/17q22 Dual Color Probe

Intended Purpose
The ZytoLight SPEC TP53/17q22 Dual Color Probe (PL156) is intended to be used for the qualitative detection of deletions involving the human TP53 gene as well as the detection of gains of chromosome 17q22 specific sequences in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC TP53/17q22 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 17p13.1* (chr17:7,495,749-7,663,022) harboring the TP53 gene region (see Fig.1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 17q22* (chr17:56,124,338-56,594,220) (see Fig.1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC TP53 Probe map; Bottom: SPEC 17q22 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the TP53 gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletion affecting the TP53 gene region, a reduced number of orange signals will be detected. Deletions affecting only parts of the TP53 gene region might result in a normal signal pattern with orange signals of reduced size. A gain of 17q involving the 17q22 region will result in an increased number of green signals. Isochromosome 17q is indicated by three green signals and one orange signal.



Example of an aberrant signal pattern: SPEC TP53/17q22 Dual Color Probe hybridized to bone marrow tissue section with deletion of the TP53 gene as indicated by one orange signal and two green signals in each nucleus.



Example of an aberrant signal pattern: SPEC TP53/17q22 Dual Color Probe hybridized to a bone marrow smear with isochromosome 17q as indicated by three green signals and one orange signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	5 (50 μl)	Z-2198-50

ZytoLight® SPEC TP53/CEN 17 Dual Color Probe

Intended Purpose
The ZytoLight SPEC TP53/CEN 17 Dual Color Probe (PL109) is intended to be used for the qualitative detection of deletions involving the human TP53 gene as well as the detection of chromosome 17 alpha satellites in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC TP53/CEN 17 Dual Color Probe is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 17p13.1* (chr17:7,495,749-7,663,022) harboring the TP53 gene region (see Fig.1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC TP53 Probe map (not to scale)




Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the TP53 gene region, two orange signals and two green signals appear.
Aberrant situation: In a cell with deletion affecting the TP53 gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the TP53 gene region might result in a normal signal pattern with orange signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: SPEC TP53/CEN 17 Dual Color Probe hybridized to bone marrow tissue section with deletion of the TP53 gene as indicated by one orange signal and two green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	5 (50 μl)	Z-2153-50
			20 (200 μl)	Z-2153-200

ZytoLight® SPEC USP6 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC USP6 Dual Color Break Apart Probe (PL107) is intended to be used for the qualitative detection of translocations involving the human USP6 gene at 17p13.2 in formalin-fixed, paraffin-embedded specimens, such as aneurysmal bone cyst (ABC) or nodular fasciitis (NF), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ABC or NF and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC USP6 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 17p13.2* (chr17:4,489,889-5,017,582) distal to the USP6 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 17p13.2* (chr17:5,087,046-5,361,104) proximal to the USP6 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC USP6 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the USP6 gene region, two green/orange fusion signals appear.
Aberrant situation: One USP6 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to aneurysmal bone cyst tissue section with polysomy of chromosome 17 but without translocation affecting the 17p13.2 locus as indicated by multiple orange/green fusion signals per nucleus.



Aneurysmal bone cyst tissue section with translocation affecting the 17p13.2 locus as indicated by one orange/green fusion (non-rearranged) signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	5 (50 μl)	Z-2151-50
			20 (200 μl)	Z-2151-200

ZytoLight® SPEC YWHAE Dual Color Break Apart Probe

Intended Purpose
The *ZytoLight* SPEC YWHAE Dual Color Break Apart Probe (PL134) is intended to be used for the qualitative detection of translocations involving the human YWHAE gene at 17p13.3 in formalin-fixed, paraffin-embedded specimens, such as endometrial stromal sarcoma (ESS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ESS and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The *ZytoLight* SPEC YWHAE Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 17p13.3* (chr17:1,339,752-1,716,668) proximal to the YWHAE breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 17p13.3* (chr17:791,171-1,124,746) distal to the YWHAE breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC YWHAE Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the YWHAE gene region, two green/orange fusion signals appear.
Aberrant situation: One YWHAE gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Endometrial stromal sarcoma tissue section with translocation affecting the YWHAE gene as indicated by one non-rearranged orange/green fusion signal, one orange, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	5 (50 μl)	Z-2175-50

ZytoLight® SPEC ERBB2/CEN 17 Dual Color Probe

Intended Purpose
The *ZytoLight* SPEC ERBB2/CEN 17 Dual Color Probe (PL8) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene as well as the detection of chromosome 17 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer and gastric/gastroesophageal junction cancer, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and gastric/gastroesophageal junction cancer and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The *ZytoLight* SPEC ERBB2/CEN 17 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 17q12-q21.1* (chr17:37,572,531-38,181,308) harboring the ERBB2 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.5 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC ERBB2 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the ERBB2 gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the ERBB2 gene region an increased number of green signals or green signal clusters will be observed.



Normal interphase cells, ERBB2 (green), CEN 17 (orange).



Breast carcinoma tissue section, ERBB2 gene cluster (green), CEN 17 (orange).

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	5 (50 μl)	Z-2015-50
			20 (200 μl)	Z-2015-200

The product is also available as a kit in two different sizes.

ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit	Tests	Order No.
Incl: Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 210 ml; Probe, 0.05 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml	5	Z-2020-5
Incl: Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 560 ml; Probe, 0.2 ml; 25x Wash Buffer A, 100 ml; DAPI/DuraTect-Solution, 0.8 ml	20	Z-2020-20

ZytoLight® CEN 17/SPEC ERBB2 Dual Color Probe

Intended Purpose
The *ZytoLight* CEN 17/SPEC ERBB2 Dual Color Probe (PL36) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene as well as the detection of chromosome 17 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer and gastric/gastroesophageal junction cancer, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and gastric/gastroesophageal junction cancer and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The *ZytoLight* CEN 17/SPEC ERBB2 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 17q12-q21.1* (chr17:37,572,531-38,181,308) harboring the ERBB2 gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC ERBB2 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the ERBB2 gene region, two green and two orange signals appear.
Aberrant situation: In cells with an amplification of the ERBB2 gene region, an increased number of orange signals or orange signal clusters will be observed.



Normal interphase cells, ERBB2 (orange), CEN 17 (green).



Breast carcinoma tissue section, ERBB2 gene cluster (orange), CEN 17 (green).

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	5 (50 μl)	Z-2077-50
			20 (200 μl)	Z-2077-200

ZytoLight® SPEC ERBB2/D17S122 Dual Color Probe

Intended Purpose
The *ZytoLight* SPEC ERBB2/D17S122 Dual Color Probe (PL148) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene as well as the detection of the D17S122 locus in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the *ZytoLight* FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The *ZytoLight* SPEC ERBB2/D17S122 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 17q12-q21.1* (chr17:37,572,531-38,181,308) harboring the ERBB2 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 17p12* (chr17:14,954,785-15,434,017) harboring the D17S122 locus (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC ERBB2 Probe map; Bottom: SPEC D17S122 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the ERBB2 gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the ERBB2 gene region an increased number of green signals or green signal clusters will be observed.



Normal interphase cells, ERBB2 (green), D17S122 (orange).



Example of an aberrant signal pattern: Breast carcinoma tissue section, ERBB2 gene cluster (green), D17S122 (orange).

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	5 (50 μl)	Z-2190-50
			20 (200 μl)	Z-2190-200

ZytoLight® SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe

Intended Purpose

The ZytoLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe (PL52) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene and the human TOP2A gene as well as the detection of chromosome 17 alpha satellites in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 17q12-q21.1* (chr17:37,572,531-38,181,308) harboring the ERBB2 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 17q21.1-q21.2* (chr17:38,323,741-38,818,030) harboring the TOP2A gene region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12.0 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ERBB2/TOP2A Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without an amplification involving the respective gene regions, two green, two orange and two blue signals appear.

Aberrant situation: In cells with an amplification of the ERBB2 gene region an increased number of green signals or green signal clusters will be observed. In cells with an amplification of the TOP2A gene region an increased number of orange signals or orange signal clusters will be observed. Deletions affecting the TOP2A gene region result in a reduced number of orange signals.



Hybridized to normal interphase cells as indicated by two green, two orange, and two blue signals per nucleus.



Example of an aberrant signal pattern: Breast cancer tissue section with two copies of chromosome 17 (blue) and TOP2A (orange) and ERBB2 gene clusters (green) in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
  	17	CE/IVD	5 (50 μl)	Z-2093-50
			20 (200 μl)	Z-2093-200

ZytoLight® SPEC COL1A1 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC COL1A1 Dual Color Break Apart Probe (PL78) is intended to be used for the qualitative detection of translocations involving the human COL1A1 gene at 17q21.33 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC COL1A1 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 17q21.33* (chr17:47,669,218-48,223,465) distal to the COL1A1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 17q21.33* (chr17:48,347,800-48,744,021) distal to the COL1A1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC COL1A1 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the COL1A1 gene region, two green/orange fusion signals appear.

Aberrant situation: One COL1A1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Example of an aberrant signal pattern: DFSP tissue section with translocation affecting the 17q21.33 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.



Example of an aberrant signal pattern: DFSP tissue section with amplification affecting the 17q21-qter and 22q10-q13.1 sequences probably due to a COL1A1-PDGFB fusion product on the ring chromosome.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17	CE/IVD	20 (200 μl)	Z-2121-200

Image kindly provided by Dr. Schildhaus, Essen, Germany.

ZytoLight® SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe

Intended Purpose

The ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe (PL73) is intended to be used for the qualitative detection of the translocation t(17;22)(q21.3;q13.1) involving the human COL1A1 and PDGFB genes in formalin-fixed, paraffin-embedded specimens, such as dermatofibrosarcoma protuberans (DFSP), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of DFSP and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 17q21.33* (chr17:47,820,343-48,744,021) harboring the COL1A1 gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 22q13.1* (chr22:38,928,973-40,267,687) harboring the PDGFB gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC COL1A1 Probe map; Bottom: SPEC PDGFB Probe map (not to scale); Right: Ideograms of chromosome 17 and 22 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.

Aberrant situation: A reciprocal translocation is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals. A signal pattern showing one fusion signal and variable numbers of separate green and orange signals can occur due to unbalanced translocations. Multiple copies of orange/green fusion signals indicate the presence of supernumerary ring chromosomes comprising low-level amplified COL1A1-PDGFB fusion sequences.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



DFSP tissue section with translocation affecting the COL1A1/PDGFB loci as indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	17, 22	CE/IVD	5 (50 μl)	Z-2116-50
			20 (200 μl)	Z-2116-200

ZytoLight® SPEC SS18 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC SS18 Dual Color Break Apart Probe (PL56) is intended to be used for the qualitative detection of translocations involving the human SS18 gene at 18q11.2 in formalin-fixed, paraffin-embedded specimens, such as synovial sarcoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of synovial sarcoma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC SS18 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 18q11.2* (chr18:23,109,942-23,466,217) proximal to the SS18 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 18q11.2* (chr18:23,772,255-24,137,169) distal to the SS18 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC SS18 Probe map (not to scale); Right: Ideogram of chromosome 18 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the SS18 gene region, two green/orange fusion signals appear.

Aberrant situation: One SS18 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Synovial sarcoma tissue section with translocation affecting the 18q11.2 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	18	CE/IVD	5 (50 μl)	Z-2097-50
			20 (200 μl)	Z-2097-200

ZytoLight® SPEC SS18/SSX1 TriCheck™ Probe

Intended Purpose

The ZytoLight SPEC SS18/SSX1 TriCheck™ Probe is designed to detect translocations involving the chromosomal region 18q11.2 harboring the SS18 (SS18, nBAF chromatin remodeling complex subunit, a.k.a. SYT) gene and the chromosomal region Xp11.23 harboring the SSX1 (SSX family member 1) gene.

Synovial sarcoma is characterized by the t(X;18) found in more than 95% of these tumors and juxtaposing the SS18 gene in 18q11.2 either next to the SSX1 or the SSX2 gene, or very rarely to the SSX4 locus.

Synovial sarcoma is one of the most common soft tissue tumors that typically occurs in the extremities of young adults with greater prevalence in males, even though, the occurrence of synovial sarcoma has also been described in a wide variety of anatomical locations and in all ages.

In combination with histopathological diagnosis, detection of SS18 rearrangements via FISH is a valuable tool to confirm the diagnosis of synovial sarcoma. Moreover, patients with SS18-SSX1 fusions were shown to have a higher risk of developing metastases compared to those with SS18-SSX2 fusions. Hence, detection of the SS18 fusion gene variant by FISH may also be of prognostic significance.

References

Amary MF, et al. (2007) Mod Pathol 20: 482-96.
Clark J, et al. (1994) Nat Genet 7: 502-8.
Kawai A, et al. (1998) N Engl J Med 338: 153-60.
Panagopoulos I, et al. (2001) Genes Chromosomes Cancer 31: 362-72.
Surace C, et al. (2004) Lab Invest 84: 1185-92.
Torres L, et al. (2008) Cancer Genet Cytogenet 187: 45-9.

Probe Description

The ZytoLight SPEC SS18/SSX1 TriCheck Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 18q11.2* (chr18:23,109,942-23,466,217) proximal to the SS18 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 18q11.2* (chr18:23,772,255-24,137,169) distal to the SS18 breakpoint region (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/μl), which target sequences mapping in Xp11.23* (chrX:48,265,856-48,792,674) proximal to the SSX1 breakpoint region (see Fig. 1).

* according to Human Genome Assembly GRCh37/hg19


Results

Normal situation: In interphases of normal female cells or female cells without SS18-SSX1 rearrangement, two green/orange fusion signals and two blue signals appear. In interphases of normal male cells or male cells without SS18-SSX1 rearrangement, two green/orange fusion signals and one blue signal appear.

Aberrant situation: An SS18-SSX1 or an SS18-SSX4 fusion is indicated by one separate orange signal co-localizing with one blue signal and one separate green signal. An SS18-SSX2 fusion is indicated by one separate green signal, one separate orange signal, and a blue signal in close proximity of the separated green signal. An SS18 translocation without involvement of SSX1, SSX2, or SSX4 is indicated by the split of one green/orange fusion signal without co-localization of the separated orange or green signal with one blue signal.



Male synovial sarcoma tissue section with SS18-SSX1 or SS18-SSX4 fusion as indicated by orange/blue fusion signals.



Female synovial sarcoma tissue section with SS18-SSX2 fusion as indicated by green/blue fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
  	18, X	CE/IVD	5 (50 μl)	Z-2184-50

ZytoLight® SPEC MALT1 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC MALT1 Dual Color Break Apart Probe (PL154) is intended to be used for the qualitative detection of translocations involving the human MALT1 gene at 18q21.32 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC MALT1 Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 18q21.31-q21.32* (chr18:55,690,725-56,330,582) proximal to the MALT1 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 18q21.32* (chr18:56,427,386-56,859,755) distal to the MALT1 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC MALT1 Probe map (not to scale); Right: Ideogram of chromosome 18 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the MALT1 gene region, two green/orange fusion signals appear.

Aberrant situation: One MALT1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Example of an aberrant signal pattern: Lymphoma tissue section with translocation of the MALT1 gene as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	18	CE/IVD	5 (50 μl)	Z-2196-50
			20 (200 μl)	Z-2196-200

ZytoLight® SPEC MALT1/IGH Dual Color Dual Fusion Probe

Intended Purpose

The ZytoLight SPEC MALT1/IGH Dual Color Dual Fusion Probe (PL279) is intended to be used for the qualitative detection of the translocation t(14;18) (q32.3;q21.3) involving the human MALT1 gene at 18q21.32 and the human IGH locus at 14q32.33 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoLight SPEC MALT1/IGH Dual Color Dual Fusion Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 18q21.31-q21.32* (18:56,021,766-56,724,408) harboring the MALT1 gene region (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC MALT1 Probe map; Bottom: SPEC IGH Probe map (not to scale); Right: Ideograms of chromosome 18 and 14 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate orange and green signals appear.

Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



Example of an aberrant signal pattern: SPEC MALT1/IGH Dual Color Dual Fusion Probe hybridized to MALT lymphoma tissue section with t(14;18) as indicated by one separate orange signal, one separate green signal and two orange/green fusion signals indicating the MALT1/IGH translocation.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	14, 18	CE/IVD	5 (50 μl)	Z-2325-50

ZytoLight® SPEC BCL2 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC BCL2 Dual Color Break Apart Probe (PL150) is intended to be used for the qualitative detection of translocations involving the human BCL2 gene at 18q21.33 in formalin-fixed, paraffin-embedded specimens, such as B-cell lymphoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of B-cell lymphoma and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC BCL2 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 18q21.33* (chr18:60,046,152-60,589,273) proximal to the BCL2 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 18q21.33-q22.1* (chr18:60,994,528-61,658,503) distal to the BCL2 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC BCL2 Probe map (not to scale); Right: Ideogram of chromosome 18 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the BCL2 gene region, two green/orange fusion signals appear.
Aberrant situation: One BCL2 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.
Neck lymph node tissue section with translocation of the BCL2 gene as indicated by two non-rearranged orange/green fusion signals, one orange and one separate green signal.

ZytoLight® SPEC BCL2/IGH Dual Color Dual Fusion Probe

Intended Purpose
The ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe (PL71) is intended to be used for the qualitative detection of the translocation t(14;18) (q32.3;q21.3) involving the human IGH and BCL2 genes in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~12 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~6 ng/μl), which target sequences mapping in 18q21.33* (chr18:60,180,078-61,507,691) harboring the BCL2 gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC BCL2 Probe map; Bottom: IGH Probe map (not to scale); Right: Ideograms of chromosome 14 and 18 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the respective gene regions, two separate green and orange signals appear.
Aberrant situation: A gene fusion is indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.
Example of an aberrant signal pattern: Bone marrow biopsy section with translocation affecting the BCL2/IGH loci as indicated by one separate orange signal, one separate green signal, and two orange/green fusion signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	14, 18	CE/IVD	5 (50 μl)	Z-2114-50
			20 (200 μl)	Z-2114-200

ZytoLight® SPEC CIC Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC CIC Dual Color Break Apart Probe (PL240) is intended to be used for the qualitative detection of translocations involving the human CIC gene at 19q13.2 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC CIC Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 19q13.2* (chr19:42,835,047-43,374,575) distal to the CIC breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 19q13.2* (chr19:41,980,301-42,751,339) proximal to the CIC breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC CIC Probe map (not to scale); Right: Ideogram of chromosome 19 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the CIC gene region, two green/orange fusion signals appear.
Aberrant situation: One CIC gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.
Example of an aberrant signal pattern: Undifferentiated roand cell sarcoma 'Ewing-like' tissue section with translocation of the 19q13.2 locus as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	19	CE/IVD	5 (50 μl)	Z-2285-50

ZytoLight® SPEC C19MC/19p13 Dual Color Probe

Intended Purpose
The ZytoLight SPEC C19MC/19p13 Dual Color Probe (PL230) is intended to be used for the qualitative detection of amplifications involving the C19MC locus as well as the detection of chromosome 19p13 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC C19MC/19p13 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 19q13.42* (chr19:54,156,239-54,768,983) harboring the C19MC locus (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 19p13.3* (chr19:658,555-1,144,465) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC C19MC Probe map; Bottom: SPEC 19p13 Probe map (not to scale); Right: Ideogram of chromosome 19 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the C19MC locus, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the C19MC locus, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two orange and two green signals per nucleus and to metaphase chromosomes of a normal cell.
Example of an aberrant signal pattern: Primitive neuroectodermal tumor tissue section with amplification of the C19MC miRNA cluster.

Labels	Chromosome	Reg. Status	Tests	Order No.
	19	CE/IVD	5 (50 μl)	Z-2274-50

Specimen kindly provided by Hannu Haapasalo, MD, PhD, Fimlab Laboratories, Finland.

ZytoLight® SPEC PTPRT/20q11 Dual Color Probe

Intended Purpose
The ZytoLight SPEC PTPRT/20q11 Dual Color Probe (PL171) is intended to be used for the detection of deletions involving the human PTPRT gene as well as the detection of chromosome 20q11 specific sequences in cytologic specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC PTPRT/20q11 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 20q12* (chr20:40,807,040-41,419,634) harboring the PTPRT gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 20q11.21* (chr20:30,301,019-30,719,110) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC PTPRT Probe map; Bottom: SPEC 20q11 Probe map (not to scale)



Results
Normal situation: In interphases of normal cells or cells without a deletion involving the PTPRT gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletion affecting the PTPRT gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the PTPRT gene region might result in a normal signal pattern with green signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.

Example of an aberrant signal pattern: Lymphocytes of a myelodysplastic syndrome showing a 20q deletion indicated by one single green and two orange signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	20	CE/IVD	5 (50 μl)	Z-2213-50

Material kindly provided from Dr. Saurabh Bhattacharya, Lal PathLabs, India.

ZytoLight® SPEC ERG Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC ERG Dual Color Break Apart Probe (PL95) is intended to be used for the qualitative detection of translocations involving the human ERG gene at 21q22.2 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC ERG Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 21q22.2* (chr21:40,078,039-40,850,582) distal to the ERG breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 21q22.13-q22.2* (chr21:39,171,790-39,733,849) proximal to the ERG breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ERG Probe map (not to scale); Right: Ideogram of chromosome 21 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the ERG gene region, two orange/green fusion signals appear.
Aberrant situation: One ERG gene region affected by a deletion resulting in the TMPRSS2-ERG fusion is indicated by the loss of one green signal. Deletions affecting only parts of the ERG gene region might result in a normal signal pattern with green signals of reduced size. A signal pattern consisting of one orange/green fusion signal, a separate green, and a separate orange signal indicates an ERG translocation without involvement of TMPRSS2 (e.g., SLC45A3-ERG).



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Example of an aberrant signal pattern: Prostate cancer tissue section with interstitial deletion of the chromosomal region 21q22.2 resulting in the TMPRSS2-ERG fusion as indicated by the loss of one green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	21	CE/IVD	20 (200 μl)	Z-2138-200

ZytoLight® SPEC ERG/TMPRSS2 TriCheck™ Probe

Intended Purpose
The ZytoLight SPEC ERG/TMPRSS2 TriCheck Probe (PL92) is intended to be used for the qualitative detection of rearrangements involving the human ERG gene at 21q22.2 and the human TMPRSS2 gene at 21q22.3 in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC ERG/TMPRSS2 TriCheck Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 21q22.2* (chr21:40,078,039-40,850,582) distal to the ERG breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 21q22.13-q22.2* (chr21:39,171,790-39,733,849) proximal to the ERG breakpoint region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/μl), which target sequences mapping in 21q22.3* (chr21:43,301,411-44,195,531) distal to the TMPRSS2 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ERG/TMPRSS2 Probe map (not to scale); Right: Ideogram of chromosome 21 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without rearrangements affecting the respective gene regions, two orange/green fusion signals and two blue signals in close proximity of the respective fusion signals appear.
Aberrant situation: One 21q22 locus affected by a 21q22.2 deletion resulting in the TMPRSS2-ERG fusion is indicated by one separate orange signal co-localizing with one blue signal, and the loss of one green signal. An ERG translocation without involvement of TMPRSS2 is indicated by a separated orange signal and a blue signal co-localizing with the separate green signal. A non-ERG translocation affecting TMPRSS2 is indicated by a separated blue signal not co-localizing with the ERG fusion signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals and two blue signals in close proximity of the respective fusion signals.

Example of an aberrant signal pattern: Prostate cancer tissue section with one 21q22 locus affected by an interstitial deletion of the chromosomal region 21q22.2 resulting in the TMPRSS2-ERG fusion as indicated by one separate orange signal co-localizing with one blue signal, and the loss of one green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	21	CE/IVD	20 (200 μl)	Z-2135-200

ZytoLight® SPEC DiGeorge/Phelan McDermid Dual Color Probe

Intended Purpose

The ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe is designed to detect deletions affecting the chromosomal regions 22q11.21 harboring the HIRA (a.k.a. TUPLE1) gene and 22q13.33 harboring the SHANK3 (a.k.a. prosap2) gene, respectively.

The 22q11.2 deletion syndrome (22q11.2DS), also known as velocardiofacial syndrome (VCF) and DiGeorge syndrome, is a genetic disorder caused by hemizygous microdeletions on chromosome 22q11.2, with population prevalence of about 1 to 4,000 births. The characteristic phenotype of 22q11.2DS includes cardiac defects, velopharyngeal insufficiency, immune deficiency due to thymic aplasia, growth restriction, and deficits in cognitive abilities.

The 22q11.2 deletion usually occurs by meiotic non-allelic homologous recombination events between low copy repeats on chromosome 22q11.2 termed LCR22. There are eight LCR22s that span the 22q11.2 region termed LCR22A through LCR22H. The majority (90%) of 22q11.2DS patients show a recurrent 3 Mb deletion between LCR22A and LCR22D harboring the HIRA gene.

The 22q13.3 deletion syndrome (Phelan-McDermid syndrome) typically results from deletions of 100 kb to 9 Mb involving the distal long arm of chromosome 22. Almost all of these deletions include the gene SHANK3 that encodes a scaffold protein in the postsynaptic densities of excitatory synapses, connecting membrane-bound receptors to the actin cytoskeleton. This syndrome is characterized by neurological deficits, which include global developmental delay, moderate to severe intellectual impairment, absent or severely delayed speech, and neonatal hypotonia.

References

Burnside RD (2015) Cytogenet Genome Res 146: 89-99.
Morrow BE, *et al.* (2018) Am J Med Genet A 176: 2070-81.
Phelan K & McDermid HE (2012) Mol Syndromol 2: 186-201.
Scambler PJ, *et al.* (1991) Genomics 10: 201-6.
Watt JL, *et al.* (1985) J Med Genet 22: 283-7.

Probe Description


The ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 22q11.21* (chr22:19,191,435-19,506,869) harboring the HIRA gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 22q13.33* (chr22:50,924,766-51,188,029) harboring the SHANK3 gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19


Results

Normal situation: In interphases of normal cells or cells without a deletion involving the respective gene regions, two green and two orange signals appear.

Aberrant situation: In a cell with deletions affecting the HIRA gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the HIRA gene region might result in a normal signal pattern with green signals of reduced size. In a cell with deletions affecting the SHANK3 gene region, a reduced number of orange signals will be observed. Deletions affecting only parts of the SHANK3 gene region might result in a normal signal pattern with orange signals of reduced size



Lymphocytes and metaphase chromosomes from a DiGeorge syndrome case showing a HIRA deletion as indicated by the loss of one green signal. *



Lymphocytes and metaphase chromosomes from a Phelan-McDermid syndrome case showing a SHANK3 deletion as indicated by the loss of one orange signal. **

Labels	Chromosome	Reg. Status	Tests	Order No.
	22	CE/IVD	5 (50 μl)	Z-2299-50

*Kindly provided by Dr. Liehr, Jena, Germany.

**Kindly provided by Dr. Kazmierczak, Bremen, Germany.





Fig. 1: Top: SPEC HIRA Probe map; Bottom: SPEC SHANK3 Probe map (not to scale)



ZytoLight® SPEC DiGeorge Triple Color Probe

Intended Purpose

The ZytoLight SPEC DiGeorge/Phelan McDermid Dual Color Probe is designed to detect deletions affecting the chromosomal regions 22q11.21 harboring the HIRA (a.k.a. TUPLE1) gene and 22q13.33 harboring the SHANK3 (a.k.a. prosap2) gene, respectively.

The 22q11.2 deletion syndrome (22q11.2DS), also known as velocardiofacial syndrome (VCF) and DiGeorge syndrome, is a genetic disorder caused by hemizygous microdeletions on chromosome 22q11.2, with population prevalence of about 1 to 4,000 births. The characteristic phenotype of 22q11.2DS includes cardiac defects, velopharyngeal insufficiency, immune deficiency due to thymic aplasia, growth restriction, and deficits in cognitive abilities. The 22q11.2 deletion usually occurs by meiotic non-allelic homologous recombination events between low copy repeats on chromosome 22q11.2 termed LCR22. There are eight LCR22s that span the 22q11.2 region termed LCR22A through LCR22H. The majority (90%) of 22q11.2DS patients show a recurrent 3 Mb deletion between LCR22A and LCR22D harboring the HIRA gene.

The 22q13.3 deletion syndrome (Phelan-McDermid syndrome) typically results from deletions of 100 kb to 9 Mb involving the distal long arm of chromosome 22. Almost all of these deletions include the gene SHANK3 that encodes a scaffold protein in the postsynaptic densities of excitatory synapses, connecting membrane-bound receptors to the actin cytoskeleton. This syndrome is characterized by neurological deficits, which include global developmental delay, moderate to severe intellectual impairment, absent or severely delayed speech, and neonatal hypotonia.

References

Ben-Shachar S, *et al.* (2008) Am J Hum Genet 82: 214-21.
Burnside RD (2015) Cytogenet Genome Res 146: 89-99.
Michaelovsky E, *et al.* (2012) BMC Med Genet 13: 122.
Morrow BE, *et al.* (2018) Am J Med Genet A 176: 2070-81.
Scambler PJ, *et al.* (1991) Genomics 10: 201-6.

Probe Description

The ZytoLight SPEC DiGeorge Triple Color Probe is composed of:

- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37 ng/μl), which target sequences mapping in 22q11.21* (chr22:19,191,435-19,932,689) harboring the HIRA gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 22q11.21* (chr22:21,096,895-21,454,102) harboring the CRKL gene region (see Fig. 1).
 - ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 22q11.21-q11.22* (chr22:21,931,816-22,587,439) harboring the MAPK1 gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: SPEC HIRA/SPEC CRKL Probe map; Bottom: SPEC MAPK1 Probe map (not to scale); Right: Ideogram of chromosome 22 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a deletion involving the respective gene regions, two blue, two orange, and two green signals appear.

Aberrant situation: In a cell with deletion affecting the HIRA gene region, a reduced number of blue signals will be detected. Deletions affecting both the HIRA and CRKL gene regions will result in a reduced number of blue and orange signals. A reduced number of green signals indicates deletion affecting the MAPK1 gene region. Deletions affecting only parts of the respective gene region might result in a normal signal pattern with signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange, two green, and two blue signals and to metaphase chromosomes of a normal cell.



Lymphocytes and metaphase chromosomes from a DiGeorge syndrome case showing a HIRA/CRKL deletion as indicated by the loss of one blue and one orange signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	22	CE/IVD	5 (50 μl)	Z-2289-50

Kindly provided by Dr. Liehr, Jena, Germany.

ZytoLight® SPEC IGL Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC IGL Dual Color Break Apart Probe (PL241) is intended to be used for the qualitative detection of translocations involving the human IGL locus at 22q11.22 in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with ZytoLight FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC IGL Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 22q11.21-q11.22* (chr22:21,931,816-22,942,402) proximal to the IGL breakpoint region (see Fig.1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 22q11.22-q11.23* (chr22:23,324,781-23,679,042) distal to the IGL breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC IGL Probe map (not to scale); Right: Ideogram of chromosome 22 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the IGL locus, two green/orange fusion signals appear.
Aberrant situation: One IGL locus affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals in each nucleus and to metaphase chromosomes of a normal cell.

Example of an aberrant signal pattern: Cell line with an IGL translocation affecting the 22q11.21-q11.23 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	22	CE/IVD	5 (50 μl)	Z-2286-50

ZytoLight® SPEC SMARCB1/22q12 Dual Color Probe

Intended Purpose
The ZytoLight SPEC SMARCB1/22q12 Dual Color Probe (PL137) is intended to be used for the qualitative detection of deletions involving the human SMARCB1 gene as well as the detection of chromosome 22q12 specific sequences in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC SMARCB1/22q12 Dual Color Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 22q11.23* (chr22:23,887,951-24,431,064) harboring the SMARCB1 gene region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 22q12.1-22q12.2* (chr22:29,340,078-29,673,440) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC SMARCB1 Probe map; Bottom: SPEC 22q12 Probe map (not to scale); Right: Ideogram of chromosome 22 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the SMARCB1 gene region, two green and two orange signals appear.
Aberrant situation: In a cell with deletion affecting the SMARCB1 gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the SMARCB1 gene region might result in a normal signal pattern with green signals of reduced size.



Hybridized to normal interphase cells as indicated by two orange and two green signals per nucleus.

Example of an aberrant signal pattern: SPEC SMARCB1/22q12 Dual Color Probe hybridized to epithelioid sarcoma tissue section with biallelic deletion of the SMARCB1 gene as indicated by missing green signals in the nuclei.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	22	CE/IVD	5 (50 μl)	Z-2178-50

ZytoLight® SPEC EWSR1 Dual Color Break Apart Probe

Intended Purpose
The ZytoLight SPEC EWSR1 Dual Color Break Apart Probe (PL55) is intended to be used for the qualitative detection of translocations involving the human EWSR1 gene at 22q12.2 in formalin-fixed, paraffin-embedded specimens, such as Ewing sarcoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of Ewing sarcoma and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC EWSR1 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 22q12.2* (chr22:29,779,841-30,179,900) distal to the EWSR1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 22q12.1-q12.2* (chr22:29,191,431-29,673,440) proximal to the EWSR1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC EWSR1 Probe map (not to scale); Right: Ideogram of chromosome 22 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the EWSR1 gene region, two green/orange fusion signals appear.
Aberrant situation: One EWSR1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Ewing sarcoma tissue section with translocation affecting the 22q12.1-q12.2 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	22	CE/IVD	5 (50 μl)	Z-2096-50
			20 (200 μl)	Z-2096-200

ZytoLight® SPEC EWSR1/FLI1 TriCheck™ Probe

Intended Purpose
The ZytoLight SPEC EWSR1/FLI1 TriCheck Probe (PL141) is intended to be used for the qualitative detection of rearrangements involving the human EWSR1 gene at 22q12.2 and the human FLI1 gene at 11q24.3 in formalin-fixed, paraffin-embedded specimens, such as Ewing sarcoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of Ewing sarcoma and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoLight SPEC EWSR1/FLI1 TriCheck Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 22q12.2* (chr22:29,779,841-30,179,900) distal to the EWSR1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 22q12.1-q12.2* (chr22:29,191,431-29,673,440) proximal to the EWSR1 breakpoint region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/μl), which target sequences mapping in 11q24.3* (chr11:128,707,454-129,346,602) distal to the FLI1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC EWSR1 Probe map; Bottom: SPEC FLI1 Probe map (not to scale); Right: Ideograms of chromosome 22 and 11 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without FLI1-EWSR1 rearrangement, two green/orange fusion signals and two blue signals are expected.
Aberrant situation: A FLI1-EWSR1 fusion is indicated by one separate orange signal co-localizing with one blue signal and one separate green signal. An EWSR1 translocation without involvement of FLI1 is indicated by the split of one green/orange fusion signal without co-localization of the separated orange signal with one blue signal.



Ewing sarcoma tissue section with FLI1-EWSR1 fusion as indicated by orange/blue fusion signals.

Ewing sarcoma tissue section with a non-FLI1 EWSR1 rearrangement as indicated by the lack of co-localization of the separated orange signal with one blue signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
  	11, 22	CE/IVD	5 (50 μl)	Z-2183-50

ZytoLight® SPEC PDGFB Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC PDGFB Dual Color Break Apart Probe (PL76) is intended to be used for the qualitative detection of translocations involving the human PDGFB gene at 22q13.1 in formalin-fixed, paraffin-embedded specimens, such as dermatofibrosarcoma protuberans (DFSP), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.


The probe is intended to be used as an aid to the differential diagnosis of DFSP and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoLight SPEC PDGFB Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 22q13.1* (chr22:39,720,415-40,267,687) distal to the PDGFB breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 22q13.1* (chr22:38,928,973-39,526,228) proximal to the PDGFB breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Fig. 1: SPEC PDGFB Probe map (not to scale)




Ideogram of chromosome 22 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the PDGFB gene region, two green/orange fusion signals appear.

Aberrant situation: One PDGFB gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Dermatofibrosarcoma protuberans tissue section with translocation affecting the 22q13.1 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	22	CE/IVD	5 (50 μl)	Z-2119-50
			20 (200 μl)	Z-2119-200

ZytoLight® SPEC TFE3 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC TFE3 Dual Color Break Apart Probe (PL66) is intended to be used for the qualitative detection of translocations involving the human TFE3 gene at Xp11.23 in formalin-fixed, paraffin-embedded specimens, such as renal cell carcinomas (RCC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of RCC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoLight SPEC TFE3 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in Xp11.23* (chrX:48,906,685-49,509,699) proximal to the TFE3 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in Xp11.23* (chrX:48,287,169-48,792,674) distal to the TFE3 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Fig. 1: SPEC TFE3 Probe map (not to scale); Right: Ideogram of chromosome X indicating the hybridization locations




Results

Normal situation: In interphases of normal female cells or cells without a translocation involving the TFE3 gene region, two green/orange fusion signals appear. In interphases of normal male cells or cells without a translocation involving the TFE3 gene region one orange/green fusion signal appears.

Aberrant situation: One TFE3 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Renal cell carcinoma section with translocation affecting the Xp11.23 locus as indicated by one non-rearranged green/orange fusion signal, one separate green signal, and one separate orange signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	X	CE/IVD	5 (50 μl)	Z-2109-50
			20 (200 μl)	Z-2109-200

ZytoLight® SPEC CRLF2 Dual Color Break Apart Probe

Intended Purpose

The ZytoLight SPEC CRLF2 Dual Color Break Apart Probe is designed to detect rearrangements involving the chromosomal regions Xp22.33 and Yp11.32 harboring the CRLF2 (cytokine receptor-like factor 2, a.k.a. CRL2, TSLPR) gene. The CRLF2 protein interacts with IL7R to form a receptor for TSLP, binding of which activates cell signaling through JAK/STAT pathways.

Approximately 7% of patients with B-cell precursor ALL (B-ALL) and more than 50% of B-ALL in children with Down syndrome harbor alterations involving the CRLF2 gene. These include the translocations t(X;14)(p22.33;q32.3) or t(Y;14)(p11.32;q32.3) which fuse the entire CRLF2 gene to the immunoglobulin heavy chain enhancer region (IGH-CRLF2).

Another common alteration is an interstitial deletion involving the pseudoautosomal region (PAR1) of the sex chromosomes upstream of CRLF2, juxtaposing the first non-coding exon of P2RY8 to the entire coding region of CRLF2 (P2RY8-CRLF2).

These rearrangements, which are often accompanied by JAK mutations, result in overexpression of CRLF2 and were shown to contribute to lymphoid transformation. Patients with CRLF2 rearrangements and JAK mutations have a poor event-free and overall survival. Moreover, the detection of CRLF2 rearrangements by FISH may help in selecting B-ALL patients eligible for therapy with inhibitors of the JAK/STAT pathway.

References

- Harvey RC, *et al.* (2010) Blood 115: 5312-21.
Mullighan CG, *et al.* (2009) Nat Genet 41: 1243-6.
Roberts KG, *et al.* (2014) N Engl J Med 371: 1005-15.
Russell LJ, *et al.* (2009) Blood 114: 2688-98.
Tasian SK, *et al.* (2012) Blood 120: 833-42.

Probe Description

- The ZytoLight SPEC CRLF2 Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in Xp22.33 (chrX:513,125-1,245,395) and Yp11.32* (chrY:463,125-1,195,395) distal to the CRLF2 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in Xp22.33 (chrX:1,497,976-1,660,328) and Yp11.32* (chrY:1,498,976-1,657,328) proximal to the CRLF2 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19






Fig. 1: Top: chromosome X SPEC CRLF2 Probe map; Bottom: chromosome Y SPEC CRLF2 Probe map (not to scale); Right: Ideograms of chromosome X and Y indicating the hybridization locations



Results

Normal situation: In interphases of normal female cells or cells without a translocation involving the CRLF2 gene region, two orange/green fusion signals appear. In interphases of normal male cells or cells without a translocation involving the CRLF2 gene region, two orange/green fusion signals appear.

Aberrant situation: A signal pattern consisting of one orange/green fusion signal, one orange signal, and a separate green signal indicates one normal CRLF2 gene region and one CRLF2 gene region affected by a translocation.



Bone marrow smear with translocation affecting the CRLF2 gene locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	X,Y	CE/IVD	5 (50 μl)	Z-2201-50
			20 (200 μl)	Z-2201-200

ZytoLight® SPEC BCOR Dual Color Break Apart Probe

Intended Purpose

The *ZytoLight* SPEC BCOR Dual Color Break Apart Probe is designed to detect rearrangements involving the chromosomal region Xp11.4 harboring the BCOR (BCL6 corepressor, a.k.a. KIAA1575) gene. In the 2020 WHO classification of soft tissue and bone tumors, BCOR-rearranged sarcoma is recognized as a distinct entity due to particular morphological, immunohistochemical, and molecular features and differing clinical outcomes compared to other undifferentiated sarcomas.

A fusion between BCOR and CCNB3 can be found in about 60% of all BCOR-rearranged sarcomas. The BCOR-CCNB3 fusion results from an X-chromosomal paracentric inversion. In vitro studies suggest that the BCOR-CCNB3 fusion protein is oncogenic and drives proliferation in these sarcomas. In addition, alternative fusion partners have been identified, including MAML3 and ZC3H7B.

BCOR-rearranged sarcoma usually occurs in bone or soft tissue of predominantly male children with a median age in the second decade of life. There are considerable overlapping morphological and immunohistochemical features with classical Ewing sarcoma, other subtypes of small round cell tumors, as well as lymphomas and carcinomas. Therefore, the evaluation of the BCOR rearrangement status by FISH may be of diagnostic relevance.

References

Antonescu CR, *et al.* (ed.) (2020) WHO Classification of Tumours Soft Tissue and Bone Tumours (5th Edition).
Pierron G, *et al.* (2012) Nat Genet 44: 461-6.
Renzi S, *et al.* (2019) J Cell Physiol 234: 7999-8007.
Sirisena UDN, *et al.* (2021) Skeletal Radiol 50: 1291-301.
Specht K, *et al.* (2016) Am J Surg Pathol 40: 433-42.


Probe Description

- The *ZytoLight* SPEC BCOR Dual Color Break Apart Probe is composed of:
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in Xp11.4* (chrX:39,262,996-39,850,787) distal to the BCOR breakpoint region.
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in Xp11.4* (chrX:39,998,508-40,345,270) proximal to the BCOR breakpoint region.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19


Results

Normal situation: In interphases of normal female cells or cells without a translocation involving the BCOR gene region, two green/orange fusion signals appear. In interphases of normal male cells or cells without a translocation involving the BCOR gene region, one green/orange fusion signal appears.

Aberrant situation: One BCOR gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Sarcoma tissue section with translocation affecting the BCOR gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal indicating the translocation.

Labels	Chromosome	Reg. Status	Tests	Order No.
	X	CE/IVD	5 (50 μl)	Z-2310-50




Fig. 1: SPEC BCOR Probe map (not to scale)



Ideogram of chromosome X indicating the hybridization locations

ZytoLight® Probes for Chromosome Enumeration - CE/IVD

The *ZytoLight*® Chromosome Enumeration Probes are designed for identification and enumeration of human chromosomes in interphase cells and as an adjunct to standard karyotyping in metaphases. These probes will produce sharp, bright signals specific for each individual chromosome.

CEN Probe Description

For most chromosomes, direct labeled *ZytoLight*® CEN™ Probes hybridizing to highly repetitive human satellite DNA sequences mainly located at the centromeric regions of chromosomes are applicable.

SPEC Probe Description

As several chromosomes share the same repetitive sequences resulting in cross-hybridization signals, they cannot be differentiated by centromere specific probes. Instead, these chromosomes can be identified by direct labeled *ZytoLight*® SPEC™ Probes hybridizing in close proximity to the respective satellite DNA sequences or to other chromosome specific loci.

Results

In a normal interphase nucleus, two signals are expected using Chromosome Enumeration Probes specific for autosomes.

Using chromosome Y specific probes will result in normal male cells in one signal and in normal female cells in no signal.


Using chromosome X specific probes will result in normal male cells in one signal and in normal female cells in two signals per nucleus.

Other signal patterns indicate numerical aberrations of the respective chromosome.

ZytoLight® SPEC Probe maps

The *ZytoLight* SPEC 13q12 Probe is composed of:


- ZyOrange (excitation 547 nm/emission 467 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 13q12.11* (chr13:20,200,365-20,892,494).
- Formamide based hybridization buffer



SPEC 13q12 Probe map (not to scale)

The SPEC 18q21 Probe, included in the *ZytoLight* SPEC 18/CEN X/Y Triple Color Probe, is composed of:


- ZyBlue (excitation 418 nm, Emission 467 nm) labeled polynucleotides (~37 ng/μl), which target sequences mapping in 18q21.31-q21.32* (chr18:55,690,725-56,455,119).



SPEC 18q21 Probe map (not to scale)

The *ZytoLight* SPEC 21q22 Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/μl), which target sequences mapping in 21q22.13-q22.2* (chr21:39,372,983-39,784,773).
- Formamide based hybridization buffer



SPEC 21q22 Probe map (not to scale)

* according to Human Genome Assembly GRCh37/hg19

CE/IVD Probes

Product Name	Alpha/Class. Sat.	Chromosome band	Reg. Status	Labels	Tests	Order No.
ZytoLight CEN 8 Probe	D8Z2	8p11.1-q11.1	CE/IVD		5/20 (50/200 µl)	Z-2004-50/-200
ZytoLight CEN 11 Probe	D11Z1	11p11.11-q11			20 (200 µl)	Z-2005-200
ZytoLight CEN 12 Probe	D12Z3	12p11.1-q11			20 (200 µl)	Z-2050-200
ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe	D18Z1	13q12.11/18p11.1/21q22.13-q22.2			5/20 (50/200 µl)	Z-2095-50/-200
ZytoLight SPEC 13/21 Dual Color Probe	-	13q12.11/21q22.13-q22.2			20 (200 µl)	Z-2164-200
ZytoLight SPEC 18/CEN X/Y Triple Color Probe	DXZ1/DYZ3	18q21.31-q21.32/Xp11.1-q11.1/Yp11.1-q11.1			20 (200 µl)	Z-2163-200
ZytoLight SPEC 21/CEN X/Yq12 Triple Color Probe	DXZ1/ III DYZ1	21q22.13-q22.2/Xp11.1-q11.1/Yq12			20 (200 µl)	Z-2180-200
ZytoLight CEN X/Yq12 Dual Color Probe	DXZ1/III DYZ1	Xp11.1-q11.1/Yq12			5/20 (50/200 µl)	Z-2016-50/-200
ZytoLight CEN X/Y Dual Color Probe	DXZ1/ DYZ3	Xp11.1-q11.1/Yp11.1-q11.1			20 (200 µl)	Z-2120-200

Related Products				
ZytoLight Aneuploidy Panel 18/X/Y and 13/21 Incl. ZytoLight SPEC 18/CEN X/Y Triple Color Probe, 0.2 ml (Z-2163-200); ZytoLight SPEC 13/21 Dual Color Probe, 0.2 ml (Z-2164-200)	CE/IVD	-	20	Z-2279-20
ZytoLight Aneuploidy Panel X/Y and 13/18/21 Incl. ZytoLight CEN X/Yq12 Dual Color Probe, 0.05 ml (Z-2016-50); ZytoLight SPEC 13/CEN 18/ SPEC 21 Triple Color Probe, 0.05 ml (Z-2095-50)			5	Z-2104-5
ZytoLight Aneuploidy Panel X/Y and 13/18/21 Incl. ZytoLight CEN X/Yq12 Dual Color Probe, 0.2 ml (Z-2016-200); ZytoLight SPEC 13/CEN 18/ SPEC 21 Triple Color Probe, 0.2 ml (Z-2095-200)			20	Z-2104-20
ZytoLight FISH-Tissue Implementation Kit Incl. Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 210 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml			5	Z-2028-5
ZytoLight FISH-Tissue Implementation Kit Incl. Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 560 ml; 25x Wash Buffer A, 100 ml; DAPI/DuraTect-Solution, 0.8 ml			20	Z-2028-20
ZytoLight FISH-Cytology Implementation Kit Incl. Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x MgCl ₂ , 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/ DuraTect-Solution, 0.8 ml			20	Z-2099-20

ZytoLight® Probes for Chromosome Enumeration - RUO

The ZytoLight® Chromosome Enumeration Probes are designed for identification and enumeration of human chromosomes in interphase cells and as an adjunct to standard karyotyping in metaphases. These probes will produce sharp, bright signals specific for each individual chromosome.

CEN Probe Description


For most chromosomes, direct labeled ZytoLight® CEN™ Probes hybridizing to highly repetitive human satellite DNA sequences mainly located at the centromeric regions of chromosomes are applicable.

SPEC Probe Description

As several chromosomes share the same repetitive sequences resulting in cross-hybridization signals, they cannot be differentiated by centromere specific probes. Instead, these chromosomes can be identified by direct labeled ZytoLight® SPEC™ Probes hybridizing in close proximity to the respective satellite DNA sequences or to other chromosome specific loci.


ZytoLight® SPEC Probe maps

- The **ZytoLight SPEC 1p12 Probe** is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 1p12* (chr1:119,537,102-119,823,147).
 - Formamide based hybridization buffer




SPEC 1p12 Probe map (not to scale)

- The **ZytoLight SPEC 13q12 Probe** is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 13q12.11* (chr13:20,200,365-20,892,494).
 - Formamide based hybridization buffer




SPEC 13q12 Probe map (not to scale)

- The **ZytoLight SPEC 2q11 Probe** is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 2q11.2* (chr2:100,132,806-100,621,725).
 - Formamide based hybridization buffer




SPEC 2q11 Probe map (not to scale)

- The **ZytoLight SPEC 21q22 Probe** is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 21q22.13-q22.2* (chr21:39,372,983-39,784,773).
 - Formamide based hybridization buffer



SPEC 21q22 Probe map (not to scale)

- The **ZytoLight SPEC 4p11 Probe** is composed of:
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 4p11* (chr4:48,329,914- 48,762,386).
 - Formamide based hybridization buffer



SPEC 4p11 Probe map (not to scale)

Results

In a normal interphase nucleus, two signals are expected using Chromosome Enumeration Probes specific for autosomes.

Using chromosome Y specific probes will result in normal male cells in one signal and in normal female cells in no signal.

Using chromosome X specific probes will result in normal male cells in one signal and in normal female cells in two signals per nucleus.

Other signal patterns indicate numerical aberrations of the respective chromosome.

* according to Human Genome Assembly GRCh37/hg19

RUO Probes

Product Name	Alpha/Class. Sat.	Chromosome band	Reg. Status	Labels	Tests	Order No.
ZytoLight SPEC 1p12 Probe	-	1p12	RUO	●	20 (200 µl)	Z-2101-200
ZytoLight SPEC 2q11 Probe	-	2q11.2		●	20 (200 µl)	Z-2049-200
ZytoLight CEN 3 Probe	D3Z1	3p11.1-q11.1		●	20 (200 µl)	Z-2001-200
ZytoLight SPEC 4p11 Probe	-	4p11		●	20 (200 µl)	Z-2083-200
ZytoLight CEN 6 Probe	D6Z1	6p11.1-q11		●	20 (200 µl)	Z-2002-200
ZytoLight CEN 7 Probe	D7Z1	7p11.1-q11.1		●	20 (200 µl)	Z-2003-200
ZytoLight CEN 9 Probe	III D9Z3	9q12		●	20 (200 µl)	Z-2067-200
ZytoLight CEN 10 Probe	D10Z1	10p11.1-q11.1		●	20 (200 µl)	Z-2079-200
ZytoLight CEN 17 Probe	D17Z1	17p11.1-q11.1		●	20 (200 µl)	Z-2006-200
ZytoLight CEN 18 Probe	D18Z1	18p11.1-q11.1		●	20 (200 µl)	Z-2007-200
ZytoLight SPEC 21q22 Probe	-	21q22.13-q22.2	RUO	●	20 (200 µl)	Z-2086-200
ZytoLight CEN X Probe	DXZ1	Xp11.1-q11.1		●	20 (200 µl)	Z-2008-200
ZytoLight CEN Yq12 Probe	III DYZ1	Yq12		●	20 (200 µl)	Z-2010-200

ZytoLight® Aneuploidy Panel 18/X/Y and 13/21

Intended Purpose

The *ZytoLight* Aneuploidy Panel 18/X/Y and 13/21 is intended to be used for the qualitative detection of human chromosome 13 and 21 specific sequences as well as human chromosome 18 specific sequences and alpha satellites of chromosome X and Y in cytologic or formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH).

The probes are intended to be used in combination with *ZytoLight* FISH Implementation Kits (Prod. No. Z-2028-5/-20, or Z-2099-20). Interpretation of the results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

Probe Description

The *ZytoLight* Aneuploidy Panel 18/X/Y and 13/21 is a set comprising two separate probes:

- *ZytoLight* SPEC 18/CEN X/Y Triple Color Probe (Prod. No. Z-2163-200)
- *ZytoLight* SPEC 13/21 Dual Color Probe (Prod. No. Z-2164-200)

The *ZytoLight* SPEC 18/CEN X/Y Triple Color Probe (PL119) is composed of:

- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~37,0 ng/µl), which target sequences mapping in 18q21.31-q21.32* (chr18:55,690,725-56,455,119) (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/µl), which target sequences mapping in Xp11.1-q11.1 specific for the alpha satellite centromeric region DXZ1 of chromosome X.
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1,5 ng/µl), which target sequences mapping in Yp11.1-q11.1 specific for the alpha satellite centromeric region DYZ3 of chromosome Y.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19

The *ZytoLight* SPEC 13/21 Dual Color Probe (PL120) is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/µl), which target sequences mapping in 13q12.11* (chr13:20,200,365-20,892,494) (see Fig. 2).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/µl), which target sequences mapping in 21q22.13-q22.2* (chr21:39,372,983-39,784,773) (see Fig. 2).
- Formamide based hybridization buffer


* according to Human Genome Assembly GRCh37/hg19

Results


Normal situation: In an interphase nucleus of a normal cell using the *ZytoLight* SPEC 13/21 Dual Color Probe, two green and two orange signals are expected.

In an interphase nucleus of a normal cell, using the *ZytoLight* SPEC 18/CEN X/Y Triple Color Probe, two blue signals are expected. Two green signals are expected in a normal female cell, or one single green and one single orange signal is expected in a normal male cell.

Aberrant situation: Other signal patterns indicate numerical aberration of the respective chromosomes.



SPEC 18/CEN X/Y Triple Color Probe hybridized to interphase nuclei of normal male cells and to chromosomes of a metaphase spread.



SPEC 13/21 Dual Color Probe hybridized to interphase cells with trisomy of chromosome 21.

ZytoLight Aneuploidy Panel 18/X/Y and 13/21

Incl. *ZytoLight* SPEC 18/CEN X/Y Triple Color Probe, 0.2 ml (Z-2163-200); *ZytoLight* SPEC 13/21 Dual Color Probe, 0.2 ml (Z-2164-200)

Reg. Status	Tests	Order No.
CE/IVD	2 x 200 µl (20 Tests)	Z-2279-20




Fig. 1: SPEC 18q21 Probe map (not to scale)



Fig. 2: Top: SPEC 13q12 Probe map; Bottom: SPEC 21q22 Probe map (not to scale)

ZytoLight® Aneuploidy Panel X/Y and 13/18/21

Intended Purpose

The *ZytoLight* Aneuploidy Panel X/Y and 13/18/21 is designed for enumeration of the chromosomes 13, 18, 21, X, and Y. Trisomies of the autosomes 13, 18, or 21 (Down Syndrome) are common genomic aberrations. Aberrant numbers of the gonosomes X and Y are resulting in disorders of sex development (DSD). Diseases such as Ulrich-Turner-Syndrome (45, X) or Triple X Syndrome (47, XXX) may cause severe developmental and metabolic disorders. The prevalence of chromosomal abnormalities detectable in the newborn including chromosomes 13, 18, 21, X, and Y, is about 0.92%.

Probe Description

The *ZytoLight* Aneuploidy Panel X/Y and 13/18/21 is a set comprising two separate probes:

- *ZytoLight* CEN X/Yq12 Dual Color Probe (Prod. No. Z-2016-50/200)
- *ZytoLight* SPEC 13/CEN 18/SPEC 21 Triple Color Probe (Prod. No. Z-2095-50/200)

The *ZytoLight* CEN X/Yq12 Dual Color Probe (PL3) is composed of:

- ZyOrange (excitation 547 nm /Emission 572 nm) labeled polynucleotides (~1,5 ng/μl), which target sequences mapping in Xp11.1-q11.1 specific for the alpha satellite centromeric region DXZ1 of chromosome X.
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in Yq12, specific in the classical satellite III centromeric region DYZ1 of chromosome Y.
- Formamide based hybridization buffer

The *ZytoLight* SPEC 13/CEN 18/SPEC 21 Triple Color Probe (PL54) is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10,0 ng/μl), which target sequences mapping in 13q12.11* (chr13:20,200,365-20,892,494) (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~12,0 ng/μl), which target sequences mapping in 18p11.1-q11.1 specific for the alpha satellite centromeric region D18Z1 of chromosome 18.
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,5 ng/μl), which target sequences mapping in 21q22.13-q22.2* (chr21:39,372,983-39,784,773) (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC 13q12 Probe map; Bottom: SPEC 21q22 Probe map (not to scale)


Results

Normal situation: In an interphase nucleus, using the *ZytoLight* CEN X/Yq12 Dual Color Probe, two orange signals are expected in a normal female cell whereas one single orange and one single green signal is expected in a normal male cell. In an interphase nucleus of a normal cell, using the *ZytoLight* SPEC 13/CEN 18/SPEC 21 Triple Color Probe, two green, two blue, and two orange signals are expected.

Aberrant situation: Other signal patterns indicate numerical aberrations of the respective chromosomes.



SPEC 13/CEN 18/ SPEC 21 Triple Color Probe hybridized to normal interphase cells.



CEN X/Yq12 Dual Color Probe hybridized to metaphase chromosomes of a normal male cell.

ZytoLight Aneuploidy Panel X/Y and 13/18/21

Incl. *ZytoLight* CEN X/Yq12 Dual Color Probe, 0.05 ml (Z-2016-50); *ZytoLight* SPEC 13/CEN 18/SPEC 21 Triple Color Probe, 0.05 ml (Z-2095-50)

ZytoLight Aneuploidy Panel X/Y and 13/18/21



Incl. *ZytoLight* CEN X/Yq12 Dual Color Probe, 0.2 ml (Z-2016-200); *ZytoLight* SPEC 13/CEN 18/SPEC 21 Triple Color Probe, 0.2 ml (Z-2095-200)

Reg. Status	Tests	Order No.
CE/IVD	2 x 50 μl (5 Tests)	Z-2104-5
	2 x 200 μl (20 Tests)	Z-2104-20

Additional *ZytoLight*® RUO Probes


Product Name	Chromosomenbande	Reg. Status	Labels	Tests	Order No.
<i>ZytoLight</i> SPEC TNFRSF14/1q25 Dual Color Probe	1p36.3	RUO	● ●	5 (50 μl)	Z-2323-50
<i>ZytoLight</i> SPEC EML4 Dual Color Break Apart Probe	2p21		● ●	5 (50 μl)	Z-2136-50
<i>ZytoLight</i> SPEC FHIT/CEN 3 Dual Color Probe	3p14.2		● ●	20 (200 μl)	Z-2062-200
<i>ZytoLight</i> SPEC TERC/CEN 3 Dual Color Probe	3q26.2		● ●	20 (200 μl)	Z-2284-200
<i>ZytoLight</i> SPEC SOX2/CEN 3 Dual Color Probe	3q26.3		● ●	20 (200 μl)	Z-2127-200
<i>ZytoLight</i> SPEC RICTOR/5q31.1 Dual Color Probe	5p13.1		● ●	20 (200 μl)	Z-2278-200
<i>ZytoLight</i> SPEC PLAG1 Dual Color Break Apart Probe	8q12		● ●	5 (50 μl)	Z-2326-50
<i>ZytoLight</i> SPEC KIF5B Dual Color Break Apart Probe	10p11.2		● ●	5 (50 μl)	Z-2131-50
<i>ZytoLight</i> SPEC BCL2L1/CEN 20 Dual Color Probe	20q11.2		● ●	20 (200 μl)	Z-2171-200
			● ●		

Z-2323-50




Top: SPEC TNFRSF14 Probe map (not to scale); Bottom: SPEC 1q25 Probe map (not to scale)

Z-2136-50




SPEC EML4 Probe map (not to scale)

Z-2062-200




SPEC FHIT Probe map (not to scale)

Z-2284-200




SPEC TERC Probe map (not to scale)

Z-2127-200




SPEC SOX2 Probe map (not to scale)

Z-2278-200




Top: SPEC RICTOR Probe map (not to scale); Bottom: SPEC 5q31.1 Probe map (not to scale)

Z-2326-50




SPEC PLAG1 Probe map (not to scale)

Z-2131-50



SPEC KIF5B Probe map (not to scale)

Z-2171-200



SPEC BCL2L1 Probe map (not to scale)

Accessories

ZytoLight® Implementation Kits

For the detection of *ZytoLight* probes

Product	Manufacturer	Tests	Reg. Status	Order No.
ZytoLight FISH-Tissue Implementation Kit				
Incl.: Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 210 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml	ZytoVision GmbH	5	CE/IVD	Z-2028-5
ZytoLight FISH-Tissue Implementation Kit				
Incl.: Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 560 ml; 25x Wash Buffer A, 100 ml; DAPI/DuraTect-Solution, 0.8 ml	ZytoVision GmbH	20	CE/IVD	Z-2028-20
ZytoLight FISH-Cytology Implementation Kit				
Incl.: Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x MgCl ₂ , 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/DuraTect-Solution, 0.8 ml	ZytoVision GmbH	20	CE/IVD	Z-2099-20

ZytoLight® Pretreatment Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
Pepsin Solution	ZytoVision GmbH	4 ml	CE/IVD	ES-0001-4
Pepsin Solution	ZytoVision GmbH	50 ml	CE/IVD	ES-0001-50
Pepsin Solution	ZytoVision GmbH	1000 ml	CE/IVD	ES-0001-1000
Cytology Pepsin Solution	ZytoVision GmbH	4 ml	CE/IVD	ES-0002-4
Cytology Pepsin Solution	ZytoVision GmbH	50 ml	CE/IVD	ES-0002-50
Heat Pretreatment Solution Citric	ZytoVision GmbH	2 x 500 ml	CE/IVD	PT-0001-1000
Formaldehyde Dilution Buffer Set Incl.: 10x MgCl ₂ , 50 ml; 10x PBS, 50 ml	ZytoVision GmbH	2 x 50 ml	CE/IVD	PT-0006-100


ZytoLight® Wash Buffers & Ancillary Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
DAPI/DuraTect™-Solution	ZytoVision GmbH	150 ng DAPI/ml, 0.8 ml	CE/IVD	MT-0007-0.8
DAPI/DuraTect™-Solution (ultra)	ZytoVision GmbH	1360 ng DAPI/ml, 0.8 ml	CE/IVD	MT-0008-0.8
Wash Buffer SSC	ZytoVision GmbH	560 ml	CE/IVD	WB-0001-560
25x Wash Buffer A	ZytoVision GmbH	50 ml	CE/IVD	WB-0002-50
20x Wash Buffer TBS	ZytoVision GmbH	50 ml	CE/IVD	WB-0005-50
Cytology Stringency Wash Buffer SSC	ZytoVision GmbH	500 ml	CE/IVD	WB-0007-500
Cytology Wash Buffer SSC	ZytoVision GmbH	500 ml	CE/IVD	WB-0008-500

RUO *ZytoLight*® Wash Buffers & Ancillary Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
20x SSC Solution	ZytoVision GmbH	50 ml	RUO	WB-0003-50

ZyBlack™ Quenching Solution




Kidney tissue section hybridized with the *ZytoLight* SPEC PTEN/CEN 10 Dual Color Probe.

Product Description

ZyBlack™ Quenching Solution is a ready-to-use solution to reduce autofluorescence on formalin-fixed paraffin-embedded specimens. It can be easily incorporated into the manual FISH protocol by applying it after the proteolytic pretreatment.

One of the major concerns of Fluorescence *in situ* Hybridization (FISH)-based diagnostic assays is the interference by autofluorescence. Several types of tissue tend to emit intense autofluorescence, including brain, liver, kidney and myocardium, making it difficult to evaluate FISH results.

ZyBlack™ Quenching Solution reduces autofluorescence without adversely affecting tissue integrity or specific fluorescence signals.



Product	Manufacturer	Amount	Reg. Status	Order No.
ZyBlack Quenching Solution	ZytoVision GmbH	8 ml		BS-0002-8

FlexISH®

FlexISH® products are designed for identification of chromosomal aberrations on formalin-fixed, paraffin-embedded (FFPE) specimens by FISH. Using the FlexISH® products gives you the flexibility to choose between a 1-day (2 h hybridization) or a 2-day (overnight hybridization) protocol by adapting the hybridization time just according to your individual needs!

 ZytoVision GmbH

Advantages of FlexISH®


- Hybridization time can be varied between 2 hours and overnight.
- With a hybridization temperature of 37°C, the FlexISH® protocol is fully compatible with routine workflows in pathology laboratories.

FlexISH®-Kit – Convenient Solution


- FlexISH® probes can be combined with the FlexISH®-Tissue Implementation Kit to obtain reliable results already within 4.5 hours.
- The FlexISH® protocol can also be incorporated into the routine workflow with overnight hybridization providing the highest flexibility.

High-Quality FISH Results with flexible Hybridization Time

The hybridization time can be varied freely between 2 hours and overnight.



Chromosome Index




Chr. Band	Product Name	Reg. Status	Order No.	Page
1q23.1	FlexISH NTRK1/NTRK3 DistinguishISH™ Probe	CE/IVD	Z-2314-50/-200	119




Chr. Band	Product Name	Reg. Status	Order No.	Page
2p23	FlexISH ALK/ROS1 DistinguishISH™ Probe	CE/IVD	Z-2203-50/-200	120




Chr. Band	Product Name	Reg. Status	Order No.	Page
3q27	FlexISH BCL2/BCL6 DistinguishISH™ Probe	CE/IVD	Z-2283-50/-200	121



Chr. Band	Product Name	Reg. Status	Order No.	Page
6q22.1	FlexISH ALK/ROS1 DistinguishISH™ Probe	CE/IVD	Z-2203-50/-200	120



Chr. Band	Product Name	Reg. Status	Order No.	Page
8q24.21	FlexISH MYC/IGH TriCheck™ Probe	CE/IVD	Z-2293-50	122



Chr. Band	Product Name	Reg. Status	Order No.	Page
10p11.2	FlexISH RET/KIF5B TriCheck™ Probe	CE/IVD	Z-2269-200	123
10q11.2	FlexISH RET/KIF5B TriCheck™ Probe	CE/IVD	Z-2269-200	123

Chr. Band	Product Name	Reg. Status	Order No.	Page
14q32.3	FlexISH MYC/IGH TriCheck™ Probe	CE/IVD	Z-2293-50	122

Chr. Band	Product Name	Reg. Status	Order No.	Page
15q25	FlexISH NTRK1/NTRK3 DistinguISH™ Probe	CE/IVD	Z-2314-50/-200	119

Chr. Band	Product Name	Reg. Status	Order No.	Page
17q12	FlexISH ERBB2/CEN 17 Dual Color Probe	CE/IVD	Z-2166-50/-200	124

Chr. Band	Product Name	Reg. Status	Order No.	Page
18q21.3	FlexISH BCL2/BCL6 DistinguISH™ Probe	CE/IVD	Z-2283-50/-200	121

Gene Index

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
ALK	CD246	FlexISH ALK/ROS1 DistinguISH™ Probe	CE/IVD	Z-2203-50/-200	120
BCL2	Bcl-2, PPP1R50	FlexISH BCL2/BCL6 DistinguISH™ Probe	CE/IVD	Z-2283-50/-200	121
BCL6	ZNF51, LAZ3	FlexISH BCL2/BCL6 DistinguISH™ Probe	CE/IVD	Z-2283-50/-200	121
ERBB2	HER2, HER-2, NEU	FlexISH ERBB2/CEN 17 Dual Color Probe	CE/IVD	Z-2166-50/-200	124
IGH	IGH@	FlexISH MYC/IGH TriCheck™ Probe	CE/IVD	Z-2293-50	122
KIF5B	KNS1	FlexISH RET/KIF5B TriCheck™ Probe	CE/IVD	Z-2269-200	123
MYC	CMYC, bHLHe39, c-Myc	FlexISH MYC/IGH TriCheck™ Probe	CE/IVD	Z-2293-50	122
NTRK1	MTC, TRK	FlexISH NTRK1/NTRK3 DistinguISH™ Probe	CE/IVD	Z-2314-50/-200	119
NTRK3	TRKC	FlexISH NTRK1/NTRK3 DistinguISH™ Probe	CE/IVD	Z-2314-50/-200	119
RET	HSCR1, CDHF12	FlexISH RET/KIF5B TriCheck™ Probe	CE/IVD	Z-2269-200	123
ROS1	MCF3, ROS	FlexISH ALK/ROS1 DistinguISH™ Probe	CE/IVD	Z-2203-50/-200	120

FlexISH® NTRK1/NTRK3 DistinguISH™ Probe

Intended Purpose
The FlexISH NTRK1/NTRK3 DistinguISH™ Probe is designed to detect rearrangements affecting the chromosomal region 1q23.1 and 15q25.3 harboring the NTRK1 (neurotrophic receptor tyrosine kinase 1, a.k.a. TRKA, TRK) and NTRK3 (neurotrophic receptor tyrosine kinase 3, a.k.a. TRKC) gene region, respectively.

The neurotrophic tyrosine receptor kinase genes (NTRK1, NTRK2, and NTRK3) encode a family of receptor tyrosine kinases that serve important roles in cell survival, proliferation, and cellular differentiation in healthy human cells. The tumor types in which NTRK gene fusions have been detected are diverse, and include, e.g., breast cancer, non-small cell lung cancer, sarcoma, melanoma, and thyroid carcinoma. The treatment of patients with NTRK fusion-positive cancers with a NTRK inhibitor, such as the FDA-approved drugs larotrectinib or entrectinib, is associated with high response rates, regardless of NTRK gene, fusion partner, and tumor type.

Hence, detection of NTRK1 and NTRK3 rearrangements by FISH may be of therapeutic significance.

References
Haller F, *et al.* (2016) J Pathol 238: 700-10.
Hsiao SJ, *et al.* (2019) J Mol Diagn 21: 553-71.
Knezevich SR, *et al.* (1998) Nat Genet 18: 184-7.
Martin-Zanca D, *et al.* (1986) Nature 319: 743-8.
Solomon JP & Hechtman JF (2019) Cancer Res 79: 3163-8.

Probe Description
The FlexISH NTRK1/NTRK3 DistinguISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 1q22-q23.1* (chr1:156,245,849-156,781,745), proximal to the NTRK1 breakpoint region, and in 15q25.3-q26.1* (chr15:88,825,346-89,475,889), distal to the NTRK3 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2,5 ng/μl), which target sequences mapping in 1q23.1* (chr1:156,854,527-157,296,918), distal to the NTRK1 breakpoint region, and in 15q25.3* (chr15:87,976,717-88,471,002), proximal to the NTRK3 breakpoint region (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~70,0 ng/μl), which target sequences mapping in 15q25.3-q26.1* (chr15:87,845,459-89,475,889) harboring the NTRK3 gene region (see Fig. 1).

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: NTRK1 Probe map; Bottom: NTRK3 Probe map (not to scale)



Ideograms of chromosome 1 and 15 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a rearrangement involving the NTRK1 or NTRK3 gene region, four green/orange fusion signals appear when using an appropriate dual bandpass filter set, and two blue signals appear when using an appropriate single bandpass filter set.
Aberrant situation: One NTRK1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal not co-localizing with blue signals. Loss of one green signal resulting in an isolated orange signal is the result of a deletion proximal to the NTRK1 breakpoint region. One NTRK3 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal, each co-localizing with a blue signal. Loss of one green signal resulting in an isolated orange signal co-localizing with a blue signal is the result of a deletion distal to the NTRK3 breakpoint region.



Cell which shows two green/orange/blue fusion signals (NTRK3) and one green/orange/blue fusion signal (NTRK1). NTRK1 rearrangement is indicated by one isolated orange signal, not co-localizing with a blue signal.

Cell which shows two green/orange fusion signals and one green/orange/blue fusion signal. NTRK3 rearrangement is indicated by one separate orange and one separate green signal, both co-localizing with blue signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
● ● ●	1, 15	CE/IVD	5 (50 μl)	Z-2314-50
			20 (200 μl)	Z-2314-200

FlexISH® ALK/ROS1 DistinguISH™ Probe

Intended Purpose

The FlexISH ALK/ROS1 DistinguISH Probe (PL161) is intended to be used for the qualitative detection of translocations involving the human ALK gene at 2p23.1-p23.2 and the human ROS1 gene at 6q22.1 in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the FlexISH-Tissue Implementation Kit (Prod. No. Z-2182-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The FlexISH ALK/ROS1 DistinguISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 2p23.1-p23.2* (chr2:29,460,144-30,095,822), proximal to the ALK breakpoint region, and in 6q22.1* (chr6:116,912,298-117,627,255), proximal to the ROS1 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2,5 ng/μl), which target sequences mapping in 2p23.2* (chr2:29,174,204-29,383,335), distal to the ALK breakpoint region, and in 6q22.1* (chr6:117,659,135-117,871,701), distal to the ROS1 breakpoint region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~70.0 ng/μl), which target sequences mapping in 6q22.1* (chr6:116,671,642-117,260,761) proximal to the ROS1 breakpoint region co-localizing with the green-labeled ROS1 polynucleotides and in 6q22.1-q22.2* (chr6:117,765,211-118,444,005) distal to the ROS1 breakpoint region co-localizing with the orange-labeled ROS1 polynucleotides (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: ALK Probe map; Bottom: ROS1 Probe map (not to scale)




Ideograms of chromosome 2 and 6 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a rearrangement involving the ALK or ROS1 gene region, four green/orange fusion signals appear when using an appropriate dual bandpass filter set, and two blue signals appear when using an appropriate single bandpass filter set.

Aberrant situation: One ALK gene region affected by a translocation is indicated by one separate green signal and one separate orange signal not co-localizing with blue signals. Loss of one green signal resulting in an isolated orange signal is the result of a deletion of 5'-ALK sequences. One ROS1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal, each co-localizing with a blue signal. Loss of one orange signal resulting in an isolated green signal co-localizing with a blue signal is the result of a deletion distal to the ROS1 breakpoint region or due to an unbalanced translocation affecting this chromosomal region.



H3122 cell line which shows two green/orange/blue fusion signals and one orange/green fusion signal. An ALK rearrangement is indicated by one separate orange and one separate green signal, both not co-localizing with blue signals.



Paraffin-embedded HCC78 cell line which shows two green/orange fusion signals and one green/orange/blue fusion signal. ROS1 rearrangement is indicated by one separate orange and one separate green signal, both co-localizing with blue signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
  	2, 6	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2203-50 Z-2203-200

FlexISH® BCL2/BCL6 DistinguISH™ Probe

Intended Purpose

The FlexISH BCL2/BCL6 DistinguISH Probe (PL238) is intended to be used for the qualitative detection of translocations involving the human BCL2 gene at 18q21.33 and the human BCL6 gene at 3q27.3 in formalin-fixed, paraffin-embedded specimens, such as B-cell lymphoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the FlexISH-Tissue Implementation Kit (Prod. No. Z-2182-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of B-cell lymphoma and therapeutic measures should not be initiated based on the test result alone.


Probe Description

The FlexISH BCL2/BCL6 DistinguISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 18q21.33* (chr18:60,046,152-60,779,138), proximal to the BCL2 breakpoint region, and in 3q27.3* (chr3:186,578,337-187,403,834), proximal to the BCL6 breakpoint region (see Fig. 1).
 - ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2,5 ng/μl), which target sequences mapping in 18q21.33-q22.1* (chr18:60,994,528-61,658,503), distal to the BCL2 breakpoint region, and in 3q27.3-q28* (chr3:187,744,962-188,411,425), distal to the BCL6 breakpoint region (see Fig. 1).
 - ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~70.0 ng/μl), which target sequences mapping in 3q27.3* (chr3:186,578,337-187,403,834) proximal to the BCL6 breakpoint region co-localizing with the green-labeled BCL6 polynucleotides and in 3q27.3-q28* (chr3:187,744,962-188,411,425) distal to the BCL6 breakpoint region co-localizing with the orange-labeled BCL6 polynucleotides (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Fig. 1: Top: BCL2 Probe map; Bottom: BCL6 Probe map (not to scale)




Ideograms of chromosome 3 and 18 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a rearrangement involving the BCL2 or BCL6 gene region, four green/orange fusion signals appear when using an appropriate dual bandpass filter set, and two blue signals appear when using an appropriate single bandpass filter set.

Aberrant situation: One 18q21.33-q22.1 locus affected by a BCL2 translocation is indicated by one separate green signal and one separate orange signal not co-localizing with blue signals. One 3q27.3-q28 locus affected by a BCL6 translocation is indicated by one separate green signal and one separate orange signal, each co-localizing with a blue signal.



Lymphoma tissue which shows two green/orange/blue fusion signals and one green/orange fusion signal. BCL2 rearrangement is indicated by one separate green and one separate orange signal, both not colocalizing with blue signals.



DLBCL tissue which shows one green/orange/blue fusion signal and one green/orange fusion signal. BCL6 rearrangement is indicated by one separate green and one separate orange signal, both co-localizing with blue signals. Additionally, one separate orange and one separate green signal indicate a further BCL2 positivity, confirming a BCL2/BCL6 co-rearrangement.

Labels	Chromosome	Reg. Status	Tests	Order No.
  	3, 18	CE/IVD	5 (50 μl) 20 (200 μl)	Z-2283-50 Z-2283-200

Specimen kindly provided by Dr. Rontogianni, Athens, Greece.

FlexISH® MYC/IGH TriCheck™ Probe

Intended Purpose

The FlexISH MYC/IGH TriCheck Probe (PL247) is intended to be used for the qualitative detection of human MYC rearrangements with and without participation of the human IGH locus in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the FlexISH-Tissue Implementation Kit (Prod. No. Z-2182-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The FlexISH MYC/IGH TriCheck Probe is composed of:


- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 8q24.21* (chr8:130,373,051-130,930,673) distal to the MYC breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2,5 ng/μl), which target sequences mapping in 8q24.21* (chr8:127,888,765-128,363,281) proximal to the MYC breakpoint region (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~70,0 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-106,995,000) harboring the IGH locus. (see Fig. 1).

* according to Human Genome Assembly GRCh37/hg19


Results

Normal situation: In interphases of normal cells or cells without MYC rearrangements, two green/orange fusion signals appear when using an appropriate dual bandpass filter set and two blue signals appear when using an appropriate single bandpass filter set. When using an appropriate triple color bandpass filter set, two green/orange fusion signals and two blue signals can be observed.

Aberrant situation: A rearrangement of the MYC gene region not involving the IGH locus is indicated by a green/orange fusion signal, separated single green and orange signals and two blue signals. A MYC rearrangement involving the IGH locus is indicated by one green/orange fusion signal, one separate blue signal and one green/blue fusion signal as well as one orange/blue fusion signal. In a cell with a rearrangement of the IGH locus not involving the MYC gene or a trisomy of chromosome 14, two green/orange fusion signals and three blue signals can be observed.



Example of an aberrant signal pattern: Non-Hodgkin lymphoma tissue section with t(8;14) as indicated by one separate green and one separate orange signal, and one additional blue signal.



Example of an aberrant signal pattern: Non-Hodgkin lymphoma tissue section with translocation of the MYC gene without IGH involvement as indicated by one separate green and one separate orange signal, without an additional blue signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
  	8, 14	CE/IVD	5 (50 μl)	Z-2293-50







Fig. 1: Top: MYC Probe map; Bottom: IGH Probe map (not to scale)



Ideograms of chromosome 8 and 14 indicating the hybridization locations

FlexISH® RET/KIF5B TriCheck™ Probe

Intended Purpose

The FlexISH RET/KIF5B TriCheck Probe (PL226) is intended to be used for the qualitative detection of rearrangements involving the human RET gene with and without participation of the human KIF5B gene in formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the FlexISH-Tissue Implementation Kit (Prod. No. Z-2182-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The FlexISH RET/KIF5B TriCheck Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2,5 ng/μl), which target sequences mapping in 10q11.21* (chr10:43,340,888-43,510,171) proximal to the RET breakpoint region (see Fig. 1).
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 10q11.21* (chr10:43,626,274-44,112,146) distal to the RET breakpoint region (see Fig. 1).
- ZyBlue (excitation 418 nm/emission 467 nm) labeled polynucleotides (~70 ng/μl), which target sequences mapping in 10p11.22* (chr10:31,640,467-33,085,804) harboring the KIF5B gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19







Fig. 1: Top: RET Probe map; Bottom: KIF5B Probe map (not to scale); Right: Ideogram of chromosome 10 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without RET-KIF5B rearrangements, two green/orange fusion signals appear when using an appropriate dual bandpass filter set, and two blue signals appear when using an appropriate single bandpass filter set. When using an appropriate triple color bandpass filter set, two green/orange fusion signals, each in close proximity to a blue signal, can be observed, indicating the normal chromosomes 10.

Aberrant situation: A KIF5B-RET fusion caused by a RET-KIF5B inversion is indicated by a green/blue fusion signal and an orange/blue fusion signal in close proximity. A RET translocation without involvement of KIF5B is indicated by an orange and a blue signal in close proximity and a separate green signal. Isolated green signals are the result of deletions proximal to the RET breakpoint region.



FlexISH RET/KIF5B TriCheck™ Probe on normal interphase cells with non-rearranged RET loci (two green/orange fusion signals), and non-rearranged KIF5B loci (two blue signals).



Example of an aberrant signal pattern: NSCLC tissue section with a KIF5B-RET inversion as indicated by one green, one separated orange, and an additional blue signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
  	10	CE/IVD	20 (200 μl)	Z-2269-200

Specimen kindly provided by Dr. Schildhaus, Essen, Germany.

FlexISH® ERBB2/CEN 17 Dual Color Probe

Intended Purpose

The FlexISH ERBB2/CEN 17 Dual Color Probe (PL122) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene as well as the detection of chromosome 17 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer and gastric/gastroesophageal junction cancer, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the FlexISH-Tissue Implementation Kit (Prod. No. Z-2182-5/-20).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and gastric/gastroesophageal junction cancer and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The FlexISH ERBB2/CEN 17 Dual Color Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/μl), which target sequences mapping in 17q12-q21.1* (chr17:37,572,531-38,181,308) harboring the ERBB2 gene region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~1.0 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19




Fig. 1: ERBB2 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without an amplification involving the ERBB2 gene region, two green and two orange signals appear.

Aberrant situation: In cells with an amplification of the ERBB2 gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized for 2 hours on a breast cancer tissue section with ERBB2 (green) amplification.



Hybridized overnight on a breast cancer tissue section with ERBB2 (green) amplification.

Labels	Chromosome	Reg. Status	Tests	Order No.
	17	CE/IVD	5 (50 μl)	Z-2166-50
			20 (200 μl)	Z-2166-200

Accessories

FlexISH® Implementation Kits

For the detection FlexISH® probes

Product	Manufacturer	Tests	Reg. Status	Order No.
FlexISH-Tissue Implementation Kit				
Incl.: Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; 5x FlexISH Wash Buffer, 150 ml; DAPI/DuraTect-Solution, 0.2 ml	ZytoVision GmbH	5	CE/IVD	Z-2182-5
FlexISH-Tissue Implementation Kit				
Incl.: Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; 5x FlexISH Wash Buffer, 500 ml; DAPI/DuraTect-Solution, 0.8 ml	ZytoVision GmbH	10	CE/IVD	Z-2182-20

FlexISH® Pretreatment Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
Pepsin Solution	ZytoVision GmbH	4 ml	CE/IVD	ES-0001-4
Pepsin Solution	ZytoVision GmbH	50 ml	CE/IVD	ES-0001-50
Pepsin Solution	ZytoVision GmbH	1000 ml	CE/IVD	ES-0001-1000
Heat Pretreatment Solution Citric	ZytoVision GmbH	2 x 500 ml	CE/IVD	PT-0001-1000

FlexISH® Wash Buffers & Ancillary Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
DAPI/DuraTect™-Solution	ZytoVision GmbH	150 ng DAPI/ml, 0.8 ml	CE/IVD	MT-0007-0.8
DAPI/DuraTect™-Solution (ultra)	ZytoVision GmbH	1360 ng DAPI/ml, 0.8 ml	CE/IVD	MT-0008-0.8
5 x FlexISH Wash Buffer	ZytoVision GmbH	500 ml	CE/IVD	WB-0010-500

ZytoMation®

The **ZytoMation® probes** combine the quality of the ZytoVision probes for Fluorescence *in situ* Hybridization (FISH) with an automated workflow. They are designed for fully automated FISH to detect genetic aberrations such as translocations and amplifications in formalin-fixed, paraffin-embedded (FFPE) tissue sections on the BOND™ fully automated systems (BOND™-III, BOND™-MAX, and BOND™ RXm (RUO)) by Leica Biosystems (ZytoVision is not affiliated or associated with Leica Biosystems).


 ZytoVision GmbH

Advantages of ZytoMation®

- + 5 h protocol on fully automated BOND™ systems
- + Ready-to-use probes
- + Reduced hands-on time

Workflow

To successfully use the ZytoMation® probes, the BOND™ FISH Kit (DS9636) is required. Prior to evaluation, the hybridized slides should be mounted using a DAPI/DuraTect™-Solution (MT-0007-0.8/MT-0008-0.8).




Chromosome Index



Chr. Band	Product Name	Reg. Status	Order No.	Page
2p23	ZytoMation ALK Dual Color Break Apart FISH Probe	CE/IVD	Z-2315-5.1ML	129



Chr. Band	Product Name	Reg. Status	Order No.	Page
3q27	ZytoMation BCL6 Dual Color Break Apart FISH Probe	CE/IVD	Z-2313-5.1ML	130




Chr. Band	Product Name	Reg. Status	Order No.	Page
6q22.1	ZytoMation ROS1 Dual Color Break Apart FISH Probe	CE/IVD	Z-2298-5.1ML	131



Chr. Band	Product Name	Reg. Status	Order No.	Page
7q31.2	ZytoMation MET/CEN 7 Dual Color FISH Probe	CE/IVD	Z-2321-5.1ML	131



Chr. Band	Product Name	Reg. Status	Order No.	Page
8q24.21	ZytoMation MYC Dual Color Break Apart FISH Probe	CE/IVD	Z-2312-5.1ML	132



Chr. Band	Product Name	Reg. Status	Order No.	Page
10q11.2	ZytoMation RET Dual Color Break Apart FISH Probe	CE/IVD	Z-2316-5.1ML	133

Chr. Band	Product Name	Reg. Status	Order No.	Page
14q32.3	ZytoMation IGH Dual Color Break Apart FISH Probe	CE/IVD	Z-2317-5.1ML	134

Chr. Band	Product Name	Reg. Status	Order No.	Page
17q12	ZytoMation ERBB2/CEN 17 Dual Color FISH Probe	CE/IVD	Z-2292-5.1ML	134

Chr. Band	Product Name	Reg. Status	Order No.	Page
18q21.3	ZytoMation BCL2 Dual Color Break Apart FISH Probe	CE/IVD	Z-2306-5.1ML	135

Gene Index

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
ALK	CD246	ZytoMation ALK Dual Color Break Apart FISH Probe	CE/IVD	Z-2315-5.1ML	129
BCL2	Bcl-2, PPP1R50	ZytoMation BCL2 Dual Color Break Apart FISH Probe	CE/IVD	Z-2306-5.1ML	135
BCL6	ZNF51, LAZ3	ZytoMation BCL6 Dual Color Break Apart FISH Probe	CE/IVD	Z-2313-5.1ML	130
ERBB2	HER2, HER-2, NEU	ZytoMation ERBB2/CEN 17 Dual Color FISH Probe	CE/IVD	Z-2292-5.1ML	134
IGH	IGH@	ZytoMation IGH Dual Color Break Apart FISH Probe	CE/IVD	Z-2317-5.1ML	134
MET	HGFR, RCCP2	ZytoMation MET/CEN 7 Dual Color FISH Probe	CE/IVD	Z-2321-5.1ML	131
MYC	CMYC, bHLHe39, c-Myc	ZytoMation MYC Dual Color Break Apart FISH Probe	CE/IVD	Z-2312-5.1ML	132
RET	HSCR1, CDHF12	ZytoMation RET Dual Color Break Apart FISH Probe	CE/IVD	Z-2316-5.1ML	133
ROS1	MCF3, ROS	ZytoMation ROS1 Dual Color Break Apart FISH Probe	CE/IVD	Z-2298-5.1ML	131

ZytoMation® ALK Dual Color Break Apart FISH Probe

Intended Purpose
The ZytoMation ALK Dual Color Break Apart FISH Probe is designed to detect rearrangements involving the chromosomal region 2p23.1-p23.2 harboring the ALK (ALK receptor tyrosine kinase, a.k.a. CD246) gene.

ALK encodes a transmembrane receptor tyrosine kinase. This gene exerts characteristic oncogenic activities through fusion to several gene partners or mutations both in hematopoietic and non-hematopoietic solid tumors.

Translocations affecting the ALK gene locus are frequently found in anaplastic large cell lymphoma (ALCL), an aggressive non-Hodgkin lymphoma arising from T-cells. The most frequent translocation t(2;5) results in a fusion with the NPM1 gene located on chromosome 5q35.

This rearrangement results in a NPM1/ALK fusion protein, which is constitutively activated through autophosphorylation, and that in turn mediates malignant cell formation by activating downstream effectors like e.g. STAT3.

Additionally, inversions affecting the ALK gene located on the short arm of chromosome 2 [inv(2)(p21p23)] have been frequently detected in non-small cell lung cancer (NSCLC) and lead to the formation of EML4-ALK fusion transcripts.

ALK kinase targeted therapies may represent a very effective therapeutic strategy in NSCLC patients carrying EML4-ALK rearrangements.

References
Ach T, *et al.* (2013) Virchows Arch 462: 65-72.
Cooper CS, *et al.* (1984) Nature 311: 29-32.
Engelman JA, *et al.* (2007) Science 316: 1039-43.
Ettl T, *et al.* (2014) Head Neck 36: 517-23.
Garcia S, *et al.* (2007) Int J Oncol 31: 49-58.
Hara T, *et al.* (1998) Lab Invest 78: 1143-53.
Lacroix L, *et al.* (2014) PLoS One 1: e84319.
Lee D, *et al.* (2015) Cancer Res Treat 47: 120-5.
Preusser M, *et al.* (2014) Histopathology 65: 684-92.
Schildhaus HU, *et al.* (2015) Clin Cancer Res 21: 907-15.

Probe Description
The ZytoMation ALK Dual Color Break Apart FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~6 ng/µl), which target sequences mapping in 2p23.1-p23.2* (chr2:29,460,144-30,095,822) proximal to the ALK breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4,0 ng/µl), which target sequences mapping in 2p23.2* (chr2:29,174,204-29,383,335) distal to the ALK breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: ALK Probe map (not to scale)




Ideogram of chromosome 2 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the ALK gene region, two green/orange fusion signals appear.
Aberrant situation: One ALK gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. An EML4-ALK inversion with deletion of 5'-ALK sequences is indicated by one or multiple isolated orange signals.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Lung cancer tissue section with translocation of the ALK gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
 	2	CE/IVD	Up to 20 (5.1 ml)	Z-2315-5.1ML

ZytoMation® BCL6 Dual Color Break Apart FISH Probe

Intended Purpose
The ZytoMation BCL6 Dual Color Break Apart FISH Probe is designed for the detection of translocations involving the chromosomal region 3q27.3 harboring the BCL6 (BCL6 transcription repressor, a.k.a. ZNF51, LAZ3) gene. The BCL6 protein acts as a transcriptional repressor that is involved in the regulation of lymphoid development and function.

Chromosomal rearrangements of the BCL6 gene region were found to occur in different types of non-Hodgkin lymphoma (NHL), including diffuse large B-cell lymphoma (DLBCL) and follicular lymphoma (FL). The most common BCL6 translocation t(3;14)(q27;q32.3) results in the IGH-BCL6 gene fusion. In addition, more than 20 partner loci have been identified including immunoglobulin (Ig) genes but also a number of non-Ig genes. As a result of these translocations, the rearranged BCL6 gene comes under the control of the promoter of the partner gene leading to deregulated expression of BCL6.

In DLBCL, the most common histologic subtype of NHL, BCL6 translocations represent one of the most frequent cytogenetic abnormality, occurring in 20% to 40% of the cases. Several studies reported a correlation of BCL6 translocation with an inferior overall survival. Moreover, DLBCL, which are positive for both BCL6 and MYC rearrangements, have been shown to have an extremely poor prognosis.

Hence, the detection of BCL6 rearrangements by FISH may help in predicting the clinical outcome in patients with NHL.

References
Akyurek N, *et al.* (2012) Cancer 118: 4173-83.
Cady FM, *et al.* (2008) J Clin Oncol 26: 4814-9.
Ohno H (2004) Histol Histopathol 19: 637-50.
Ohno H (2006) J Clin Exp Hematop 46: 43-53.

Probe Description
The ZytoMation BCL6 Dual Color Break Apart FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~6.0 ng/μl), which target sequences mapping in 3q27.3* (chr3:186,737,897-187,403,834) proximal to the BCL6 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2,5 ng/μl), which target sequences mapping in 3q27.3-q28* (chr3:187,744,962-188,097,195) distal to the BCL6 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: BCL6 Probe map (not to scale)



Ideogram of chromosome 3 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the BCL6 gene region, two green/orange fusion signals appear.
Aberrant situation: One BCL6 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Lymphoma tissue section with translocation of the BCL6 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	3	CE/IVD	Up to 20 (5.1 ml)	Z-2313-5.1ML

ZytoMation® ROS1 Dual Color Break Apart FISH Probe

Intended Purpose
The ZytoMation ROS1 Dual Color Break Apart FISH Probe (PL251) is intended to be used for the qualitative detection of translocations involving the human ROS1 gene at 6q22.1 in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC), by fluorescence *in situ* hybridization (FISH).

The probe is intended to be used in combination with the Bond FISH Kit (DS9636) on the automated Bond-MAX or Bond III system by Leica Biosystems.

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and therapeutic measures should not be initiated based on the test result alone

Probe Description
The ZytoMation ROS1 Dual Color Break Apart FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~6,0 ng/μl), which target sequences mapping in 6q22.1* (chr6:116,912,298-117,627,255) proximal to the ROS1 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2,5 ng/μl), which target sequences mapping in 6q22.1* (chr6:116,912,298-117,627,255) distal to the ROS1 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: ROS1 Probe map (not to scale); Right: Ideogram of chromosome 6 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the ROS1 gene region, two green/orange fusion signals appear.
Aberrant situation: One ROS1 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated green signals are the result of deletions distal to the ROS1 breakpoint region or are due to unbalanced translocations affecting this chromosomal region.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Lung cancer tissue section with translocation of the ROS1 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	6	CE/IVD	Up to 20 (5.1 ml)	Z-2298-5.1ML

ZytoMation® MET/CEN 7 Dual Color FISH Probe

Intended Purpose
The ZytoMation MET/CEN 7 Dual Color FISH Probe is designed for the detection of MET gene amplifications found in a variety of human tumors.

The MET gene (a.k.a. c-Met) is located in the chromosomal region 7q31.2 and encodes a transmembrane tyrosine kinase receptor for the hepatocyte growth factor (HGF). HGF and MET play an important role in angiogenesis and tumor growth.

Activation or upregulation of MET was found in a number of carcinomas including lung, breast, colorectal, prostate, and gastric carcinomas as well as in gliomas, melanomas and some sarcomas. MET overexpression is known as a negative prognostic indicator in patients with various carcinomas, multiple myeloma, or glioma. Therefore, several inhibitors of the HGF/MET signaling pathway are being studied and developed as potent therapies to inhibit angiogenesis and tumor growth.

In addition, it was shown that MET amplification leads to resistance to gefitinib or erlotinib in lung cancer by driving ERBB3-dependent activation of the PI3K pathway.

References
Ach T, *et al.* (2013) Virchows Arch 462: 65-72.
Cooper CS, *et al.* (1984) Nature 311: 29-32.
Engelman JA, *et al.* (2007) Science 316: 1039-43.
Ettl T, *et al.* (2014) Head Neck 36: 517-23.
Garcia S, *et al.* (2007) Int J Oncol 31: 49-58.
Hara T, *et al.* (1998) Lab Invest 78: 1143-53.
Lacroix L, *et al.* (2014) PLoS One 1: e84319.
Lee D, *et al.* (2015) Cancer Res Treat 47: 120-5.
Preusser M, *et al.* (2014) Histopathology 65: 684-92.
Schildhaus HU, *et al.* (2015) Clin Cancer Res 21: 907-15.

Probe Description
The ZytoMation MET/CEN 7 Dual Color FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~5 ng/μl), which target sequences mapping in 7q31.2* (chr7:115,925,700-116,718,699) harboring the MET region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~0,2 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19




Fig. 1: MET Probe map (not to scale); Right: Ideogram of chromosome 7 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the MET gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the MET gene region an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two green and two orange signals in each nucleus.

Lung adenocarcinoma tissue section with amplification of the MET gene locus as indicated by multiple copies of the green signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	7	CE/IVD	Up to 20 (5.1 ml)	Z-2321-5.1ML

ZytoMation® MYC Dual Color Break Apart FISH Probe

Intended Purpose
The ZytoMation MYC Dual Color Break Apart FISH Probe is designed to detect translocations involving the chromosomal region 8q24.21 harboring the MYC gene.

The MYC proto-oncogene (MYC proto-oncogene, bHLH transcription factor, a.k.a. CMYC) encodes a transcription factor essential for cell growth and proliferation and is broadly implicated in tumorigenesis.

Translocations involving the MYC gene are considered cytogenetic hallmarks for Burkitt lymphoma but are also found in other types of lymphomas. The most frequent translocation involving the MYC gene region is t(8;14)(q24.21;q32.3) juxtaposing the MYC gene in 8q24.21 next to the IGH (immunoglobulin heavy chain locus) gene in 14q32.33.

Further translocations affecting the MYC gene are t(8;22)(q24.21;q11.2) and t(2;8)(p11.2;q24.21), both of which involve one of the two immunoglobulin light chain loci. All three translocations bring the MYC gene under the control of a regulatory element from one of the immunoglobulin loci resulting in constitutive overexpression of MYC.

References
Boerma EG, *et al.* (2009) Leukemia 23: 225-34.
Dalla-Favera R, *et al.* (1982) Proc Natl Acad Sci U S A 79: 6497-501.
Haralambieva E, *et al.* (2004) Genes Chromosomes Cancer 40: 10-8.
Veronese ML, *et al.* (1995) Blood 85: 2132-8.

Probe Description
The ZytoMation MYC Dual Color Break Apart FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~6.0 ng/µl), which target sequences mapping in 8q24.21* (chr8:130,373,051-130,930,673) distal to the MYC breakpoint region (see Fig.1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2.5 ng/µl), which target sequences mapping in 8q24.21* (chr8:127,888,765-128,363,281) proximal to the MYC breakpoint region (see Fig. 1).
- Hybridisierungsbuffer auf Basis von Formamid

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: MYC Probe map (not to scale)



Results
Normal situation: In interphases of normal cells or cells without a translocation involving the MYC gene region, two green/orange fusion signals appear.
Aberrant situation: One MYC gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Burkitt lymphoma tissue section with translocation of the MYC gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	8	CE/IVD	Up to 20 (5.1 ml)	Z-2312-5.1ML

ZytoMation® RET Dual Color Break Apart FISH Probe

Intended Purpose
The ZytoMation RET Dual Color Break Apart FISH Probe is designed to detect translocations involving the chromosomal region 10q11.21 harboring the RET (retproto-oncogene) gene. RET encodes a tyrosine kinase (TK) receptor. Translocations involving RET were first described in papillary thyroid carcinoma (PTC) where somatic rearrangements result in the fusion of its TK catalytic domain with an N-terminal dimerization domain encoded by various fusion partner genes.

In addition, recurrent inversions [inv(10)(p11.2q11.2)] fusing the coiled-coil domains of the kinesin family member 5B (KIF5B) gene to the RET kinase domain have been detected in lung adenocarcinoma.

The resulting KIF5B-RET fusion protein can form homodimers through the coiled-coil domains of KIF5B, causing an aberrant activation of the TK of RET, a mechanism known from KIF5B-ALK fusions which is also found in lung adenocarcinoma.

RET translocations are responsible for 1-2% of non-squamous NSCLCs. Similarly to ALK and ROS1, they are more characteristic for young non-smokers and females. This category of cancers is known to be responsive to treatment with RET tyrosine kinase inhibitors.

References
Gautschi O, *et al.* (2013) J Thorac Oncol 8: e43-4.
Imyanitov EN, *et al.* (2021) Crit Rev Oncol Hematol 157: 103194.
Ju YS, *et al.* (2012) Genome Res 22: 436-45.
Kohn T, *et al.* (2012) Nat Med 18: 375-7.
Lee SE, *et al.* (2015) Mod Pathol 28: 468-79.
Nikiforov YE (2002) Endocr Pathol 13: 3-16.
Takahashi M, *et al.* (1985) Cell 42: 581-8.
Takeuchi K, *et al.* (2012) Nat Med 18: 378-81.

Probe Description
The ZytoMation RET Dual Color Break Apart FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~6.0 ng/µl), which target sequences mapping in 10q11.21* (chr10:43,626,274-43,902,346) distal to the RET breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.0 ng/µl), which target sequences mapping in 10q11.21* (chr10:43,340,888-43,510.171) proximal to the RET breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: RET Probe map (not to scale)



Ideogram of chromosome 10 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the RET gene region, two green/orange fusion signals appear.
Aberrant situation: One RET gene region affected by a translocation is indicated by one separate green signal and one separate orange signal. Isolated green signals are the result of deletions proximal to the RET breakpoint region. Inversions involving genes in close proximity such as CCDC6 might not be observable due to just subtle signal separation.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Lung adenocarcinoma tissue section with rearrangement of the RET gene as indicated by isolated green signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
	10	CE/IVD	Up to 20 (5.1 ml)	Z-2316-5.1ML

ZytoMation® IGH Dual Color Break Apart FISH Probe

Intended Purpose
The ZytoMation IGH Dual Color Break Apart FISH Probe is designed to detect translocations involving the chromosomal region 14q32.33 harboring the IGH locus. Rearrangements involving the IGH (immunoglobulin heavy locus, a.k.a. IGH@) gene locus are considered to be cytogenetic hallmarks for non-Hodgkin lymphoma (NHL). NHLs represent 50% of all hematological malignancies. IGH locus rearrangements have been identified in about 50% of NHLs and are associated with specific subtypes of NHLs.

Translocation t(11;14)(q13.3;q32.3) can be found in about 95% of mantle cell lymphoma (MCL), t(14;18)(q32.3;q21.3) in 80% of follicular lymphoma (FL), t(3;14)(q27;q32.3) in diffuse large B-cell lymphoma (DLBCL), and t(8;14)(q24.21;q32.3) in Burkitt lymphoma. In all of these translocations an oncogene located near the breakpoint of the translocation partner is activated by juxtaposing to IGH regulatory sequences.

Rearrangements involving 14q32.33 have unique biological characteristics and correlate with clinical, morphological, and immunophenotypic features. Fluorescence *in situ* Hybridization is a helpful tool for the diagnosis, selecting treatment, and giving prognostic information.

References
Bernicot I, *et al.* (2007) Cytogenet Genome Res 118: 345-52.
Hohne S, *et al.* (2012) Pathol Res Pract 208: 510-7.
Lu S, *et al.* (2004) Cancer Genet and Cytogenet 152: 141-5.
Nishida K, *et al.* (1997) Blood 90: 526-34.
Quintero-Rivera F, *et al.* (2009) Cancer Genet and Cytogenet 190: 33-9.

Probe Description
The ZytoMation IGH Dual Color Break Apart FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~6.0 ng/µl), which target sequences mapping in 14q32.33* (chr14:106,690,778-107,268,412) distal to the IGH breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.0 ng/µl), which target sequences mapping in 14q32.33* (chr14:105,296,741-105,909,611) proximal to the IGH breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: IGH Probe map (not to scale); Right: Ideogram of chromosome 14 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the IGH locus, two green/orange fusion signals appear.
Aberrant situation: One IGH locus affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Burkitt lymphoma tissue section with translocation of the IGH gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	14	CE/IVD	Up to 20 (5.1 ml)	Z-2317-5.1ML

ZytoMation® ERBB2/CEN 17 Dual Color FISH Probe

Intended Purpose
The ZytoMation ERBB2/CEN 17 Dual Color FISH Probe (PL246) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene as well as the detection of chromosome 17 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer and gastric/gastroesophageal junction cancer, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the Bond FISH Kit (DS9636) on the automated Bond-MAX or Bond III system by Leica Biosystems.

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and gastric/gastroesophageal junction cancer and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoMation ERBB2/CEN 17 Dual Color FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/µl), which target sequences mapping in 17q12-q21.1* (chr17:37,572,531-38,181,308) harboring the ERBB2 gene region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~0.2 ng/µl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: ERBB2 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the ERBB2 gene region, two green signals and two orange signals appear.
Aberrant situation: In cells with an amplification of the ERBB2 gene region an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two green and two orange signals in each nucleus.

Breast cancer tissue section with amplification of the ERBB2 gene locus as indicated by multiple copies of the green signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
	17	CE/IVD	Up to 20 (5.1 ml)	Z-2292-5.1ML

ZytoMation® BCL2 Dual Color Break Apart FISH Probe

Intended Purpose
The ZytoMation BCL2 Dual Color Break Apart FISH Probe (PL260) is intended to be used for the qualitative detection of translocations involving the human BCL2 gene at 18q21.33 in formalin-fixed, paraffin-embedded specimens, such as B-cell lymphoma, by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the Bond FISH Kit (DS9636) on the automated Bond-MAX or Bond III system by Leica Biosystems.

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of B-cell lymphoma and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoMation BCL2 Dual Color Break Apart FISH Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~6.0 ng/µl), which target sequences mapping in 18q21.33* (chr18:60,046,152-60,589,273) proximal to the BCL2 breakpoint region (see Fig. 1).
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~2.5 ng/µl), which target sequences mapping in 18q21.33-q22.1* (chr18:60,994,528-61,658,503) distal to the BCL2 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: BCL2 Probe map (not to scale); Right: Ideogram of chromosome 18 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the BCL2 gene region, two green/orange fusion signals appear.
Aberrant situation: One BCL2 gene region affected by a translocation is indicated by one separate green signal and one separate orange signal.



Hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.

Follicular lymphoma tissue section with translocation of the BCL2 gene as indicated by one non-rearranged orange/green fusion signal, one orange and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
	18	CE/IVD	Up to 20 (5.1 ml)	Z-2306-5.1ML

ZytoDot®

Reliable and Simple Detection of Genomic Alterations using Light Microscopy!

The **ZytoDot® products** are designed for the detection of aneuploidies and gene amplifications by Chromogenic *in situ* Hybridization (CISH) in formalin-fixed, paraffin-embedded (FFPE) tissue sections.

ZytoVision GmbH


Advantages of CISH

- + Simultaneous observation of tissue morphology and CISH signals
- + Storage of slides at room temperature - CISH signals are permanent
- + No costly fluorescent microscope needed

ZytoDot® Kits – Convenient Solutions

For making CISH analysis reliable and user-friendly, all ZytoDot® CISH probes can be combined with the ZytoDot® CISH Implementation Kit (C-3018-40) which includes all necessary pretreatment solutions, wash buffers, antibodies, chromogenic substrates, counterstaining solution, mounting solution and a detailed protocol to perform successful CISH experiments.

Method



The ZytoDot® system uses Digoxigenin (DIG)-labeled probes [1], which are, after blocking [2], detected using a Mouse-anti-DIG antibody [3]. This antibody is detected by a polymerized HRP-Goat-anti-Mouse antibody [4]. The enzymatic reaction of DAB [5] leads to the formation of strong permanent brown signals that can be visualized by light microscopy using a 40x objective.

Protocol Overview



ZytoDot® 2C™

2-Color CISH for the Detection of Genomic Alterations

The **ZytoDot® 2C™ products** are designed for the simultaneous detection of two different genomic targets by Chromogenic *in situ* Hybridization (CISH) in formalin-fixed, paraffin-embedded (FFPE) tissue sections. This two color system is especially useful for the differentiation of aneuploidies from gene amplifications, and the detection of deletions and translocations.

ZytoVision GmbH


Advantages of ZytoDot® 2C™

- + Simultaneous observation of tissue morphology and CISH signals at 40x using light microscopy
- + Two targets detected simultaneously
- + High contrasting distinct red and green signals
- + Standardized and complete kits
- + No costly fluorescent microscope needed

ZytoDot® 2C™ Kits – Standardized Solutions


For making CISH analysis reliable and user-friendly, any ZytoDot® 2C™ probe can be combined with the ZytoDot® 2C™ CISH Implementation Kit resulting in target specific kit solutions. The ZytoDot® 2C™ CISH Implementation Kit includes all necessary pretreatment solutions, wash buffers, antibodies, chromogenic substrates, counterstaining and mounting solution, and a detailed protocol.

Method



The ZytoDot® 2C™ system uses DIG-and DNP-labeled probe cocktails targeting different genomic sections [1] which are detected using a Mouse-anti-DIG/Rabbit-anti-DNP cocktail [2]. These antibodies are detected by a unique cocktail of polymerized HRP-Goat-anti-Mouse/AP-Goat-anti-Rabbit antibodies [3]. The enzymatic reaction of AP-Red [4] and HRP-Green [5] leads to the formation of strong permanent red respectively green signals that can be visualized by light microscopy using a 40x objective.


Protocol Overview



Chromosome Index




Chr. Band	Product Name	Reg. Status	Order No.	Page
1p36.3	ZytoDot 2C Glioma 1p/19q Probe Set	CE/IVD	C-3076-10/-40	142
	ZytoDot 2C SPEC 1p36/1q25 Probe	CE/IVD	C-3036-100/-400	143
1q25.3	ZytoDot 2C Glioma 1p/19q Probe Set	CE/IVD	C-3076-10/-40	142
	ZytoDot 2C SPEC 1p36/1q25 Probe	CE/IVD	C-3036-100/-400	143



Chr. Band	Product Name	Reg. Status	Order No.	Page
2p24	ZytoDot SPEC MYCN Probe	CE/IVD	C-3029-400	144
2p23	ZytoDot 2C SPEC ALK Break Apart Probe	CE/IVD	C-3055-100/-400	144




Chr. Band	Product Name	Reg. Status	Order No.	Page
3q27	ZytoDot 2C SPEC BCL6 Break Apart Probe	CE/IVD	C-3074-100	145



Chr. Band	Product Name	Reg. Status	Order No.	Page
6q22.1	ZytoDot 2C SPEC ROS1 Break Apart Probe	CE/IVD	C-3063-100/-400	145




Chr. Band	Product Name	Reg. Status	Order No.	Page
7p11.2	ZytoDot SPEC EGFR Probe	RUO	C-3007-400	155
	ZytoDot 2C SPEC EGFR/CEN 7 Probe	CE/IVD	C-3033-100/-400	146
7p11.1-q11.1	ZytoDot CEN 7 Probe	RUO	C-3008-400	156
7q31.2	ZytoDot 2C SPEC MET/CEN 7 Probe	CE/IVD	C-3057-400	146




Chr. Band	Product Name	Reg. Status	Order No.	Page
8p11.1-q11.1	ZytoDot CEN 8 Probe	RUO	C-3016-400	156
8q24.21	ZytoDot SPEC MYC Probe	RUO	C-3013-400	155
	ZytoDot 2C SPEC MYC Break Apart Probe	CE/IVD	C-3066-400	147




Chr. Band	Product Name	Reg. Status	Order No.	Page
9p21	ZytoDot 2C SPEC CDKN2A/CEN 9 Probe	CE/IVD	C-3067-400	147




Chr. Band	Product Name	Reg. Status	Order No.	Page
10q11.2	ZytoDot 2C SPEC RET Break Apart Probe	CE/IVD	C-3064-100	148
10q23.3	ZytoDot 2C SPEC PTEN/CEN 10 Probe	RUO	C-3053-400	155
10q26.1	ZytoDot 2C SPEC FGFR2/CEN 10 Probe	CE/IVD	C-3056-400	148




Chr. Band	Product Name	Reg. Status	Order No.	Page
11q13.3	ZytoDot 2C SPEC CCND1 Break Apart Probe	RUO	C-3075-100	149




Chr. Band	Product Name	Reg. Status	Order No.	Page
12q13.3	ZytoDot 2C SPEC DDIT3 Break Apart Probe	CE/IVD	C-3047-100	149
12q14	ZytoDot 2C SPEC CDK4/CEN 12 Probe	CE/IVD	C-3062-400	150
12q15	ZytoDot SPEC MDM2 Probe	CE/IVD	C-3012-400	150
	ZytoDot 2C SPEC MDM2/CEN 12 Probe	CE/IVD	C-3049-100/-400	151




Chr. Band	Product Name	Reg. Status	Order No.	Page
13q14.1	ZytoDot 2C SPEC FOXO1 Break Apart Probe	RUO	C-3065-100	155




Chr. Band	Product Name	Reg. Status	Order No.	Page
14q32.3	ZytoDot 2C SPEC IGH Break Apart Probe	CE/IVD	C-3071-100	139



Chr. Band	Product Name	Reg. Status	Order No.	Page
17p13	ZytoDot 2C SPEC USP6 Break Apart Probe	CE/IVD	C-3077-100	152
17q12	ZytoDot SPEC ERBB2 Probe	CE/IVD	C-3001-400	152
	ZytoDot SPEC ERBB2 Probe Kit	CE/IVD	C-3003-40	152
	ZytoDot 2C SPEC ERBB2/CEN 17 Probe	CE/IVD	C-3032-100/-400	153
	ZytoDot 2C SPEC ERBB2/CEN 17 Probe Kit	CE/IVD	C-3022-10/-40	153



Chr. Band	Product Name	Reg. Status	Order No.	Page
18q11.2	ZytoDot 2C SPEC SS18 Break Apart Probe	CE/IVD	C-3046-100	153
18q21.3	ZytoDot 2C SPEC BCL2 Break Apart Probe	CE/IVD	C-3073-100	154
	ZytoDot 2C SPEC MALT1 Break Apart Probe	RUO	C-3072-100	154



Chr. Band	Product Name	Reg. Status	Order No.	Page
19p13.3	ZytoDot 2C Glioma 1p/19q Probe Set	CE/IVD	C-3076-10/-40	142
	ZytoDot 2C SPEC 19q13/19p13 Probe	CE/IVD	C-3037-100/-400	143
19q13.3	ZytoDot 2C Glioma 1p/19q Probe Set	CE/IVD	C-3076-10/-40	142
	ZytoDot 2C SPEC 19q13/19p13 Probe	CE/IVD	C-3037-100/-400	143



Chr. Band	Product Name	Reg. Status	Order No.	Page
22q12.2	ZytoDot 2C SPEC EWSR1 Break Apart Probe	CE/IVD	C-3043-100	154



Chr. Band	Product Name	Reg. Status	Order No.	Page
Xp11.1-q11.1	ZytoDot CEN X Probe	RUO	C-3025-400	156
	ZytoDot 2C CEN X/Y Probe	RUO	C-3048-400	156



Chr. Band	Product Name	Reg. Status	Order No.	Page
Yp11.1-q11.1	ZytoDot 2C CEN X/Y Probe	RUO	C-3048-400	156
Yq12	ZytoDot CEN Yq12 Probe	RUO	C-3020-400	156

Gene Index

HUGO Name	Synonym	Product Name	Reg. Status	Order No.	Page
ALK	CD246	ZytoDot 2C SPEC ALK Break Apart Probe	CE/IVD	C-3055-100/-400	144
BCL2	Bcl-2, PPP1R50	ZytoDot 2C SPEC BCL2 Break Apart Probe	CE/IVD	C-3073-100	154
BCL6	ZNF51, LAZ3	ZytoDot 2C SPEC BCL6 Break Apart Probe	CE/IVD	C-3074-100	145
CCND1	BCL1, PRAD1	ZytoDot 2C SPEC CCND1 Break Apart Probe	RUO	C-3075-100	149
CDK4	PSK-J3	ZytoDot 2C SPEC CDK4/CEN 12 Probe	CE/IVD	C-3062-400	150
CDKN2A	p16, ARF, INK4	ZytoDot 2C SPEC CDKN2A/CEN 9 Probe	CE/IVD	C-3067-400	147
DDIT3	CHOP, GADD153	ZytoDot 2C SPEC DDIT3 Break Apart Probe	CE/IVD	C-3047-100	149
EGFR	HER1, ERBB1	ZytoDot SPEC EGFR Probe	RUO	C-3007-400	152
		ZytoDot 2C SPEC EGFR/CEN 7 Probe	CE/IVD	C-3033-100/-400	138
ERBB2	HER2, HER-2, NEU	ZytoDot SPEC ERBB2 Probe	CE/IVD	C-3001-400	152
		ZytoDot SPEC ERBB2 Probe Kit	CE/IVD	C-3003-40	152
		ZytoDot 2C SPEC ERBB2/CEN 17 Probe	CE/IVD	C-3032-100/-400	153
		ZytoDot 2C SPEC ERBB2/CEN 17 Probe Kit	CE/IVD	C-3022-10/-40	153
EWSR1	EWS	ZytoDot 2C SPEC EWSR1 Break Apart Probe	CE/IVD	C-3043-100	154
FGFR2	BEK, CD332	ZytoDot 2C SPEC FGFR2/CEN 10 Probe	CE/IVD	C-3056-400	148
FOXO1	FKHR, FKH1	ZytoDot 2C SPEC FOXO1 Break Apart Probe	RUO	C-3065-100	155
IGH	IGH@	ZytoDot 2C SPEC IGH Break Apart Probe	CE/IVD	C-3071-100	151
MALT1	MLT	ZytoDot 2C SPEC MALT1 Break Apart Probe	RUO	C-3072-100	155
MDM2	HDM2	ZytoDot SPEC MDM2 Probe	CE/IVD	C-3012-400	150
		ZytoDot 2C SPEC MDM2/CEN 12 Probe	CE/IVD	C-3049-100/-400	151
MET	HGFR, RCCP2	ZytoDot 2C SPEC MET/CEN 7 Probe	CE/IVD	C-3057-400	146
MYC	CMYC, bHLHe39, c-Myc	ZytoDot SPEC MYC Probe	RUO	C-3013-400	155
		ZytoDot 2C SPEC MYC Break Apart Probe	CE/IVD	C-3066-400	147
MYCN	NMYC, N-myc	ZytoDot SPEC MYCN Probe	CE/IVD	C-3029-400	144
PTEN	MMAC1, TEP1	ZytoDot 2C SPEC PTEN/CEN 10 Probe	RUO	C-3053-400	155
RET	HSCR1, CDHF12	ZytoDot 2C SPEC RET Break Apart Probe	CE/IVD	C-3064-100	148
ROS1	MCF3, ROS	ZytoDot 2C SPEC ROS1 Break Apart Probe	CE/IVD	C-3063-100/-400	145
SS18	SYT, SSXT	ZytoDot 2C SPEC SS18 Break Apart Probe	CE/IVD	C-3046-100	153
USP6	Tre-2, TRE17	ZytoDot 2C SPEC USP6 Break Apart Probe	CE/IVD	C-3077-100	152

ZytoDot® 2C Glioma 1p/19q Probe Set

Intended Purpose

The ZytoDot 2C Glioma 1p/19q Probe Set is intended to be used for the qualitative detection of deletions involving the human chromosomal region 1p36.31 as well as deletions involving the human chromosomal region 19q13.32-q13.33 in formalin-fixed, paraffin-embedded specimens, such as glioma, by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of glioma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoDot 2C Glioma 1p/19q Probe Set is a set comprising two separate probes:

- ZytoDot 2C SPEC 1p36/1q25 Probe (Prod. No. C-3036-100/-400)
- ZytoDot 2C SPEC 19q13/19p13 Probe (Prod. No. C-3037-100/-400)

The ZytoDot 2C SPEC 1p36/1q25 Probe (PD21) is composed of:

- Digoxigenin-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 1q25.3* (chr1:184,562,510-184,752,938) (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 1p36.31* (chr1:5,808,946-6,176,336) (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19




Fig. 1: Top: SPEC 1p36 Probe map; Bottom: SPEC 1q25 Probe map (not to scale)

The ZytoDot 2C SPEC 19q13/19p13 Probe (PD22) is composed of:

- Digoxigenin-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 19p13.3* (chr19:815,938-962,244) (see Fig. 2).
- Dinitrophenyl-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 19q13.32-q13.33* (chr19:47,857,776-48,339,398) (see Fig. 2).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 2: Top: SPEC 19p13 Probe map; Bottom: SPEC 19q13 Probe map (not to scale)

Results

Normal situation: Using the SPEC 1p36/1q25 Probe or the SPEC 19q13/19p13 Probe in a normal interphase nucleus, two red and two green signals are expected.

Aberrant situation: In a cell with deletions affecting the 1p36 or 19q13 locus, one or no copy of the red signal will be observed.



SPEC 1p36/1q25 Probe hybridized to glioma tissue section with 1p36 deletion as indicated by one red signal in each nucleus.

SPEC 19q13/19p13 Dual Color Probe hybridized to glioma tissue section with 19q13 deletion as indicated by one red signal in each nucleus.

Labels	Contents	Reg. Status	Tests	Order No.
DNP/DIG	ZytoDot 2C SPEC 1p36/1q25 Probe, 0.1 ml; ZytoDot 2C SPEC 19q13/19p13 Probe, 0.1 ml	CE/IVD	10	C-3076-10
	ZytoDot 2C SPEC 1p36/1q25 Probe, 0.4 ml; ZytoDot 2C SPEC 19q13/19p13 Probe, 0.4 ml		40	C-3076-40

Images kindly provided by Prof. W. Müller, University Leipzig, Germany.

ZytoDot® 2C SPEC 1p36/1q25 Probe

Intended Purpose

The ZytoDot 2C SPEC 1p36/1q25 Probe (PD21) is intended to be used for the qualitative detection of deletions involving the human chromosomal region 1p36.31 as well as chromosome 1q25.3 specific sequences in formalin-fixed, paraffin-embedded specimens, such as glioma, by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of glioma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoDot 2C SPEC 1p36/1q25 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 1q25.3* (chr1:184,562,510-184,752,938) (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 1p36.31* (chr1:5,808,946-6,176,336) (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC 1p36 Probe map; Bottom: SPEC 1q25 Probe map (not to scale); Right: Ideogram of chromosome 1 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without deletion involving the 1p36 locus, two distinct dot-shaped red and two distinct dot-shaped green signals appear.

Aberrant situation: In a cell with deletion affecting the 1p36 locus, a reduced number of red signals will be observed. Deletions affecting only parts of the 1p36 locus might result in a normal signal pattern with red signals of reduced size.



Hybridized to glioma tissue section with 1p36 deletion as indicated by one red signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
DNP/DIG	1	CE/IVD	10 (100 μl)	C-3036-100
			40 (400 μl)	C-3036-400

Image kindly provided by Prof. W. Müller, University Leipzig, Germany.

ZytoDot® 2C SPEC 19q13/19p13 Probe

Intended Purpose

The ZytoDot 2C SPEC 19q13/19p13 Probe (PD22) is intended to be used for the qualitative detection of deletions involving the human chromosomal region 19q13.32-q13.33 as well as chromosome 19p13.3 specific sequences in formalin-fixed, paraffin-embedded specimens, such as glioma, by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of glioma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoDot 2C SPEC 19q13/19p13 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 19p13.3* (chr19:815,938-962,244) (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1.7 ng/μl), which target sequences mapping in 19q13.32-q13.33* (chr19:47,857,776-48,339,398) (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: Top: SPEC 19p13 Probe map; Bottom: SPEC 19q13 Probe map (not to scale); Right: Ideogram of chromosome 19 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without deletion involving the 19q13 locus, two distinct dot-shaped red and two distinct dot-shaped green signals appear.

Aberrant situation: In a cell with deletion affecting the 19q13 locus, a reduced number of red signals will be observed. Deletions affecting only parts of the 19q13 locus might result in a normal signal pattern with red signals of reduced size.



Hybridized to glioma tissue section with 19q13 deletion as indicated by one red signal in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
DNP/DIG	19	CE/IVD	10 (100 μl)	C-3037-100
			40 (400 μl)	C-3037-400

Image kindly provided by Prof. W. Müller, University Leipzig, Germany.

ZytoDot® SPEC MYCN Probe

Intended Purpose

The ZytoDot SPEC MYCN Probe (PD17) is intended to be used for the qualitative detection of amplifications involving the human MYCN gene in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot CISH Implementation Kit (Prod. No. C-3018-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoDot SPEC MYCN Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1,8 ng/μl), which target sequences mapping in 2p24.3* (chr2:15,846,046-16,213,717) harboring the MYCN gene region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC MYCN Probe map (not to scale); Right: Ideogram of chromosome 2 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without an amplification involving the MYCN gene region, two distinct dot-shaped brown signals appear.

Aberrant situation: In cells with an amplification of the MYCN gene region or polyploidy of chromosome 2, an increased number of brown signals or brown signal clusters will be observed.



Normal nuclei each with two MYCN signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG	2	CE/IVD	40 (400 μl)	C-3029-400

ZytoDot® 2C SPEC ALK Break Apart Probe

Intended Purpose

The ZytoDot 2C SPEC ALK Break Apart Probe (PD35) is intended to be used for the qualitative detection of translocations involving the human ALK gene at 2p23.2 in formalin-fixed, paraffin-embedded specimens, such as non-small cell lung cancer (NSCLC), by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of NSCLC and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoDot 2C SPEC ALK Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 2p23.2* (chr2:29,460,144-29,681,581) proximal to the ALK breakpoint region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~0,75 ng/μl), which target sequences mapping in 2p23.2* (chr2:29,174,204-29,383,335) distal to the ALK breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ALK Probe map (not to scale); Right: Ideogram of chromosome 2 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the ALK gene region, two red/green fusion signals appear.

Aberrant situation: One ALK gene region affected by a translocation is indicated by one separate green signal and one separate red signal. EML4-ALK inversion with deletion of 5'-ALK sequences is indicated by one or multiple isolated red signals.



Lung carcinoma tissue section with translocation affecting the 2p23.2 locus as indicated by one red/green fusion (non-rearranged) signal, one red signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	2	CE/IVD	10 (100 μl)	C-3055-100
			40 (400 μl)	C-3055-400

ZytoDot® 2C SPEC BCL6 Break Apart Probe

Intended Purpose

The ZytoDot 2C SPEC BCL6 Break Apart Probe (PD54) is intended to be used for the qualitative detection of translocations involving the human BCL6 gene at 3q27.3 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoDot 2C SPEC BCL6 Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 3q27.3* (chr3:187,028,236-187,403,834) proximal to the BCL6 breakpoint region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~0,75 ng/μl), which target sequences mapping in 3q27.3-q28* (chr3:187,744,962-188,097,195) distal to the BCL6 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC BCL6 Probe map (not to scale); Right: Ideogram of chromosome 3 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the BCL6 gene region, two red/green fusion signals appear.

Aberrant situation: One BCL6 gene region affected by a translocation is indicated by one separate green signal and one separate red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	3	CE/IVD	10 (100 μl)	C-3074-100

ZytoDot® 2C SPEC ROS1 Break Apart Probe

Intended Purpose

The ZytoDot 2C SPEC ROS1 Break Apart Probe (PD43) is intended to be used for the qualitative detection of translocations involving the human ROS1 gene at 6q22.1 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoDot 2C SPEC ROS1 Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 6q22.1* (chr6:117,448,964-117,627,255) proximal to the ROS1 breakpoint region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~0,75 ng/μl), which target sequences mapping in 6q22.1* (chr6:117,659,135-117,871,701) distal to the ROS1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC ROS1 Probe map (not to scale); Right: Ideogram of chromosome 6 indicating the hybridization locations


Results

Normal situation: In interphases of normal cells or cells without a translocation involving the ROS1 gene region, two red/green fusion signals appear.

Aberrant situation: One ROS1 gene region affected by a translocation is indicated by one separate green signal and one separate red signal. Isolated green signals are the result of deletions distal to the ROS1 breakpoint region or are due to unbalanced translocations affecting this chromosomal region.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.



Example of an aberrant signal pattern: Lung cancer tissue section with rearrangement of the ROS1 gene as indicated by isolated green signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	6	CE/IVD	10 (100 μl)	C-3063-100
			40 (400 μl)	C-3063-400

ZytoDot® 2C SPEC EGFR/CEN 7 Probe

Intended Purpose
The ZytoDot 2C SPEC EGFR/CEN 7 Probe (PD18) is intended to be used for the qualitative detection of amplifications involving the human EGFR gene as well as the detection of chromosome 7 alpha satellites in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC EGFR/CEN 7 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1,1 ng/μl), which target sequences mapping in 7p11.2* (chr7:55,082,262-55,278,647) harboring the EGFR region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1,1 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC EGFR Probe map (not to scale); Right: Ideogram of chromosome 7 indicating the hybridization locations

Results
Normal situation: In Interphases of normal cells or cells without an amplification involving the EGFR gene region, two distinct dot-shaped green and two distinct dot-shaped red signals appear.
Aberrant situation: In cells with an amplification of the EGFR gene region, an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	7	CE/IVD	10 (100 μl)	C-3033-100
			40 (400 μl)	C-3033-400

ZytoDot® 2C SPEC MET/CEN 7 Probe

Intended Purpose
The ZytoDot 2C SPEC MET/CEN 7 Probe (PD37) is intended to be used for the qualitative detection of amplifications involving the human MET gene as well as the detection of chromosome 7 alpha satellites in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC MET/CEN 7 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 7q31.2* (chr7:116,298,989-116,718,699) harboring the MET region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19




Fig. 1: SPEC MET Probe map (not to scale); Right: Ideogram of chromosome 7 indicating the hybridization locations

Results
Normal situation: In Interphases of normal cells or cells without an amplification involving the MET gene region, two distinct dot-shaped green and two distinct dot-shaped red signals appear.
Aberrant situation: In cells with an amplification of the MET gene region, an increased number of green signals or green signal clusters will be observed.



Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	7	CE/IVD	40 (400 μl)	C-3057-400

ZytoDot® 2C SPEC MYC Break Apart Probe

Intended Purpose
The ZytoDot 2C SPEC MYC Break Apart Probe (PD46) is intended to be used for the qualitative detection of translocations involving the human MYC gene at 8q24.21 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC MYC Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 8q24.21* (chr8:130,373,051-130,847,951) distal to the MYC breakpoint region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~0,75 ng/μl), which target sequences mapping in 8q24.21* (chr8:127,888,765-128,363,281) proximal to the MYC breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC MYC Probe map (not to scale); Right: Ideogram of chromosome 8 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the MYC gene region, two red/green fusion signals appear.
Aberrant situation: One MYC gene region affected by a translocation is indicated by one separate distinct dot-shaped green signal and one separate distinct dot-shaped red signal.



Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	8	CE/IVD	40 (400 μl)	C-3066-400

ZytoDot® 2C SPEC CDKN2A/CEN 9 Probe

Intended Purpose
The ZytoDot 2C SPEC CDKN2A/CEN 9 Probe is designed for the detection of CDKN2A deletions frequently observed in most tumor cell lines as well as in primary human malignancies.

The CDKN2A gene, often referred to as p16 or INK4a/ARF, is located in the chromosomal region 9p21.3. Using alternative first exons and an alternative reading frame, the gene encodes for two distinct tumor suppressor proteins p16INK4a and p14ARF, both involved in cell cycle regulation. CDKN2A has been identified as a major susceptibility gene for melanoma.

The tumor suppressor gene CDKN2A is inactivated by homozygous deletions with high frequency in a variety of human primary tumors e.g. bladder and renal cell carcinoma, prostate and ovarian adenocarcinoma, non-small cell lung cancer, sarcoma, glioma, mesothelioma, and melanoma. Furthermore, deletion of the CDKN2A gene is found in up to 80% of T-cell acute lymphoblastic leukemia cases and is associated with poor prognosis and relapse of the disease.

References
Arif Q & Husain AN (2015) Arch Pathol Lab Med 139: 978-80.
Cowan JM *et al.* (1988) J Natl Cancer Inst 80: 1159-64.
Holley T, *et al.* (2012) PLoS One 7: e50586.
Hussussian CJ, *et al.* (1994) Nat Genet 8: 15-21.
Kamb A, *et al.* (1994) Science 264: 436-40.
Nobori T, *et al.* (1994) Nature 368: 753-6.
Quelle DE, *et al.* (1995) Cell 83: 993-1000.
Rocco JW & Sidransky D (2001) Exp Cell Res 264: 42-55.
Schopmeyer K, *et al.* (1999) Neoplasia 1: 128-37.
Schwarz S, *et al.* (2008) Cytometry A 73: 305-11.
Sharpless NE (2005) Mutat Res 576: 22-38.

Probe Description
The ZytoDot 2C SPEC CDKN2A/CEN 9 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 9p21.3* (chr9:21,742,629-22,056,853) harboring the CDKN2A gene region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in chromosomal region 9q12 specific for the classical satellite III region D9Z3 of chromosome 9.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC CDKN2A Probe map (not to scale); Right: Ideogram of chromosome 9 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a deletion involving the CDKN2A gene region, two distinct dot-shaped green and two distinct dot-shaped red signals appear.
Aberrant situation: In a cell with deletion affecting the CDKN2A gene region, a reduced number of green signals will be observed. Deletions affecting only parts of the CDKN2A gene region might result in a normal signal pattern with green signals of reduced size.



Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	9	CE/IVD	40 (400 μl)	C-3067-400

ZytoDot® 2C SPEC RET Break Apart Probe

Intended Purpose
The ZytoDot 2C SPEC RET Break Apart Probe (PD44) is intended to be used for the qualitative detection of translocations involving the human RET gene at 10q11.21 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC RET Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.5 ng/μl), which target sequences mapping in 10q11.21* (chr10:43,687,278-43,856,587) distal to the RET breakpoint region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~0.75 ng/μl), which target sequences mapping in 10q11.21* (chr10:43,340,888-43,510,171) proximal to the RET breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC RET Probe map (not to scale); Right: Ideogram of chromosome 10 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the RET gene region, two red/green fusion signals appear.
Aberrant situation: One RET gene region affected by a translocation is indicated by one separate green signal and one separate red signal. Isolated green signals are the result of deletions proximal to the RET breakpoint region.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.

ZytoDot® 2C SPEC FGFR2/CEN 10 Probe

Intended Purpose
The ZytoDot 2C SPEC FGFR2/CEN 10 Probe (PD36) is intended to be used for the qualitative detection of amplifications involving the human FGFR2 gene as well as the detection of chromosome 10 alpha satellites in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC FGFR2/CEN 10 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 10q26.12-10q26.13* (chr10:123,080,085-123,492,398) harboring the FGFR2 region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 10p11.1-q11.1 specific for the alpha satellite centromeric region D10Z1 of chromosome 10.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC FGFR2 Probe map (not to scale); Right: Ideogram of Chromosome 10 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the FGFR2 gene region, two distinct dot-shaped green and two distinct dot-shaped red signals appear.
Aberrant situation: In cells with an amplification of the FGFR2 gene region, an increased number of green signals or green signal clusters will be observed.



Example of an aberrant signal pattern: Breast carcinoma tissue section with FGFR2 (green) amplification.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	10	CE/IVD	10 (100 μl)	C-3064-100

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	10	CE/IVD	40 (400 μl)	C-3056-400

ZytoDot® 2C SPEC CCND1 Break Apart Probe

Intended Purpose
The ZytoDot 2C SPEC CCND1 Break Apart Probe (PD55) is intended to be used for the qualitative detection of translocations involving the human CCND1 gene at 11q13.3 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

Probe Description
The ZytoDot 2C SPEC CCND1 Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 11q13.3* (chr11:68,522,105-68,705,283) proximal to the CCND1 breakpoint region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~0.75 ng/μl), which target sequences mapping in 11q13.3* (chr11:69,453,301-70,031,240) distal to the CCND1 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC CCND1 Probe map (not to scale); Right: Ideogram of chromosome 11 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the CCND1 gene region, two red/green fusion signals appear.
Aberrant situation: One CCND1 gene region affected by a translocation is indicated by one separate green signal and one separate red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.



Example of an aberrant signal pattern: Mantle cell lymphoma tissue section with translocation affecting the 11q13.3 locus as indicated by one non-rearranged red/green fusion signal, one red signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	11	RUO	10 (100 μl)	C-3075-100

ZytoDot® 2C SPEC DDIT3 Break Apart Probe

Intended Purpose
The ZytoDot 2C SPEC DDIT3 Break Apart Probe (PD27) is intended to be used for the qualitative detection of translocations involving the human DDIT3 gene at 12q13.3 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC DDIT3 Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 12q13.3-q14.1* (chr12:58,024,366-58,486,511) distal to the DDIT3 breakpoint region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~0.75 ng/μl), which target sequences mapping in 12q13.3* (chr12:57,386,302-57,865,800) proximal to the DDIT3 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC DDIT3 Probe map (not to scale); Right: Ideogram of chromosome 12 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the DDIT3 gene region, two red/green fusion signals appear.
Aberrant situation: One DDIT3 gene region affected by a translocation is indicated by one separate distinct dot-shaped green signal and one separate distinct dot-shaped red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.



Example of an aberrant signal pattern: Myxoid liposarcoma tissue section with translocation affecting the 12q13.3-q14.1 locus as indicated by one non-rearranged red/green fusion signal, one red signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	12	CE/IVD	10 (100 μl)	C-3047-100

ZytoDot® 2C SPEC CDK4/CEN 12 Probe


Intended Purpose
The ZytoDot 2C SPEC CDK4/CEN 12 Probe (PD42) is intended to be used for the qualitative detection of amplifications involving the human CDK4 gene as well as the detection of chromosome 12 alpha satellites in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description


- The ZytoDot 2C SPEC CDK4/CEN 12 Probe is composed of:
- Digoxigenin-labeled polynucleotides (~1,1 ng/μl), which target sequences mapping in 12q13.3-q14.1* (chr12:58,004,553-58,313,271) harboring the CDK4 region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~1,1 ng/μl), which target sequences mapping in 12p11.1-q11 specific for the alpha satellite centromeric region D12Z3 of chromosome 12.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without an amplification involving the CDK4 gene region, two distinct dot-shaped green and two distinct dot-shaped red signals appear.

Aberrant situation: In cells with an amplification of the CDK4 gene region, an increased number of green signals or green signal clusters will be observed.



Example of an aberrant signal pattern: Liposarcoma tissue section with CDK4 amplification as indicated by large green clusters.

ZytoDot® SPEC MDM2 Probe


Intended Purpose
The ZytoDot SPEC MDM2 Probe (PD9) is intended to be used for the qualitative detection of amplifications involving the human MDM2 gene in formalin-fixed, paraffin-embedded specimens, such as atypical lipomatous tumor/well-differentiated liposarcoma (ALT/WDLPS) and dedifferentiated liposarcoma (DDLPS), by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot CISH Implementation Kit (Prod. No. C-3018-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ALT/WDLPS and DDLPS and therapeutic measures should not be initiated based on the test result alone.

Probe Description


- The ZytoDot MDM2 Probe is composed of:
- Digoxigenin-labeled polynucleotides (~1,8 ng/μl), which target sequences mapping in 12q15* (chr12:69,190,708-69,430,185) (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without an amplification involving the MDM2 gene region, two distinct dot-shaped brown signals appear.

Aberrant situation: In cells with an amplification of the MDM2 gene region or aneuploidy of chromosome 12, an increased number of the brown signal or brown signal clusters will be observed.



Normal nuclei each with two MDM2 signals.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	12	CE/IVD	40 (400 μl)	C-3062-400

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG	12	CE/IVD	40 (400 μl)	C-3012-400

ZytoDot® 2C SPEC MDM2/CEN 12 Probe


Intended Purpose
The ZytoDot 2C SPEC MDM2/CEN 12 Probe (PD29) is intended to be used for the qualitative detection of amplifications involving the human MDM2 gene as well as the detection of chromosome 12 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as atypical lipomatous tumor/well-differentiated liposarcoma (ALT/WDLPS) and dedifferentiated liposarcoma (DDLPS), by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of ALT/WDLPS and DDLPS and therapeutic measures should not be initiated based on the test result alone.

Probe Description


- The ZytoDot 2C SPEC MDM2/CEN 12 Probe is composed of:
- Digoxigenin-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 12q15* (chr12:69,190,708-69,430,185) harboring the MDM2 region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 12p11.1-q11 specific for the alpha satellite centromeric region D12Z3 of chromosome 12.
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19



Results

Normal situation: In interphases of normal cells or cells without an amplification involving the MDM2 gene region, two distinct dot-shaped green and two distinct dot-shaped red signals appear.

Aberrant situation: In cells with an amplification of the MDM2 gene region, an increased number of green signals or green signal clusters will be observed.



Liposarcoma tissue section with MDM2 amplification as indicated by large green clusters.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	12	CE/IVD	10 (100 μl)	C-3049-100
			40 (400 μl)	C-3049-400

ZytoDot® 2C SPEC IGH Break Apart Probe


Intended Purpose
The ZytoDot 2C SPEC IGH Break Apart Probe (PD51) is intended to be used for the qualitative detection of translocations involving the human IGH locus at 14q32.33 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description


- The ZytoDot 2C SPEC IGH Break Apart Probe is composed of:
- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 14q32.33* (chr14:106,690,778-106,883,535) distal to the IGH breakpoint region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~0,75 ng/μl), which target sequences mapping in 14q32.33* (chr14:105,462,169-105,909,611) proximal to the IGH breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Results

Normal situation: In interphases of normal cells or cells without a translocation involving the IGH gene region, two red/green fusion signals appear.

Aberrant situation: One IGH gene region affected by a translocation is indicated by one separate distinct dot-shaped green signal and one separate distinct dot-shaped red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.



Example of an aberrant signal pattern: Burkitt lymphoma tissue section with translocation affecting the 14q32.33 locus as indicated by one red/green fusion (non-rearranged) signal, one red signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	14	CE/IVD	10 (100 μl)	C-3071-100

ZytoDot® 2C SPEC USP6 Break Apart Probe

Intended Purpose
The ZytoDot 2C SPEC USP6 Break Apart Probe (PD56) is intended to be used for the qualitative detection of translocations involving the human USP6 gene at 17p13.2 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC USP6 Break Apart Probe is composed of:

- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 17p13.2* (chr17:4,825,753-5,017,582) distal to the USP6 breakpoint region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~0.75 ng/μl), which target sequences mapping in 17p13.2* (chr17:5,087,046-5,361,104) proximal to the USP6 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC USP6 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the USP6 gene region, two red/green fusion signals appear.
Aberrant situation: One USP6 gene region affected by a translocation is indicated by one separate green signal and one separate red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	17	CE/IVD	10 (100 μl)	C-3077-100

ZytoDot® 2C SPEC ERBB2 Probe

Intended Purpose
The ZytoDot SPEC ERBB2 Probe (PD1) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene in formalin-fixed, paraffin-embedded specimens, such as breast cancer, by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot CISH Implementation Kit (Prod. No. C-3018-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot ERBB2 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1.8 ng/μl), which target sequences mapping in 17q12* (chr17:37,725,661-37,973,541) harboring the ERBB2 gene region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC ERBB2 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the ERBB2 gene region, two distinct dot-shaped brown signals appear.
Aberrant situation: In cells with an amplification of the ERBB2 gene region or polysomy of chromosome 17, an increased number of the brown signal or brown signal clusters will be observed.



Normal nuclei each with two ERBB2 signals.



Breast cancer tissue section with ERBB2 amplification.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG	17	CE/IVD	40 (400 μl)	C-3001-400

The product is also available as a kit.

ZytoDot® SPEC ERBB2 Probe Kit	Tests	Order No.
Incl.: Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; Probe, 0.4 ml; Wash Buffer SSC, 560 ml; PBS/Tween, good for 2000 ml; Blocking Solution, 4 ml; Mouse-anti-DIG, 4 ml; Anti-Mouse-HRP-Polymer, 4 ml; DAB Solution A, 0.3 ml; DAB Solution B, 10 ml; Mayer's Hematoxylin Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	40	C-3003-40

ZytoDot® 2C SPEC ERBB2/CEN 17 Probe

Intended Purpose
The ZytoDot 2C SPEC ERBB2/CEN 17 Probe (PD12) is intended to be used for the qualitative detection of amplifications involving the human ERBB2 gene as well as the detection of chromosome 17 alpha satellites in formalin-fixed, paraffin-embedded specimens, such as breast cancer, by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of breast cancer and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC ERBB2/CEN 17 Probe is composed of:

- Digoxigenin-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 17q12* (chr17:37,725,661-37,973,541) harboring the ERBB2 region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~1.1 ng/μl), which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19






Fig. 1: SPEC ERBB2 Probe map (not to scale); Right: Ideogram of chromosome 17 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without an amplification involving the ERBB2 gene region, two distinct dot-shaped green and two distinct dot-shaped red signals appear.
Aberrant situation: In cells with an amplification of the ERBB2 gene region, an increased number of green signals or green signal clusters will be observed.



Hybridized to normal interphase cells as indicated by two red and two green signals per nucleus.



Breast cancer tissue section with ERBB2 amplification as indicated by multiple green signals in each nucleus.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	17	CE/IVD	10 (100 μl)	C-3032-100
			40 (400 μl)	C-3032-400

The product is also available as a kit in two different sizes.

ZytoDot® 2C SPEC ERBB2/CEN 17 Probe Kit	Tests	Order No.
Incl.: Heat Pretreatment Solution EDTA, 150 ml; Pepsin Solution, 1 ml; Probe, 0.1 ml; Wash Buffer SSC, 210 ml; 20x Wash Buffer TBS, 50 ml; Anti-DIG/DNP-Mix, 1 ml; HRP/AP-Polymer-Mix, 1 ml; AP-Red Solution A, 0.1 ml; AP-Red Solution B, 4 ml; HRP-Green Solution A, 0.2 ml; HRP-Green Solution B, 4 ml; Nuclear Blue Solution, 4 ml; Mounting Solution (alcoholic), 1 ml	10	C-3022-10
Incl.: Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; Probe, 0.4 ml; Wash Buffer SSC, 560 ml; 20x Wash Buffer TBS, 2x 50 ml; Anti-DIG/DNP-Mix, 4 ml; HRP/AP-Polymer-Mix, 4 ml; AP-Red Solution A, 0.4 ml; AP-Red Solution B, 15 ml; HRP-Green Solution A, 0.8 ml; HRP-Green Solution B, 15 ml; Nuclear Blue Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	40	C-3022-40

ZytoDot® 2C SPEC SS18 Break Apart Probe

Intended Purpose
The ZytoDot 2C SPEC SS18 Break Apart Probe (PD26) is intended to be used for the qualitative detection of translocations involving the human SS18 gene at 18q11.2 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description
The ZytoDot 2C SPEC SS18 Break Apart Probe is composed of:


- Digoxigenin-labeled polynucleotides (~0.5 ng/μl), which target sequences mapping in 18q11.2* (chr18:23,109,942-23,262,464) proximal to the SS18 breakpoint region (see Fig. 1).
- Dinitrophenyl-labeled polynucleotides (~0.75 ng/μl), which target sequences mapping in 18q11.2* (chr18:23,772,255-24,137,169) distal to the SS18 breakpoint region (see Fig. 1).
- Formamide based hybridization buffer

* according to Human Genome Assembly GRCh37/hg19




Fig. 1: SPEC SS18 Probe map (not to scale); Right: Ideogram of chromosome 18 indicating the hybridization locations

Results
Normal situation: In interphases of normal cells or cells without a translocation involving the SS18 gene region, two red/green fusion signals appear.
Aberrant situation: One SS18 gene region affected by a translocation is indicated by one separate distinct dot-shaped green signal and one separate distinct dot-shaped red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.



Example of an aberrant signal pattern: Synovial sarcoma tissue section with translocation affecting the 18q11.2 locus as indicated by one non-rearranged red/green fusion signal, one red signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	18	CE/IVD	10 (100 μl)	C-3046-100

ZytoDot® 2C SPEC BCL2 Break Apart Probe

Intended Purpose

The ZytoDot 2C SPEC BCL2 Break Apart Probe (PD53) is intended to be used for the qualitative detection of translocations involving the human BCL2 gene at 18q21.33 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoDot 2C SPEC BCL2 Break Apart Probe is composed of:
- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 18q21.33* (chr18:60,415,418-60,589,273) proximal to the BCL2 breakpoint region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~0.75 ng/μl), which target sequences mapping in 18q21.33* (chr18:60,994,528-61,507,691) distal to the BCL2 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19





Fig. 1: SPEC BCL2 Probe map (not to scale); Right: Ideogram of chromosome 18 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the BCL2 gene region, two red/green fusion signals appear.

Aberrant situation: One BCL2 gene region affected by a translocation is indicated by one separate distinct dot-shaped green signal and one separate distinct dot-shaped red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.

Example of an aberrant signal pattern: Lymphoma tissue section with translocation affecting the 18q21.33 locus as indicated by one red/green fusion (non-rearranged) signal, one red signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	18	CE/IVD	10 (100 μl)	C-3073-100

ZytoDot® 2C SPEC EWSR1 Break Apart Probe

Intended Purpose

The ZytoDot 2C SPEC EWSR1 Break Apart Probe (PD24) is intended to be used for the qualitative detection of translocations involving the human EWSR1 gene at 22q12.2 in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoDot 2C CISH Implementation Kit (Prod. No. C-3044-10/-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

- The ZytoDot 2C SPEC EWSR1 Break Apart Probe is composed of:
- Digoxigenin-labeled polynucleotides (~0.50 ng/μl), which target sequences mapping in 22q12.2* (chr22:29,779,841-30,057,928) distal to the EWSR1 breakpoint region (see Fig. 1).
 - Dinitrophenyl-labeled polynucleotides (~0.75 ng/μl), which target sequences mapping in 22q12.1-22q12.2* (chr22:29,413,831-29,673,440) proximal to the EWSR1 breakpoint region (see Fig. 1).
 - Formamide based hybridization buffer
- * according to Human Genome Assembly GRCh37/hg19




Fig. 1: SPEC EWSR1 Probe map (not to scale); Right: Ideogram of chromosome 22 indicating the hybridization locations

Results

Normal situation: In interphases of normal cells or cells without a translocation involving the EWSR1 gene region, two red/green fusion signals appear.

Aberrant situation: One EWSR1 gene region affected by a translocation is indicated by one separate distinct dot-shaped green signal and one separate distinct dot-shaped red signal.



Hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.


Example of an aberrant signal pattern: Ewing sarcoma tissue section with translocation affecting the 22q12.1-q12.2 locus as indicated by one non-rearranged red/green fusion signal, one red signal, and one separate green signal.

Labels	Chromosome	Reg. Status	Tests	Order No.
DIG/DNP	22	CE/IVD	10 (100 μl)	C-3043-100

ZytoDot® RUO Probes


Product Name	Chr.band	Reg. Status	Labels	Tests	Order No.
ZytoDot SPEC EGFR Probe	7p11.2	RUO	DIG	40 (400 μl)	C-3007-400
ZytoDot SPEC MYC Probe	8q24.21		DIG	40 (400 μl)	C-3013-400
ZytoDot 2C SPEC PTEN/CEN 10 Probe	10q23.3		DIG/DNP	40 (400 μl)	C-3053-400
ZytoDot 2C SPEC FOXO1 Break Apart Probe	13q14.1		DIG/DNP	10 (100 μl)	C-3065-100
ZytoDot 2C SPEC MALT1 Break Apart Probe	18q21.3		DIG/DNP	10 (100 μl)	C-3072-100

C-3007-400




SPEC EGFR Probe map (not to scale)

C-3013-400




SPEC MYC Probe map (not to scale)

C-3053-400




SPEC PTEN Probe map (not to scale)

C-3065-100



SPEC FOXO1 Probe map (not to scale)

C-3072-100



SPEC MALT1 Probe map (not to scale)

ZytoDot® Probes for Chromosome Enumeration - RUO

The ZytoDot® Chromosome Enumeration Probes are designed for identification and enumeration of human chromosomes in interphase cells and as an adjunct to standard karyotyping in metaphases. These probes will produce sharp, bright signals specific for each individual chromosome.

CEN Probe Description


For most chromosomes, direct labeled ZytoDot® CEN™ Probes hybridizing to highly repetitive human satellite DNA sequences mainly located at the centromeric regions of chromosomes are applicable.

Results

In a normal interphase nucleus, two signals are expected using Chromosome Enumeration Probes specific for autosomes.

Using chromosome Y specific probes will result in normal male cells in one signal and in normal female cells in no signal.

Using chromosome X specific probes will result in normal male cells in one signal and in normal female cells in two signals per nucleus. Other signal patterns indicate numerical aberrations of the respective chromosome.



CEN X/Y Probe hybridized on normal male interphase cells as indicated by one red (chromosome X) and one green (chromosome Y) signal per nucleus.

CEN X/Y Probe hybridized on normal female interphase cells as indicated by two red (chromosome X) signals per nucleus.

RUO Probes

Product Name	Alpha/Class. Sat.	Chr.band	Reg. Status	Labels	Tests	Order No.
ZytoDot CEN 7 Probe	D7Z1	7p11.1-q11.1	RUO	DIG	40 (400 µl)	C-3008-400
ZytoDot CEN 8 Probe	D8Z2	8p11.1-q11.1		DIG	40 (400 µl)	C-3016-400
ZytoDot CEN X Probe	DXZ1	Xp11.1-q11.1		DIG	40 (400 µl)	C-3025-400
ZytoDot CEN Yq12 Probe	III DYZ1	Yq12		DIG	40 (400 µl)	C-3020-400
ZytoDot 2C CEN X/Y Probe	DXZ1/DYZ3	Xp11.1-q11.1/Yp11.1-q11.1		DNP/DIG	40 (400 µl)	C-3048-400

Accessories

ZytoDot® Implementation Kits

For the detection of ZytoDot® probes

Product	Manufacturer	Tests	Reg. Status	Order No.
ZytoDot CISH Implementation Kit For the detection of Digoxigenin-labeled ZytoDot Probes Incl.: Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 560 ml; PBS/Tween, ausreichend für 2000 ml; Blocking Solution, 4 ml; Mouse-anti-DIG, 4 ml; Anti-Mouse-HRP-Polymer, 4 ml; DAB Solution A, 0.3 ml; DAB Solution B, 10 ml; Mayer’s Hematoxylin Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	ZytoVision GmbH	40	CE/IVD	C-3018-40
ZytoDot 2C CISH Implementation Kit For the detection of Digoxigenin/Dinitrophenyl-labeled ZytoDot 2C Probes Incl.: Heat Pretreatment Solution EDTA, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 210 ml; 20x Wash Buffer TBS, 50 ml; Anti-DIG/DNP-Mix, 1 ml; HRP/AP-Polymer-Mix, 1 ml; AP-Red Solution A, 0.2 ml; AP-Red Solution B, 4 ml; HRP-Green Solution A, 0.2 ml; HRP-Green Solution B, 4 ml; Nuclear Blue Solution, 4 ml; Mounting Solution (alcoholic), 1 ml	ZytoVision GmbH	10	CE/IVD	C-3044-10
ZytoDot 2C CISH Implementation Kit For the detection of Digoxigenin/Dinitrophenyl-labeled ZytoDot 2C Probes Incl.: Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 560 ml; 20x Wash Buffer TBS, 2x 50 ml; Anti-DIG/DNP-Mix, 4 ml; HRP/AP-Polymer-Mix, 4 ml; AP-Red Solution A, 0.5 ml; AP-Red Solution B, 15 ml; HRP-Green Solution A, 0.8 ml; HRP-Green Solution B, 15 ml; Nuclear Blue Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	ZytoVision GmbH	40	CE/IVD	C-3044-40

ZytoDot® Pretreatment Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
Pepsin Solution	ZytoVision GmbH	4 ml	CE/IVD	ES-0001-4
Pepsin Solution	ZytoVision GmbH	50 ml	CE/IVD	ES-0001-50
Pepsin Solution	ZytoVision GmbH	1000 ml	CE/IVD	ES-0001-1000
Heat Pretreatment Solution EDTA	ZytoVision GmbH	500 ml	CE/IVD	PT-0002-500

ZytoDot® Wash Buffers & Ancillary Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
Mouse-anti-DIG	ZytoVision GmbH	4 ml	CE/IVD	AB-0001-4
Anti-Mouse-HRP-Polymer	ZytoVision GmbH	4 ml	CE/IVD	AB-0002-4
HRP/AP-Polymer-Mix	ZytoVision GmbH	4 ml	CE/IVD	AB-0013-4
Anti-DIG/DNP-Mix	ZytoVision GmbH	4 ml	CE/IVD	AB-0014-4
Blocking Solution	ZytoVision GmbH	4 ml	CE/IVD	BS-0001-4
DAB Solution Set Incl.: DAB Solution A, 0.3 ml; DAB Solution B, 10 ml; good for 10 ml DAB Solution	ZytoVision GmbH	-	CE/IVD	C-3015-100
ZytoDot AP-Red Solution Set Incl.: AP-Red Solution A, 0.5 ml; AP-Red Solution B, 15 ml; good for 15 ml AP-Red Solution	ZytoVision GmbH	-	CE/IVD	C-3038-100
ZytoDot HRP-Green Solution Set Incl.: HRP-Green Solution A, 0.8 ml; HRP-Green Solution B, 15 ml; good for 15 ml HRP-Green Solution	ZytoVision GmbH	-	CE/IVD	C-3039-100
Mayer’s Hematoxylin Solution	ZytoVision GmbH	20 ml	CE/IVD	CS-0001-20
Nuclear Blue Solution	ZytoVision GmbH	20 ml	CE/IVD	CS-0002-20
Mounting Solution (alcoholic)	ZytoVision GmbH	4 ml	CE/IVD	MT-0004-4
Wash Buffer SSC	ZytoVision GmbH	560 ml	CE/IVD	WB-0001-560
PBS/Tween	ZytoVision GmbH	good for 1000 ml	CE/IVD	WB-0004-1000
20x Wash Buffer TBS	ZytoVision GmbH	50 ml	CE/IVD	WB-0005-50

RUO ZytoDot® Accessories

Product	Manufacturer	Amount	Reg. Status	Order No.
ERBB2 Control Slide Set	ZytoVision GmbH	2 pcs.	RUO	E-4007-2

ZytoFast®

The **ZytoFast® products** are designed for outstandingly fast detection and determination of lymphocyte clonality by detecting IGK and IGL light chain RNA by Chromogenic *in situ* Hybridization (CISH) in formalin-fixed, paraffin-embedded (FFPE) specimens.

ZytoVision GmbH


ZytoFast®
Outstandingly fast CISH

- Optimized protocols and faster tissue penetration due to short oligonucleotide probes of the ZytoFast® system, make the ZytoFast® CISH procedure outstandingly fast.
- Results can be achieved within just 5 hours, hands-on time is about 3 hours!

Advantages of CISH


- + Simultaneous observation of tissue morphology and CISH signals
- + No risk of false positives due to mispriming or contamination as with PCR
- + Easy method comparable to IHC
- + No costly equipment needed
- + Ability to test archival specimens

Method



The ZytoFast® system uses oligonucleotide probes tagged with Biotin and Digoxigenin [1], which are detected using HRP-conjugated antibodies and AP-conjugated streptavidin targeting the tags [2]. The enzymatic reaction of chromogenic substrates [3], e.g. BCIP/NBT and AEC, leads to the formation of strong color precipitates that can be visualized by light microscopy.

Protocol Overview



ZytoFast® PLUS

The **ZytoFast® PLUS products** are designed for outstandingly fast detection and discrimination of human pathogen viruses, e.g. HPV and EBV and the determination of lymphocyte clonality by detecting IGK and IGL light chain RNA by Chromogenic *in situ* Hybridization (CISH) in formalin-fixed, paraffin-embedded (FFPE) specimens. The signal intensity of ZytoFast® probes is increased even more when using the ZytoFast® PLUS Implementation Kits.

ZytoVision GmbH


ZytoFast® PLUS
Outstandingly fast and sensitive CISH

- Depending on the time required for dewaxing and pre-treatment of tissue sections, ZytoFast® PLUS protocols can be performed within approx. 4 hours!
- Thus, due to optimized protocols, the ZytoFast® PLUS method takes only slightly more time compared to ZytoFast® protocols while being much more sensitive!

ZytoFast® PLUS
Flexibility that meets your needs

- Several ZytoFast® PLUS CISH Implementation Kits using different enzyme/substrate combinations can be combined with any separately available Digoxigenin-labeled ZytoFast® probe to meet your preferences concerning the detection chemistry and counterstaining.
- Each ZytoFast® PLUS CISH Implementation Kit includes a detailed protocol and all necessary reagents for versatile use in DNA as well as RNA *in situ* hybridizations.

Method



The ZytoFast® PLUS system uses Digoxigenin-labeled probes [1], which are detected using primary antibodies [2]. These antibodies are detected by polymerized enzyme-conjugated secondary antibodies [3]. The enzymatic reaction of chromogenic substrates [4], e.g. DAB, leads to the formation of strong color precipitates that can be visualized by light microscopy.

Protocol Overview



Virus Index

Virus	Product	Label	Amount	Reg. Status	Order No.	Page
HPV	ZytoFast HPV type 6/11 Probe	DIG	400 µl	CE/IVD	T-1055-400	161
	ZytoFast HPV type 16/18 Probe	DIG	400 µl	CE/IVD	T-1056-400	161
	ZytoFast HPV type 31/33 Probe	DIG	400 µl	CE/IVD	T-1057-400	161
	ZytoFast HPV High-Risk (HR) Types Probe (specific for HPV type 16/18/31/33/35/39/45/51/52/56/58/59/66/68/82)	DIG	400 µl	CE/IVD	T-1140-400	162
	ZytoFast HPV Screening Probe (specific for HPV type 6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82)	DIG	400 µl	CE/IVD	T-1144-400	162
EBV	ZytoFast EBV Probe	DIG	400 µl	CE/IVD	T-1114-400	162
CMV	ZytoFast CMV Probe	DIG	400 µl	RUO	T-1113-400	163

mRNA Index

mRNA	Product	Label	Amount	Reg. Status	Order No.	Page
Ig-kappa	ZytoFast human Ig-kappa Probe	DIG	400 µl	CE/IVD	T-1115-400	164
	ZytoFast human Ig-kappa/Ig-lambda Probe	DIG/Biotin	400 µl	CE/IVD	T-1017-400	164
	ZytoFast human Ig-kappa/Ig-lambda CISH Kit	DIG/Biotin	40 tests	CE/IVD	T-1005-40	164
	ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit	DIG/Biotin	40 tests	CE/IVD	T-1105-40	164
Ig-lambda	ZytoFast human Ig-lambda Probe	DIG	400 µl	CE/IVD	T-1116-400	164
	ZytoFast human Ig-kappa/Ig-lambda Probe	DIG/Biotin	400 µl	CE/IVD	T-1017-400	164
	ZytoFast human Ig-kappa/Ig-lambda CISH Kit	DIG/Biotin	40 tests	CE/IVD	T-1005-40	164
	ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit	DIG/Biotin	40 tests	CE/IVD	T-1105-40	164

ZytoFast® HPV-CISH System

Intended Purpose

The ZytoFast HPV type 6/11 Probe (PF25) is intended to be used for the qualitative detection of human papillomavirus (HPV) type 6/11 DNA in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH).

The ZytoFast HPV type 16/18 Probe (PF26) is intended to be used for the qualitative detection of human papillomavirus (HPV) type 16/18 DNA in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH).

The ZytoFast HPV type 31/33 Probe (PF27) is intended to be used for the qualitative detection of human papillomavirus (HPV) type 31/33 DNA in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH).

The probes are intended to be used in combination with the ZytoFast PLUS CISH Implementation Kit HRP-DAB (Prod. No. T-1063-40).

The products are intended for professional use only. All tests using these products should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probes are intended to be used as an aid to the differential diagnosis of various cancers and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoFast HPV type 6/11 Probe is composed of:

- Digoxigenin-labeled oligonucleotides (~0.6 ng/µl) specific for HPV type 6/11, which target DNA sequences encoding for HPV 6/11 proteins E6, E7, and/or L1.
- The probe also targets the respective RNA sequences of E6, E7, and/or L1 proteins, which are expressed during some stages of infection.
- Formamide based hybridization buffer

The ZytoFast HPV type 16/18 Probe is composed of:


- Digoxigenin-labeled oligonucleotides (~0.6 ng/µl) specific for HPV type 16/18, which target DNA sequences encoding for the HPV 16/18 proteins E6, E7, and/or L1.
- The probe also targets the respective RNA sequences of E6, E7, and/or L1 proteins, which are expressed during some stages of infection.
- Formamide based hybridization buffer

The ZytoFast HPV type 31/33 Probe is composed of:

- Digoxigenin-labeled oligonucleotides (~0.6 ng/µl) specific for HPV type 31/33, which target DNA sequences encoding for HPV 31/33 proteins E6, E7, and/or L1.
- The probe also targets the respective RNA sequences of E6, E7, and/or L1 proteins, which are expressed during some stages of infection.
- Formamide based hybridization buffer

Results

A positive reactivity for HPV DNA in epithelial cells is indicated by a distinctly stained nucleus. Due to the detection of HPV DNA as well as E6, E7, and/or L1 RNAs, depending on the infection stage, cytoplasmic staining might be observed additionally. Colored precipitates, which can be clearly distinguished from the background, will be dark brown when using DAB as substrate for the detection.



Example of an HPV positive signal pattern: HPV infected cervix tissue hybridized with the ZytoFast HPV type 6/11 Probe, detected with the ZytoFast PLUS CISH Implementation Kit HRP-DAB.

Product Name	Labels	Reg. Status	Tests	Order No.
ZytoFast HPV type 6/11 Probe	DIG	CE/IVD	40 (400 µl)	T-1055-400
ZytoFast HPV type 16/18 Probe	DIG	CE/IVD	40 (400 µl)	T-1056-400
ZytoFast HPV type 31/33 Probe	DIG	CE/IVD	40 (400 µl)	T-1057-400




Fig. 1: Schematic representation of the HPV genome with E and L open reading frames. Genomic regions targeted by ZytoFast HPV specific oligonucleotides are indicated in blue.

ZytoFast® HPV High-Risk (HR) Types Probe

Intended Purpose
The ZytoFast HPV High-Risk (HR) Types Probe (PF46) is intended to be used for the qualitative detection of human papillomavirus (HPV) type 16/18/31/33/35/39/45/51/52/56/58/59/66/68/82 DNA in formalin-fixed, paraffin-embedded specimens, such as cervical carcinoma and oropharyngeal squamous cell cancer (OSCC), by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoFast PLUS CISH Implementation Kit HRP-DAB (Prod. No. T-1063-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of cervical carcinoma and OSCC and therapeutic measures should not be initiated based on the test result alone.

- Probe Description**
The ZytoFast HPV High-Risk (HR) Types Probe is composed of:
- Digoxigenin-labeled oligonucleotides (~ 2.2 ng/µl) specific for HPV type 16/18/31/33/35/39/45/51/52/56/58/59/66/68/82, which target DNA sequences encoding for HPV 16/18/31/33/ 35/39/45/51/52/56/ 58/59/66/68/82 proteins E6, E7, and/or L1.
 - The probe also targets the respective RNA sequences of E6, E7, and/or L1 proteins, which are expressed during some stages of infection.
 - Formamide based hybridization buffer





Fig. 1: Schematic representation of the HPV genome with E and L open reading frames. Genomic regions targeted by ZytoFast HPV specific oligonucleotides are indicated in blue.

Results
A positive reactivity for HPV DNA in epithelial cells is indicated by a distinctly stained nucleus. Due to the detection of HPV DNA as well as E6, E7, and/or L1 RNAs, depending on the infection stage, cytoplasmic staining might be observed additionally. Colored precipitates, which can be clearly distinguished from the background, will be dark brown when using DAB as substrate for the detection.



HPV infected cervix tissue hybridized with the ZytoFast HPV High-Risk (HR) Types Probe, detected with the ZytoFast PLUS CISH Implementation Kit HRP-DAB.

Labels	Reg. Status	Tests	Order No.
DIG	CE/IVD	40 (400 µl)	T-1140-400

ZytoFast® HPV Screening Probe

Intended Purpose
The ZytoFast HPV Screening Probe is designed for the detection and discrimination of human papillomavirus (HPV) DNA in paraffin-embedded tissue sections or cell samples.

At least 50 percent of sexually active men and women acquire some form of genital HPV infection at some point in their lives.

Most of the approx. 30 identified genital HPV types, predominantly types 6 and 11, are called “low-risk” types, and may cause mild Pap test abnormalities or genital warts.

Until now, approximately 10–15 HPV types are associated with lesions that can progress to cancer. Among those are the HPV types 16/18/31/33/35/39/45/ 51/52/56/58/59/66/68/82.

These cancer-associated HPV types are designated as high-risk HPV (hr-HPV) types. The infection with the HPV hr-types can lead to development of cancer of the cervix, vulva, vagina, anus, or penis. The majority of malignant cervical carcinomas (approx. 70%) occur as a result of infections with HPV types 16 or 18.

References
Cothran MM & White JP (2002) Health Care Women Int 23: 306-19.
Cubie HA & Norval M (1988) J Virol Methods 20: 239-49.
Faulkner-Jones BE, et al. (1993) J Virol Methods 41: 277-96.
Francis IM, et al. (2013) Sultan Qaboos Univ Med J 13: 527-33.
Grandmeier N, et al. (2011) Dermatology 223: 293-300.
Kaspersen MD, et al. (2011) PLoS One 6: e18095.
Mirasoli M, et al. (2009) Anal Bioanal Chem 394: 981-7.
Montag M, et al. (2011) Arch Gynecol Obstet 284: 999-1005.
Poljak M & Kocjan BJ (2010) Expert Rev Anti Infect Ther 8: 1139-62.
Reinholz M, et al. (2013) Arch Dermatol Res 305: 723-32.
Walboomers JMM, et al. (1999) J Pathol 189: 12-9.

- Probe Description**
The ZytoFast HPV Screening Probe is composed of:
- Digoxigenin-labeled oligonucleotides (~2.6 ng/µl) specific for HPV type 6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82, which target DNA sequences encoding for HPV 6/11/16/18/31/33/35/39/45/ 51/52/56/58/59/66/68/82 proteins E6, E7, and/or L1.
 - The probe also targets the respective RNA sequences of E6, E7, and/or L1 proteins, which are expressed during some stages of infection.
 - Formamide based hybridization buffer





Fig. 1: Schematic representation of the HPV genome with E and L open reading frames. Genomic regions targeted by ZytoFast HPV specific oligonucleotides are indicated in blue.

Results
A positive reactivity for HPV DNA in epithelial cells is indicated by a distinctly stained nucleus. Due to the detection of HPV DNA as well as E6, E7, and/or L1 RNAs, depending on the infection stage, cytoplasmic staining might be observed additionally. Depending on the detection chemistry that is used, colored precipitates, which can be clearly distinguished from the background, will be dark brown when using DAB as substrate for the detection.



HPV infected cells with signals from integrated and episomal HPV hybridized with the ZytoFast HPV Screening Probe, detected with the ZytoFast PLUS CISH Implementation Kit HRP-DAB.

Labels	Reg. Status	Tests	Order No.
DIG	CE/IVD	40 (400 µl)	T-1144-400

ZytoFast® EBV Probe

Intended Purpose
The ZytoFast EBV Probe (PF29) is intended to be used for the qualitative detection of human Epstein-Barr virus (EBV) EBER RNA in formalin-fixed, paraffin-embedded specimens, such as diffuse large B-cell lymphomas (DLBCL) or Hodgkin lymphomas, by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoFast PLUS CISH Implementation Kit HRP-DAB (Prod. No. T-1063-40).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The probe is intended to be used as an aid to the differential diagnosis of DLBCL or Hodgkin lymphomas and therapeutic action should not be initiated on the basis of the test result alone.

- Probe Description**
The ZytoFast EBV Probe is composed of:
- Digoxigenin-labeled oligonucleotides (~ 0.2 ng/µl), which target mRNA sequences encoding EBER-1 and EBER-2 regions.





Fig. 1: Schematic representation of the EBV genome with the EBER-1 and EBER-2 encoding region indicated in blue. U1-U5 indicate unique nucleotide sequences, hatched boxes represent terminal and internal repeats.

Results
A positive reactivity for Epstein-Barr-virus (EBV) EBER RNA in the target cells is indicated by a distinctly stained nucleus. Colored precipitates, which can be clearly distinguished from the background, will be dark brown when using DAB for detection.



EBV infected tonsil tissue hybridized with ZytoFast EBV Probe, detected with ZytoFast PLUS CISH Implementation Kit HRP-DAB.

Labels	Reg. Status	Tests	Order No.
DIG	CE/IVD	40 (400 µl)	T-1114-400

ZytoFast® CMV Probe

Intended Purpose
The ZytoFast CMV Probe (PF28) is intended to be used for the qualitative detection of human cytomegalovirus (CMV (a.k.a. human herpesvirus-5, HHV-5)) DNA in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH).

The ZytoFast CMV Probe is intended to be used in combination with one of the ZytoFast PLUS CISH Implementation Kits, either the ZytoFast PLUS CISH Implementation Kit AP-NBT/BCIP (Prod. No. T-1061-40), the ZytoFast PLUS CISH Implementation Kit HRP-DAB (Prod. No. T-1063-40), or the ZytoFast PLUS CISH Implementation Kit AP-Permanent Red (Prod. No. T-1151-40).

- Probe Description**
The ZytoFast CMV Probe is composed of:
- Digoxigenin-labeled oligonucleotides (~ 0.2 ng/µl), which target DNA sequences encoding the β2.7 gene, the most abundantly transcribed early CMV gene.
 - Formamide based hybridization buffer





Fig. 1: Schematic representation of the CMV genome with the β2.7 encoding region indicated in blue. UL and US indicate unique nucleotide sequences, hatched boxes represent terminal and internal repeats.

Results
Due to the detection of CMV DNA as well as of the abundantly transcribed β2.7 RNA, a positive reactivity for cytomegalovirus (CMV) in the target cells is indicated by a cytoplasmic and/or nuclear staining pattern. Depending on the detection chemistry that is used, colored precipitates, which can be clearly distinguished from the background, will be dark violet-blue when using NBT/BCIP as substrate, or strong red when using Permanent Red.



CISH analysis of paraffin-embedded adrenal gland tissue using the ZytoFast CMV Probe, detected with ZytoFast PLUS CISH Implementation Kit AP-NBT/BCIP.

Labels	Reg. Status	Tests	Order No.
DIG	RUO	40 (400 µl)	T-1113-400

ZytoFast® Ig-kappa/Ig-lambda CISH System

Intended Purpose

The ZytoFast Ig-kappa/Ig-lambda-CISH System is intended to be used for the qualitative detection of human Ig-kappa (κ) and Ig-lambda (λ) light chain mRNA in formalin-fixed, paraffin-embedded specimens, such as multiple myeloma, by chromogenic *in situ* hybridization (CISH).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The product is intended to be used as an aid to the differential diagnosis of multiple myeloma and therapeutic measures should not be initiated based on the test result alone.

Probe Description

The ZytoFast human Ig-kappa Probe (PF30) is composed of:


- Digoxigenin-labeled oligonucleotides (~ 0.2 ng/μl), which target mRNA sequences encoding Ig-kappa light chain constant regions.

The ZytoFast human Ig-lambda Probe (PF31) is composed of:


- Digoxigenin-labeled oligonucleotides (~ 0.2 ng/μl), which target mRNA sequences encoding Ig-lambda light chain constant regions.

The ZytoFast human Ig-kappa/Ig-lambda Probe (PF22) is composed of:


- Digoxigenin-labeled oligonucleotides (~ 1 ng/μl), which target mRNA sequences encoding Ig-kappa light chain constant regions.
- Biotin-labeled oligonucleotides (~ 1ng/μl), which target mRNA sequences encoding Ig-lambda light chain constant regions.



CISH analysis of a paraffin-embedded bone marrow biopsy specimen using the ZytoFast human Ig-kappa/Ig-lambda CISH Kit.



CISH analysis of a paraffin-embedded tonsil tissue using the ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit.



Multiple myeloma tissue with B-cells expressing Ig-kappa hybridized with ZytoFast human Ig-kappa Probe, detected with ZytoFast PLUS CISH Implementation Kit HRP-DAB.




Fig. 1: Basic immunoglobulin structure indicating the heavy chains (black), I (blue) and k (red) lights chains. The light chain constant regions (C) whose encoding mRNA sequences are targeted by ZytoFast Ig-lambda and Ig-kappa probes are indicated in dark blue and red respectively, the variable regions (V) in light blue and red.

Results

A positive reactivity in the target cells is indicated by cytoplasmic staining. Depending on the detection chemistry that is used, colored precipitates, which can be clearly distinguished from the background, will be dark violet-blue when using NBT/BCIP as substrate, strong red when using AEC, dark brown when using DAB, green when using HRP-Green, or strong red when using Permanent Red.

Using the ZytoFast human Ig-kappa Probe, B-cells expressing antibodies with κ light chains will result in cytoplasmic staining whereas IGL expressing B-cells are not stained.

Using the ZytoFast human Ig-lambda Probe, B-cells expressing antibodies with λ light chains will result in cytoplasmic staining whereas IGL expressing B-cells are not stained.

Using the ZytoFast human Ig-kappa/Ig-lambda CISH Kit, B-cells expressing antibodies with κ light chains will result in a red cytoplasmic staining and simultaneously IGL expressing B-cells will result in a dark violet-blue cytoplasmic staining.

Using the ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit, B-cells expressing antibodies with κ light chains will result in a green cytoplasmic staining and simultaneously IGL expressing B-cells will result in permanent red cytoplasmic staining.

Digoxigenin-labeled Probes

Product Name	Reg. Status	Tests	Order No.
ZytoFast human Ig-kappa Probe	CE/IVD	40 (400 μl)	T-1115-400
ZytoFast human Ig-lambda Probe	CE/IVD	40 (400 μl)	T-1116-400

Biotin-/Digoxigenin-labeled Probes and Probe Kits

Product Name	Reg. Status	Tests	Order No.
ZytoFast human Ig-kappa/Ig-lambda Probe	CE/IVD	40 (400 μl)	T-1017-400
ZytoFast human Ig-kappa/Ig-lambda CISH Kit	CE/IVD	40	T-1005-400
Incl.: Ig-kappa/Ig-lambda Probe (DIG/Biotin-labeled), 0.4 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 2x 50 ml; Anti-Biotin/DIG-Mix, 4 ml; AEC Solution, 4 ml; NBT/BCIP Solution, 4 ml; Nuclear Green Solution, 20 ml; Mounting Solution (aqueous), 4 ml			
ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit	CE/IVD	40	T-1105-400
Incl: Ig-kappa/Ig-lambda Probe (DIG/Biotin-labeled), 0.4 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 2x 50 ml; Anti-Biotin/DIG-Mix, 4 ml; HRP-Green-Solution A, 0.8 ml; HRP-Green-Solution B, 15 ml; Permanent Red Solution A, 0.25 ml; Permanent Red Solution B, 15 ml; Nuclear Blue Solution, 20 ml; Mounting Solution (alcoholic), 4 ml			


Accessories

ZytoFast® PLUS Implementation Kits for Use in Diagnostic Procedures

For the detection of Digoxigenin-labeled ZytoFast® probes

Product	Manufacturer	Tests	Reg. Status	Order No.
ZytoFast PLUS CISH Implementation Kit HRP-DAB Incl.: Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 4x 50 ml; Mouse-anti-DIG, 4 ml; Anti-Mouse-HRP-Polymer, 4 ml; DAB Solution A, 0.3 ml; DAB Solution B, 10 ml; Nuclear Blue Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	ZytoVision GmbH	40	CE/IVD	T-1063-40

Peroxidase (HRP)



ZytoFast® Pretreatment Reagents

Product	Manufacturer	Amount	Reg. Status	Order No.
Pepsin Solution	ZytoVision GmbH	4 ml	CE/IVD	ES-0001-4
Pepsin Solution	ZytoVision GmbH	50 ml	CE/IVD	ES-0001-50
Pepsin Solution	ZytoVision GmbH	1000 ml	CE/IVD	ES-0001-1000
Heat Pretreatment Solution EDTA	ZytoVision GmbH	500 ml	CE/IVD	PT-0002-500

ZytoFast® Wash Buffers & Ancillary Reagents


Product	Manufacturer	Amount	Reg. Status	Order No.
Mouse-anti-DIG	ZytoVision GmbH	4 ml	CE/IVD	AB-0001-4
Anti-Mouse-HRP-Polymer	ZytoVision GmbH	4 ml	CE/IVD	AB-0002-4
Anti-Biotin/DIG-Mix	ZytoVision GmbH	4 ml	CE/IVD	AB-0015-4
DAB Solution Set Incl.: DAB Solution A, 0.3 ml; DAB Solution B, 10 ml; good for 10 ml DAB Solution	ZytoVision GmbH	-	CE/IVD	C-3015-100
ZytoDot HRP-Green Solution Set Incl: HRP-Green Solution A, 0.8 ml; HRP-Green Solution B, 15 ml; good for 15 ml HRP-Green Solution	ZytoVision GmbH	-	CE/IVD	C-3039-100
Mayer's Hematoxylin Solution	ZytoVision GmbH	20 ml	CE/IVD	CS-0001-20
Nuclear Blue Solution	ZytoVision GmbH	20 ml	CE/IVD	CS-0002-20
Nuclear Green Solution	ZytoVision GmbH	20 ml	CE/IVD	CS-0004-20
Mounting Solution (alcoholic)	ZytoVision GmbH	4 ml	CE/IVD	MT-0004-4
NBT/BCIP Solution	ZytoVision GmbH	4 ml	CE/IVD	SB-0004-4
AEC Solution	ZytoVision GmbH	4 ml	CE/IVD	SB-0005-4
20x Wash Buffer TBS	ZytoVision GmbH	50 ml	CE/IVD	WB-0005-50

ZytoFast® PLUS Implementation Kits for Research Use Only


For the detection of Digoxigenin-labeled ZytoFast® probes

Product	Manufacturer	Tests	Reg. Status	Order No.
ZytoFast PLUS CISH Implementation Kit AP-NBT/BCIP Incl.: Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 4x 50 ml; Rabbit-anti-DIG, 4 ml; Anti-Rabbit-AP-Polymer, 4 ml; NBT/BCIP Solution, 4 ml; Nuclear Red Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	ZytoVision GmbH	40	RUO	T-1061-40
ZytoFast PLUS CISH Implementation Kit AP-Permanent Red Incl: Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 4x50 ml; Rabbit-anti-DIG, 4 ml; Anti-Rabbit-AP-Polymer, 4 ml; Permanent Red Solution A, 0.25 ml; Permanent Red Solution B, 15 ml; Mayer’s Hematoxylin Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	ZytoVision GmbH	40	RUO	T-1151-40

Alkaline Phosphatase (AP)



T-1061-40



T-1151-40

ZytoFast® Wash Buffers & Ancillary Reagents for Research Use Only

Product	Manufacturer	Amount	Reg. Status	Order No.
Rabbit-anti-DIG	ZytoVision GmbH	4 ml	RUO	AB-0011-4
Anti-Rabbit-AP-Polymer	ZytoVision GmbH	4 ml	RUO	AB-0012-4
Nuclear Red Solution	ZytoVision GmbH	20 ml	RUO	CS-0003-20

ZytoFast® Control Probes

Product	Labels	Manufacturer	Tests	Reg. Status	Order No.
ZytoFast DNA (+) Control Probe	DIG	ZytoVision GmbH	10 (100 µl)	RUO	T-1053-100
ZytoFast DNA (-) Control Probe	DIG	ZytoVision GmbH	10 (100 µl)	RUO	T-1054-100
ZytoFast 28S rRNA (+) Control Probe	DIG	ZytoVision GmbH	10 (100 µl)	RUO	T-1120-100
ZytoFast RNA (-) Control Probe	DIG	ZytoVision GmbH	10 (100 µl)	RUO	T-1119-100

Arrays

DNA arrays: modern technologies for reliable results

Our diverse product portfolio is fully aligned with the current standards of molecular pathology diagnostics and research. The sensitive VisionArray® system enables the reliable detection of various pathogens in infected tissues.



DNA Microarray Chips

168

VisionArray®

The **VisionArray® products** are designed for the qualitative detection of specific DNA sequences by DNA/DNA hybridization on immobilized catcher molecules which are arranged on a glass chip. All capture sequences and positive controls are set up on the VisionArray® Chips as duplicates.

 ZytoVision GmbH

Advantages of VisionArray®


-  Quick & easy 1 hour protocol
-  Automated evaluation using a VisionArray® Analyzer Software – simple visualization & quick analysis in just a few minutes

Sample Collection


For the detection of DNA sequences with the VisionArray® system, the following raw material can be used for DNA extraction; depending on the VisionArray® Chip used:
Formalin-fixed, paraffin-embedded (FFPE) tissue or cell samples

Method


Step 1: Amplification and Labeling in a PCR
The DNA is extracted from, e.g., FFPE samples and is used as a template for PCR. Biotinylated primers are used to amplify and label different sections of the target sequences. The human HLA-DQA1 gene is also amplified and serves as a PCR positive control and as a genomic control.



Step 2: Hybridization on the Glass Chip
After amplification, the biotinylated sequences hybridize to complementary DNA capture sequences on the glass chip.



Step 3: Detection and Visualization
Specifically bound and biotinylated sequences are visualized by secondary marking with a streptavidin-peroxidase conjugate and a staining with tetra-methylbenzidine. After color development, evaluation is performed using a VisionArray® Analyzer Software.



VisionArray® HPV Chip 1.0


Intended Purpose
The VisionArray® HPV Chip 1.0 is intended to be used for the qualitative detection and genotyping of PCR-amplificates of 41 clinically relevant human papil-loma virus (HPV) genotypes that have been produced with the help of the VisionArray® HPV PreCise Master Mix (Prod. No. ES-0007-50) from formalin-fixed, paraffin-embedded specimens, such as cervical carcinoma or head and neck squamous cell carcinoma. The chip is intended to be used in combination with a VisionArray® Software.

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

The product is intended to be used as an aid to the differential diagnosis of cervical carcinoma or head and neck squamous cell carcinoma and therapeutic measures should not be initiated based on the test result alone.

Chip Description
The components of the product are the chip as well as the VisionArray® HPV Chip File 1.0.

Positioning of the capture sequences on the chip:



-  High Risk
-  Probably High Risk
-  Low Risk
-  Guide Dots (GD)/Positive Control (+)



*HPV 55 is classified by now as subtype of HPV 44, but is still labeled HPV 55 for consistency reasons.

Product	Manufacturer	Amount	Reg. Status	Order No.
VisionArray® HPV Chip 1.0	ZytoVision GmbH	10 Arrays (10 Tests)	CE/IVD	VA-0001-10

VisionArray® Myco Chip 2.0


Intended Purpose
The VisionArray® MYCO Chip 2.0 is intended to be used with a VisionArray® Analysis Package for the qualitative detection and identification of PCR amplicates of the genera Mycobacterium, Mycobacteroides, Mycolicibacillus, Mycolicibacter, and Mycolicibacterium as well as several clinically relevant mycobacterial species that have been produced with the help of the VisionArray® MYCO PreCise Master Mix 2.0.

The mycobacterial genera comprise more than 140 species, which, for the purpose of diagnosis and treatment, have been grouped into three categories: M. tuberculosis complex (MTC), M. leprae, and non-tuberculous mycobacteria (NTM). The majority of the Mycobacterium species belongs to the NTM group and can be found in different environments. Many of these bacteria cause life-threatening infections in humans and in recent years, the mortality and morbidity associated with NTMs has increased especially in immunocompromised patients worldwide. Treatment of NTMs is specific to each species and therefore a clear distinction between the present species is of extreme importance.

Reliable and rapid molecular diagnostics are the basis of an adequate therapy that is given by the VisionArray® MYCO Chip 2.0.

References
Griffith DE, *et al.* (2007) Am J Respir Crit Core Med 175: 367-416.
Gupta RS, *et al.* (2018) Front Microbiol 9: 67.
Oren A & Carrity GM (2019) Int J Syst Evol Microbiol 69: 597-9.
Perez-Martinez I, *et al.* (2013) BMC Res Notes 6: 531.
Simons S, *et al.* (2011) Emerg Infect Dis 17: 343-9.
Tortoli E (2009) Clin Microbiol Infect 15: 906-10.

Chip Description
The VisionArray® MYCO Chip 2.0 is designed to detect several clinically relevant mycobacterial species. All capture sequences and the positive control are set up on the Chip as duplicates and the guide dots as triplicates. The signals are visible on the Chip as dark blue areas. The automated evaluation of the results is performed by a VisionArray® Software.




M. tuberculosis (MTC) complex	
1	M. tuberculosis complex (ITS Region)
2	M. tuberculosis complex (IS6110 Region)
Nontuberculous Mycobacteria (NTM)	
3	M. abscessus
4	M. avium / M. intracellulare complex
5	M. chelonae
6	M. chimaera
7	M. fortuitum
8	M. genavense
9	M. goodii
10	M. haemophilum
11	M. kansasii
12	M. malmoense
13	M. marinum / M. ulcerans
14	M. scrofulaceum / M. parascrofulaceum
15	M. simiae
16	M. smegmatis
17	M. szulgai
18	M. xenopi

Product	Manufacturer	Amount	Reg. Status	Order No.
VisionArray® MYCO Chip 2.0	ZytoVision GmbH	10 Arrays (10 Tests)	CE/IVD	VA-0005-10

VisionArray® FUNGI Chip 1.0

Intended Purpose
The VisionArray® FUNGI Chip 1.0 is intended to be used for the qualitative detection and genotyping of PCR-amplicates of 30 clinically relevant fungi genotypes that have been produced with the help of the VisionArray® FUNGI PreCise Master Mix 1.0 (Prod. No. ES-0009-50) from formalin-fixed, paraffin-embedded specimens. The chip is intended to be used in combination with a VisionArray® Software.

Chip Description
The components of the product are the chip as well as the VisionArray® FUNGI Chip File 1.0. Positioning of the capture sequences on the chip:



Aspergillus	
1	Aspergillus flavus
2	Aspergillus fumigatus
3	Aspergillus nidulans / quadrilineatus
4	Aspergillus niger
5	Aspergillus terreus
6	Aspergillus versicolor
Candida	
7	Candida albicans
8	Candida auris
9	Candida dubliniensis
10	Candida parapsilosis
11	Candida tropicalis
12	Clavispora lusitanae
13	Kluyveromyces marxianus
14	Meyerozyma guilliermondii
15	Nakaseomyces glabratus
16	Pichia fermentans
17	Pichia kudriavzevii
18	Pichia norvegensis
19	Wickerhamomyces anomalus
Other	
20	Cryptococcus neoformans
21	Fusarium spp.
22	Lichtheimia corymbifera
23	Mucor spp.
24	Paecilomyces variotii
25	Pneumocystis jirovecii
26	Purpureocillium lilacinum
27	Rhizomucor pusillus
28	Rhizopus spp.
29	Scedosporium spp.
30	Trichophyton/Microsporum

Product	Manufacturer	Amount	Reg. Status	Order No.
VisionArray® FUNGI Chip 1.0	ZytoVision GmbH	10 Arrays (10 Tests)	RUO	VA-0006-10

VisionArray® PCR and Detection

VisionArray® Detection Kit

For hybridization and detection of PCR products on VisionArray® Chips

Product	Manufacturer	Amount	Reg. Status	Order No.
VisionArray® Detection Kit Incl.: Hybridization Solution, 1 ml; Detection Solution, 5 ml; Blue Spot Solution, 5 ml; 100x Wash Buffer, 250 ml	ZytoVision GmbH	1 Kit (50 Tests)	CE/IVD	VK-0003-50

VisionArray® PCR Reagents

For contamination-free amplification and biotinylation of target sequences with a high quality heat stable Taq polymerase

Product	Manufacturer	Amount	Reg. Status	Order No.
VisionArray® HPV PreCise Master Mix Incl.: VisionArray® HPV Primer; dNTP/dUTP Solution; PreCise Taq DNA Polymerase; PCR-Buffer; MgCl ₂ ; Uracil-DNA Glycosylase	ZytoVision GmbH	1 Kit (50 Tests)	CE/IVD	ES-0007-50
VisionArray® MYCO PreCise Master Mix 2.0 Incl.: VisionArray® MYCO Primer; dNTP/dUTP Solution; PreCise Taq DNA Polymerase; PCR-Buffer; MgCl ₂ ; Uracil-DNA Glycosylase	ZytoVision GmbH	1 Kit (50 Tests)	CE/IVD	ES-0008-50


VisionArray® PCR Reagents for Research Use Only

For contamination-free amplification and biotinylation of target sequences with a high quality heat stable Taq polymerase

Product	Manufacturer	Amount	Reg. Status	Order No.
VisionArray® FUNGI PreCise Master Mix 1.0 Incl.: VisionArray® FUNGI Primer; dNTP/dUTP Solution; PreCise Taq DNA Polymerase; PCR-Buffer; MgCl ₂ ; Uracil-DNA Glycosylase	ZytoVision GmbH	1 Kit (50 Tests)	CE/IVD	ES-0009-50

VisionArray® SingleScan Software


- Simple visualization and quick analysis of the VisionArray® Chip data
- Analysis of a Chip and the report of the results can be achieved in just a few minutes
- Program navigation is easy and intuitive for the user



Product	Reg. Status	Order No.
VisionArray® SingleScan Software	CE/IVD	E-4301-1

VisionArray® MultiScan Software

- Simple visualization and quick analysis of up to 6 VisionArray® Chips simultaneously
- All available VisionArray® Chips can be automatically detected by the software offering maximum flexibility
- Analysis of the Chips and the report of the results can be achieved in just a few minutes



Product	Reg. Status	Order No.
VisionArray® MultiScan Software	CE/IVD	E-4302-1



Publisher: ZYTOMICS

Contact Details: www.zytomics.com/distribution



Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx