

Biological Evaluation of MILLIPLEX® MAP Human and Mouse Th17 Panels




Robert Keith*, Alex Munmun Banerjee, Tim Warmke, and Qiang Xiao

R&D Dept., EMD Millipore Corporation, St. Charles, MO.
*Corresponding Author: Robert.Keith@merckgroup.com

Introduction


Th17 cells are recently discovered T helper cells that play important roles in the establishment and maximization of the capabilities of the immune system. We recently developed MILLIPLEX® MAP Human and Mouse Th17 multiplex assays for simultaneous measurement of the cytokines that are secreted from Th17 cells and/or regulate Th17 cell differentiation and activation. In this study, we used the two Th17 panels to analyze the cytokine secretion in human and mouse PBMCs stimulated with lipopolysaccharide (LPS), concanavalin A (Con A) or Phytohemagglutinin (PHA). We also examined the effect of Adiponectin pretreatment on the LPS- or Con A-induced PBMC response. We examined cytokine levels in human plasma from sepsis, rheumatoid arthritis and lupus patients and examined *in vivo* mouse plasma cytokine response to LPS treatment. Our results demonstrated that LPS, Con A or PHA treatments induced significantly increased secretion of many cytokines including Th17 cell specific cytokines from human and mouse PBMCs (e.g. IL-6, TNF α , IL-17A, IL-17F and IL-22). The diseased human and the *in vivo* LPS-treated mice plasma samples also expressed increased levels of Th17 cell specific plasma cytokines. Adiponectin may regulate certain LPS-induced or Con A-induced cytokine secretion from PBMCs. In conclusion, the Th17 panels will be useful tools for Th17 cell-related cytokine profiling in various biological samples.



LPS, a component of Gram-negative bacterial cell walls, stimulates the innate immune response via TLR (Toll-like receptor) 4. Con A and PHA are lectins which bind to cell membrane carbohydrates and induce cell agglutination and T-cell mitosis and differentiation. Adiponectin, also called GBP-28, apM1, AdipoQ, and Acrp30, is a novel adipose tissue-specific protein that has structural homology to collagen VIII and X and Complement factor C1q, and circulates in human plasma at high levels. There has been a reported adiponectin attenuation of the TNF α response *in vivo* which indicates an anti-inflammatory role for the protein (Xu et al., 2003).

Methods

Luminex® 200 system. This is a compact unit consisting of an analyzer, a computer, and software (Luminex Corporation, Austin, TX).
Microspheres. Magnetic microsphere beads were purchased from Luminex Corp. Each set of beads is distinguished by different ratios of two internal dyes yielding a unique fluorescent signature to each bead set. Capture antibodies were covalently coupled to the carboxylate-modified magnetic microsphere beads.



Sample Preparation. PBMCs (Bioreclamation) were thawed, washed and resuspended in RPMI Media (Gibco) containing 10% FCS and 1% Penicillin/Streptomycin at 1 million cells per mL. The PBMCs were aliquoted at one million cells per well and incubated overnight at 37°C. The next day either LPS, Con A or PHA (Sigma) were added for final concentrations of 10, 20 or 200 ng/mL, respectively. For Adiponectin pretreatment studies, the Adiponectin (Sigma) was added to select wells to 10 ng/mL final concentration, 5 min prior to the addition of LPS, Con A or PHA. Conditioned media was collected at various time points, centrifuged at 15,000 RPM for 10 min, and the cell-free supernatants were stored at -80°C prior to assaying in Multiplex format. Human disease plasma and *in vivo* LPS-challenged mouse plasma samples were obtained from Bioreclamation.

Immunoassay Procedure. The multiplex assay was performed in a 96-well plate. The detailed protocol is as follows: wet the plate with 150 μ L assay buffer for 10 min and decant. Add 25 μ L standards or samples, 25 μ L beads, 25 μ L assay buffer in sample wells or 25 μ L matrix in standard wells and incubate at O/N at 4°C. Wash the beads two times then add 25 μ L biotinylated detection Ab cocktail and incubate at RT for 1 hour. Add 25 μ L Streptavidin-Phycoerythrin and further incubate at RT for 30 min. Lastly, wash beads two times, add 100 μ L sheath fluid and read on Luminex instrumentation.

MILLIPLEX is a registered trademark and EMD Millipore and the M logo are trademarks of Merck KGaA, Darmstadt, Germany. Luminex is a trademark of Luminex Corp, Austin, TX.
© 2011 EMD Millipore Corporation. All rights reserved. Lit. No. PS321002EN00_09/11

Results: Human Th17 Multiplex Panel

Table 1. Stimulation of Human PBMC

Human PBMCs were treated with LPS, Con A or PHA at 37°C for 48 hr and cell-free samples were collected and assayed as described. Approximate cytokine responses are indicated as (+) 40 to 100 pg/mL, (++) 101 to 1000 pg/mL or (+++) >1000 pg/mL. Unstimulated samples were 0-40 pg/mL for most analytes.

Cytokine Response	LPS	Con A	PHA
IL-1 β	+++	++	++
IL-2	+	+++	+++
IL-4		++	+
IL-5		+	+
IL-6	+++	+++	+++
IL-9		++	++
IL-10	++	++	++
IL-13		+++	++
IL-17A		++	++
IL-17F		+++	++
IL-21		+	+
IL-23	+	+	
IL-27	++		
IL-31	++		
IFN γ		+++	+++
GM-CSF	++	+++	+++
MIP-3 α	++	++	++
TNF α	+++	+++	+++
TNF β		++	++

Table 2. Stimulation of Human PBMC: Representative Study

Human PBMCs were treated with LPS or Con A alone or after 5 min pretreatment with Adiponectin then incubated for 48 hr at 37°C and cell-free samples were collected and assayed as described. Values are pg/mL.

Sample	IL-1 β	IL-2	IL-4	IL-5	IL-6	IL-9	IL-10	IL-13	IL-17A	IL-17F	IL-21	IL-23	IL-27	IL-31	IFN γ	GM-CSF	MIP-3 α	TNF α	TNF β
Unstimulated	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
LPS	19	37	8	2	10	7	23	4	0	38	0	2	0	5	5	25	150	4	80
Con A	61	37	10	11	20	42	107	8	6	187	10	5	11	33	23	47	366	5	108
PHA	13	16	7	2	6	3	12	0	0	0	0	0	0	0	0	150	4	88	0
Adiponectin	15	20	10	2	18	1	18	0	0	150	0	0	0	0	0	176	0	75	0
LPS + Adiponectin	19	37	8	2	10	7	23	4	0	38	0	2	0	5	5	25	150	4	80
Con A + Adiponectin	61	37	10	11	20	42	107	8	6	187	10	5	11	33	23	47	366	5	108
PHA + Adiponectin	13	16	7	2	6	3	12	0	0	0	0	0	0	0	0	150	4	88	0

Table 3. Human Disease Samples

Human disease plasma samples were obtained from Bioreclamation and assayed in the Human Th17 Panel as described. Values are the sample range in pg/mL for 16 sepsis, 4 rheumatoid arthritis (RA) and 4 lupus samples.

Cytokine Response	Sepsis (n=16)	RA (n=4)	Lupus (n=4)
IL-1 β	0-129	0-143	0-30
IL-2	0-160	0-514	0-108
IL-4	0-735	0-3	ND
IL-5	0-88	0-5	ND
IL-6	31-703	0-373	0-11
IL-9	0-185	0-143	ND
IL-10	0-160(3)	0-2	0-8
IL-12p70	0-126	0-185	0-23
IL-13	12-300	0-11	0-1
IL-15	9-80	0-89	0-13
IL-17A	0-111	0-1	ND
IL-17E/IL-25	0-12479	0-18634	0-951
IL-17F	15-1038	0-1237	0-79
IL-21	0-21	0-34	0-3
IL-22	37-2042	0-729	ND
IL-23	0-13183	0-12472	0-1685
IL-27	2333-17179	505-8777	429-804
IL-28A	200-2974	0-1123	1-313
IL-31	0-4614	0-3013	0-17
IL-33	0-1202	0-320	0-10
IFN γ	0-445	0-30	ND
GM-CSF	61-3262	0-1493	0-330
MIP-3 α	48-1207	0-127	2-48
TNF α	20-4355	4-72	2-16
TNF β	0-386	0-474	0-13

Results: Mouse Th17 Multiplex Panel

Table 4. Stimulation of Mouse PBMC

Mouse PBMCs were treated with LPS, Con A or PHA at 37°C for 48 hr, then cell-free samples were collected and assayed as described. Approximate cytokine responses are indicated as (+) 1 to 20 pg/mL, (++) 21 to 500 pg/mL or (+++) >500 pg/mL. Unstimulated samples were not detectable for most samples except IL-6 at approximately 30 pg/mL.

Cytokine Response	LPS	Con A	PHA
IL-1 β	+		
IL-2		+	+
IL-4		+	
IL-5		+	
IL-6	+++	++	++
IL-10	+		
IL-17A		++	+
IL-17F		++	+
IL-22		+	
MIP-3 α	+		
TNF α	++		

Table 5. *in vivo* LPS Effect on Mouse Plasma Cytokine Levels


8-12 week old CD-1 Mice were injected with 1 mg/kg E.Coli 055:B5 LPS (Sigma-Aldrich) suspension in saline, IP, disodium EDTA plasma was collected at various time points and assayed as described. N=8, pg/mL (SEM)

hr	IL-1 β	IL-6	IL-10	IL-12p70	IL-13	IL-15
0	3 (3)	31 (22)	53 (16)	44 (17)	198 (55)	57 (5)
1	11 (3)	7648 (352)	3053 (591)	106 (38)	363 (58)	75 (17)
2	102 (46)	7129 (871)	1855 (516)	80 (12)	339 (30)	48 (5)
4	60 (13)	5542 (1017)	1192 (429)	77 (10)	108 (20)	44 (6)
6	20 (6)	2616 (921)	442 (93)	114 (29)	0 (0)	124 (48)

hr	IL-17A	IL-17E	IL-17F	IL-22	IL-23	IL-27
0	45 (26)	3914 (1456)	40 (8)	158 (34)	460 (120)	1383 (494)
1	94 (14)	6049 (694)	15 (6)	190 (96)	962 (190)	2093 (535)
2	102 (11)	6392 (476)	82 (37)	923 (163)	680 (77)	2183 (433)
4	113 (21)	5044 (617)	379 (69)	1278 (250)	612 (70)	1961 (384)
6	132 (22)	3590 (510)	365 (151)	661 (84)	538 (48)	2805 (992)

hr	IL-28B	CD40L	GM-CSF	IFN γ	MIP-3a	TNFa
0	137 (15)	157 (24)	45 (20)	20 (12)	96 (19)	20 (6)
1	171 (38)	167 (48)	70 (12)	27 (10)	157 (73)	1871 (489)
2	239 (59)	131 (8)	188 (52)	38 (9)	3307 (758)	1017 (376)
4	199 (17)	110 (5)	57 (13)	1052 (370)	8360 (5979)	164 (42)
6	228 (45)	201 (26)	39 (12)	592 (246)	352 (61)	70 (20)

Figures 1 & 2. *in vivo* LPS-effect on Mouse Plasma Cytokine levels. Data from Table 3



Summary

- LPS, Con A or PHA treatment induced secretion of cytokines from Human or Mouse PBMC including Th17 cell specific cytokines.
- Human sepsis, RA and lupus plasma samples contain significant levels of many cytokines including Th17 cell specific cytokines.
- Mouse *in vivo* LPS-stimulated plasma samples contain significant levels of many cytokines including Th17 cell specific cytokines.
- Adiponectin may regulate cytokine secretion from PBMC.
- MILLIPLEX Th17 panels are powerful tools for cytokine profiling in biological samples.

References:

- Kom T, et al. IL-17 and Th17 Cells. Annu Rev Immunol 27:485-517 (2009).
- Xu A, et al. The fat-derived hormone adiponectin alleviates alcoholic and nonalcoholic liver diseases in mice. J Clin Invest 112:91-100 (2003).