
# Adam, Rouilly SERVING MEDICAL EDUCATION WORLDWIDE



NEW

AR405/AR405-B


DIGITAL EYE AND EAR EXAMINATION TRAINER SET



The ideal training set for eye and ear examination



# AR405 DIGITAL EYE AND EAR EXAMINATION TRAINER SET, WHITE AR405-B DIGITAL EYE AND EAR EXAMINATION TRAINER SET, BLACK



Adam, Rouilly's new **AR402 DIGITAL EAR EXAMINATION TRAINER** and popular **AR403 DIGITAL EYE RETINOPATHY TRAINER** are now available together in an **educational set** ideal for any clinical skills setting.

The AR402 DIGITAL EAR EXAMINATION TRAINER features 48 ear conditions and the AR403 DIGITAL EYE RETINOPATHY TRAINER features 36 eye conditions (see rear of this leaflet) totalling 84 ear and eye conditions that form the ideal basis for both introductory and advanced ear and eye examination, as well as training in ophthalmoscope and otoscope use.

Both models feature high resolution displays on which each condition is represented, as well as life-size proportions and realistic features to greatly enhance any training scenario.



AR405 Digital Eye and Ear Examination Trainer Set, White

# **Features**



Cost saving advantage as a two model set



Simple to set up and use



48 ear conditions



36 eye conditions



High resolution digital displays



Easy to use, digital control for eye and ear conditions



Examination cover to hide displays of condition numbers



Battery or worldwide mains power compatible



Sleep mode to conserve power








Simple digital control for eye and ear conditions





Cover to hide condition LED indicators under examination conditions

# **Includes**

- Instruction manuals describing details of each condition
- Low voltage mains adaptors with worldwide plug fixings
- Removable Shoulder base for Ear Model
- Rigid carrying cases

Optional Extras
ARR00001 OTOSCOPE
ARR00002 OPHTHALMOSCOPE



Worldwide mains power



Battery Power (batteries not included)

# **Conditions and Diseases of the Retina** Diabetic Retinopathy (According to National Screening Programmes England, Europe and WHO)

- Background Diabetic Retinopathy 1.
- Maculopathy 2.
- 3 Pre-Proliferative Diabetic Retinopathy
- 4. Proliferative Diabetic Retinopathy
- New Vessels Disc 5
- 6.
- Photocoagulation 7
- Ungradable

### **Important and Common Retinal Conditions**

- 9 Normal
- 10. Glaucoma
- Papilloedema 11
- 12. Optic Atrophy
- Age-related Macular Degeneration (Dry) 13.
- 14. Hypertensive Retinopahty
- Central Retinal Vein Occlusion 15.
- Central Retinal Artery Occlusion 16.
- 17. Drusen
- 18. Retinitis Pigmentosa
- Myelinated Nerve Fibres 19.
- 20 High Myopia
- 21 Branch Retinal Vein Occlusion

## **Important and Less Common Retinal Conditions**

- Pre-Retinal Vein Occlusion 22
- Multiple Retinal Haemorrhages
- Retinal Detachment 24
- Angioid Streaks
- Benign Disc Neavus 26.
- 27. Malignant Melanoma
- 28. Macular Hamorrhage
- 29. Choroidal Naevus
- 30. Macular Scar (Toxoplasma)
- 31. Cytomegalovirus Retinitis
- 32. Lipaemia Retinalis
- 33 Medusa Head
- 34. Myopic Crescent - Normal Choroidal Vessels
- Sub Hyaloid Haemorrhage Resolving 35
- 36. Macular Burn

### Conditions and Diseases of the Ear

- 1 Normal I
- 2 Normal II
- 3. Ear Wax (Cerumen)
- Swimmer's Osteoma 4.
- 5. Fungal Ear 1
- Fungal Ear 2 6.
- Acute Viral Ear
- Acute Secretory Otitis Media I 8.
- 9. Resolving Secretory Otitis Media
- 10. Acute Secretory Otitis Media II
- 11. Acute Suppurative Otitis Media III
- 12. Perforation following an Acute Suppurative Otitis Media (ASOM)
- Childhood Glue Far 13.
- Glue Ear in a Child with a Dermoid Cyst in the Fardrum
- Adult Glue Ear
- A Standard Ventilation Tube in the Membrane 16.
- 17. Infected Mini Grommet with Otitis Externa Secondary to the Mucus Discharge

- 18 Permayent Tube in Place
- 19 Large Perforation of the Tympanic Membrane
- 20. A Posterior Perforation of the Tympanic
- 21. Two Small Traumatic Perforations following a Blow to the Ear
- 22. Subtotal Perforation of the Tympanic Membrane
- Perforation with Tympanosclerosis 23.
- 24. Grommet Scar Healed
- 25. Tympanosclerosis of Tympanic Membrane
- 26. Posterior Retraction
- 27. Retraction onto Long Process of the Incus
- Retraction with Loss of Long Process of the Incus and Keratin Trail 28
- 29. Retraction with loss of Long Process of the Incus
- Posterior Retraction Pocket onto Jugular Bulb and 30. with Middle Ear Fluid
- 31. Retraction with Early Keratin Build Up
- 32. Childhood Attic Retraction
- Deep Attic Retraction

- 34 Attic retraction Accumulation Keratin -Underlying Cholesteatoma
- Extensive Accumulation with Cholesteatoma in Middle Ear 35.
- 36. Wet Cholesteatoma
- Clean Dry Reconstructed Mastoid Cavity 37.
- 38. Old Style Mastoid Cavity with Residual Cholesteatoma
- 39. Mastoid Cavity with Fistula into Lateral Semicircular Canal
- 40. Congenital Cholesteatoma
- Large Congenital Cholesteatoma 41.
- 42. Ear Canal Cholesteatoma I
- Ear Canal Cholesteatoma II 43.
- 44. Keratosis Obturans
- Glomus Tympanicum Tumours 45.
- Glomus Jugulare Tumour
- 47. Foreign bodies
- 48. Aural Polyp



# 2 Year Guarantee

All products manufactured by Adam, Rouilly are covered by our full 2 Year Guarantee. This guarantee applies to models which have been used correctly and covers durability and functionality.



B16 L1

