



MILESTONE


H E L P I N G
P A T I E N T S

FormSAFE: still formalin... no fumes!



Multilayered liquid container

Stratified compound system: floating shield
+10% Neutral Buffered Formalin



Floating barrier →
Formalin

- Immiscible in formalin
- Non Toxic
- “Inert” for biopsies

60 ML CONTAINER W 30 ML FORMALIN
SAFETY FOR **OR** AND **AP**
STANDARD PROCEDURE AND EASY TO USE
EASY SCREW CAP
OPTIMAL FIXATION, NO CHANGE OF FIXATIVE

MILESTONE



HELPING
PATIENTS