

HMGB1

Multi-tasker of the

innate immune system.

IMMUNOLOGY / CYTOKINES



At the center of HMGB1 research.

Over the years, HMGB1 has increasingly attracted the attention of academic as well as clinical researchers and has become one of the most intriguing molecules within the complete arsenal of proteins of the innate immune system.


HMGB1 is a nuclear protein that can bind and bend DNA, but it can also be released to the extracellular environment where it exerts additional functions.

The unraveling of the post-translational modifications has led to a better understanding of the mechanism of its translocation and its function within the immune system.

As of today it is known that acetylation of HMGB1 is necessary, so that it can be translocated from the nucleus to the cytosol, that the fully reduced HMGB1 isoform has chemotactic effects and causes the migration of immune cells to the danger zone and that the disulfide HMGB1 isoform acts as a pro-inflammatory cytokine inducing production of further cytokines.

It will now also allow more in-depth studies into the role of HMGB1 in numerous pathologies, including

- Respiratory diseases such as COVID-19, ARDS, pneumonia or else
- Autoimmune diseases
- Sepsis
- Cancer
- Diabetes
- Cardial infarction
- Stroke
- ... any inflammatory process



Different HMGB1 redox isoforms determine their location and biological function.

TECAN AT THE CENTER OF HMGB1 RESEARCH

Through its collaboration with HMGB1 experts from academia and commercial partners, Tecan is at the forefront of new developments. We at Tecan offer the most complete range of products - including a highly sensitive ELISA for the quantitative measurement of HMGB1 - which have been widely used and are cited in many publications.

(If you wish to receive a complete list of the referenced articles, please contact us at ibl@tecan.com)



HMGB PRODUCTS

Immunoassay	Catalog#	Determ.	Assay range	Incubation time	Sample type
HMGB1 Express ELISA	30164033	96 well	2.5 – 80 ng/mL	1st: 37 °C, 2h, 2nd: RT, 1h, 3rd: RT, 20 min	Serum, plasma (EDTA, Citrate)

For Research use only

Antibodies*	Catalog#	Quantity	Intended use
Anti-HMGB1 Rabbit IgG PoAb	ST326052219	50 µg	WB (1-2 µg/mL), IHC ²⁴ , immunofluorescence ²⁵ and immunoprecipitation ²⁶
Anti-HMGB1 Chicken IgY PoAb	ST326052226	50 µg	WB ²⁷ (1-2 µg/mL)
Anti-HMGB1 Chicken IgY Neutralising PoAb	ST326052233	1 mg	WB and neutralization experiments ²⁸⁻³¹ (2 mg/kg/mouse)
Anti-HMGB1,2 Mouse IgG1 MoAb	ST326052240	50 µg	WB (1-2 µg/mL)
Anti-HMGB1 [DPH1.1] Mouse IgG1 MoAb	REHM901	50 µg	WB ³² (1 µg/mL), immunofluorescence ³² , IHC, blocking experiments in cell migration assay ³² and blocking recruitment of inflammatory cells to sites of necrosis and infection in vitro and in vivo (220 µg/mouse) ³²
	REHM902	250 µg	
	REHM903	1 mg	

Protein isoforms and related proteins*	Catalog#	Quantity	Intended use
Fully reduced HMGB1, LPS-free	REHM114	500 µg	To study HMGB1 cell migratory effects in vitro ³³ and in vivo
	REHM115	100 µg	
	REHM116	50 µg	
Disulfide HMGB1, LPS-free	REHM120	500 µg	To study HMGB1 induced cytokine effects in proinflammatory processes ³⁵
	REHM121	100 µg	
	REHM122	50 µg	
BoxA from HMGB1, LPS-free	REHM012	100 µg	To study HMGB1-RAGE interaction by using BoxA as an antagonist for HMGB1 ^{33,34}
	REHM013	500 µg	
	REHM014	2 mg	
BoxB from HMGB1, LPS-free	REHM052	100 µg	To study HMGB1 induced cytokine activity by measuring BoxB as a read out ³⁶
	REHM051	250 µg	
	REHM050	1 mg	
HMGB2, LPS-free	REHM151	500 µg	To study HMGB2 induced fibroblast migration ³⁷
	REHM152	100 µg	
	REHM153	50 µg	

* Distributed by IBL International

Manufactured by

IBL International GmbH
Flughafenstrasse 52a
22335 Hamburg
Germany

Phone: +49 (0)40-53 28 91-0
Fax: +49 (0)40-53 28 91-11
Email: IBL@Tecan.com
www.tecan.com/immunoassays

.....
Australia +61 3 9647 4100 **Austria** +43 62 46 89 330 **Belgium** +32 15 42 13 19 **China** +86 21 220 63 206 **France** +33 4 72 76 04 80 **Germany** +49 79 51 94 170
Italy +39 02 92 44 790 **Japan** +81 44 556 73 11 **Netherlands** +31 18 34 48 17 4 **Nordic** +46 8 750 39 40 **Singapore** +65 644 41 886 **Spain** +34 93 93 595 25 31
Switzerland +41 44 922 89 22 **UK** +44 118 9300 300 **USA** +1 919 361 5200 **Other countries** +41 44 922 81 11
.....

Tecan Group Ltd. makes every effort to include accurate and up-to-date information; however, it is possible that omissions or errors might have occurred. Tecan Group Ltd. cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided. Changes can be made at any time without notice. All mentioned trademarks are protected by law. In general, the trademarks and designs referenced herein are trademarks, or registered trademarks, of Tecan Group Ltd., Maennedorf, Switzerland. A complete list may be found at www.tecan.com/trademarks. Product names and company names that are not contained in the list but are noted herein may be the trademarks of their respective owners. For technical details and detailed procedures of the specifications provided please contact your Tecan representative. This may contain reference to applications and products which are not available in all markets. Please check with your local sales representative.

© 2021 Tecan Trading AG, Switzerland, all rights reserved. For disclaimer and trademarks please visit www.tecan.com
.....

www.tecan.com

 **TECAN.**



Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx