





THE ALL NEW TissueSAFE plus

The best is now better

1. Easier operation



New intuitive 7" touch display with dedicated Software, icon driven and user-friendly. Before selecting the vacuum program, user can press the icon "i" to list the suitable program for each tissue type.

2. ID data tracking




Enter specimen and/or patient case ID through barcode reader.

3. Specimen safety



To start a vacuum program, user must scan the bag ID. TissueSAFE plus is a closed system for vacuum bags.

4. Time/Temperature tracking



Data Logger Card: Specimen ID and temperature monitoring for documented and standardized transportation (valid only for TissueSAFE plus with integrated datalogger card version).

THE ALL NEW TissueSAFE plus

The best is now better

5. Full documentation



Built-in printer, for protocol summary, delivers an adhesive label.



Anatomical Pathology: additional details are displayed for both delivered and received Cases ID.



USB ports: for full documentation, allow to export all saved runs and events log of interest.

6. Accounts configuration




Create more users protected by username and password.

7. Unit maintenance



Replacement interval of Hepa and Charcoal activated filters controlled by software.



MILESTONE



MILESTONE Srl - Via Fatebenefratelli, 1/5 - 24010 Sorisole (BG) - Italy
 Tel: +39 035 4128264 - Fax: +39 035 575498
 www.milestonemedsrl.com - email: medical@milestonesrl.com

UNI EN ISO 9001: 2008 / UNI EN ISO 13485:2012 CERTIFIED

MILESTONE MEDICAL TECHNOLOGIES, INC.
 6475 Technology Avenue, Suite F, Kalamazoo, MI 49009 - USA
 Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949
 www.milestonemed.com - email: info@milestonemed.com

CAT369EN-001