

Features to support your specific requirements

LIFE GUARD
INDICATORS

Logistics:

Shipping and handling at **room temperature**, storing at cool temperature (ideally together with blood bags).



Application:

In warm or cool environment possible, activated or non-activated for **optimized working conditions**.



Activation:

After cooling below 6°C **touchless and effortless** activation with UV-C light by mobile or stationary devices in less than 5 seconds.



Safe Handling:

The patented "Life Guard" insulated surface rules out any temperature influence from outside. **Reliable protection** against unintended contact with the hand.



Unique color code:

Intuitive colors support the **right decisions** in the operating process:
non-activated – yellow
activated – dark blue
temperature exceeded – orange



Quality assurance:

Maximum safety through functional principle. Activation only possible after sufficient cooling of the blood bag at 6°C and below.



Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx

on point medicals GmbH

- Founded 2008 in Klagenfurt/Austria
- Supplier of worldwide leading temperature indicator technology
- Patents pending: EP2010/007009, EP2010/007010
- Partner of first-rate R&D institutions in Europe and the US
- Financing by private investors

LIFE GUARD
INDICATORS

- Please contact:
Alexander Höfinger
Managing Director

on point medicals GmbH
Phone: +43 (463) 410355-100
E-Mail: alexander.hoefinger@onpoint.at
Web: www.life-guard.at

LIFE GUARD
INDICATORS

"Life Guard"



Monitoring System for Blood Temperature



The multiple benefits of single blood bag temperature monitoring

Stop the wastage of blood

Closing the circuit of the blood bag supply chain makes it possible to reuse blood which otherwise is discarded.



Increase patient safety

Permanent monitoring of blood temperature from storage to distribution and back ensures an uncompromised level of quality. This removes the guesswork of what might have happened to your blood outside of the bloodbank.



Full compliance with regulatory standards

By using the "Life Guard" system it is possible to prove that your handling of blood has been compliant with all relevant regulatory standards at all times.



The way it works in a typical blood bag life cycle

