

Ordering information

Instruments


• CAPILLARYS 3 OCTA	1245
• CAPILLARYS 3 TERA	1246
• CAPILLARYS 3 TERA TLA	1248
• CAPILLARYS 3 TERA Loader	1249
• CAPILLARYS 3 TERA Conveyor	1250

Consumables

• CAPILLARYS 3 CAPICLEAN	2060
• CAPILLARYS 3 Wash Solution	2062
• CAPILLARYS 3 Reagent Cups (24x14 units)	2582
• CAPILLARYS 3 Bins for used Reagent Cups	2581

CAPILLARYS 3 | Beyond Separation


CAPILLARYS 3
OCTA



CAPILLARYS 3
TERA



CAPILLARYS 3
TERA TLA



CAPILLARYS 3
TERA MC3



Sebia | RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008 Evry cedex | France | Tel.: +33 1 69 89 80 80 | E-mail: sebia@sebia.com | www.sebia.com

Benelux | Tel.: +32 2 702 64 64 (BE) | Tel.: +31 20 610 8099 (NL) | E-mail: sebia.benelux@sebia.com • **Brasil** | Tel.: +55 11 3849 0148 | E-mail: brasil@sebia.com
Canada | Tel.: +1 (888) 725-3013 | Fax: +1 (438) 380-3780 | E-mail: canada@sebia.com • **China** | Tel.: +86 (21) 6360 1436 | Fax: +86 (21) 6361 2011 | E-mail: china@sebia.com • **Germany & Austria** | Tel.: +49 (0) 661 33081 | Fax: +49 (0) 661 31881 | E-mail: sebia@sebia.de • **India** | Tel.: +91 124 4262100 | E-mail: sebiaindia@sebia.com • **Italy** | Tel.: +39 055 24851 | Fax: +39 055 0982083 | E-mail: info@sebia.it • **Japan** | Tel.: (03) 5830-0139 | Fax: (03) 5830-0140
Portugal | Tel.: +351 211 224 653 | Fax: +351 211 305 079 | E-mail: portugal@sebia.com • **Spain** | Tel.: +34 93 208 15 52 | Fax: +34 93 458 55 86 | E-mail: sebia@sebia.es • **Singapore** | Tel.: +65 6223 2273 | E-mail: singapore@sebia.com • **Switzerland** | Tel.: +41 44 787 88 10 | E-mail: contact.ch@sebia.com
UK | Tel.: +44 (0)1276 600636 | Fax: +44 (0)1276 38827 | E-mail: sales@sebia.co.uk • **USA** | Tel.: +971 (0)4 878 69 77 • **USA** | Tel.: +1 (770) 446-3707 - Fax: +1 (770) 446-8511 | E-mail: info@sebia-usa.com

©2021 Sebia - Non-contractual illustration - 10/2021 - Ref: C3BIEN - Sebia reserves the right to modify at any time and without notice the information contained in this document intended for healthcare professionals. This document contains information intended for wide distribution and may therefore contain product details or information that is not available or valid in your country. Carefully read the instructions in the reagent package inserts and instrument manuals.

sebia
The new language of life

Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx

Your journey towards full lab automation with Sebia

Flexibility of solutions

Stand alone, workcell or Total Laboratory Automation, which one is the best for you today?

Service & support

Do you have a question? Just call our Sebia Knowledge Center

Traceability & security

Be assured that everything is under control, with full traceability for samples and reagents, capped tube management and mixing for whole blood

Result interpretation software

Powerful Phoresis Core software makes your patient management easy

Quality of results

Have confidence when reporting your patients' results with the best-in-class separation assay

Optimized workflow & productivity

High throughput autonomous instruments with less technician intervention – your patients benefit from faster generation of results

Multiparametric instruments with large range of tests*:

- Hb A1c (capillary blood & venous blood)
- Serum & Urine Protein Electrophoresis (SPE & UPE)
- Serum & Urine Immunotyping
- CDT/CDT_{IFCC}
- Hemoglobin Electrophoresis

Incomparable separation for a unique quality of result




CAPILLARYS
OCTA/TERA 3


Multiparametric stand alone instruments with 8/12 capillaries in parallel. A proven analytical excellence meets high throughput

CAPILLARYS
TERA MC 3

Workcell configuration with up to 3 CAPILLARYS 3 TERAs and a tube loader. Unbeatable autonomy: 700 samples can be processed without technician intervention



CAPILLARYS
TERA TLA 3



Track connectivity with point-in-space instrument. The analytical excellence is now fully automated