


# PRODUCT CATALOGUE





## **Content**

Bovine	4
Ovine	7
Porcine	8
Equine	9
Avian	10





### **Improved standards worldwide**

Boehringer Ingelheim Svanova contributes to the improvement of animal health and food safety on a global scale by delivering diagnostic solutions to prevent, monitor and eradicate animal diseases. Under the trade names, SVANOVIR® and SVANODIP® our tests have been successfully applied in eradication and control programmes for more than 20 years. In 2011, Svanova became a part of Boehringer Ingelheim Vetmedica, further strengthening our activities and services worldwide.



# Bovine

## IBR

### Infectious Bovine Rhinotracheitis Antibody Test

An indirect ELISA for the detection of specific antibodies to Infectious Bovine Rhinotracheitis in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® IBR-Ab	10-2100-02	strips	2	96	88
SVANOVIR® IBR-Ab	10-2100-10	plate	10	480	440
SVANOVIR® IBR-Ab (screen)	10-2100-50	strips	10	960	920

## BVDV

### Bovine Virus Diarrhoea Virus Antibody Test

An indirect ELISA for the detection of specific antibodies to Bovine Virus Diarrhoea Virus in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BVDV-Ab	10-2200-02	strips	2	96	88
SVANOVIR® BVDV-Ab	10-2200-10	plate	10	480	440
SVANOVIR® BVDV-Ab (screen)	10-2200-50	strips	10	960	920

### Bovine Virus Diarrhoea Virus p80 Antibody Test

A blocking ELISA for the detection of specific antibodies to Bovine Virus Diarrhoea p80 (NS3) in bovine, caprine and ovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BVDV p80-Ab	10-2210-02	strips	2	192	184

## BLV

### Bovine Leukaemia Virus gp51 Antibody Test

An indirect ELISA for the detection of specific antibodies to Bovine Leukaemia Virus gp51 in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BLV gp51-Ab	10-2351-02	strips	2	96	88
SVANOVIR® BLV gp51-Ab	10-2351-10	plate	10	480	440
SVANOVIR® BLV gp51-Ab (screen)	10-2351-50	strips	10	960	920

## BCV

### Bovine Coronavirus Antibody Test

An indirect ELISA for the detection of specific antibodies to Bovine Coronavirus in bovine serum and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BCV-Ab	10-2400-02	strips	2	96	88

## BRSV

### Bovine Respiratory Syncytial Virus Antibody Test

An indirect ELISA for the detection of specific antibodies to Bovine Respiratory Syncytial Virus in bovine serum and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BRSV-Ab	10-2500-02	strips	2	96	88

Tests - number of tests. Samples - number of samples, wells for kit controls excluded.

## PIV3

### **Bovine Parainfluenza Virus Type 3 Antibody Test**

An indirect ELISA for the detection of specific antibodies to Bovine Parainfluenza Virus Type 3 in bovine serum and milk samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® PIV3-Ab	10-2600-02	strips	2	96	88

## BRUCELLA

### **Brucella I-ELISA Antibody Test**

An indirect ELISA for the detection of specific antibodies to *Brucella* in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab I-ELISA	10-2700-10	strips	10	960	920

### **Brucella C-ELISA Antibody Test**

A competitive ELISA for the detection of specific antibodies to *Brucella* in bovine, ovine, caprine and porcine serum and plasma samples; which is capable of discriminating between naturally infected and S19 vaccinated cattle.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	10-2701-02	strips	2	192	176
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	10-2701-10	plate	10	960	880

## FMDV

### **Foot and Mouth Disease Virus 3ABC Antibody Test – Ruminant**

An indirect ELISA for the detection of specific antibodies to Foot and Mouth Disease Virus 3ABC in bovine, ovine and caprine serum samples; which is capable of discrimination between naturally infected and vaccinated animals.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® FMDV 3ABC-Ab ruminant	10-2801-10	plate	10	480	440

### **IAEA Standard Bovine FMDV Positive Control Sera**

A panel of reference sera for FMDV assays (lyophilised).

Product name	Article No.	Quantity
IAEA standard FMDV type A	33-2800-01	0.5 ml
IAEA standard FMDV type O	33-2800-02	0.5 ml
IAEA standard FMDV type Asia1	33-2800-03	0.5 ml
IAEA standard FMDV type SAT1	33-2800-04	0.5 ml
IAEA standard FMDV type SAT2	33-2800-05	0.5 ml
IAEA standard FMDV type SAT3	33-2800-06	0.5 ml

### **Foot and Mouth Disease Virus Antigen Test**

A penside test for the detection of Foot and Mouth Disease Virus in swab and tissue samples.

Product name	Article No.	No of tests
SVANODIP® FMDV-Ag	10-4100-13	20 test devices

### **Foot and Mouth Disease Virus Tissue Sample Extraction Kit**

A tissue sample extraction kit designed to extract the antigen from tissues.

Product name	Article No.	No of tests
SVANODIP® FMDV-Ag Extraction	36-4100-11	20 extraction vials

## RPV

### Rinderpest Virus Antigen Test

A penside test for the detection of Rinderpest virus in eye swab samples.

Product name	Article No.	No of tests
SVANODIP® RPV-Ag	10-4300-13	50 test devices

## OSTERTAGIA

### Ostertagia ostertagi Antibody Test

An indirect ELISA for the detection of specific antibodies to *Ostertagia ostertagi* in bovine bulk tank milk.

Product name	Article no.	Format	Plates	Tests	Samples
SVANOVIR® <i>O. ostertagi</i> -Ab	10-2940-02	strips	2	192	184

## NEOSPORA

### Neospora caninum iscom ELISA Antibody Test

An indirect ELISA for the detection of specific antibodies to *Neospora caninum* in bovine serum and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Neospora</i> -Ab	10-2950-02	strips	2	192	184

## ANAPLASMA

### Anaplasma marginale Antibody Test

An indirect ELISA for the detection of specific antibodies to *Anaplasma marginale* in bovine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>A. marginale</i> -Ab	10-2960-02	strips	2	192	184

## BABESIA

### Babesia bigemina Antibody Test

An indirect ELISA for the detection of specific antibodies to *Babesia bigemina* in bovine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>B. bigemina</i> -Ab	10-2970-02	strips	2	192	184

Tests - number of tests. Samples - number of samples, wells for kit controls excluded.

# OVINE / CAPRINE

## BDV

### **Border Disease Virus Antibody Test**

An indirect ELISA for the detection of antibodies to Border Disease Virus in ovine and caprine serum and milk samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BDV-Ab	10-8100-02	strips	2	96	88

## BRUCELLA

### **Brucella I-ELISA Antibody Test (small ruminants)**

An indirect ELISA for the detection of specific antibodies to *Brucella* in ovine and caprine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab I-ELISA o/c	10-8200-02	strips	2	192	184

### **Brucella C-ELISA Antibody Test**

An competitive ELISA for the detection of specific antibodies to *Brucella* in ovine, caprine, bovine and porcine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	10-2701-02	strips	2	192	176
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	10-2701-10	plate	10	960	880

## FMDV

### **Foot and Mouth Disease Virus 3ABC Antibody Test – Ruminant**

An indirect ELISA for the detection of specific antibodies to Foot and Mouth Disease Virus 3ABC in bovine, ovine and caprine serum samples; which is capable of discrimination between naturally infected and vaccinated animals.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® FMDV 3ABC-Ab ruminant	10-2801-10	plate	10	480	440



# PORCINE

## PRV

### **Pseudorabies/Aujeszky Disease Virus gE Antibody Test**

A blocking ELISA for the detection of specific antibodies to Pseudorabies/Aujeszky Disease Virus gE in porcine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® PRV gE-Ab	10-7161-02	strips	2	192	184
SVANOVIR® PRV gE-Ab	10-7161-10	plate	10	960	920

### **Pseudorabies/Aujeszky Disease Virus gB Antibody Test**

A blocking ELISA for the detection of specific antibodies to Pseudorabies/Aujeszky Disease Virus gB in porcine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® PRV gB-Ab	10-7262-02	strips	2	192	184
SVANOVIR® PRV gB-Ab	10-7262-10	plate	10	960	920

## ASFV

### **African Swine Fever Antibody Test**

An indirect ELISA for the detection of specific antibodies to African Swine Fever Virus in porcine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® ASFV-Ab confirmation	10-7300-02	strips	2	96	88
SVANOVIR® ASFV-Ab combi*	10-7300-50	plate+strips	9+1	912	872

\*Includes 9 screening plates and 1 confirmation plate

## TGEV/PRCV

### **Transmissible Gastroenteritis Virus/Porcine Respiratory Coronavirus Differentiating test**

A blocking ELISA for the detection and/or differentiation of specific antibodies to Transmissible Gastroenteritis and Porcine Respiratory Coronavirus in porcine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® TGEV/PRCV-Ab	10-7500-02	strips	2	96	88

## BRUCELLA

### **Brucella C-ELISA Antibody Test**

A competitive ELISA for the detection of specific antibodies in porcine, bovine, ovine and caprine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® Brucella-Ab C-ELISA	10-2701-02	strips	2	192	176
SVANOVIR® Brucella-Ab C-ELISA	10-2701-10	plate	10	960	880

## SALMONELLA

### **Porcine Salmonella Antibody Test**

An indirect ELISA for the detection of specific antibodies to *Salmonella* in pig serum and meat juice samples.

Product name	Article No.	Format	Plates	Tests	Samples
VETSIGN™ Salmonella-Ab	10-7601-02	strips	2	192	184
VETSIGN™ Salmonella-Ab	10-7601-10	plate	10	960	920

Tests - number of tests. Samples - number of samples, wells for kit controls excluded.

## EHV

### Equine Herpesvirus type 1 and 4 Discriminating test

An indirect ELISA for the detection of and/or discrimination of specific antibodies to Equine Herpesvirus type 1 and type 4 in equine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® EHV1/EHV4-Ab	10-3100-02	strips	2	64	52

# AVIAN

## APV

### **Avian Pneumovirus Antibody Test**

A blocking ELISA for the detection of specific antibodies to Avian Pneumovirus in serum and egg yolk from various avian species.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® APV-Ab	10-1100-02	strips	2	192	184

## MG

### **Mycoplasma gallisepticum Antibody Test**

A blocking ELISA for the detection of specific antibodies to *Mycoplasma gallisepticum* in serum and egg yolk from various avian species.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® MG-Ab	10-1300-02	strips	2	192	184

## MS

### **Mycoplasma synoviae Antibody Test**

A blocking ELISA for the detection of specific antibodies to *Mycoplasma synoviae* in serum and egg yolk from various avian species.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® MS-Ab	10-1400-02	strips	2	192	184

Tests - number of tests. Samples - number of samples, wells for kit controls excluded.



### **We test our tests**

The quality of our processes, products and services are continuously optimized and improved in accordance with our customer's needs and the high standards of veterinary diagnostic laboratories. All products are developed, evaluated and manufactured according to ISO 9001 to ensure consistency and reliability of our tests. With a combination of long experience and professionalism, Boehringer Ingelheim Svanova offers customers an ongoing, personal support. Our competence is part of our product.

[www.svanova.com](http://www.svanova.com)

svanova 

Distributed by Abacus ALS

AUS 1800 222 287 | NZ 0800 222 170

[info@abacus-als.com](mailto:info@abacus-als.com) | [www.abacus-als.com](http://www.abacus-als.com)

