

Novel Genitourinary (GU) Markers



NKX3.1 (EP356)

NKX3.1 is a transcription factor that is specific for prostate cancer. It is useful in distinguishing prostate from bladder and other non-prostate tissue in instances of metastasis. This nuclear marker has higher sensitivity and specificity than the already established PSA for prostate carcinoma.

Description	Cat. No.
0.1 ml concentrate	441R-14
0.5 ml concentrate	441R-15
1 ml concentrate	441R-16
1 ml predilute	441R-17
7 ml predilute	441R-18
5 Positive Control Slides	4415

ERG (EP111)

ERG is an important novel marker for the identification of vascular neoplasms due its strong and specific nuclear expression in endothelial cells. ERG is strongly expressed in Kaposi sarcoma, which is usually associated with HHV-8, as well as other vascular tumors such as hemangioendothelioma and

Description	Cat. No.
0.1 ml concentrate	434R-14
0.5 ml concentrate	434R-15
1 ml concentrate	434R-16
1 ml predilute	434R-17
7 ml predilute	434R-18
5 Positive Control Slides	4345

angiosarcoma. ERG has shown to be a valuable addition to an endothelial panel that includes Factor VIII, CD31, CD34, and D2-40. ERG is also helpful in identifying prostate adenocarcinoma. With its nuclear staining properties, ERG is able to differentiate prostate adenocarcinoma from low grade PIN and hyperplastic cells. ERG is a more specific marker than the established P504s for prostate cancer and should be used in a panel with PSA, P501s, P504s, and NKX3.1.

SALL4 (6E3)

SALL4 (6E3) is used in the identification of carcinomas of the gastrointestinal tract. It is seen to demonstrate high sensitivity for tumor cells in intratubular germ cell neoplasia, seminomas/dysgerminomas, embroyonal carcinomas, and yolk sac tumor. Anti-SALL4 also stains teratomas and mononucleated trophoblastic cells in choriocarcinomas.


Description	Cat. No.
0.1 ml concentrate	385M-14
0.5 ml concentrate	385M-15
1 ml concentrate	385M-16
1 ml predilute	385M-17
7 ml predilute	385M-18
5 Positive Control Slides	385S



GLUT3 (polyclonal)

GLUT3 is valuable in the identification of germ cell tumors and differentiation non-germ cell tumors and neoplasms. In IHC studies conducted at Cell Marque using IHC, GLUT3 demonstrated positive expression in testis and spermatozoa but had negative staining on brain tissue. GLUT3 is commonly used in panels with SALL4 and Oct-4.


Description	Cat. No.
0.1 ml concentrate	413A-14
0.5 ml concentrate	413A-15
1 ml concentrate	413A-16
1 ml predilute	413A-17
7 ml predilute	413A-18
5 Positive Control Slides	413S



Oct-4 (MRQ-10)

Oct-4 is a nuclear transcription factor that maintains and regulates pluripotency in embryonic stem and germ cells. It has a high sensitivity and specificity for seminoma/dysgerminoma, embryonal carcinoma, and the germ cell component of gonadoblastoma.


Description	Cat. No.
0.1 ml concentrate	309M-14
0.5 ml concentrate	309M-15
1 ml concentrate	309M-16
1 ml predilute	309M-17
7 ml predilute	309M-18
5 Positive Control Slides	309S



p57^{Kip2} (Kp10)

p57 is useful in the discrimination of complete hydatidiform mole (CHM) from partial hydatiform mole (PHM) and hydropic abortion.


Description	Cat. No.
0.1 ml concentrate	457M-94
0.5 ml concentrate	457M-95
1 ml concentrate	457M-96
1 ml predilute	457M-97
7 ml predilute	457M-98
5 Positive Control Slides	457S



Smoothelin (R4A)

Smoothelin is useful in distinguishing bladder muscualris mucosae (MM) from muscularis propria (MP) muscle bundles as it is exclusively observed in MP. Smoothlin's staining pattern of MP (positive) and MM (negative) makes IHC an attractive diagnostic tool for the sometimes difficult task of staging bladder urothelial carcinoma.

Description	Cat. No.
0.1 ml concentrate	377M-14
0.5 ml concentrate	377M-15
1 ml concentrate	377M-16
1 ml predilute	377M-17
7 ml predilute	377M-18
5 Positive Control Slides	377S
5 Positive Control Sildes	3//5



Glypican-3 (1G12)

Glypican-3 (GPC3) is a glycosylphospatidyl inositol-anchored membrane protein, which may also be found in a secreted form. In the context of testicular germ cell tumors, GPC3 expression is up-regulated in certain histologic subtypes, specifically yolk sac tumors and choriocarcinoma.

Description	Cat. No.
0.1 ml concentrate	261M-94
0.5 ml concentrate	261M-95
1 ml concentrate	261M-96
1 ml predilute	261M-97
7 ml predilute	261M-98
5 Positive Control Slides	2615





Rev. 0.0

©2016 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH and Cell Marque are trademarks of Sigma-Aldrich Co. LLC.

