CryoStore™ Freezing Bags

The OriGen CryoStore™ bag is indicated for freezing blood components. It has the traceability, durability, and safety that you need.

Traceability - Heat Welded Lots

A permanent lot number is heat welded into the tail of each bag for easy traceability.

Aliquots - Matched Bag and Tube Printing

Each bag has a donor tube with printed segments of approximately 100uL and a matched number on the bag. This can aid in keeping aliquot samples with the bag.

Durability - Better Material, Proven Performance


The hermetically sealed spike ports and film used in every CryoStore™ bag are made from a rugged EVA. Our customers have had great success with the strength of the bag over time.

Safety - Needle Free Injection Port

Standard tubing sets now have luer-actuated needle-free injection ports. These swabable ports maintain sterility of the sample through numerous actuations (white paper available) and can eliminate the possibility of needle sticks.

Ease of Use - Designed In

The donor tube is intentionally placed in the corner of the bag for easier de-bubbling and sealing. OriGen can provide your choice of several stock or non-stock tube sets (see reverse). You can order prefabricated multi-bag sets to reduce processing time and the possibility of contamination in the lab.







Technical Information

Product Code	Width, (cm)*	Length, (cm)*	Freeze volume min-max, ml
CS 25	9.0	7.3	10-25
CS 50	7.6	14.5	10-30
CS 250	12.7	15.3	30-70
CS 500	12.7	22.8	55-100
CS 750	12.7	25.3	80-190
CS 1000	15.0	28.5	120-275
CS 2000	23.0	33.5	180-450

^{*} Length and width are the overall bag footprint or cassette inside dimensions

The OriGen CryoStore $\ensuremath{^{\text{TM}}}$ bag is indicated for cryogenically freezing blood components. Carefully read Instructions for Use before using this product. Sterile, single use only, made in the USA.



Product Code	Description	Tube Set Description	Stock/ Non-Stock
CS25F	CryoStore 25 EVA Freezing Bag, 10 - 25 ml, "F" Tube Set, case of 24	Female luer with caps	Stock
CS25N	CryoStore 25 EVA Freezing Bag, 10 - 25 ml, "N" Tube Set, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS25NS	CryoStore 25 EVA Freezing Bag, 10 - 25 ml, "N" Tube Set, double aseptically packed, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS25NSW	CryoStore 25 EVA Freezing Bag, 10 - 25 ml, "N" Tube Set, double aseptically packed in OW, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS25XN	CryoStore 25 EVA Freezing Bag, 10 - 25 ml, "N" Tube Set, 8" donor tube, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port	Non-Stock
CS25Y4	CryoStore 25 EVA Freezing Bag, 10 - 25 ml, "Y4" Tube Set, case of 24	One spike, one male luer with cap and 50 cm SCD Tube with a needle-free injection port and pinch clamp	Non-Stock
CS50	CryoStore 50 EVA freezing Bag, 10 — 30 ml, "Standard" Tube Set, case of 24	Two female luers with caps	Stock
CS50M3	CryoStore 50 EVA Freezing Bag, 10 - 30 ml, "M3" Tube Set, case of 24	One male luer with cap, two female luers with caps, a needle-free injection port and a pinch clamp	Non-Stock
CS50N	CryoStore 50 EVA Freezing Bag, 10 - 30 ml, "N" Tube Set, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS50NS	CryoStore 50 EVA Freezing Bag, 10 - 30 ml, "N" Tube Set, double aseptically packed, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS250	CryoStore 250 EVA Freezing Bag, 30 - 70 ml, "Standard" Tube Set, case of 24	Two female luers with caps	Stock
CS250N	CryoStore 250 EVA Freezing Bag, 30 - 70 ml, "N" Tube Set, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS250NS	CryoStore 250 EVA Freezing Bag, 30 - 70 ml, "N" Tube Set, double aseptically packed, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS250NSW	CryoStore 250 EVA Freezing Bag, 30 - 70 ml, "N" Tube Set, double aseptically packed in OW, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS250S	CryoStore 250 EVA Freezing Bag, 30 - 70 ml, "Standard" Tube Set, double aseptically packed, one bag, case of 12	Two female luers with caps	Non-Stock
CS500	CryoStore 500 EVA Freezing Bag, 55 - 100 ml, "Standard" Tube Set, case of 24	Two female luers with caps	Stock
CS500N	CryoStore 500 EVA Freezing Bag, 55 - 100 ml, "N" Tube Set, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS500NS	CryoStore 500 EVA Freezing Bag, 55 - 100 ml, "N" Tube Set, double aseptically packed, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS500NSW	CryoStore 500 EVA Freezing Bag, 55 - 100 ml, "N" Tube Set, double aseptically packed in OW, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS500V2	CryoStore 500 EVA 2-Bag set with spikes, 55 - 100 ml, 2 sets per pouch, case of 12 sets	One nylon spike with cover, one male luer with cap, a needle-free injection port and pinch clamps.	Non-Stock
CS750N	CryoStore 750 EVA Freezing Bag, 80 - 190 ml, "N" Tube Set, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS750NS	CryoStore 750 EVA Freezing Bag, 80 - 190 ml, "N" Tube Set, double aseptically packed, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS750NSW	CryoStore 750 EVA Freezing Bag, 80 - 190 ml, "N" Tube Set, double aseptically packed in OW, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS750	CryoStore 750 EVA Freezing Bag, 80 - 190 ml, "Standard" Tube Set, case of 24	Two female luers with caps	Stock
CS1000N	CryoStore 1000 EVA Freezing Bag, 120 - 275 ml, "N" Tube Set, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port with roller clamps	Stock
CS1000NS	CryoStore 1000 EVA Freezing Bag, 120 - 275 ml, "N" Tube Set, double aseptically packed, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port with roller clamps	Non-Stock
CS2000N	CryoStore 2000 EVA Freezing Bag, 180 - 450 ml, "N" Tube Set, case of 24	One male luer with cap, two female luers with caps, and a needle-free injection port	Stock
CS2000NS	CryoStore 2000 EVA Freezing Bag, 180 - 450 ml, "N" Tube Set, double aseptically packed, case of 12	One male luer with cap, two female luers with caps, and a needle-free injection port	Non-Stock
	, , , , , , , , , , , , , , , , , , , ,		

